VIASURE

Real Time PCR Detection Kit

Cryptosporidium, Giardia & E. histolytica
for ELITe InGenius®

CE IVD

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Content

1.	Intended use.....	5
2.	Summary and Explanation	5
3.	Principle of the procedure	5
4.	Reagents provided	6
5.	Reagents and equipment to be supplied by the user	7
6.	Transport and storage conditions.....	7
7.	Precautions for users	7
8.	Test procedure	9
8.1.	Specimen collection, transport and storage	9
8.2.	Sample Preparation and DNA extraction.....	9
8.3.	Extraction control	10
8.4.	Lyophilized Positive Control	10
8.5.	Lyophilized Reaction-Mix tube.....	10
8.6.	PCR Condition	11
8.6.1.	PCR Test program for VIASURE <i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Real Time PCR Detection Kit for ELITE InGenius®.....	11
8.6.2.	ELITE InGenius® Instrument set up.....	12
9.	Result interpretation	13
10.	Limitations of the test	15
11.	Quality control	16
12.	Performance characteristics.....	17
12.1.	Clinical sensitivity and specificity.....	17
12.2.	Analytical sensitivity	17
12.3.	Analytical specificity	18
12.4.	Analytical reactivity	19
1.	Usó previsto.....	20

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ENGLISH

1. Intended use

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® is an automated in-vitro diagnostic test designed for qualitative detection of *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in human faecal samples from patients with signs and symptoms of gastrointestinal infection. This test is intended for use as an aid in the diagnosis of the *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in combination with clinical and epidemiological risk factors. The assay is validated in association with ELITE InGenius®, an automated and integrated system for extraction, amplification, and detection of nucleic acids. DNA is extracted from faecal specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica*.

2. Summary and Explanation

Enteric protozoa continue to be the most commonly encountered parasitic diseases affecting millions of people each year and causing significant morbidity and mortality worldwide. Among them, *Cryptosporidium* spp., *Entamoeba histolytica*, and *Giardia lamblia* are the major etiological agents.

Cryptosporidium can be acquired through person-to-person, contact with companion and farm animals, ingestion of contaminated food, drinking water and recreational water. *Cryptosporidium* can be present without symptoms or cause diarrhoea and abdominal discomfort with weight loss and malabsorption. *Giardia lamblia* infection occurs by the faecal-oral route transmission, either through direct contact or by ingestion of contaminated food or water. Diarrhoea is the major symptom of the disease, but the clinical manifestations vary; from non-symptomatic infections to severe diarrhoea, abdominal cramps, bloating and flatulence often accompanied by nausea and weight loss. *Entamoeba histolytica* disease is spread primarily by food or water contaminated with cysts but may also be transmitted from person to person by faecal-oral route. Clinical features of amoebiasis range from asymptomatic colonization (90% of the cases) to amoebic dysentery (stomach pain, nausea or vomiting and severe diarrhoea with bloody and slimy stools, fever) and invasive extraintestinal disease, most commonly in the form of liver abscesses.

Microscopy has been considered to be the gold standard method for diagnosis of these parasites. However, is time-consuming, not very sensitive and unable to distinguish the invasive pathogenic *E. histolytica* from the commensal parasite *E. dispar*. Real-time PCR assay is less labour-intensive and has higher sensitivity and specificity, making it an attractive alternative.

3. Principle of the procedure

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® is designed for the qualitative detection of DNA from *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* in faecal specimens. After DNA isolation, the identification of these pathogens is performed by the amplification of a conserved region of the 18S rRNA gene for *Cryptosporidium*, *Giardia lamblia* and *Entamoeba histolytica* respectively, using specific primers and a fluorescent-labelled probe.

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence is measured on ELITE InGenius®.

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer and polymerase) in a stabilized format, as well as an Extraction Control to monitor PCR inhibition. *Cryptosporidium* DNA targets are amplified and detected in Cy5 channel, *Giardia lamblia* DNA targets are amplified and detected in FAM channel, *Entamoeba histolytica* DNA targets are amplified and detected in ROX channel and the extraction control (EC) in HEX, VIC or JOE channel.

Target	ELITE InGenius® Channel	Compatible Fluorophore	Gene
<i>Cryptosporidium</i>	Ch 5	Cy5	18S rRNA
<i>Giardia lamblia</i>	Ch 1	FAM	18S rRNA
<i>Entamoeba histolytica</i>	Ch 4	ROX	18S rRNA
Extraction Control (EC)	Ch 2	HEX, VIC or JOE *	-

Table 1. Target, channel and genes.

4. Reagents provided

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® includes the following materials and reagents detailed in Table 2:

Reagent/Material	Description	Colour	Amount
<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
<i>Cryptosporidium</i> , <i>Giardia</i> & <i>E. histolytica</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	2 vials x 1 mL

Table 2. Reagents and materials provided in VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® with Cat. N°. VS-2340KGE196TE.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius®.

- Automated sample-to-result solution dedicated to Molecular Diagnostics: ELITE InGenius® Instrument (ELITechGroup S.p.A., ref. INT030).
- Extraction cartridges «ELITE InGenius® SP 200» (ELITechGroup S.p.A., ref. INT032SP200),
- Consumables for extraction and amplification «ELITE InGenius® SP 200 Consumable Set» (ELITechGroup S.p.A, ref. INT032CS),
- Amplification cartridges «ELITE InGenius® PCR Cassette» (ELITechGroup S.p.A, ref. INT035PCR),
- Tips «300 µL Filter Tips Axygen» (Axygen BioScience Inc., CA, USA, ref. TF-350-L-R-S),
- Boxes «ELITE InGenius® Waste Box» (ELITechGroup S.p.A, ref. F2102-000).
- Sarstedt 2.0 mL tube skirted screw-cap (Sarstedt Ref. 72.694.005).
- Sarstedt 0.5 mL tube skirted screw-cap (Sarstedt Ref. 72.730.005).
- 50 mL conical tube
- Specimen Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Microcentrifuge.
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Laminar airflow hood.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C.
- Store the kit at 2-40°C until expiration date stated in the label, except the rehydrated Reaction-Mix tube. Once the *Cryptosporidium*, *Giardia* & *E. histolytica* Reaction-Mix tube has been reconstituted, it may be kept at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, we recommend storing at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).
- The rehydrated Reaction-Mix should be recapped and stored inside the pouch with desiccant material in order to protect from light and humidity.
- Once the positive control has been re-suspended, store it at -20°C. We recommend separating it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR, human genetics and *in vitro* diagnostic procedures (including training on ELITE InGenius®).

- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use. Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius®, ELITE InGenius® SP 200 extraction cartridge, any additional reagents required for testing, and ELITE InGenius® are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents and cartridges.
- Please make sure that ELITE InGenius® SP 200 extraction cartridge and ELITE InGenius® PCR Cassette are assembled correctly (see ELITE InGenius® Ref. INT030 Instructions for Use).
- To avoid contamination of the environment by amplicons, do not remove the cap of the ELITE InGenius® PCR Cassette after use. The cap of the ELITE InGenius® PCR Cassette are designed to prevent contamination.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the

environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.

- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.
- Consult the ELITE InGenius® User's Manual and 'ELITE InGenius® Quick Reference Guide' for additional warnings, precautions and procedures before using VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius®.
- Please make sure that the software of ELITE InGenius® System is the last updated version.

8. Test procedure

8.1. Specimen collection, transport and storage

The VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® has been validated on faecal specimens. Another different types of samples must be validated by the user.

Collection, storage, and transport of specimens should be maintained per the conditions validated by the user. Overall, faecal specimens should be collected and labelled appropriately in clean containers. After collection, specimens should be placed in a biohazard bag and should be transported and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at Room Temperature (RT) for up to 2 hours, or at 2 to 8°C for up to 24 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 24 hours), we recommend shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 24 hours or frozen at -20°C or ideally at -70°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The clinical specimens must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC guidelines (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>), the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94), and García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

8.2. Sample Preparation and DNA extraction

Faecal samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend using fresh samples.

For longer storage, the samples must be frozen between -20°C and -70 °C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise faecal sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations in the instructions for use of the extraction kit used, ELITE InGenius® and 'ELITE InGenius® Quick Reference Guide'. Note that some other samples may require pre-processing. For faecal samples, use 'Stool sample: pre-treatment for bacterial and parasitic target with ELITE InGenius®' protocol. Application-specific extraction preparation procedures should be developed and validated by the user.

8.3. Extraction control

Reconstitute the extraction control (green vial) by adding 500 µL of Water RNase/DNase free (white vial) and vortex thoroughly. Transfer the extraction control rehydrated to the 0.5 mL tube skirted screw-cap (Ref. 72.730.005).

8.4. Lyophilized Positive Control

Cryptosporidium, *Giardia* & *E. histolytica* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Cryptosporidium*, *Giardia* & *E. histolytica* Positive Control (red vial) by adding 950 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles

8.5. Lyophilized Reaction-Mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *Cryptosporidium*, *Giardia* & *E. histolytica* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

8.6. PCR Condition

Note: Please, refer to the ELITE InGenius® User's Manual and 'ELITE InGenius® Quick Reference Guide' for detailed instructions.

8.6.1. PCR Test program for VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius®

The product VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® can be used with the ELITE InGenius® system in order to perform:

- A. Integrated run (Extract + PCR)
- B. Amplification run (PCR only)
- C. Extraction run (Extraction only)

To setup an integrated run (A), carry out the following steps as per the Graphical user Interface (GUI).

- 1) Power on the instrument and login with your credentials.
- 2) Select 'Perform Run' from the 'Home' screen.

Note: Verify that the Amplification Controls (Positive Control and Negative Control) have been run and have not expired (status) along with the Amplification Reagent lot to be used. If you do not have valid and approved amplification controls, they should be run as outlined in the following sections.

- 3) Verify the 'Extraction Input Volume' is 200 µL and the 'Extracted Elute Volume' is 100 µL.
- 4) For each 'Track' scan or key the sample barcodes.
- 5) Select assay protocol in the 'Assay' column: VIASURE_KGE_00. Make sure the selected protocol is 'Extract + PCR'.
- 6) Select the sample loading position 'Sample position' column, and then select 'Extraction Tube'.
- 7) Load and verify all required consumables and reagents.
 - Load the Extraction Control and the *Cryptosporidium*, *Giardia* & *E. histolytica* Reaction-Mix tube into the selected 'Inventory Block' following the instructions on the interface.
 - Load and check the Tip Racks in the 'Inventory Area' selected by following the GUI instruction.
 - Load the 'PCR Cassettes', the 'ELITE InGenius® SP 200' extraction cartridges, all the required consumables and the samples to be extracted, following the GUI instruction
- 8) Close door and start run.
- 9) View and approve results.

All the parameters needed for the session are included in the Assay protocol available on the instrument and are automatically recalled when the Assay protocol is selected.

A summary of program settings is shown below:

- I. Assay protocol name: VIASURE_KGE_00

II. Extraction

- a. Extraction Cassette Name: ELITE InGenius® SP 200.
- b. Extraction Input Volume: 200 µL
- c. Extraction Elution Volume: 100 µL

III. PCR

- a. Total Reaction Volume: 20 µL
- b. PCR Template Volume: 5 µL
- c. PCR protocol conditions:

Cycles	Step	Time (s)	Temperature
1	Initial denaturation	120	95°C
45	Denaturation	10	95°C
	Annealing/Extension (Data collection*)	50	60°C

Table 3. PCR Protocol.

- d. Read Fluorescence on the 'Annealing/Extension (Data collection*)' step. Fluorogenic data should be collected during the extension step (*) through the *Cryptosporidium*, *Giardia lamblia*, *E. histolytica* and Extraction Control channels.

PCR only mode (B) implies all the steps mentioned above considering the following variations:

- Step 3: Even if no extraction will be carried out, ensure that the 'Extraction Input Volume' is 200 µL and the 'Extraction Elute Volume' is 100 µL.
- Step 6: Select 'PCR only' and set sample position as 'Elution tube' (bottom row).

Extraction only mode (C) implies all the steps mentioned above considering the following variations:

- Step 3: Ensure that the 'Extraction Input Volume' is 200 µL and the 'Extraction Elute Volume' is 100 µL.
- Step 6: Select 'Extraction only' and either 'Primary or Extraction tube'.

8.6.2. ELITE InGenius® Instrument set up

For instruction on how to configure the ELITE InGenius® System, including Setting up the instrument, Running an assay, Loading the cartridge rack and sample rack on ELITE InGenius®, review and approval of results, results output, finishing the run and maintenance, and Shutting down the instrument, refer to the ELITE InGenius® User's Manual and 'ELITE InGenius® Quick Reference Guide'.

9. Result interpretation

For a detailed description on how to analyse and approve the data, refer to the ELITE InGenius® User's Manual.

All the result of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Extraction Control (EC) signal to verify the extraction procedure and/or correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions. Using the following tables 4 and 5 read and analyse the results.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control tube and the presence of signal for *Cryptosporidium*, *Giardia lamblia* and *E. histolytica* in the positive control tube. For a valid diagnostic test run, the following control conditions should be met:

Controls	<i>Cryptosporidium</i> ¹	<i>Giardia lamblia</i> ¹	<i>E. histolytica</i> ¹	Extraction Control	Interpretation of Controls
Positive Control (PC)	Positive	Positive	Positive	Positive ²	Passed
Negative Control (NC)	Negative or no signal	Negative or no signal	Negative or no signal	n.a.	Passed

Table 4. Expected Performance of Controls. no signal = no amplification curve

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control tube for any target channel), all results are reported as 'failed' and retesting is required.

2 The Extraction Control (EC) should show an amplification signal (Positive value in positive control wells (PC)). Differences can be observed in the values of Ct in the extraction controls between the positive control and the clinical samples, due to the extraction process.

Assessment of clinical samples test results should be performed after the positive and negative controls have been validated and approved. It should be verified that amplification Positive Control and Negative Control were run with the lot of amplification reagent to be used. Both controls results will expire after 15 days.

If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of individual patient sample results, use the following criteria for read and analyse the results:

<i>Cryptosporidium</i>	<i>Giardia lamblia</i>	<i>E. histolytica</i>	Extraction Control	Interpretation for individual patients' samples	
Positive	Positive	Positive	Positive or no signal ¹	Passed	<i>Cryptosporidium, Giardia lamblia and E. histolytica</i> Positive
Negative	Negative	Negative	Positive	Passed	<i>Cryptosporidium, Giardia lamblia and E. histolytica</i> Negative
Positive	Negative	Negative	Positive or no signal ¹	Passed	<i>Cryptosporidium</i> Positive, <i>Giardia lamblia</i> and <i>E. histolytica</i> Negative
Positive	Positive	Negative	Positive or no signal ¹	Passed	<i>Cryptosporidium</i> and <i>Giardia lamblia</i> Positive, and <i>E. histolytica</i> Negative
Positive	Negative	Positive	Positive or no signal ¹	Passed	<i>Cryptosporidium</i> and <i>E. histolytica</i> Positive, and <i>Giardia lamblia</i> Negative
Negative	Positive	Negative	Positive or no signal ¹	Passed	<i>Giardia lamblia</i> Positive, <i>Cryptosporidium</i> and <i>E. histolytica</i> Negative
Negative	Positive	Positive	Positive or no signal ¹	Passed	<i>Giardia lamblia</i> and <i>E. histolytica</i> Positive, <i>Cryptosporidium</i> Negative
Negative	Negative	Positive	Positive or no signal ¹	Passed	<i>E. histolytica</i> Positive, <i>Cryptosporidium</i> and <i>Giardia lamblia</i> Negative
Negative no signal	Negative or no signal	Negative or no signal	Negative or no signal ²	Failed	Test Failure – Repeat Testing ² Reported as 'Invalid' in the graphical user interface (GUI) in the ELITE InGenius® System.


Table 5. Interpretation of patient sample results. Ct values. no signal = no amplification curve.

1 The Extraction Control (EC) shows or not an amplification signal (Positive or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids. Differences can be observed in the values of Ct in the extraction controls between the controls and the clinical samples, due to the extraction process.

2 In the case of *Cryptosporidium*, *Giardia lamblia* and *E. histolytica* target genes negative, EC must show an amplification signal with. If there is an absence of signal or Ct value Negative of the Extraction Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the ELITE InGenius® User's Manual and 'ELITE InGenius® Quick Reference Guide' to verify the correct performance of extraction process, each RT-qPCR step and review the parameters, and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 1. Correct run of negative and positive control run on the ELITE InGenius®.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from human faecal samples.
- The performance evaluation characteristics have been evaluated following the PCR protocol named: VIASURE_KGE_00. In the case of using a protocol other than the one mentioned, the test must be validated by the end user.
- For good test performance, the lyophilized product should be at the bottom of the tube and not adhered to the top area of the tube or the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- The quality of the test depends on the quality of the sample: nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica*, either the high number of DNA template copies which contains each *Cryptosporidium*, *Giardia lamblia* & *Entamoeba histolytica* Positive Control vial, samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Cryptosporidium*, *Giardia lamblia* & *Entamoeba histolytica* Positive Control, which contains high copies of template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in an established order and in a separate laboratory area.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimen collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the DNA during sample shipping/storage and/or processing.
 - A bacterial load in the specimen below the limit of detection for the assay.
 - The presence of RT-qPCR inhibitors or other types of interfering substances. The impact of antibiotics used to prevent *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.

If in doubt, refer to section 9 to check the correct interpretation of the results.

- A negative EC signal does not preclude the presence of *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* DNA in a clinical specimen.
- A positive test result does not necessarily indicate the presence of viable pathogens and does not imply that these pathogens are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of target bacterial sequences (*18S rRNA gene*).
- Negative results do not preclude *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak pathogen levels during infections caused by *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogens.
- If diagnostic tests for other gastrointestinal diseases are negative and the patient's clinical presentation and epidemiological information suggest that *Cryptosporidium*, *Giardia lamblia* and/or *Entamoeba histolytica* infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Fluorescence values may vary due to multiple factors such as: PCR equipment, extraction system, type of sample, previous treatment of the sample, etc... among others.

11. Quality control

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Extraction Control (EC) in each specimen tube confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® was tested using clinical stool samples from patients with suspected gastrointestinal infection.

In order to determine the clinical diagnostic accuracy, different evaluations have been conducted through collaboration with national entities. A summary of the sites, sample type and workflow are included in the following table. The results were as follows:

	Site	Sample type	Workflow	Target
1	Hospital Clínico Universitario Lozano Blesa (Zaragoza, Spain)	Stool swabs	ELITE InGenius® SP 200 + ELITE InGenius® System	<i>Cryptosporidium</i>
				<i>Giardia</i>
				<i>E. histolytica</i>

Table 6. Site, sample type, workflow, and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	Overall agreement	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	Ziehl-Neelsen stain and microscopy visualization + RIDA®GENE Parasitic Stool Panel II (r-Biopharm)	<i>Cryptosporidium</i>	0.99 (0.95-0.99)	29	99	1*	0	1 (0.88-1)	0.99 (0.94-1)	0.96 (0.8-0.99)	1 (0.95-1)
		<i>Giardia</i>	1 (0.97-1)	32	97	0	0	1 (0.89-1)	1 (0.96-1)	1 (0.86-1)	1 (0.95-1)
		<i>E. histolytica</i>	1 (0.97-1)	20	109	0	0	1 (0.83-1)	1 (0.96-1)	1 (0.79-1)	1 (0.95-1)
		<i>Cryptosporidium</i> with $CT_{cut-off}$	1 (0.97-1)	29	100	0	0	1 (0.88-1)	1 (0.96-1)	1 (0.85-1)	1 (0.95-1)

Table 7. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV), Predictive Negative Values (NPV) for VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius®.

* Discrepant samples were tested using the reference molecular method RIDA®GENE Bacterial Stool panel (r-Bpahrm). Results indicates that there were 1 false positive for *Cryptosporidium*. Using a cut-off 37.5 or 37 (only in Cy5 channel), all false positive values disappear.

Results show high agreement to detect *Cryptosporidium*, *Giardia lamblia* and *E. histolytica* using VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® in stool samples.

12.2. Analytical sensitivity

VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® has a detection limit of 40 copies per reaction for *Cryptosporidium*, 25 copies per reaction for *Giardia lamblia* and 0.6 cells per reaction for *E. histolytica*.

12.3. Analytical specificity

The specificity of the *Cryptosporidium*, *Giardia lamblia* and *E. histolytica* assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter lari</i>	-	Enteropathogenic <i>Escherichia coli</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter fetus</i>	-	<i>Klebsiella oxytoca</i>	-
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter coli</i>	-	<i>Listeria monocytogenes</i>	-
<i>Helicobacter heilmannii</i>	-	<i>Campylobacter jejuni subsp. jejuni</i>	-	<i>Candida albicans</i>	-
<i>Shigella flexneri</i>	-	<i>Campylobacter upsaliensis</i>	-	<i>Arcobacter butzleri</i>	-
<i>Shigella dysenteriae</i>	-	<i>Proteus vulgaris</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Salmonella typhi</i>	-	<i>Aeromonas hydrophila subsp. Hydrophila</i>	-	<i>Enterococcus faecalis</i>	-
<i>Salmonella paratyphi A</i>	-	<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3 y O:9	-
<i>Salmonella paratyphi B</i>	-	<i>Staphylococcus aureus subsp. Aureus</i>	-	Enteroinvasive <i>Escherichia coli</i>	-
<i>Salmonella typhimurium</i>	-	<i>Serratia liquefaciens</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella bongori</i>	-	<i>Vibrio parahaemolyticus</i>	-	Adenovirus 1-5, 8, 15, 31, 40 y 41	-
<i>Salmonella enteritidis</i>	-	<i>Clostridium difficile</i>	-	Rotavirus A	-
<i>Salmonella enterica subsp. Entérica</i>	-	<i>Clostridium perfringens</i>	-	Norovirus Genotypes I y II	-
<i>Salmonella pullorum</i>	-	Enterotoxigenic <i>Escherichia coli</i>	-	Astrovirus Genotype 1-8	-
<i>Salmonella gallinarum</i>	-	Enterohemorrhagic <i>Escherichia coli</i>	-	Ecovirus 30	-
Enterovirus 68 and 71	-	<i>Aeromonas caviae</i>	-	Sapovirus	-
<i>Dientamoeba fragilis</i>	-	<i>Blastocystis hominis</i>	-	<i>Entamoeba dispar</i>	-
<i>Campylobacter hyointestinalis</i>	-	Coxsackievirus A24, A9 y B3	-		

Table 8. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® for *Giardia lamblia* was evaluated against DNA extracted from *Giardia intestinalis* culture (synonym: *Giardia lamblia*, *Giardia duodenalis*), showing positive results.

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® for *Entamoeba histolytica* was evaluated against DNA extracted from *Entamoeba histolytica*, showing positive results.

The reactivity of VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit for ELITE InGenius® for *Cryptosporidium* was evaluated against DNA extracted from *Cryptosporidium parvum*, showing positive results.

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.