

Instructions for use

MRSA/SA - ELITe Positive Control

plasmid DNA control for qualitative assay



REF M800356

UDI 08033891486549



CHANGE HISTORY

Rev.	Notice of change	Date (dd/mm/yy)
07	Update of the paragraph: "Other product required"	07/11/25
06	Extended use of the product in association with «ELITe BeGenius®» instrument (REF INT040) Update of the paragraph "Symbols" with the symbol "Consult instructions for use" New graphics and content setting of the IFU	27/03/25
05	The number of tubes and the volume of Positive Control has been modified from 4 x 65 µL to 2 x 160 µL.	05/07/18
00–04	New product development and subsequent changes	—

NOTE

The revision of this IFU is also compatible with the previous versions of the kit

TABLE OF CONTENT

1 INTENDED USE 4

2 PRODUCT DESCRIPTION 4

3 MATERIALS PROVIDED IN THE PRODUCT 4

4 MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT..... 4

5 OTHER PRODUCTS REQUIRED..... 5

6 WARNINGS AND PRECAUTIONS 5

7 PROCEDURE..... 6

8 REFERENCES..... 6

9 SYMBOLS 7

1 INTENDED USE

The product **MRSA/SA - ELITE Positive Control** is an *in vitro* diagnostic medical device intended to be used by healthcare professionals as DNA positive control in nucleic acids Real-Time PCR assay for the detection of the DNA of *Staphylococcus aureus* (SA) and methicillin-resistant *Staphylococcus aureus* (MRSA, including the LGA251 strain) in association with **MRSA/SA ELITE MGB® Kit** and **ELITE InGenius®**, **ELITE BeGenius®** and **ABI 7500 Fast Dx Real-Time PCR Instruments**.

2 PRODUCT DESCRIPTION

The product supplies the **MRSA/SA Positive Control** and the **LGA251/SA Positive Control**, plasmid DNAs at known titre in stabilizing solutions based on Tris-HCl and EDTA, aliquoted into **two ready-to-use test tubes** each.

One plasmid solution contains a region of a *Staphylococcus aureus*-specific gene and a region of the *mecA* gene; the other plasmid solution contains the same region of a *Staphylococcus aureus*-specific gene and a region of *mecC* gene. Both *mecA* and *mecC* genes (*mecA*_{LGA251}) confer resistance to methicillin and to other beta-lactam antibiotics.

The detection of target DNAs, using **MRSA/SA ELITE MGB Kit** product in association with **ELITE InGenius** and **ELITE BeGenius** instruments, attests the system ability to detect the SA and MRSA DNA, including LGA251 strain, and consequently the verification of the system (product batch and instrument).

The product contains sufficient reagents for **10 separate sessions** on **ELITE InGenius** and **ELITE BeGenius** (5 sessions each tube), and **24 separate analytic sessions** in association with other systems, by using 10 µL per reaction.

NOTE

The plasmid DNAs concentration in copies / mL was determined through absorbance measurement by spectrophotometer. There are no WHO approved standards for the target genomic DNAs.

3 MATERIALS PROVIDED IN THE PRODUCT

Table 1

Component	Description	Quantity	Classification of Hazards
MRSA/SA Positive Control Ref. M800356	plasmid DNAs solution in tube with RED cap	2 x 160 µL	-
LGA251/SA Positive Control Ref. M800356	plasmid DNAs solution in tube with BLACK cap	2 x 160 µL	-

4 MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT

- Laminar airflow hood.
- Disposable powderless nitrile gloves or similar material.
- Vortex mixer.
- Bench microcentrifuge (~13,000 RPM).
- Micropipettes and sterile tips with aerosol filter or sterile positive displacement tips (volume range: 0.5-1000 µL).
- 2.0 mL sterile screw capped tubes (Sarstedt, Germany, ref. 72.694.005).
- 0.5 mL sterile screw capped tubes (Sarstedt, Germany, ref. 72.730.005).
- Molecular biology grade water.

5 OTHER PRODUCTS REQUIRED

The reagents for Real-Time amplification reaction and the consumables **are not** included in this product.

To perform the assay the following products are required:

Table 2

Instruments and softwares	Products and reagents
ELITe InGenius (ELITechGroup S.p.A., EG SpA, ref. INT030) ELITe InGenius Software version 1.3.0.19 (or later) MRSA-SA ELITe_PC_200_100 or MRSA-SA ELITe_PC_200_50 , Assay Protocols with parameters for Positive Control analysis.	MRSA/SA ELITe MGB Kit (EG SpA, ref. M800351) ELITe InGenius and ELITe BeGenius Consumables (see ELITe InGenius and ELITe BeGenius Instruction for use)
ELITe BeGenius (EG SpA, ref. INT040) ELITe BeGenius Software version 2.3.0 (or later) MRSA-SA ELITe_Be_PC_200_100 or MRSA-SA ELITe_Be_PC_200_50 , Assay Protocols with parameters for Positive Control analysis.	
7500 Fast Dx Real-Time PCR Instrument (ThermoFisher Scientific, ref. 4406985) NucliSENS® easyMAG (bioMérieux SA, Ref. 200111)	MRSA/SA ELITe MGB Kit (EG SpA, ref. M800351) MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL (Life Technologies, ref. 4346906)

6 WARNINGS AND PRECAUTIONS

This product is designed for in vitro use only.

6.1 Warnings and general precautions

- Handle and dispose of all reagents and all materials used to carry out the assay as if they were infectious. Avoid direct contact with the reagents. Avoid splashing or spraying. Waste must be handled and disposed of in compliance with adequate safety standards. Disposable combustible material must be incinerated. Liquid waste containing acids or bases must be neutralized before disposal.
- Wear suitable protective clothes and gloves and protect eyes and face.
- Never pipette solutions by mouth.
- Do not eat, drink, smoke or apply cosmetic products in the work areas.
- Carefully wash hands after handling samples and reagents.
- Dispose of leftover reagents and waste in compliance with the regulations in force.
- Carefully read all the instructions provided before running the assay.
- While running the assay, follow the product instructions provided.
- Do not use the product after the indicated expiry date.
- Only use the reagents provided with the product and those recommended by the manufacturer.
- Do not use reagents from different batches.
- Do not use reagents from other manufacturers.

6.2 Warnings and precautions for molecular biology

Molecular biology procedures, require qualified and trained staff to avoid the risk of erroneous results, especially due to sample nucleic acids degradation or sample contamination by PCR products.

When amplification session is manually setup, it is necessary to have available separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplification products. Never introduce an amplification product in the area designated for extraction / preparation of amplification reactions.

Laboratory coats, gloves and tools dedicated to work session setup are needed.

The reagents must be handled under a laminar airflow hood. The pipettes used to handle the reagents must be exclusively used for this purpose. The pipettes must be of the positive displacement type or be used with aerosol filter tips. The tips used must be sterile, free from DNases and RNases, free from DNA and RNA.

The PCR Cassette must be handled carefully and never opened to prevent PCR product diffusion and carryover contamination.

6.3 Warnings and precautions specific for the components

Table 3

Component	Storage temperature	Use from first opening	Freeze / thaw cycles	On board stability (ELiTe InGenius and ELiTe BeGenius)
MRSA/SA Positive Control	-20°C or below	one month	up to twelve	up to four separate sessions* of three hours each
LGA251/SA Positive Control	-20°C or below	one month	up to twelve	up to four separate sessions* of three hours each

*with intermediate freezing

7 PROCEDURE

The product **MRSA/SA - ELiTe Positive Control** must be used in association with the product **MRSA/SA ELiTe MGB Kit**.

The components **MRSA/SA Positive Control** and **LGA251/SA Positive Control** are ready to use: a volume of **10 µL** each is directly added to two separate reaction mixtures (**MRSA/SA PCR Mix**, component of **MRSA/SA ELiTe MGB Kit**) by the instrument ELiTe InGenius or ELiTe BeGenius, or manually when other instruments are used.

Before use, take and thaw the **MRSA/SA Positive Control** and **LGA251/SA Positive Control** tubes at room temperature (+16 / +26°C) for 30 minutes. Mix gently, spin down the content for 5 seconds and keep them on ice or in a cool block.

The complete assay procedure is described in detail in the instructions for use of the product **MRSA/SA ELiTe MGB Kit**.

The performance characteristics and procedure limitations of the complete assay are described in detail in the instructions for use of the product **MRSA/SA ELiTe MGB Kit**.

NOTE

The results of the **MRSA/SA - ELiTe Positive Control** will be stored by the ELiTe InGenius and ELiTe BeGenius instruments and used to set up the Control Charts monitoring the amplification step performances. For each batch of the product **MRSA/SA ELiTe MGB Kit**, the amplification of Positive Control is required. The stored results of the Positive Control amplification will expire **after 15 days**.

8 REFERENCES

Garcia-Alvarez, L. et al. Methicillin-resistant *Staphylococcus aureus* with a novel *mecA* homologue in human and bovine populations in the UK and Denmark: a descriptive study. *Lancet Infect Dis* 2011, 11:595-603

Ito T. et al. Guidelines for reporting novel *mecA* gene homologues. *Antimicrob Agents Chemother.* 2012 October; 56(10): 4997-4999.

9 SYMBOLS



Catalogue Number.



Upper limit of temperature.



Batch code.



Use by (last day of month).



in vitro diagnostic medical device.



Fulfilling the requirements of the European Directive 98/79/EC for *in vitro* diagnostic medical device.



Unique Device Identification



Contains sufficient for "N" tests.



Consult instructions for use.



Contents.



Manufacturer.



ELITechGroup S.p.A.
C.so Svizzera, 185, 10149 Torino ITALY
Tel. +39-011 976 191
Fax +39-011 936 76 11
E. mail: emd.support@elitechgroup.com
WEB site: www.elitechgroup.com