

Instructions for use

EBV ELITe MGB[®] Kit

reagents for DNA Real-Time PCR



REF RTS020PLD

UDI 08033891483579

CE IVD
0123

CHANGE HISTORY

Rev.	Notice of change	Date (dd/mm/yyyy)
24-R	Expansion of the use of the product in association with the ELITE InGenius® and ELITE BeGenius® instruments and Cerebrospinal Fluid (CSF) matrix. Update of the paragraph "Notice To Purchaser: Limited License". Update to "Model 2 version 4.0.2" of Assay protocols for the instruments ELITE InGenius® (ref INT030) and ELITE BeGenius® (ref INT040).	03/03/2026
23-R	Expansion of the use of the product in association with MyGenius PRO instrument (REF INT050) with whole blood matrix.	08/10/2025
22-R	Update for inclusion of the extraction kit ELITE InGenius SP 1000 (ref INT033SP1000) and related data in association with plasma and ELITE InGenius with "1000 in 100" protocol. Upgrade of the analytical and diagnostic performances in PERFORMANCE CHARACTERISTICS paragraph (LoD, Linear measuring range, Diagnostic Sensitivity and Specificity) Update of the paragraph "Symbols" with the symbol "Consult instructions for use"	04/03/2025
21-R	Update for compliance with the Regulation (EU) 2017/746 on in vitro diagnostic medical devices (IVDR) requirements. Upgrade of the analytical and diagnostic performances in PERFORMANCE CHARACTERISTICS paragraph Update of the Intended use. New graphics and content setting of the IFU.	29/08/2024
20	Description of the IC cut off value in the section: Diagnostic Specificity	22/12/2022
19	Modification of ULoQ value with whole blood samples in association with ELITE InGenius instrument	01/04/2022
18	Expanded use on the automated and integrated instrument ELITE BeGenius	30/09/2021
17	Addition of analytical data (LoD, linearity, ...) obtained in association with Roche cobas z 480 analyzer	07/05/2020
00 — 16	new product development and succeeding changes	-

TABLE OF CONTENT

1 INTENDED USE	4
2 ASSAY PRINCIPLE	4
3 PRODUCT DESCRIPTION	4
4 MATERIALS PROVIDED IN THE PRODUCT	5
5 MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT.....	5
6 OTHER PRODUCTS REQUIRED.....	5
7 WARNINGS AND PRECAUTIONS	7
8 SPECIMENS AND CONTROLS for ELITe InGenius, ELITe BeGenius AND MyGenius PRO.....	8
9 ELITe InGenius PROCEDURE.....	12
10 ELITe BeGenius PROCEDURE	17
11 MyGenius PRO PROCEDURE.....	23
12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO.....	28
13 SPECIMENS AND CONTROLS FOR ABI 7500 Fast Dx Real-Time PCR Instrument.....	41
14 ABI 7500 Fast Dx Real-Time PCR Instrument PROCEDURE	42
15 PERFORMANCE CHARACTERISTICS WITH ABI 7500 Fast Dx Real-Time PCR Instrument.....	47
16 REFERENCES.....	50
17 PROCEDURE LIMITATIONS	50
18 TROUBLESHOOTING	51
19 SYMBOLS	57
20 NOTICE TO THE USERS.....	58
21 NOTICE TO PURCHASER: LIMITED LICENSE	58
Appendix A QUICK START GUIDE.....	59
Appendix B QUICK START GUIDE.....	64
Appendix C QUICK START GUIDE.....	67

1 INTENDED USE

The product **EBV ELITe MGB® Kit** is an *in vitro* diagnostic medical device intended to be used by healthcare professionals as quantitative nucleic acids Real-Time PCR assay for the detection and quantification of the **DNA of Epstein-Barr virus (EBV)** extracted from clinical specimens.

The assay is validated in association with the **ELITe InGenius®** and **ELITe BeGenius®** instruments, automated and integrated systems for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid.

The assay is also validated in association with **MyGenius PRO®** (registration name ELIVERSE®) instrument, automated and integrated system for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA.

The assay is also validated in association with the **ELITe GALAXY**, automatic extraction and PCR set-up system and **7500 Fast Dx Real-Time PCR Instrument**, Real-Time PCR platform, using human specimens of whole blood collected in EDTA.

The product is intended for use as an aid in the diagnosis and monitoring of EBV infections in patients suspected of having or undergoing monitoring of EBV infections.

The results must be interpreted in combination with all relevant clinical observation and laboratory outcomes.

2 ASSAY PRINCIPLE

The assay is a quantitative Real-Time PCR detecting EBV DNA isolated from specimens and amplified using the assay reagent **EBV Q - PCR Mix**, that contains primers and probes with ELITe MGB® technology.

The ELITe MGB probes are activated when hybridize with the related PCR products. **ELITe InGenius**, **ELITe BeGenius** and **MyGenius PRO** instruments monitor fluorescence increase and calculate the threshold cycles (Ct) and the melting temperatures (Tm). The EBV quantity is calculated based on a stored calibration curve.

7500 Fast Dx Real-Time PCR Instrument measures and records the increase of fluorescence emission. The subsequent data processing allows the detection and quantification of EBV in the primary specimen.

In the ELITe MGB probes the fluorophores are quenched in the random-coiled, single-stranded state of probe. The fluorophores are active in the probe / amplicon duplex as the quencher is spatially separated from the fluorophore. Note the fluorophore is not cleaved during PCR and can be utilized for dissociation analysis and melting temperature calculation.

3 PRODUCT DESCRIPTION

The **EBV ELITe MGB Kit** provides the assay reagent **EBV Q-PCR Mix**, an optimized and stabilized PCR mixture that contains the specific primers and probes for:

- a region of the **EBNA-1 gene** of EBV, detected in Channel **EBV**; the probe is stabilized by MGB, quenched by the Eclipse Dark Quencher®, and labelled by FAM dye.
- Internal Control, specific for **promoter and 5' UTR region** of the **human beta-globin gene**, detected in Channel **IC**; the probe is stabilized by MGB, quenched by the Eclipse Dark Quencher and labelled by AquaPhluor® 525 (AP525) dye.

The **EBV Q-PCR Mix** also contains buffer, magnesium chloride, triphosphates nucleotides, AP593 fluorophore, (used instead of ROX or Cy5) as passive reference for fluorescence normalization), the enzyme Uracil N-glycosidase (UNG) to inactivate contamination by the amplification product and the "hot start" DNA Polymerase enzyme.

The product **EBV ELITe MGB Kit** contains sufficient reagents for **96 tests** on **ELITe InGenius**, **ELITe BeGenius** and **MyGenius PRO** with **20 µL** used per reaction.

The product **EBV ELITe MGB Kit** contains sufficient reagents for **100 tests on other systems**, with **20 µL** used per reaction.

NOTE

A conversion factor is required to express the results of the quantitative analysis in International Units of EBV of "1st WHO International Standard for Human Epstein Barr Virus for Nucleic Acid Amplification Techniques" (NIBSC ref. 09/260, United Kingdom).

4 MATERIALS PROVIDED IN THE PRODUCT**Table 1**

Component	Description	Quantity	Classification of hazards
EBV Q-PCR Mix ref. RTS020PLD	Mixture of reagents for Real-Time PCR in tube with NATURAL cap	4 x 540 µL	-

5 MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT

- Laminar airflow hood.
- Disposable nitrile powder-free gloves or similar material.
- Vortex mixer.
- Bench centrifuge (~5,000 RPM).
- Bench microcentrifuge (~13,000 RPM).
- Micropipettes and sterile tips with aerosol filter or sterile positive displacement tips (volume range: 0.5-1000 µL).
- 2.0 mL sterile screw capped tubes (Sarstedt, Germany, ref. 72.694.005).
- 0.5 mL sterile screw capped tubes (Sarstedt, Germany, ref. 72.730.005)
- Molecular biology grade water.

6 OTHER PRODUCTS REQUIRED

The reagents for the extraction of sample DNA, the extraction and inhibition internal control, the amplification positive and negative controls, the DNA standards and the consumables **are not** provided with this product.

For the extraction of nucleic acids, Real-Time PCR and result interpretation of samples, the following products are required:

Table 2

Instruments and softwares	Products and reagents
<p>ELITE InGenius (ELITechGroup S.p.A., EG SpA, ref. INT030) ELITE InGenius Software version 1.3.0.19 (or later) EBV ELITE_STD, Assay Protocol with parameters for Calibrators analysis EBV ELITE_PC, Assay Protocol with parameters for Positive Control analysis EBV ELITE_NC, Assay Protocol with parameters for Negative Control analysis EBV ELITE_STD_1000_100, Assay Protocol with parameters for Calibrators analysis (only for plasma specimens, sample volume 1000 µL) EBV ELITE_PC_1000_100, Assay Protocol with parameters for Positive Control analysis (only for plasma specimens, sample volume 1000 µL) EBV ELITE_NC_1000_100, Assay Protocol with parameters for Negative Control analysis (only for plasma specimens, sample volume 1000 µL) EBV ELITE_WB_200_100 Assay Protocol with parameters for whole blood specimen analysis EBV ELITE_PL_200_100 Assay Protocols with parameters for plasma specimen analysis EBV ELITE_PL_1000_100 Assay Protocol with parameters for plasma specimen analysis EBV ELITE_CSF_200_100 Assay Protocols with parameters for cerebrospinal fluid specimen analysis</p>	<p>EBV ELITE Standard (EG SpA, ref. STD020PLD) EBV - ELITE Positive Control (EG SpA, ref. CTR020PLD) CPE - Internal Control (EG SpA, ref. CTRCPE) ELITE InGenius SP 200 (EG SpA, ref. INT032SP200) ELITE InGenius SP1000 (EG SpA, ref. INT033SP1000) ELITE InGenius and ELITE BeGenius Consumables (see ELITE InGenius and ELITE BeGenius Instruction for use)</p>
<p>ELITE BeGenius (EG SpA ref. INT040) ELITE BeGenius Software version 2.3.0 (or later) EBV ELITE_Be_STD, Assay Protocol with parameters for Calibrators analysis EBV ELITE_Be_PC, Assay Protocol with parameters for Positive Control analysis EBV ELITE_Be_NC, Assay Protocol with parameters for Negative Control analysis EBV ELITE_Be_WB_200_100, Assay Protocol with parameters for whole blood specimen analysis EBV ELITE_Be_PL_200_100, Assay Protocol with parameters for plasma specimen analysis EBV ELITE_Be_CSF_200_100, Assay Protocol with parameters for cerebrospinal fluid specimen analysis</p>	
<p>MyGenius PRO (EG SpA ref: INT050). MyGenius PRO Software version BB-04 (or later) EBV ELITE_My_STD, Assay Protocol with parameters for Calibrators analysis EBV ELITE_My_PC, Assay Protocol with parameters for Positive Control analysis EBV ELITE_My_NC, Assay Protocol with parameters for Negative Control analysis EBV ELITE_My_WB_IU_200_100, Assay Protocol with parameters for whole blood specimen analysis and result in IU/mL EBV ELITE_My_WB_cmL_200_100, Assay Protocol with parameters for whole blood specimen analysis and result in copies/mL</p>	<p>EBV ELITE Standard (EG SpA, ref. STD020PLD) EBV - ELITE Positive Control (EG SpA, ref. CTR020PLD) Negative Control (EG SpA, ref. CTRNEG) Internal Control Maxi (EG SpA, ref. ICMAXI) MyGenius PRO Consumables (see MyGenius PRO Instruction for use)</p>

Table 2 (continued)

Instruments and softwares	Products and reagents
<p>7500 Fast Dx Real-Time PCR Instrument (ThermoFisher Scientific, ref. 4406985)</p> <p>ELITE GALAXY (EG SpA, ref. INT020) with software version 1.3.1 (or later). Extraction Protocol for ELITE GALAXY, xNA Extraction (Universal)</p>	<p>EBV ELITE Standard (EG SpA, ref. STD020PLD)</p> <p>EBV - ELITE Positive Control (EG SpA, ref. CTR020PLD)</p> <p>CPE - Internal Control (EG SpA, ref. CTRCPE)</p> <p>ELITE GALAXY 300 Extraction Kit (EG SpA, ref. INT021EX)</p> <p>ELITE GALAXY Consumables (see ELITE GALAXY Instruction for use)</p> <p>MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL (Life Technologies, ref. 4346906), microplates with 0.1 mL wells and adhesive sealing sheets for real time amplification.</p>

7 WARNINGS AND PRECAUTIONS

This product is designed for in-vitro use only.

7.1 General warnings and precautions

Handle and dispose of all biological samples as if they were infectious. Avoid direct contact with biological samples. Avoid splashing or spraying. Tubes, tips and other materials that come into contact with the biological samples must be treated for at least 30 minutes with 3% sodium hypochlorite (bleach) or autoclaved for one hour at 121°C before disposal.

Handle and dispose of all reagents and all materials used to carry out the assay as if they were infectious. Avoid direct contact with the reagents. Avoid splashing or spraying. Waste must be handled and disposed of in compliance with adequate safety standards. Disposable combustible material must be incinerated. Liquid waste containing acids or bases must be neutralized before disposal. Do not allow extraction reagents to contact sodium hypochlorite (bleach).

- Wear suitable protective clothes and gloves and protect eyes and face,
- Never pipette solutions by mouth,
- Do not eat, drink, smoke or apply cosmetic products in the work areas,
- Carefully wash hands after handling samples and reagents,
- Dispose of leftover reagents and waste in compliance with the regulations in force,
- Carefully read all the instructions provided before running the assay,
- While running the assay, follow the product instructions provided,
- Do not use the product after the indicated expiry date,
- Only use reagents provided with the product and those recommended by the manufacturer,
- Do not use reagents from different batches,
- Do not use reagents from other manufacturers.

7.2 Warnings and precautions for molecular biology

Molecular biology procedures require qualified and trained staff to avoid the risk of erroneous results, especially due to sample nucleic acid degradation or sample contamination by PCR products.

When amplification session is manually setup, it is necessary to have available separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplification products. Never introduce an amplification product in the area designated for extraction / preparation of amplification reactions.

When amplification session is manually setup, it is necessary to have available lab coats, gloves and tools which are exclusively used for the extraction / preparation of the amplification reactions and for the amplification / detection of amplification products.

Never transfer lab coats, gloves or tools from the area designated for the amplification / detection of amplification products to the area designated for the extraction / preparation of the amplification reactions.

Laboratory coats, gloves and tools dedicated to work session setup are needed.

The samples must be suitable and, if possible, dedicated for this type of analysis. Samples must be handled under a laminar airflow hood. Pipettes used to handle samples must be exclusively used for this specific purpose. The pipettes must be of the positive displacement type or be used with aerosol filter tips. The tips used must be sterile, free from DNases and RNases, and free from DNA and RNA.

The reagents must be handled under a laminar airflow hood. The pipettes used to handle the reagents must be exclusively used for this purpose. The pipettes must be of the positive displacement type or be used with aerosol filter tips. The tips used must be sterile, free from DNases and RNases and free from DNA and RNA.

The extraction products must be handled to prevent dispersion into the environment and to avoid contamination of the instrument's working area.

The PCR Cassette must be handled carefully and never opened to prevent PCR product diffusion and carryover contamination.

7.3 Warnings and precautions specific for the components

Table 3

Component	Storage temperature	Use from first opening	Freeze / Thaw cycles	On board stability (ELITE InGenius and ELITE BeGenius)	On board stability (MyGenius PRO)
EBV Q-PCR Mix	-20 °C or below (protected from light)	one month	up to five	up to five separate* sessions of three hours each or up to 7 consecutive hours (2 sessions of 3 hours each and the time needed to start a third session)	up to 7 consecutive hours or up to 3 consecutive hours for five times*

* with intermediate freezing.

8 SPECIMENS AND CONTROLS for ELITE InGenius, ELITE BeGenius AND MyGenius PRO

8.1 Specimens and Assay Protocols

This product is intended for use on the **ELITE InGenius**, **ELITE BeGenius** and **MyGenius PRO** with the relative validated clinical specimens identified and handled according to laboratory guidelines, and collected, transported, and stored under the following conditions:

Specimen	Collection requirements	Transport/Storage conditions			
		+16 / +26 °C (room temperature)	+2 / +8 °C	-20 ± 10 °C	-70 ± 15 °C
Whole blood	EDTA	≤ 1 d	≤ 3 d	≤ 30 d	≤ 30 d
Plasma	EDTA	≤ 1 d	≤ 3 d	≤ 30 d	≤ 30 d
CSF	Minimal patient Blood contamination	≤ 1 d	≤ 2 d	≤ 30 d	≤ 30 d

EDTA, Ethylenediaminetetraacetic acid; CSF, cerebrospinal fluid, d, day.

Even if longer storage periods at -70 °C are possible, as extensively reported by scientific literature, their application should be evaluated internally by the end-users of this product.

It is recommended to divide the specimens into aliquots before freezing to prevent repeated freeze/thaw cycles. When using frozen samples, thaw the samples just before the extraction to avoid possible nucleic acid degradation.

To perform samples testing on the **ELITE InGenius** and the **ELITE BeGenius**, the following Assay Protocols must be used. These IVD protocols were specifically validated with ELITE MGB Kits and the **ELITE InGenius** or **ELITE BeGenius** with the indicated matrices.

Table 4 Assay Protocols for EBV ELITE MGB Kit and ELITE InGenius and ELITE BeGenius

Specimen	Instrument	Assay Protocol name	Report	Characteristics
Whole blood in EDTA	ELITE InGenius	EBV ELITE_WB_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Sonication: NO Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
	ELITE BeGenius	EBV ELITE_Be_WB_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
Plasma in EDTA	ELITE InGenius	EBV ELITE_PL_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Sonication: NO Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
	ELITE InGenius	EBV ELITE_PL_1000_100	copies/mL or IU/mL	Extraction Input Volume: 1000 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Sonication: NO Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
	ELITE BeGenius	EBV ELITE_Be_PL_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL

Table 4 Assay Protocols for EBV ELITE MGB Kit and ELITE InGenius and ELITE BeGenius (continued)

Specimen	Instrument	Assay Protocol name	Report	Characteristics
CSF	ELITE InGenius	EBV ELITE_CSF_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Sonication: NO Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
	ELITE BeGenius	EBV ELITE_Be_CSF_200_100	copies/mL or IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL

IU, international units

To perform samples testing on the **MyGenius PRO**, the following Assay Protocols must be used. These IVD protocols were specifically validated with ELITE MGB Kits and the **MyGenius PRO** with the indicated matrices.

Table 5 Assay Protocols for EBV ELITE MGB Kit and MyGenius PRO

Specimen	Instrument	Assay Protocol name	Report	Characteristics
Whole blood in EDTA	MyGenius PRO	EBV ELITE_My_WB_IU_200_100	IU/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Sonication: NO Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
Whole blood in EDTA	MyGenius PRO	EBV ELITE_My_WB_cmL_200_100	copies/mL	Extraction Input Volume: 200 µL Extraction Elution Volume: 100 µL Internal Control: 10 µL Dilution factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL

NOTE

Verify if the primary tube and the volume of the sample are compatible with ELITE InGenius, ELITE BeGenius or MyGenius PRO following the Instruction for use of the instruments.

The volume of the sample in a primary tube varies according to the type of the tube loaded. Refer to the instructions for use of the extraction kit for more information on how to set up and perform the extraction procedure.

If required, 200 µL of sample must be transferred into an Extraction tube (for ELITE InGenius) or 2 mL Sarstedt Tube (for ELITE BeGenius), and 260 µL of sample into a 2 mL Sarstedt Tube (for MyGenius PRO).

NOTE

Pipetting samples to the **Extraction tube** or to the **2 mL Sarstedt Tube** might **generate contamination**. Use the appropriate pipettes and follow all recommendations reported in the 7 WARNINGS AND PRECAUTIONS page 7.

Purified nucleic acids can be left at room temperature for 16 hours and stored at -20 °C or below for no longer than one month.

Refer to “Potentially Interfering Substances” in the [12 PERFORMANCE CHARACTERISTICS WITH ELITE InGenius, ELITE BeGenius and MyGenius PRO page 28](#) section to check data concerning interfering substances.

NOTE

Do not use samples collected in heparin, which is a known reverse transcription and PCR inhibitor.

8.2 PCR calibrators and controls

The Calibration curve must be generated and approved for each lot of PCR reagent.

- For the calibration curve, use the four levels of the product **EBV ELITE Standard** (not provided with this kit) with the **EBV ELITE_STD**, **EBV ELITE_STD_1000_100**, **EBV ELITE_Be_STD** or **EBV ELITE_My_STD** Assay Protocols.

NOTE

The concentration of Q – PCR Standards are expressed in copies / reaction (10^5 copies / rxn, 10^4 copies / rxn, 10^3 copies / rxn, 10^2 copies / rxn). Refer to “Standard Curve Uncertainty” in the [12 PERFORMANCE CHARACTERISTICS WITH ELITE page 28](#) section.

PCR control results must be generated and approved for each lot of PCR reagent.

- For the Positive Control, use the product **EBV - ELITE Positive Control** (not provided with this kit) with the **EBV ELITE_PC**, **EBV ELITE_PC_1000_100**, **EBV ELITE_Be_PC** or **EBV ELITE_My_PC** Assay Protocols.
- For the Negative Control, use molecular biology grade water (not provided with this kit) with the **EBV ELITE_NC**, **EBV ELITE_NC_1000_100**, **EBV ELITE_Be_NC** or use **Negative Control** (not provided with this kit) with the **EBV ELITE_My_NC** Assay Protocols.

NOTE

The **ELITE InGenius**, **ELITE BeGenius** and **MyGenius PRO** allow generation and storage of the calibration curve and PCR control validation for each lot of PCR reagent.

Calibration curves expire after **60 days**, at which time it is necessary to re-run the calibration.

PCR control results expire after **15 days**, at which time it is necessary to re-run the positive and negative controls.

The Calibrators and PCR controls must be re-run if any of the following events occur:

- a new lot of reagents is used,
- results of quality control analysis (see following paragraph) are out of specification,
- any major maintenance or service is performed on the **ELITE InGenius**, **ELITE BeGenius** and **MyGenius PRO** instruments.

8.3 Quality controls

Verification of the extraction and PCR procedure is recommended. Archived samples or certified reference material may be used. External controls should be used in accordance with local, state, and federal accrediting organizations, as applicable.

9 ELITe InGenius PROCEDURE

The procedure to use the **EBV ELITe MGB Kit** with the **ELITe InGenius** consists of three steps:

Table 6

STEP 1	Verification of the system readiness	
STEP 2	Session setup	A) Sample run (Extract + PCR)
		B) Eluted sample run (PCR Only)
		C) Calibration run (PCR Only)
		D) Positive Control and Negative Control run (PCR Only)
STEP 3	Review and approval of results	1) Validation of Calibration curve
		2) Validation of Positive Control and Negative Control results
		3) Validation of sample results
		4) Sample result reporting

9.1 STEP 1 – Verification of the system readiness

Before starting the session:

- switch on the **ELITe InGenius** and login “**CLOSED**” mode,
- in the “Calibration” menu on the Home page, verify the Calibrators (**Q – PCR Standard**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid Calibrators are available for the **PCR Mix** lot, perform calibration as described in the following sections,
- in the “Controls” menu on the Home page, verify the PCR Controls (**Positive Control, Negative Control**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid PCR Controls are available for the **PCR Mix** lot, run the PCR Controls as described in the following sections,
- select the type of run, following the instructions on the Graphical User Interface (GUI) for the session setup and use the Assay Protocols provided by EG SpA. (see “Specimens and Controls”).

If the Assay Protocol of interest is not loaded in the system, contact your local ELITechGroup Customer Service.

9.2 STEP 2 – Session Setup

The **EBV ELITe MGB Kit** can be used on **ELITe InGenius** to perform:

- Sample run (Extract + PCR),
- Eluted sample run (PCR Only),
- Calibration run (PCR Only),
- Positive Control and Negative Control run (PCR Only).

All required parameters are included in the Assay Protocols available on the instrument and are loaded automatically when the Assay Protocol is selected.

NOTE

The **ELITe InGenius** can be connected to the “Laboratory Information System” (LIS) which enables downloading the session information. Refer to the instrument manual for more details.

Before to setup a run:

Thaw the needed **PCR Mix** tubes at room temperature for 30 minutes. Each tube is sufficient for **24 tests**. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.

NOTE

Protect the **PCR Mix** from light while thawing because this reagent is photosensitive.

To set up one of the four types of run follow the steps below while referring to the GUI:

	A. Sample run (Extract + PCR)	B. Eluted sample run (PCR Only)
1	<p>Identify samples and, if needed, thaw at room temperature, mix gently, spin down the contents for 5 seconds and keep on ice or cool block. If required, transfer 200 or 1000 µL of sample in an Extraction tube previously labelled.</p> <p>Thaw the needed CPE tubes at room temperature for 30 minutes. Mix gently, spin down the contents for 5 seconds and keep on ice or cool block. Each tube is sufficient for 12 extractions.</p>	<p>Thaw the Elution tube containing the extracted nucleic acids at room temperature. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.</p>
2	Select " Perform Run " from the "Home" screen.	Select " Perform Run " from the "Home" screen.
3	Ensure the "Extraction Input Volume" is 200 or 1000 µL and the Extracted Elute Volume is 100 µL.	Ensure the "Extraction Input Volume" is 200 or 1000 µL and the Extracted Elute Volume is 100 µL.
4	For each sample, assign a Track and enter the "SampleID" (SID) by typing or by scanning the sample barcode.	For each sample, assign a Track and enter the "SampleID" (SID) by typing or by scanning the sample barcode.
5	Select the Assay Protocol in the "Assay" column (see "Specimens and Controls").	Select the Assay Protocol in the "Assay" column (see "Specimens and Controls").
6	Ensure the "Protocol" displayed is: "Extract + PCR".	Select "PCR Only" in the "Protocol" column.
7	Select the sample loading position as "Primary tube" or "Extraction Tube" in the "Sample Position" column. Ensure the " Dilution factor " is "1".	Ensure the sample loading position in the "Sample Position" column is "Elution Tube (bottom row)". Ensure the " Dilution factor " is "1".
8	Click "Next" to continue.	Click "Next" to continue.
9	Load CPE and the PCR Mix on the "Inventory Block" referring to the "Load List" and enter CPE and PCR Mix lot number, expiry date and number of reactions for each tube.	Load the PCR Mix on the "Inventory Block" referring to the "Load List" and enter PCR Mix lot number, expiry date and number of reactions for each tube.
10	Click "Next" to continue.	Click "Next" to continue.
11	Verify the tips in the "Tip Racks" in the "Inventory Area" and replace Tip Racks if necessary.	Verify the tips in the "Tip Racks" in the "Inventory Area" and replace Tip Racks if necessary.
12	Click "Next" to continue.	Click "Next" to continue.
13	Load PCR Cassette, ELITe InGenius SP 200 or SP1000 extraction cartridges, and all required consumables and samples to be extracted.	Load PCR Cassette and Elution tubes with samples extracted.
14	Click "Next" to continue.	Click "Next" to continue.
15	Close the instrument door.	Close the instrument door.
16	Press "Start".	Press "Start".

	C. Calibration run (PCR Only)	D. Positive Control and Negative Control run (PCR Only)
1	Thaw the needed Q-PCR Standard tubes (Cal1: Q-PCR Standard 10 ² , Cal2: Q-PCR Standard 10 ³ , Cal3: Q-PCR Standard 10 ⁴ , Cal4: Q-PCR Standard 10 ⁵) at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.	Thaw Positive Control tubes at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block. Prepare the Negative Control by transferring at least 50 µL of molecular biology grade water to an "Elution tube", provided with the ELiTe InGenius SP 200 Consumable Set.
2	Select "Perform Run" from the "Home" screen.	Select "Perform Run" from the "Home" screen.
3	Ensure the "Extraction Input Volume" is 200 or 1000 µL and the "Extracted Elute Volume" is 100 µL.	Ensure "Extraction Input Volume" is 200 or 1000 µL and "Extracted Elute Volume" is 100 µL.
4	For the Q-PCR Standard, assign the "Track", select the Assay Protocol (see "Specimen and Controls") in the "Assay" column and enter the reagent lot number and expiry date.	Select the Assay Protocol in the "Assay" column (see "Specimens and Controls"). Enter the lot number and expiry date of the Positive Control and of the molecular biology grade water.
5	Ensure "PCR Only" is selected in the "Protocol" column.	Ensure "PCR Only" is selected in the "Protocol" column.
6	Ensure the sample loading position in "Sample Position" column is "Elution Tube (bottom row)".	Ensure the sample loading position in the "Sample Position" column is "Elution Tube (bottom row)".
7	Load the PCR Mix on the "Inventory Block" referring to the Load List and enter the PCR Mix lot number, expiry date and number of reactions for each tube.	Load the PCR Mix on the "Inventory Block" referring to the "Load List" and enter the PCR Mix lot number, expiry date and number of reactions for each tube.
8	Click "Next" to continue.	Click "Next" to continue.
9	Verify the tips in the "Tip Rack(s)" in the "Inventory Area" and replace Tip Rack(s) if necessary.	Verify the tips in the "Tip Rack(s)" in the "Inventory Area" and replace Tip Rack(s) if necessary.
10	Click "Next" to continue.	Click "Next" to continue.
11	Load the PCR Cassette and the Q-PCR Standard tubes.	Load PCR Cassette, Positive Control and Negative Control.
12	Click "Next" to continue.	Click "Next" to continue.
13	Close the instrument door.	Close the instrument door.
14	Press "Start"	Press "Start".

When the session is finished, the **ELiTe InGenius** allows users to view, approve, store the results, print and save the report.

NOTE

At the end of the run the remaining Extracted Sample in the **Elution tube** must be removed from the instrument, capped, identified and stored at -20 ± 10 °C for no longer than one month. Avoid spilling of the Extracted Sample.

NOTE

At the end of the run the **PCR Mix** can be removed from the instrument, capped and stored at -20 °C or below or can be kept on board in the refrigerated block for up to 7 hours (2 sessions of 3 hours each and the time needed to start a third session); mix gently and spin down the content for 5 seconds before starting the next session.

NOTE

At the end of the run the remaining **Q - PCR Standard** can be removed from the instrument, capped, and stored at -20 °C or below. Avoid spilling the Q - PCR Standard.

NOTE

The **Q - PCR Standard** can be used for 4 separate sessions of 2 hours each.

NOTE

At the end of the run the remaining **Positive Control** can be removed from the instrument, capped and stored at -20 °C or below. Avoid the spilling of the Positive Control. The remaining **Negative Control** must be discarded.

NOTE

The **Positive Control** can be used for 4 separate sessions of 3 hours each.

NOTE

At the end of the run the **PCR Cassette** and the other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

9.3 STEP 3 - Review and approval of results

The **ELITE InGenius** monitors target and internal control fluorescence signals for each reaction and automatically applies the Assay Protocol parameters to generate PCR curves which are then interpreted into results.

At the end of the run, the “Results Display” screen is automatically shown. In this screen the results and the run information are shown. From this screen results can be approved, and reports printed or saved (“Sample Report” or “Track Report”). Refer to the instrument manual for more details.

NOTE

The **ELITE InGenius** can be connected to the “Laboratory Information System” (LIS) which enables uploading the session results to the laboratory data center. Refer to the instrument manual for more details.

The **ELITE InGenius** generates results with the **EBV ELITE MGB Kit** through the following procedure:

1. Validation of Calibration curve,
2. Validation of Positive Control and Negative Control results,
3. Validation of sample results,
4. Sample result reporting.

9.3.1 Validation of Calibration curve

The **ELITE InGenius software** interprets the PCR results for the target of the Calibrator reactions with the **EBV ELITE STD** or **EBV ELITE STD_1000_100** Assay Protocol parameters. The resulting Ct versus concentration produces the Calibration curve.

The Calibration curves, specific for the PCR reagent lot, are recorded in the database (Calibration). They can be viewed and approved by “Administrator” or “Analyst” users, following the GUI instructions.

The Calibration curve expires **after 60 days**.

NOTE

If the Calibration curve does not meet the acceptance criteria, the “Failed” message is shown on the “Calibration” screen. In this case, the results cannot be approved, and the Calibrator amplification reactions must be repeated. In addition, if samples were included in the run, these are not quantified and must also be repeated to generate quantitative results.

9.3.2 Validation of amplification Positive Control and Negative Control results

The **ELITe InGenius Software** interprets the PCR results for the target of the Positive Control and Negative Control reactions with the **EBV ELITe_PC** or **EBV ELITe_PC_1000_100** and **EBV ELITe_NC** or **EBV ELITe_NC_1000_100** Assay Protocols parameters. The resulting Ct values are converted to concentration and used to verify the system (reagents lot and instrument).

The Positive Control and Negative Control results, specific for the PCR reagent lot, are recorded in the database (Controls). They can be viewed and approved by "Administrator" or "Analyst" users following the GUI instructions.

The Positive Control and Negative Control results expire **after 15 days**.

The **ELITe InGenius Software** processes the Positive Control and Negative Control results and generates Control Charts. The approval of the Positive Control is based on the evaluation of the obtained logarithmic quantity that should be within the expected logarithmic quantity range (PC Chart). This ensures the system performance is within the acceptance criteria. The second chart (L-J chart) is dedicated exclusively to monitoring the Positive Control trend over time. Refer to the instrument manual for more details.

NOTE

If the Positive Control or Negative Control result does not meet the acceptance criteria, the "Failed" message is displayed on the "Controls" screen. In this case, the results cannot be approved, and the Positive Control or Negative Control runs must be repeated.

NOTE

If the Positive Control or Negative Control result is not valid and samples were included in the same run, the samples can be approved but their results are not validated. In this case, the failed Control(s) and samples must all be repeated.

9.3.3 Validation of Sample results

The **ELITe InGenius Software** interprets the PCR results for the target (Channel **EBV**) and the Internal Control (Channel **IC**) with the **EBV ELITe_WB_200_100**, **EBV ELITe_PL_200_100**, **EBV ELITe_PL_1000_100** or **EBV ELITe_CSF_200_100** Assay Protocol parameters. The resulting target Ct values are converted to concentration.

Results are shown in "Results Display" screen.

The sample results can be approved when the three conditions in the table below are true.

1) Calibration Curve	Status
EBV Q-PCR Standard	APPROVED
2) Positive Control	Status
EBV Positive Control	APPROVED
3) Negative Control	Status
EBV Negative Control	APPROVED

The sample results are automatically interpreted by the **ELITe InGenius Software** using Assay Protocol parameters.

The possible result messages are listed in the table below.

For each sample the system reports a combination of the following messages specifying if the pathogen DNAs are either detected or not detected.

Result of sample run	Interpretation
EBV:DNA Detected, quantity equal to XXXcopies/mL or IU/mL	EBV DNA was detected in the sample within the assay measurement range, its concentration is shown.

EBV:DNA Detected, quantity below LLoQcopies/mL or IU/mL	EBV DNA was detected in the sample, its concentration is below the assay Lower Limit of Quantification
EBV:DNA Detected, quantity beyond ULoQcopies/mL or IU/mL	EBV DNA was detected in the sample, its concentration is above the assay Upper Limit of Quantification
EBV:DNA Not detected or below LoDcopies/mL or IU/mL	EBV DNA was not detected in the sample. The sample is negative for EBV DNA, or its concentration is below the assay Limit of Detection.
Invalid-Retest Sample	Not valid assay result , caused by Internal Control failure (due to e.g., incorrect extraction or inhibitors carry-over). The test should be repeated.

Samples reported as “Invalid-Retest Sample”: in this case, the Internal Control DNA was not efficiently detected, which could be due to problems in sample collection, extraction or PCR steps (e. g. incorrect sampling, degradation or loss of DNA during the extraction, or inhibitors in the eluate), which may cause incorrect results.

If sufficient eluate volume remains, the eluate can be retested (as is or diluted) by an amplification run in “PCR Only” mode. If the second result is invalid, the sample must be retested starting from extraction of a new sample using “Extract + PCR” mode (see [18 TROUBLESHOOTING page 51](#)).

Samples reported as “EBV:DNA Not Detected or below “LoD” copies/mL or IU/mL” are suitable for analysis but EBV was not detected. In this case the sample may be either negative for EBV DNA or the EBV DNA is present at a concentration below the Limit of Detection of the assay (see [12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO page 28](#)).

EBV DNA positive samples at a concentration below the Limit of Detection (and Lower Limit of Quantification) of the assay, if detected, are reported as “EBV:DNA Detected, quantity below “LLoQ” copies/mL or IU/mL” (see [12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO page 28](#)).

EBV DNA positive samples within the Linear Measuring Range are detected and are reported as “EBV:DNA Detected, quantity equal to “XXX” copies/mL or IU/mL” (see [12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO page 28](#)).

EBV DNA positive samples that are above the Upper Limit of Quantification are reported as “EBV:DNA Detected, quantity beyond “ULoQ” copies/mL or IU/mL” (see [12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO page 28](#)), and they are not suitable for quantification. If needed the sample may be diluted before extraction or PCR and retested to yield results within the Linear Measuring Range of the assay.

NOTE

The results obtained with this assay must be interpreted in combination with all relevant clinical observation and laboratory outcomes.

The sample results are stored in the database and, if valid, can be approved (Results Display) by “Administrator” or “Analyst” users, following the GUI instruction. From the “Results Display” window it is possible to print and save the Sample run results as “Sample Report” and “Track Report”.

9.3.4 Sample result reporting

The sample results are stored in the database and reports can be exported as “Sample Report” and “Track Report”.

The “Sample Report” shows the results details by selected sample (SID).

The “Track Report” shows the results details by selected Track.

The "Sample Report" and "Track Report" can be printed and signed by authorized personnel.

10 ELITe BeGenius PROCEDURE

The procedure to use the **EBV ELITe MGB Kit** with the **ELITe BeGenius** consists of three steps:

Table 7

STEP 1	Verification of the system readiness	
STEP 2	Session setup	A) Sample run (Extract + PCR)
		B) Eluted sample run (PCR Only)
		C) Calibration run (PCR Only)
		D) Positive Control and Negative Control run (PCR Only)
STEP 3	Review and approval of results	1) Validation of Calibration curve
		2) Validation of Positive Control and Negative Control results
		3) Validation of sample results
		4) Sample result reporting

10.1 STEP 1 - Verification of the system readiness

Before starting the session:

- switch on the **ELITe BeGenius** and login “**CLOSED**” mode,
- in the “Calibrations” menu on the Home page, verify the Calibrators (**Q - PCR Standard**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid Calibrators are available for the **PCR Mix** lot, perform calibration as described in the following sections,
- in the “Controls” menu on the Home page, verify the PCR Controls (**Positive Control, Negative Control**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid PCR Controls are available for the **PCR Mix** lot, run the PCR Controls as described in the following sections,
- select the type of run, following the instructions on the Graphical User Interface (GUI) for the session setup and use the Assay Protocols provided by EG SpA (see “Specimens and Controls”).

If the Assay Protocol of interest is not loaded in the system, contact your local ELITechGroup Customer Service.

10.2 STEP 2 – Session Setup

The **EBV ELITe MGB Kit** can be used on the **ELITe BeGenius** to perform:

- A. Sample run (Extract + PCR),
- B. Eluted sample run (PCR Only),
- C. Calibration run (PCR Only),
- D. Positive Control and Negative Control run (PCR Only).

All the required parameters are included in the Assay Protocols available on the instrument and are loaded automatically when the Assay Protocol is selected.

NOTE

The **ELITe BeGenius** can be connected to the “Laboratory Information System” (LIS) which enables downloading the session information. Refer to the instrument manual for more details.

Before to setup a run:

Thaw the needed **PCR Mix** tubes at room temperature for 30 minutes. Each tube is sufficient for **24 tests**. Mix gently then spin down the contents for 5 seconds and keep on ice or cool block.

NOTE

Protect the **PCR Mix** from light while thawing because this reagent is photosensitive.

To set up one of the four types of run follow the steps below while referring to the GUI:

	A. Sample run (Extract + PCR)	B. Eluted sample run (PCR Only)
1	<p>Identify samples and, if needed, thaw at room temperature, mix gently, spin down the contents for 5 seconds and keep on ice or cool block. If required, transfer 200 µL of sample in a 2 mL Sarstedt tube previously labelled.</p> <p>Thaw the needed CPE tubes at room temperature for 30 minutes. Mix gently, spin down the contents for 5 seconds and keep on ice or cool block. Each tube is sufficient for 12 extractions.</p>	<p>Thaw the Elution tube containing the extracted nucleic acids at room temperature. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.</p>
2	Select “Perform Run” from the “Home” screen.	Select “Perform Run” from the “Home” screen.
3	Remove all the Racks from the “Cooler Unit” and place them on the preparation table.	Remove the “Racks” from “Lane 1, 2 and 3” (L1, L2, L3) of the “Cooler Unit” and place them on the preparation table.
4	Select the “Run mode”: “Extract + PCR”.	Select the “Run mode”: “PCR Only”.
5	Load the samples into the “Sample Rack”. (Note: when secondary tubes “2 mL Tubes” are loaded, use the blue adaptors for the “Sample Rack”).	Load the samples into the “Elution Rack”.
6	Insert the “Sample Rack” into the “Cooler Unit” starting from the “Lane 5” (L5). If needed, insert the “Sample ID” (SID) for each “Position” used. (If secondary tubes are loaded, flag “2 mL Tube”. If secondary tubes are not barcoded, type manually the “Sample ID”).	Insert the “Elution Rack” into the “Cooler Unit” starting from “Lane 3” (L3) If needed, for each “Position” enter the “Sample ID”, the “Sample matrix”, the “Extraction kit” and the “Extracted eluate vol.” (eluate volume).
7	Click “Next” to continue.	Click “Next” to continue.
8	Ensure the “Extraction Input Volume” is 200 µL and the “Extracted Elute Volume” is 100 µL.	Not applicable
9	Select the Assay Protocol in the “Assay” column (see “Specimens and Controls”).	Select the Assay Protocol in the “Assay” column (see “Specimens and Controls”).
10	Click “Next” to continue.	Click “Next” to continue.
11	When more than 12 samples are processed, repeat the procedure from point 6.	When more than 12 samples are processed, repeat the procedure from point 6.
12	Load the “Elution tubes” into the “Elution Rack” (Elution tubes can be labelled with barcode to improve traceability).	Not applicable
13	Insert the “Elution Rack” into the “Cooler Unit” starting from “Lane 3” (L3). When more than 12 samples are processed, repeat using “Lane 2” (L2).	Not applicable
14	Click “Next” to continue.	Not applicable
15	Load CPE and the PCR Mix into the “Reagent/Elution Rack”.	Load the PCR Mix into “Reagent/Elution Rack”.
16	Insert the “Reagent/Elution Rack” into the “Cooler Unit” in “Lane 2” (L2) if available or in “Lane 1” (L1). If needed, for each PCR Mix and / or CPE enter the “S/N” (serial number), the “Lot No.” (lot number), the “Exp. Date” (expiry date) and the “T/R” (number of reactions).	Insert the “Reagent/Elution Rack” into the “Cooler Unit” in “Lane 2” (L2) if available or in “Lane 1” (L1). If needed, for each PCR Mix enter the “S/N” (serial number), the “Lot No.” (lot number), the “Exp. Date” (expiry date) and the “T/R” (number of reactions).
17	Click “Next” to continue	Click “Next” to continue.
18	Verify the tips in the “Tip Rack (s)” in the “Inventory Area” and replace Tip Rack(s) if necessary.	Verify the tips in the “Tip Rack (s)” in the “Inventory Area” and replace Tip Rack(s) if necessary.

	A. Sample run (Extract + PCR)	B. Eluted sample run (PCR Only)
19	Click "Next" to continue.	Click "Next" to continue.
20	Load the "PCR Rack" with "PCR Cassette" in the Inventory Area.	Load the "PCR Rack" with "PCR Cassette" in the Inventory Area.
21	Click "Next" to continue.	Click "Next" to continue.
22	Load the "Extraction Rack" with the "ELITE InGenius SP 200" extraction cartridges and required extraction consumables.	Not applicable
23	Close the instrument door.	Close the instrument door.
24	Press "Start".	Press "Start".

	C. Calibration run (PCR Only)	D. Positive Control and Negative Control run (PCR Only)
1	Thaw the needed Q-PCR Standard tubes (Cal1: Q-PCR Standard 10 ² , Cal2: Q-PCR Standard 10 ³ , Cal3: Q-PCR Standard 10 ⁴ , Cal4: Q-PCR Standard 10 ⁵) at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.	Thaw Positive Control tubes at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block. Prepare the Negative Control by transferring at least 50 µL of molecular biology grade water to an "Elution tube", provided with the ELITE InGenius SP 200 Consumable Set.
2	Select "Perform Run" from the "Home" screen.	Select "Perform Run" from the "Home" screen.
3	Remove the "Racks" from "Lane 1, 2 and 3" (L1, L2, L3) from the "Cooler Unit" and place them on the preparation table.	Remove the "Racks" from "Lane 1, 2 and 3" (L1, L2, L3) from the "Cooler Unit" and place them on the preparation table.
4	Select the "Run mode: PCR Only".	Select the "Run mode": "PCR Only".
5	Load the Q-PCR Standard tubes into the "Elution Rack".	Load the Positive Control and Negative Control tubes into the "Elution Rack".
6	Insert the "Elution Rack" into the "Cooler Unit" starting from the "Lane 3" (L3). If needed, for each "Position" enter the "Reagent name" and the "S/N" (serial number), the "Lot No." (lot number), the "Exp. Date" (expiry date) and the "T/R" (number of reactions).	Insert the "Elution Rack" into the "Cooler Unit" starting from the "Lane 3" (L3). If needed, for each "Position" enter the "Reagent name" and the "S/N" (serial number), the "Lot No." (lot number), the "Exp. Date" (expiry date) and the "T/R" (number of reactions).
7	Click "Next" to continue.	Click "Next" to continue.
8	Select the Assay Protocol in the "Assay" column (see "Specimens and Controls").	Select the Assay Protocol in the "Assay" column (see "Specimens and Controls").
9	Click "Next" to continue.	Click "Next" to continue.
10	Load the PCR Mix into "Reagent/Elution Rack".	Load the PCR Mix into "Reagent/Elution Rack".
11	Insert the "Reagent/Elution Rack" into the "Cooler Unit" in "Lane 2" (L2) If needed, for each PCR Mix enter the "S/N" (serial number), the "Lot No." (lot number), the "Exp. Date" (expiry date) and the "T/R" (number of reactions).	Insert the "Reagent/Elution Rack" into the "Cooler Unit" in "Lane 2" (L2). If needed, for each PCR Mix enter the "S/N" (serial number), the "Lot No." (lot number), the "Exp. Date" (expiry date) and the "T/R" (number of reactions).
12	Click "Next" to continue.	Click "Next" to continue.
13	Verify the tips in the "Tip Racks" in the "Inventory Area" and replace Tip Rack(s) if necessary.	Verify the tips in the "Tip Racks" in the "Inventory Area" and replace Tip Rack(s) if necessary.
14	Click "Next" to continue.	Click "Next" to continue.
15	Load the "PCR Rack" with " PCR Cassette " in the Inventory Area.	Load the "PCR Rack" with " PCR Cassette " in the Inventory Area.
16	Click "Next" to continue.	Click "Next" to continue.
17	Close the instrument door.	Close the instrument door.
18	Press "Start".	Press "Start".

When the session is finished, the **ELITE BeGenius** allows users to view, approve, store the results, print and save the report.

NOTE

At the end of the run, the remaining Extracted Sample in the **Elution tube** must be removed from the instrument, capped, identified, and stored at -20 ± 10 °C for no longer than one month. Avoid the spilling of the Extracted Sample.

NOTE

At the end of the run the **PCR Mix** can be removed from the instrument, capped and stored at -20 °C or below or can be kept on board in the refrigerated block for up to 7 hours (2 sessions of 3 hours each and the time needed to start a third session); mix gently and spin down the content for 5 seconds before starting the next session.

NOTE

At the end of the run, the remaining **Q - PCR Standard** can be removed from the instrument, capped and stored at -20 °C or below. Avoid spilling the Q - PCR Standard.

NOTE

The **Q- PCR Standard** can be used for 4 separate sessions of 2 hours each.

NOTE

At the end of the run, the remaining **Positive Control** can be removed from the instrument, capped and stored at -20 °C or below. Avoid the spilling of the **Positive Control**. The remaining **Negative Control** must be discarded.

NOTE

The **Positive Control** can be used for 4 separate sessions of 3 hours each.

NOTE

At the end of the run, the **PCR Cassette** and the other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

10.3 STEP 3 - Review and approval of results

The **ELITe BeGenius** monitors target and internal control fluorescence signals for each reaction and automatically applies the Assay Protocol parameters to generate PCR curves which are then interpreted into results.

At the end of the run, the “Results Display” screen is automatically shown. In this screen the results and the run information are shown. From this screen results can be approved, and reports printed or saved (“Sample Report” or “Track Report”). Refer to the instrument manual for more details.

NOTE

The **ELITe BeGenius** can be connected to the “Laboratory Information System” (LIS) which enables uploading the session results to the laboratory data center. Refer to the instrument manual for more details.

The **ELITe BeGenius** generates the results with the **EBV ELITe MGB Kit** through the following procedure:

1. Validation of Calibration curve,
2. Validation of Positive Control and Negative Control results,
3. Validation of sample results,
4. Sample result reporting.

NOTE

Please, refer to the same paragraph of the **ELITe InGenius Procedure** for the details.

11 MyGenius PRO PROCEDURE

The procedure to use the **EBV ELITE MGB Kit** with the **MyGenius PRO** consists of three steps:

Table 8

STEP 1	Verification of the system readiness	
STEP 2	Session setup	A) Sample testing session (Extract + PCR)
		B) Calibration testing session (PCR Only)
		C) Positive Control and Negative Control testing session (PCR Only)
STEP 3	Review and approval of results	1) Validation of Calibration curve
		2) Validation of Positive Control and Negative Control results
		3) Validation of sample results
		4) Sample result reporting

STEP 1 - Verification of the system readiness

Before starting the testing session:

- switch on MyGenius PRO and log in; the unit will start in STAND-BY mode.
- in the “Calibrations” menu on the Home Screen, verify the Calibrators (**Q - PCR Standard**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid Calibrators are available for the **PCR Mix** lot, perform calibration as described in the following sections,
- in the “Controls” menu on the Home Screen, verify the PCR Controls (**Positive Control, Negative Control**) are approved and valid (Status) for the **PCR Mix** lot to be used. If no valid PCR Controls are available for the **PCR Mix** lot, run the PCR Controls as described in the following sections,
- follow the instructions on the Graphical User Interface (GUI) for the session setup and use the Assay Protocols provided by EG SpA (see “Specimens and Controls”).

If the Assay Protocol of interest is not loaded in the system, contact your local ELITechGroup Customer Service.

STEP 2 – Workflow Setup on the system

The **EBV ELITE MGB Kit** can be used on the **MyGenius PRO** to perform:

- A. Sample testing session (Extract + PCR),
- B. Calibration testing session (PCR Only),
- C. Positive Control and Negative Control testing session (PCR Only).

All the required parameters are included in the Assay Protocols available on the instrument and are loaded automatically when the Assay Protocol is selected.

NOTE

The **MyGenius PRO** can be connected to the “Laboratory Information System” (LIS) which enables downloading the information for the workflow. Refer to the instrument manual for more details.

Before to setup a testing session:

Thaw the needed **PCR Mix** tubes at room temperature for 30 minutes. Each tube is sufficient for **24 tests**. Mix gently then spin down the contents for 5 seconds and keep on ice or cool block.

NOTE

Protect the **PCR Mix** from light while thawing because this reagent is photosensitive.

To set up one of the three types of testing session follow the steps below while referring to the GUI:

	A. Sample testing session (Extract + PCR)	B. Calibration testing session (PCR Only)	C. Positive Control and Negative Control testing session (PCR Only)
1	Identify samples and, if necessary, thaw at room temperature, mix gently, spin down the contents for 5 seconds and keep on ice or cool block. If the primary tube is compatible with MyGenius PRO racks, insert tube with samples into the racks; if this is not possible transfer 260 µL of sample into a previously labelled secondary tube (2mL). Thaw the needed Internal Control Maxi IC MAXI tubes at room temperature for 30 minutes. Mix gently, spin down the contents for 5 seconds. Each tube is sufficient for 72 tests.	Thaw the needed Q-PCR Standard tubes (Cal1: Q-PCR Standard 102, Cal2: Q-PCR Standard 103, Cal3: Q-PCR Standard 104, Cal4: Q-PCR Standard 105) at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block.	Thaw Positive Control tubes at room temperature for 30 minutes. Mix gently, then spin down the contents for 5 seconds and keep on ice or cool block. Spin down the Negative Control tubes for 5 seconds.
2	If the instrument is connected to the LIS, the sample barcodes are automatically read and recognized by the MyGenius PRO when placed in the Sample Area. If the instrument is not connected to the LIS, manually assign the samples using the ASSIGN TEST or ASSIGN SAMPLE functions on the "Sample List" page of the GUI, before loading them into the instrument. These functions also allow assigning the correct test protocol to each sample.	Load the PCR Mix and the Q-PCR Standard tubes into the refrigerated Reagent Carousel.	Load the PCR Mix and the Positive and Negative Controls tubes into the refrigerated Reagent Carousel.
3	In STAND-BY mode, check that there is a sufficient quantity of consumables for completing the test and that solid and liquid waste levels are adequate for the instrument to operate. If needed, load the required consumables into the appropriate drawers and empty the waste box and liquid tank. For detailed loading procedures, refer to the instrument's user manual.	In STAND-BY mode, check that there is a sufficient quantity of consumables for completing the test and that solid and liquid waste levels are adequate for the instrument to operate. If needed, load the required consumables into the appropriate drawers and empty the waste box and liquid tank. For detailed loading procedures, refer to the instrument's user manual.	In STAND-BY mode, check that there is a sufficient quantity of consumables for completing the test and that solid and liquid waste levels are adequate for the instrument to operate. If needed, load the required consumables into the appropriate drawers and empty the waste box and liquid tank. For detailed loading procedures, refer to the instrument's user manual.
4	Press Start in Home page and wait for the instrument to go into status 'Preparation' and then 'Operation'.	Press Start in Home page and wait for the instrument to go into status 'Preparation' and then 'Operation'.	Press Start in Home page and wait for the instrument to go into status 'Preparation' and then 'Operation'.
5	Load the PCR Mix in the Refrigerated Reagent Carousel.	Select "Calibration" button from the "Home" screen	Select "Controls" button from the "Home" screen
6	Book the opening of the door of the Auto Sampler Area and Load the Samples and IC MAXI.	Select in the Calibration page the Q-PCR Standard tube loaded in step 2 and then press "Order".	Select in the Control page the Positive Control and Negative Control loaded in step 2 and then press "Order".
7	Close the Auto Sampler door and the testing session starts.		

At the end of the testing session, the remaining **Q - PCR Standard** can be removed from the instrument, capped and stored at -20 °C or below. Avoid spilling the Q - PCR Standard.

NOTE

At the end of the testing session the **PCR Mix** can be removed from the instrument, capped and stored at -20 °C or below or can be kept on board in the Carousel reagent for up to 7 hours.

As soon as a result is obtained, the **MyGenius PRO** allows users to view, approve, store the results and save the final report.

NOTE

NOTE

The **Q- PCR Standard** can be used for 4 calibrations leaving them on board the instrument for a maximum of 2 hours for each calibration.

NOTE

At the end of the testing session, the remaining **Positive Control** can be removed from the instrument, capped and stored at -20 °C or below. Avoid the spilling of the **Positive Control**. The remaining **Negative Control** must be discarded.

NOTE

The **Positive Control** can be used for 4 calibrations leaving them on board the instrument for a maximum of 3 hours for each calibration.

NOTE

When necessary or when the instrument requires it, remove from the drawer of the instrument the **PCR Cas- sette** and the other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

NOTE

MyGenius PRO can be connected to the "Laboratory Information System" (LIS) which enables downloading the session information. Refer to the instrument manual for more details.

STEP 3 -Review and approval of results

The **MyGenius PRO** monitors target and internal control fluorescence signals for each reaction and automatically applies the Assay Protocol parameters to generate PCR curves which are then interpreted into results.

When the analysis of each sample is complete, the results can be viewed in the 'Results' screen, showing the results and information about the samples. From this screen it is possible to approve the results and save reports. See the instrument manual for more details.

The **MyGenius PRO** generates the results with the **EBV ELITe MGB Kit** through the following procedure:

1. Validation of Calibration curve,
2. Validation of Positive Control and Negative Control results,
3. Validation of sample results,
4. Sample result reporting.

Validation of Calibration curve

The **MyGenius PRO software** interprets the PCR results for the target of the Calibrator reactions with the **EBV ELITe_My_STD** Assay Protocol parameters. The resulting Ct versus concentration produces the Calibration curve.

The Calibration curves, specific for the PCR reagent lot, are recorded in the database (Calibrations). They can be viewed and approved by the user, selecting the calibration of interest and "View chart" section.

The Calibration curve expires **after 60 days**.

NOTE

If the Calibration curve does not meet the acceptance criteria, the “Error” message is shown on the “Calibrations” screen. In this case, the results cannot be approved, and the Calibrator amplification reactions must be repeated. In addition, if samples were included in the testing session, these are not quantified: in this case, after calibration repetition, all results must be interpreted and approved (see section dedicated to samples results).

Validation of amplification Positive Control and Negative Control results

The **MyGenius PRO Software** interprets the PCR results for the target of the Positive Control and for the target and the Internal Control (Channel **IC**) of the Negative Control reactions with the **EBV ELITE_My_PC** and **EBV ELITE_My_NC** Assay Protocols parameters. The resulting Ct values are converted to concentration and used to verify the system (reagents lot and instrument).

The Positive Control and Negative Control results, specific for the PCR reagent lot, are recorded in the database (Controls). They can be viewed and approved by the user selecting the control of interest and “View chart” section.

The Positive Control and Negative Control results expire **after 15 days**.

The **MyGenius PRO Software** processes the Positive Control and Negative Control results and generates Control Charts. The results are analyzed by the software to ensure the system performances are within the acceptance criteria, shown in the Control Chart plots. Refer to the instrument manual for more details.

NOTE

If the Positive Control or Negative Control result does not meet the acceptance criteria, the “Error” message is displayed on the “Controls” screen. In this case, the results cannot be approved, and the Positive Control or Negative Control must be repeated.

NOTE

If the Positive Control or Negative Control result is not valid and samples were included in the same testing session, the samples cannot be approved. In this case, after controls repetition, all results must be interpreted and approved (see section dedicated to samples results).

Validation of Sample results

The **MyGenius PRO Software** interprets the PCR results for the target (Channel **EBV**) and the Internal Control (Channel **IC**) with the **EBV ELITE_My_WB_IU_200_100** and **EBV ELITE_My_WB_cmL_200_100** Assay Protocol parameters. The resulting target Ct values are converted to concentration.

Results are shown in “Results” section.

The sample results can be approved when the three conditions in the table below are true.

1) Calibration Curve	Status
EBV Q-PCR Standard	APPROVED
2) Positive Control	Status
EBV Positive Control	APPROVED
3) Negative Control	Status
EBV Negative Control	APPROVED

The sample results are automatically interpreted by the **MyGenius PRO Software** using Assay Protocol parameters if results of related calibration and controls are approved. If valid calibrations and controls are not present and approved when sample results are obtained, samples cannot be interpreted and approved: in this case, calibrations and controls must be approved and then, each sample must be interpreted in the “sample list” section and approved in the “Results” section using dedicated buttons.

The possible result messages are listed in the table below.

For each sample the system reports a combination of the following messages specifying if the pathogen DNAs are either detected or not detected.

Result of sample run	Interpretation
EBV:DNA DETECTED, QUANTITY BEYOND ULoQ copies/mL OR IU/mL	EBV DNA was detected in the sample within the assay measurement range, its concentration is shown.
EBV:DDNA NOT DETECTED OR BELOW LoDcopies/mL OR IU/mL	EBV DNA was detected in the sample, its concentration is below the assay Lower Limit of Quantification
EBV:DNA DETECTED, QUANTITY BEYOND ULoQ copies/mL OR IU/mL	EBV DNA was detected in the sample, its concentration is above the assay Upper Limit of Quantification
EBV:DNA NOT DETECTED OR BELOW LoDcopies/mL OR IU/mL	EBV DNA was not detected in the sample. The sample is negative for EBV DNA, or its concentration is below the assay Limit of Detection.
INVALID - CHANGE REAGENTS	Not valid assay result , caused by Internal Control failure (due to e.g., incorrect extraction or inhibitors carry-over). The test should be repeated.

Samples reported as "INVALID-CHANGE REAGENTS": in this case, the Internal Control DNA was not efficiently detected, which could be due to problems in sample collection, extraction or PCR steps (e. g. incorrect sampling, degradation or loss of DNA during the extraction, or inhibitors in the eluate), which may cause incorrect results. In this case the sample must be retested starting from extraction of a new sample.

Samples reported as "EBV:DNA NOT DETECTED OR BELOW "LOD" COPIES/ML OR IU/ML" are suitable for analysis, but EBV was not detected. In this case the sample may be either negative for EBV DNA or the EBV DNA is present at a concentration below the Limit of Detection of the assay (see [12 12 PERFORMANCE CHARACTERISTICS WITH page 28](#)).

EBV DNA positive samples at a concentration below the Limit of Detection (and Lower Limit of Quantification) of the assay, if detected, are reported as "EBV:DNA DETECTED, QUANTITY BELOW "LLOQ" COPIES/ML OR IU/ML" (see [12 12 PERFORMANCE CHARACTERISTICS WITH page 28](#)).

EBV DNA positive samples within the Linear Measuring Range are detected and are reported as "EBV:DNA DETECTED, QUANTITY EQUAL TO "XXX" COPIES/ML OR IU/ML" (see [12 12 PERFORMANCE CHARACTERISTICS WITH page 28](#)).

EBV DNA positive samples that are above the Upper Limit of Quantification are reported as "EBV: DNA DETECTED, QUANTITY BEYOND "ULOQ" COPIES/ML OR IU/ML" (see [12 12 PERFORMANCE CHARACTERISTICS WITH page 28](#)), and they are not suitable for quantification. If needed, the sample may be diluted before extraction and retested to yield results within the Linear Measuring Range of the assay.

NOTE

The results obtained with this assay must be interpreted in combination with all relevant clinical observation and laboratory outcomes.

The sample results are stored in the database and, if valid, can be approved by selecting and approving them in "Results" section by the user. From the "Results" window it is possible to print and save the Sample results as "Summary Report" and "Details Report" by selecting the samples of interest and pressing the dedicated button.

Sample result reporting

The sample results are stored in the database and reports can be exported as "Summary Report" and "Details Report".

The "Details Report" shows the detailed results for each selected sample (SID).

The "Summary Report" shows the interpretation results for all the selected samples (SID).

The "Details Report" and "Summary Report" can be printed and signed by authorized personnel.

12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO

12.1 Limit of Detection (LoD)

The Limit of Detection (LoD) of the assay in association with EDTA whole blood, EDTA plasma and CSF matrices was determined on **ELITe InGenius** or **ELITe BeGenius** instruments by testing a panel of EBV negative matrix spiked with reference material of EBV (1st WHO International Standard for Epstein-Barr Virus EBV for Nucleic Acid Amplification Techniques, NIBSC code 09/260, United Kingdom). Probit regression analysis was performed on the results, and the LoD estimated as the concentration corresponding to 95% probability of a positive call.

The results for all matrices are reported in the following tables.

Table 9 Limit of Detection (IU / mL)

Sample volume	Matrix	LoD	95% confidence range	
			lower limit	upper limit
200 µL	EDTA whole blood	104 IU / mL	75 IU / mL	175 IU / mL
	EDTA plasma	124 IU / mL	77 IU / mL	290 IU / mL
	CSF	176 IU / mL	112 IU / mL	346 IU / mL
1000 µL	EDTA plasma	18 IU / mL	13 IU / mL	28 IU / mL

The analytical sensitivity as copies / mL for each matrix is calculated by applying the specific conversion factor reported at paragraph [12.11 Conversion factor to International Units page 37](#)

The analytical sensitivity as copies / mL is reported below.

Table 10 Limit of Detection (copies / mL)

Sample volume	Matrix	LoD	95% confidence range	
			lower limit	upper limit
200 µL	EDTA whole blood	36 copies / mL	26 copies / mL	60 copies / mL
	EDTA plasma	65 copies / mL	41 copies / mL	153 copies / mL
	CSF	80 copies / mL	51 copies / mL	157 copies / mL
1000 µL	EDTA plasma	11 copies / mL	8 copies / mL	17 copies / mL

Sample volume 200 µL: the calculated LoD value was verified for each matrix by testing on ELITe InGenius and ELITe BeGenius a pool of each matrix spiked with EBV certified reference material at the claimed concentration.

The results obtained confirmed the claimed concentration for the target of EBV ELITe MGB Kit on both ELITe InGenius and ELITe BeGenius for each matrix.

Sample volume 1000 µL: the calculated LoD value was verified for Plasma matrix by testing on ELITe InGenius a pool of matrix spiked with EBV certified reference material at the claimed concentration. The results obtained confirmed the claimed concentration for the target of EBV ELITe MGB Kit on ELITe InGenius.



The LoD value calculated on ELITe InGenius and Whole Blood matrix is valid also for MyGenius PRO.

12.2 Linear measuring range and Limits of quantification

The linear measuring range of the assay was determined in association with **ELITe InGenius** and **ELITe BeGenius** using a panel of EBV reference material (Zeptomatrix in association with EDTA Plasma and 1st WHO International Standard for Epstein-Barr Virus (EBV) in association with EDTA Whole Blood and CSF) in EBV DNA - negative matrix.

The results for each matrix are reported in the following paragraphs.

Whole Blood:



The linear measuring range as copies / mL is calculated by applying the specific conversion factor reported at paragraph [12.11 Conversion factor to International Units page 37](#)



The final results are summarized in the following table.

Table 11 Linear measuring range for whole blood samples and ELITE InGenius and ELITE BeGenius

Unit	Lower limit	Upper limit
IU / mL	104	1,000,001
copies / mL	36	344,828

The Linear measuring range determined for ELITE InGenius and Whole Blood matrix is valid also for MyGenius PRO.

Plasma:



The linear measuring range as copies / mL for plasma (Sample Volume 200 µL) is calculated by applying the specific conversion factor reported at paragraph [12.11 Conversion factor to International Units page 37](#)

The final results are summarized in the following table.

Table 12 Linear measuring range for plasma samples (200 µL) and ELITE InGenius and ELITE BeGenius


Unit	Lower limit	Upper limit
IU / mL	124	10,000,000
copies / mL	65	5,263,158

The linear measuring range as copies / mL for Plasma (Sample Volume 1000 µL) is calculated by applying the specific conversion factor reported at paragraph 14.9 Conversion to International Units page 55.

The final results are summarized in the following table:

Table 13 Linear measuring range for plasma samples (1000 µL) and ELITE InGenius

Unit	Lower limit	Upper limit
IU / mL	99	1,000,000
copies / mL	62	625,000

CSF:

The linear measuring range as copies / mL for CSF is calculated by applying the specific conversion factor reported at paragraph 12.11 12.10 Conversion factor to International Units page 34 page 37

The final results are summarized in the following table.

Table 14 Linear measuring range for CSF samples (200 µL) and ELITE InGenius and ELITE BeGenius

Unit	Lower limit	Upper limit
IU / mL	176	1,000,000
copies / mL	80	454,545

12.3 Standard Curve Uncertainty

The Uncertainty value of the Standard curve was calculated by combining the random errors (SD) of all level quantifications and multiplying for the Coverage factor $k = 2$ (Expanded Combined Uncertainty) and is equal to 0.2949 Log copies / reaction.

Table 15

Standard curve levels	Theoretical	Measured	SD	Expanded Combined Uncertainty
	Log c/rxn	Log c/rxn		
EBV Q - PCR Standard 10 ⁵	5.0000	5.0032	0.0636	0.2949
EBV Q - PCR Standard 10 ⁴	4.0000	3.9860	0.0739	
EBV Q - PCR Standard 10 ³	3.0000	3.0184	0.0654	
EBV Q - PCR Standard 10 ²	2.0000	1.9924	0.0892	

12.4 Inclusivity: Efficiency of detection and quantification efficiency on different genotypes

The inclusivity of the assay, as efficiency of detection for different genotypes of EBV was evaluated by *in silico* analysis of the sequences available in nucleotide databases. The analysis showed sequence conservation and absence of significant mutations. So, efficient detection for the different strains or isolates is expected.

12.5 Potential interfering markers: Cross-reactivity

The Potential cross-reactivity of unintended organisms that may be found in clinical specimens was evaluated by *in silico* analysis. The analysis showed no significant homology with other unintended organisms (viruses, bacteria, protozoa and fungi). Therefore, no cross-reactivity should be expected.

The absence of cross-reactivity with potential interfering organisms was also verified through the analysis of the panel of unintended organisms (ATCC, Zeptomatrix) at high titre.

The results are reported in the following table.

Table 16 Potentially interfering organism

Sample	Pos. / Rep.	Outcome
CMV	0/5	No cross-reactivity
HHV6	0/5	No cross-reactivity
HHV7	0/5	No cross-reactivity
HHV8	0/5	No cross-reactivity
HSV1	0/5	No cross-reactivity
HSV2	0/5	No cross-reactivity
VZV	0/5	No cross-reactivity

All potentially interfering markers tested showed no cross-reactivity for the EBV target amplification using the EBV ELITe MGB Kit

12.6 Potential interfering markers: Inhibition

The potential inhibition of unintended organisms that may be found in clinical specimens was evaluated for the assay through the analysis of a panel of unintended organisms in EBV positive samples, from different providers (ATCC, Zeptomatrix).

The results are reported in the following table.

Table 17 Potentially interfering organism

Sample	Pos. / Rep.	Outcome
CMV	5/5	No interference
HHV6	5/5	No interference
HHV7	5/5	No interference
HHV8	5/5	No interference
HSV1	5/5	No interference
HSV2	5/5	No interference
VZV	5/5	No interference

All potential interfering organisms tested showed no interference of the EBV target detection and quantification using the EBV ELITe MGB Kit

12.7 Potential interfering substances: Cross-reactivity

The cross-reactivity by potentially interfering substances (endogenous and exogenous) that might be found in clinical specimens was evaluated for the assay by analysis of a panel of substances at relevant concentration.

The results are reported in the following table.

Table 18 Whole Blood

Samples	EBV Pos. / Rep.	Outcome
Bilirubin	0 / 3	No cross-reactivity
EDTA	0 / 3	No cross-reactivity
Ganciclovir	0 / 3	No cross-reactivity
Azithromycin	0 / 3	No cross-reactivity
Ribavirin	0 / 3	No cross-reactivity
Abacavir Sulfate	0 / 3	No cross-reactivity
Cidofovir	0 / 3	No cross-reactivity
Cyclosporine A	0 / 3	No cross-reactivity
Heparin	0 / 3	Invalid samples
Vancomycin	0 / 5	No cross-reactivity
Ciprofloxacin	0 / 5	No cross-reactivity
Cefpodoxime	0 / 5	No cross-reactivity

The test showed that most of the substances do not cross-react with the EBV target detection and quantification using the EBV ELITe MGB Kit. Heparin was confirmed to be capable of inhibiting the amplification of EBV and of the Internal Control

Table 19 Plasma

Samples	EBV Pos. / Rep.	Outcome
Bilirubin	0 / 5	No cross-reactivity
EDTA	0 / 5	No cross-reactivity
Ganciclovir	0 / 5	No cross-reactivity
Azithromycin	0 / 5	No cross-reactivity
Ribavirin	0 / 5	No cross-reactivity
Abacavir Sulfate	0 / 5	No cross-reactivity
Cidofovir	0 / 5	No cross-reactivity
Cyclosporine A	0 / 5	No cross-reactivity
Heparin	0 / 5	Invalid samples
Vancomycin	0 / 5	No cross-reactivity
Ciprofloxacin	0 / 5	No cross-reactivity
Cefpodoxime	0 / 5	No cross-reactivity

The test showed that most of the substances do not cross-react with the EBV target detection and quantification using the EBV ELITe MGB Kit. Heparin was confirmed to be capable of inhibiting the amplification of EBV and of the Internal Control

12.8 Potential interfering substances: Inhibition

The potential inhibition of interfering substances (endogenous and exogenous) that might be found in clinical specimens was evaluated for the assay by analysis of a panel of substances at relevant concentration in EBV positive samples.

The results, for each matrix, are reported in the following paragraphs.

Table 20 Whole Blood

Samples	EBV Pos. / Rep	Outcome
Bilirubin	3 / 3	No interference
EDTA	3 / 3	No interference
Ganciclovir	3 / 3	No interference
Azithromycin	3 / 3	No interference
Ribavirin	3 / 3	No interference
Abacavir Sulfate	3 / 3	No interference
Cidofovir	3 / 3	No interference
Cyclosporine A	3 / 3	No interference
Heparin	3 / 3	No interference
Vancomycin	5 / 5	No interference

Table 20 Whole Blood (continued)

Samples	EBV Pos. / Rep	Outcome
Ciprofloxacin	5 / 5	No interference
Cefpodoxime	5 / 5	No interference

The test showed that all the substances do not interfere with the EBV target detection and quantification using the EBV ELITe MGB Kit.

Table 21 Plasma

Sample	EBV Pos. / Rep.	Outcome
Bilirubin	5 / 5	No interference
EDTA	5 / 5	No interference
Ganciclovir	5 / 5	No interference
Azithromycin	5 / 5	No interference
Ribavirin	5 / 5	No interference
Abacavir Sulfate	5 / 5	No interference
Cidofovir	5 / 5	No interference
Cyclosporine A	5 / 5	No interference
Heparin	0 / 5 invalid samples	Interference
Vancomycin	5 / 5	No interference
Ciprofloxacin	5 / 5	No interference
Cefpodoxime	5 / 5	No interference

The test showed that all the substances, with the exception for heparin, do not interfere with the EBV target detection and quantification using the EBV ELITe MGB Kit.

Table 22 CSF

Sample	EBV Pos. / Rep.	Outcome
Glucose	5 / 5	No interference
Lactate	5 / 5	No interference
Albumin	5 / 5	No interference
Immunoglobulin	5 / 5	No interference
Hemoglobin	5 / 5	No interference
White blood cells	5 / 5	No interference
Whole Blood	5 / 5	No interference
Ethanol	5 / 5	No interference
Cefpodoxime	5 / 5	No interference
Ampicillin	5 / 5	No interference

Table 22 CSF (continued)

Sample	EBV Pos. / Rep.	Outcome
Azithromycin	5 / 5	No interference
Vancomycin	5 / 5	No interference
Gentamicin	5 / 5	No interference
Meropenem	5 / 5	No interference
Ganciclovir	5 / 5	No interference
Acyclovir	5 / 5	No interference
Fluconazole	5 / 5	No interference
Cyclosporine	5 / 5	No interference

The test showed that all the substances do not interfere with the EBV target detection and quantification using the EBV ELITe MGB Kit.

12.9 Repeatability

The Intra-Session and Inter-Session Repeatability of the assay was evaluated on ELITe InGenius and ELITe BeGenius by analysis of a panel of Whole Blood samples collected in EDTA, including one negative sample and two samples spiked by EBV certified reference material (1st WHO EBV International Standard, NIBSC)

An example of Intra-Session Repeatability (on one day) results is shown in the tables below.

Table 23 Intra – Session Repeatability on ELITe InGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	% Agreement
Negative	8	-	-	-	100%
3 x LoD	8	35.68	0.57	1.60	100%
10 x LoD	8	34.22	0.24	0.70	100%

Table 24 Intra – Session Repeatability on ELITe BeGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	% Agreement
Negative	8	-	-	-	100%
3 x LoD	8	37.32	0.49	1.30	100%
10 x LoD	8	35.97	0.43	1.19	100%

An example of Inter-Session Repeatability (on two days) results is shown in the tables below.

Table 25 Inter – Session Repeatability on ELITe BeGenius

Sample	EBV ELITe MGB Kit - Days 1-2				
	N	Mean EBV Ct	SD EBV Ct	% CV EBV Ct	% Agreement
Negative	14	-	-	-	100%
3 x LoD	16	37.29	0.67	1.79	100%
10 x LoD	16	35.82	0.67	1.86	100%

Table 26 Inter – Session Repeatability on ELITe InGenius

Sample	EBV ELITe MGB Kit - Days 1-2				
	N	Mean EBV Ct	SD EBV Ct	% CV EBV Ct	% Agreement
Negative	14	-	-	-	100%
3 x LoD	16	35.72	0.53	1.48	100%
10 x LoD	16	34.39	0.37	1.07	100%

In the Repeatability test, the EBV ELITe MGB Kit detected all the samples as expected and showed a maximum variability of target Ct values as %CV lower than 5%.

12.10 Reproducibility

The Reproducibility of the assay was evaluated on ELITe BeGenius and ELITe InGenius by analysis of a panel of Whole blood samples collected in EDTA, including one negative sample and two samples spiked by reference material (1st WHO International Standard, NIBSC).

A summary of Inter-Instrument Reproducibility (on two instruments) is shown in the tables below.

Table 27 Inter–Instrument Reproducibility on ELITe InGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	% Agreement
Negative	8	-	-	-	100%
3 x LoD	8	35.78	0.44	1.24	100%
10 x LoD	8	30.38	0.36	1.17	100%

Table 28 Inter–Instrument Reproducibility on ELITe BeGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	% Agreement
Negative	7	-	-	-	100%
3 x LoD	8	36.79	0.86	2.32	100%
10 x LoD	8	35.15	0.65	1.84	100%

A summary of Inter-batch Reproducibility (on two lots) is shown in the tables below:

Table 29 Inter–Batch Reproducibility on ELITe InGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	%Agreement
Negative	8	-	-	-	100%
3 x LoD	8	35.91	0.38	1.06	100%
10 x LoD	8	34.48	0.15	0.43	100%

Table 30 Inter–Batch Reproducibility on ELITe BeGenius

Sample	EBV ELITe MGB Kit				
	N	Mean Ct	SD	% CV	%Agreement
Negative	7	-	-	-	100%
3 x LoD	8	37.45	0.65	1.72	100%
10 x LoD	8	35.57	0.42	1.18	100%

In the Reproducibility test, the EBV ELITe MGB Kit detected all the samples as expected and showed a maximum variability of target Ct values as %CV lower than 5%.

12.11 Conversion factor to International Units

The conversion factor to report the quantitative results in International Units / mL starting from copies / mL, was calculated, for each matrix, using the certified calibrated reference material “1st WHO International Standard for Epstein-Barr Virus (EBV) for Nucleic Acid Amplification Techniques” (NIBSC).

The results for each matrix are summarized in the following table


Table 31 Conversion factor to International Units

Sample volume	Matrix	Fc (IU / copies)
200 µL	whole blood	2.9
	plasma	1.9
	CSF	2.2
1000 µL	plasma	1.6


The Conversion factor to International Units calculated for ELITe InGenius and Whole Blood matrix is valid also for MyGenius PRO.

Sample volume 200 µL: The Conversion factor, to report the quantitative results in International Units / mL starting from copies / mL, was verified on **ELITE InGenius** and **ELITE BeGenius** instruments using the certified calibrated reference material (1st WHO International Standard, NIBSC). The results obtained were analysed by orthogonal and linear regression in order to calculate the correlation.

The results, for each matrix, are reported in the following paragraphs.


Whole Blood

The Orthogonal Regression analysis generated an intercept equal to -0.307 (95% CI: - 0.4419; – -0.1725) and a slope equal to 1.020 (95% CI: 0.9879; - 1.0512). The linear regression analysis generated an R2 of 0.991.

Plasma:

The Orthogonal Regression analysis generated an intercept equal to -0.176 (95% CI: - 0.496 – 0.112) and a slope equal to 1.044 (95% CI: 0.981 - 1.121). The linear regression analysis generated an R2 of 0.984.

CSF:




The Orthogonal Regression analysis generated an intercept equal to -0.285 (95% (-0.4623; -0.1067) and a slope equal to 1.050 (95% CI: 1.0076; 1.0925). The linear regression analysis generated an R2 of 0.984.

12.12 MyGenius PRO: method correlation

The correlation analysis of different methods was evaluated on MyGenius PRO by analysing EBV samples from patients whose viral load was within the measuring range of the reference method (ELITe InGenius). The results obtained with the MyGenius PRO and the reference method (ELITe InGenius) were analysed by Orthogonal and Linear Regression.

The correlation study was performed at one site on 136 positive clinical samples of whole blood collected in EDTA certified positive for EBV DNA or spiked with reference material using the ELITe InGenius as comparator.

The results are summed up in the following figure.



The Orthogonal Regression analysis generated an intercept equal to -0.0767 (95% CI: -0.2726; 0.1190) and a slope equal to 1.000 (95% IC: 0.9537; 1.0465). The linear regression analysis generated an R2 of 0.930.

12.13 Diagnostic Specificity: confirmation of negative samples

The Diagnostic Specificity of the assay, as confirmation of negative samples, was evaluated in association with **ELITE InGenius** or **ELITE BeGenius**, analysing negative clinical samples of whole blood collected in EDTA, plasma collected in EDTA and CSF.

In case of extraction volume of 200 µL, as EBV products showed equivalent analytical performances on both **ELITE BeGenius** and **ELITE InGenius**, the diagnostic performances of the assay performed on the two instruments are also considered equivalent. Therefore, the Diagnostic Specificity of the assay obtained in association with ELITE InGenius is also applicable to ELITE BeGenius and vice versa.

On the other hand, in case of an extraction volume of 1000 µL, the showed results are just related to ELITE InGenius instrument.

As **MyGenius PRO** showed equivalent analytical performances to ELITE InGenius in association to EDTA whole blood matrix, the diagnostic performances of the assay performed on the two instruments and EDTA whole blood matrix are also considered equivalent. Therefore, the Diagnostic Specificity of the assay obtained in association with ELITE InGenius is also applicable to MyGenius PRO.

The results are summed up in the following table.

Table 32 Diagnostic Specificity

Samples	Sample volume	N	Positive	Negative	% Diagnostic Specificity
Whole blood collected in EDTA and negative for EBV DNA	200 µL	110	3	107	97.3
Plasma collected in EDTA and negative for EBV DNA	200 µL	113	1	112	98.3
	1000 µL	62	2	60	
CSF negative for EBV DNA	200 µL	77	0	77	100

The IC Ct cut-off value is set at 35 for whole blood samples collected in EDTA when tested with ELITE InGenius, ELITE BeGenius and MyGenius PRO.

The IC Ct cut-off value is set at 35 for plasma samples collected in EDTA and CSF when tested with ELITE InGenius and ELITE BeGenius.

12.14 Diagnostic Sensitivity: confirmation of positive samples

The Diagnostic Sensitivity of the assay, as confirmation of positive clinical samples, was evaluated in association with **ELITE InGenius** or **ELITE BeGenius** analysing positive or spiked clinical samples of whole blood collected in EDTA, plasma collected in EDTA and CSF.

In case of extraction volume of 200 µL, as EBV products showed equivalent analytical performances on both **ELITE BeGenius** and **ELITE InGenius**, the diagnostic performances of the assay performed on the two instruments are also considered equivalent. Therefore, the Diagnostic Sensitivity of the assay obtained in association with ELITE InGenius is also applicable to ELITE BeGenius and vice versa.

On the other hand, in case of an extraction volume of 1000 µL the showed results are just related to ELITE InGenius instrument.

As MyGenius PRO showed equivalent analytical performances to ELITE InGenius in association to EDTA whole blood matrix, the diagnostic performances of the assay performed on the two instruments and EDTA whole blood matrix are also considered equivalent. Therefore, the Diagnostic Sensitivity of the assay obtained in association with ELITE InGenius is also applicable to MyGenius PRO.

The results are summed up in the following table.

Table 33 Diagnostic Sensitivity

Samples	Sample volume	N	Positive	Negative	% Diagnostic Sensitivity
Whole blood collected in EDTA and positive for EBV DNA	200 µL	133	131	2	98.5
Plasma collected in EDTA and positive for EBV DNA	200 µL	12	12	0	98.9
Plasma collected in EDTA and spiked for EBV DNA	200 µL	46	46	0	
	1000 µL	30	29	1	
CSF positive for EBV DNA	200 µL	4	3	1	98.8
CSF spiked for EBV DNA		80	80	0	

NOTE

The complete data and results of the tests carried out to evaluate the product performance characteristics with matrices and instruments are recorded in the Product Technical File for the "EBV ELITE MGB® Kit", FTP020PLD.

13 SPECIMENS AND CONTROLS FOR ABI 7500 Fast Dx Real-Time PCR Instrument

13.1 Specimens

The following specimens and nucleic acid extraction methods are validated for use with the **EBV ELITE MGB Kit** using the ABI 7500 Fast Dx Real-Time PCR Instrument.

Table 34

Specimen type	Kit/Method	Protocol	Input volume (µL)	Elution volume (µL)	Primary tube minimum volume (µL)	Special instruction
Whole blood	ELITE GALAXY	xNA Extraction (Universal)	300	200	400-650	Add 10 µL/ sample of CPE to the IC + Carrier solution

13.2 Interfering substances

The DNA extracted from the sample must not contain heparin, haemoglobin, dextran, Ficoll®, ethanol or 2-propanol in order to prevent inhibition problems and the possibility of frequent invalid results.

High quantity of human genomic DNA in the DNA extracted from the sample may inhibit the amplification reaction.

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

Do not use samples collected in heparin, which is a known reverse transcription and PCR inhibitor

13.3 Amplification controls

It is mandatory to validate each amplification session with a Negative Control reaction and a Positive Control reaction.

For the Negative Control, use molecular biology grade water (not provided with this kit) added to the reaction in place of the DNA extracted from the sample.

For the Positive Control, use the **EBV - ELITe Positive Control** product or the **EBV - ELITe Standard** product.

13.4 Quality controls

Verification of the extraction and PCR procedure is recommended. Archived samples or certified reference material may be used. External controls should be used in accordance with local, state, and federal accrediting organizations, as applicable.

14 ABI 7500 Fast Dx Real-Time PCR Instrument PROCEDURE

14.1 Setting of the real time amplification session


Before starting the session, refer to the instrument documentation to:

- switch on the instrument, switch on the computer, open the dedicated software and an "absolute quantification" session and set "Run mode: Fast 7500";
- set (Detector Manager) the "detector" for the EBV probe with the "reporter" = "FAM" and the "quencher" = "none" (non fluorescent) and label it "EBV";
- set (Detector Manager) the "detector" for the internal control probe with the "reporter" = "VIC" (AP525 is analogous to VIC) and the "quencher" = "none" (non fluorescent) and label it "IC";
- for each well in use in the microplate, set (Well Inspector) the "detector" (fluorescence to be measured), the "passive reference" = "Cy5" (AP593 is used instead of Cy5 for normalisation of fluorescence levels), and the type of reaction (sample, negative control, positive control or known quantity standard).

NOTE

To quantify the DNA in the starting sample, include a series of reactions using the **Q-PCR Standards** (10^5 copies, 10^4 copies, 10^3 copies, 10^2 copies) to obtain the **Standard curve**.

See below, by way of example, how to set up the quantitative analysis of 12 samples.



Legend: S1 -S12: Samples to be analysed; **NC:** Negative Control of amplification;

10²: 10^2 standard copies; **10³:** 10^3 standard copies; **10⁴:** 10^4 standard copies; **10⁵:** 10^5 standard copies.

Refer to the instrument documentation to set the **thermal cycling** parameters (Instrument > Thermal Cycler Protocol > Thermal Profile):

- add a **20 second extension at 72 °C** (Add Step);

NOTE

Note: Fluorescence acquisition must be set during the 60°C hybridization step (Instrument > Thermal Cycler Protocol > Settings > Data Collection).

- modify thermal cycling temperatures and times as indicated in the table "**Thermal cycle**";
- set the number cycles at **45**
- set the sample volume at **30 µL**
- optional: add a dissociation stage (Add Dissociation Stage), set the starting temperature at **40 °C** and the ending temperature at **80 °C**.

Table 35 Thermal cycle

Stage	Temperatures	Timing
Decontamination	50 °C	2 min.
Initial denaturation	94 °C	2 min.
Amplification and detection (45 cycles)	94 °C	10 sec.
	60 °C (fluorescence acquisition)	30 sec.
	72 °C	20 sec.
Dissociation (optional)	95 °C	15 sec.
	40 °C	1 min.
	80 °C	15 sec.
	60 °C	15 sec.

14.2 Real-time PCR session set-up

(Performed by the **ELITe GALAXY** instrument)

To perform the PCR session set up:

- thaw the **Q-PCR Mix** tubes required for the session (each tube is sufficient for **25 reactions**)
- thaw the **Positive Control** (qualitative analysis: detection of extracted DNA) or the **Q - PCR Standard** (quantitative analysis: quantification of extracted DNA) tubes
- mix gently the reagents and spin down the contents for 5 seconds
- prepare the **Negative Control** (not provided) as per the instruction of use of the instrument
- prepare a **Q-PCR microplate**. Handle it with powderless gloves and do not damage the wells

NOTE

To prepare the PCR on the **ELITe GALAXY**, load the elution microplate containing the extracted DNA samples, the reagents and the **Q-PCR microplate** as indicated in the instrument user manual and follow the steps on the GUI.

The instrument automatically performs the PCR set-up dispensing in each well of the **Q-PCR microplate**:

- **20 µL** of **Q-PCR Mix**
- **20 µL** of **extracted DNA / Q-PCR Standard / Controls**

NOTE

If not all the Q—PCR Mix is used, store the remaining volume in the dark at -20°C for no longer than one month. Freeze and thaw the Q—PCR Mix for a maximum of **5 TIMES**.

After the PCR set-up performed by the instrument:

- seal the **Q-PCR microplate** with an optical seal

- transfer the **Q-PCR microplate** onto the **7500 Fast Dx Real-Time PCR Instrument** and start the PCR. Save the run file with a unique and recognizable name (e.g. "year-month-day-TARGET-EGSpA")

NOTE

At the end of the PCR the **Q-PCR microplate** must be discarded following all governmental and environmental regulations. In order to avoid spilling the PCR products, the **optical seal must not be removed from the Q-PCR microplate**.

14.3 General settings for analysis of results

Before starting the analysis, refer to the instrument documentation to:

- manually adjust the calculation range for the **Baseline** (fluorescence background level) from cycle 6 to cycle 15 (Results > Amplification plot > delta Rn vs Cycle);

NOTE

The FAM fluorescence of the EBV probe in a sample with a high concentration of EBV DNA may begin to increase before the cycle 15. In this case, lower the **Baseline** calculation range to the cycle at which the FAM fluorescence of the sample begins to increase (Results > Component).

- manually set the thresholds for the detectors:

set the FAM detector "EBV" **threshold** at **0.2**;

set the VIC detector "IC" **threshold** at **0.1**.

The PCR cycle at which a sample's fluorescence level reaches the **threshold** value determines the **threshold cycle (Ct)** for that sample.

The instrument software automatically analyses the fluorescence levels in the controls, standards and sample reactions and calculates Ct values.

14.4 Qualitative analysis of results

The EBV **Ct** value of the **Positive Control** is used to validate the PCR. The PCR run is valid when results are as described in the following table:

Table 36

Positive Control reaction Detector FAM "EBV "	Assay result	Amplification / Detection
Ct ≤ 25	POSITIVE	CORRECT

If the result of the **Positive Control** is **Ct > 25** or **Ct Undetermined** for Detector FAM "EBV", the session is not valid and must be repeated starting from the PCR step. This may indicate an issue during the PCR setup, the PCR or the detection step (e.g., incorrect dispensation or degradation of the Q-PCR Mix or positive control, incorrect placement of the positive control, incorrect thermal cycle settings), which could lead to incorrect results.

NOTE

When the product is used for the quantification of EBV DNA, the **Q - PCR Standard** reactions were set up instead of the **Positive Control** reaction. In this case, validate the amplification and the detection by referring to the amplification reaction of **Q - PCR Standard 10⁵** (**Ct ≤ 25**).

The EBV Ct value of the **Negative Control** is used to validate the PCR. The PCR run is valid when results are as described in the following table:

Table 37

Negative control reaction Detector FAM "EBV"	Assay result	Amplification / Detection
Ct Undetermined	NEGATIVE	CORRECT

If the result of the **Negative control** amplification reaction is different from **Ct Undetermined** for Detector FAM "EBV", the session is not valid and must be repeated starting from the PCR step. This may indicate a contamination which could lead to incorrect and false positives results.

The **Ct** value of EBV in each sample is used to detect the target DNA while the **Ct** value of the internal control is used to validate the extraction, PCR, and detection.

NOTE

Verify by the amplification plot (Results > Amplification plot > delta Rn vs Cycle) that the **Ct** of each sample was determined by a fast and regular increase in fluorescence and not by peaks or an increase in background signal (irregular or high background).

Possible sample results (Results > Report) are described in the following table:

Table 38

Sample reaction		Sample suitability	Assay sample result	EBV DNA
Detector FAM "EBV"	Detector VIC "IC"			
Ct Undetermined	Ct > 35 or Ct Undetermined	unsuitable	invalid	-
	Ct ≤ 35	suitable	valid, negative	NOT DETECTED
Ct Determined	Ct > 35 or Ct Undetermined	suitable	valid, positive	DETECTED
	Ct ≤ 35	suitable	valid, positive	DETECTED

A sample result of **Ct Undetermined** for EBV and **Ct > 35** or **Ct Undetermined** for the internal control is invalid and indicates an issue during nucleic acid extraction or PCR (e.g., degradation of sample DNA, loss of DNA during extraction, presence of inhibitors in the DNA, inefficient or absent amplification), which may lead to incorrect results. The sample is not suitable for the analysis and the assay needs to be repeated starting from nucleic acid extraction of a new sample.

A sample result of **Ct Undetermined** for EBV and **Ct ≤ 35** for the internal control is a valid result and indicates that EBV DNA was not detected in the sample. The sample may contain no EBV DNA or it contains EBV DNA at a concentration lower than the detection limit of the product (see [15 Performance Characteristics page 47](#)). A sample result of **Ct Determined (Ct ≤ 45)** for EBV and **Ct > 35**, **Ct Undetermined**, or **Ct ≤ 35** for the IC is a valid result and indicates that EBV DNA was detected in the sample

NOTE

In case of Ct Determined for EBV and Ct > 35 or Undetermined for the IC, the PCR efficiency of the IC may have been impacted by competition with the high PCR efficiency of the EBV DNA. In this case the sample is suitable, and the positive result is valid.

NOTE

The results obtained with this assay must be interpreted in combination with all relevant clinical observation and laboratory outcomes.

14.5 Quantitative analysis of the results

In the amplification reactions of the four **Q - PCR standards**, the EBV **Ct** values are used to calculate the **Standard Curve** (Results > Standard Curve) for the amplification session and to validate the amplification and the detection as described in the following table:

Table 39

Standard Curve Detector FAM " EBV "	Acceptability range	Amplification / Detection
Correlation coefficient (R2)	$0.990 \leq R2 \leq 1.000$	CORRECT

If the **Correlation coefficient (R2)** value does not fall within the limits, the session is not valid and must be repeated starting from the PCR step. This may indicate an issue during the PCR or the detection step (e. g., incorrect dispensation or degradation of the Q-PCR Mix or of the standards, incorrect placement of the standards, incorrect thermal cycle settings or cross-contamination), which could lead to incorrect results.

Table 40

Sample result for detector FAM "EBV"	EBV copies per reaction
Quantity > 1×10^6	MORE THAN 1×10^6
$1 \times 10^1 \leq \text{Quantity} \leq 1 \times 10^6$	= Quantity
Quantity < 1×10^1	LESS THAN 10

The results (**Quantity**) of each sample (Results > Report) are used to calculate the copies of EBV present in the specimen used in the extraction (**Nc**) according to this formula:

Table 41

$$Nc = \frac{Ve \times \text{Quantity}}{Vc \times Va \times Ep}$$

where:

Ve is the total volume in μL of the extracted DNA sample (elution volume)

Quantity is the **copies/reaction** of the sample calculated by the instrument software (PCR result)

Vc is the volume of the specimen used for nucleic acid extraction (input volume) expressed in the required unit of measurement

Va is the volume in μL of the extracted DNA sample (eluate) used in the PCR

Ep is the efficiency of the procedure (extraction and PCR) **expressed as a decimal**

To convert the sample quantity from copies/mL to IU/mL, multiply the copies/mL value by the **conversion factor (Fc)**. The Fc was calculated using calibrated certified reference material ("1st WHO International Standard for Epstein Barr virus for Nucleic Acid Amplification Techniques", NIBSC) (See [15 Performance Characteristics page 47](#)).

For convenience, the following are simplified formulas in which $Ve/(Vc \times Va \times Ep)$ and its conversion to IU/mL have been calculated.

Table 42

Matrix	Nucleic acid extraction method	$Ve/(Vc \times Va \times Ep)$	Formula to quantify Nc (copies/mL)	Fc (IU/copy)	Formula to quantify Nc (IU/mL)
Whole blood	ELITe GALAXY	35	$35 \times \text{Quantity}$	0.82	$31.2 \times \text{Quantity}$

15 PERFORMANCE CHARACTERISTICS WITH ABI 7500 Fast Dx Real-Time PCR Instrument

15.1 Limit of detection (LoD)

The Limit of Detection (LoD) of the assay in association to Whole Blood was tested on the ELITe GALAXY and ABI 7500 instruments, using a panel of EBV negative matrix spiked with certified reference material (1st WHO International Standard for Epstein Barr Virus for Nucleic Acid Amplification Technique, NIBSC, United Kingdom). Probit regression analysis was performed on the results, and the LoD estimated as the concentration corresponding to 95% probability of a positive call.

Table 43 Limit of Detection with ELITe GALAXY (IU / mL)

Matrix	95% positivity	95% confidence range	
		lower limit	upper limit
whole blood	99 IU / mL	57 IU / mL	376 IU / mL

The LoD as copies / mL for each matrix is calculated by applying the specific conversion factor reported at paragraph [15.3 Conversion to International Units page 48](#).

The analytical sensitivity as copies / mL is reported below


Table 44 Limit of Detection with ELITe GALAXY (copies / mL)

Matrix	95% positivity	95% confidence range	
		lower limit	upper limit
whole blood	121 copies/ mL	70 copies / mL	459 copies / mL

15.2 Linear measuring range

The linear measuring range of the assay was determined in association with Whole Blood on **ELITe GALAXY** and **ABI 7500** using a panel of dilution of a plasmid DNA containing the amplification product.

The results on Whole Blood are reported in the following graph



The linear measuring range as copies / mL is calculated by applying the specific conversion factor reported at paragraph [15.3 Conversion to International Units page 48](#).

The final results are summarized in the following table.

Table 45 Linear measuring range for whole blood samples

Unit of Measure	lower limit	upper limit
copies / reaction	10	1,000,000

15.3 Conversion to International Units

The conversion factor to report the quantitative results in International Units / mL starting from copies / mL, was calculated, for Whole blood collected in EDTA, using the certified calibrated reference material "1st WHO International Standard for Epstein Barr virus for Nucleic Acid Amplification Techniques"(NIBSC).

The results are summarized in the following table.

Table 46 Conversion factor to International Units with ABI 7500 and Whole blood

Instrument	Fc (IU / copies)
ELITe GALAXY	0.82

15.4 Inclusivity: Efficiency of detection and quantification efficiency on different genotypes

The inclusivity of the assay, as efficiency of detection for different genotypes of EBV was evaluated by *in silico* analysis of the sequences available in nucleotide databases. The analysis showed sequence conservation and absence of significant mutations. So, efficient detection for the different strains or isolates is expected.

15.5 Potential interfering organisms: Cross-reactivity

The Potential cross-reactivity of unintended organisms that may be found in clinical specimens was evaluated by *in silico* analysis. The analysis showed no significant homology with other unintended organisms (viruses, prokaryotes, fungi, phages, invertebrates and human). Therefore, no cross-reactivity should be expected.

15.6 Repeatability

The Within-run Repeatability of the assay was evaluated on ABI 7500 by analysis of a panel of samples spiked by a plasmid DNA containing the amplification product and one negative sample.

An example of Within-run Repeatability results is shown in the table below.

Table 47 Within —run Repeatability on ABI 7500

Sample	EBV				
	N	Mean Ct	SD	% CV	% Agreement
10,000,000	9	15.16	0.11	0.70	100%
1,000,000	9	18.52	0.14	0.75	100%
100,000	9	21.93	0.21	0.94	100%
10,000	9	25.21	0.09	0.35	100%
1,000	9	28.50	0.19	0.67	100%
100	9	31.88	0.23	0.72	100%
10	9	35.42	0.70	1.97	100%

In the Repeatability test, the EBV ELITe MGB Kit was detected for all the samples as expected and demonstrated a maximum variability of target Ct values with %CV less than 5%.

15.7 Reproducibility

The Reproducibility of the assay was evaluated on ABI 7500 by analysis of a panel of samples spiked by a plasmid DNA containing the amplification product and one negative sample.

A summary of Inter-batch Reproducibility (on five lots) is shown in the table below:

Table 48 Inter-batch Reproducibility on ABI 7500

Sample	EBV				
	N	Mean Ct	SD	% CV	% Agreement
50,000 copies / rxn	15	23.45	0.3793	1.6172	100%
5,000 copies / rxn	15	26.73	0.4545	1.7001	100%
500 copies / rxn	15	29.98	0.5477	1.8270	100%
Negative	15	-	-	-	100%

In the Reproducibility test, the EBV ELITe MGB Kit was detected for all the samples as expected and demonstrated a maximum variability of target Ct values with %CV less than 5%.

15.8 Diagnostic sensitivity: detection and quantification efficiency on different genotypes / subtypes

The diagnostic sensitivity of the assay, as detection and quantification efficiency on different genotypes / subtypes, was evaluated by comparison of sequences with nucleotide databases.

The analysis of the regions chosen for the hybridisation of the primers and of the fluorescent probe in the alignment of the sequences available in the database for the EBNA-1 gene of EBV showed conservation and absence of significant mutations.

15.9 Diagnostic sensitivity: confirmation of positive samples

The diagnostic sensitivity of the assay, as confirmation of positive clinical samples, was evaluated by analyzing, in association with **ELITe GALAXY** and ABI 7500 Fast Dx Real-Time PCR Instrument, samples certified positive for the target.

The results are summed up in the following table

Table 49

Samples	N	positive	negative	% Diagnostic Sensitivity
Whole blood collected in EDTA and positive for EBV DNA	54	54	0	100

NOTE

The complete data and results of the tests carried out to evaluate the product performance characteristics with matrices and instruments are recorded in the Product Technical File "EBV ELITe MGB Kit", FTP 020PLD.

15.10 Diagnostic specificity: confirmation of negative samples

The diagnostic specificity of the assay, as confirmation of negative clinical samples, was evaluated by analyzing, in association with **ELITe GALAXY** and ABI 7500 Fast Dx Real-Time PCR Instrument, samples certified negative for the target.

The results are summed up in the following table.

Table 50 Diagnostic specificity

Samples	N	positive	negative	% Diagnostic Specificity
Whole blood collected in EDTA and negative for EBV DNA	52	2	50	96.1

16 REFERENCES

- S. W. Aberle et al. (2002) *J Clin Virology* 25: S79 - S85
- C. N. Kotton et al. (2025) *Transplantation* 109: 1066-1110
- E. A. Lukhtanov et al. (2007) *Nucleic Acids Res.* 35: e30

17 PROCEDURE LIMITATIONS

Use this product only with the following clinical samples:

- whole blood collected in EDTA (all the instruments),
- plasma collected in EDTA (ELITe InGenius, ELITe BeGenius only).
- cerebrospinal fluid (CSF) (ELITe InGenius, ELITe BeGenius only).

Currently there are no data available concerning product performance with other clinical samples such as: suspensions of leukocytes, suspensions of granulocytes.

Plasma collected in EDTA shall be obtained from whole blood stored at room temperature or +2 / +8 °C for no longer than 24 hours.

Do not use DNA extracted from heparinized samples with this product: heparin inhibits the amplification reaction of nucleic acids and causes invalid results.

Do not use extracted DNA that is contaminated with haemoglobin, dextran, Ficoll®, ethanol or 2-propanol with this product: these substances inhibit the amplification reaction of nucleic acids and may cause invalid results.

Do not use with this product extracted DNA containing high quantity of human genomic DNA that may inhibit the amplification reaction of nucleic acids.

There is no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic, or immunosuppressant drugs.

The results obtained with this product depend on proper identification, collection, transport storage and processing of the samples. To avoid incorrect results, it is therefore necessary to take care during these steps and to carefully follow the instructions for use provided with the product.

Due to its high analytical sensitivity, the Real Time PCR method is susceptible to contamination from positive samples, positive controls, and PCR products. Such cross-contamination may cause false positive results. The product format is designed to reduce this, but only good laboratory practices and compliance with these Instructions for Use can prevent it.

This product must be handled by qualified personnel trained in the processing of potentially infective biological samples and chemical preparations classified as dangerous to prevent accidents with potentially serious consequences for the user and other people.

This product requires the use of personal protective equipment and areas that are suitable for the processing of potentially infective biological samples and chemical preparations classified as dangerous to prevent accidents with potentially serious consequences for the user and other people.

This product requires the use of personal protective equipment and instruments dedicated to work session setup to avoid false positive results.

To avoid incorrect results, this product must be handled by professional personnel, qualified and trained in molecular biology techniques such as extraction, PCR and detection of nucleic acids.

Due to inherent differences between technologies, it is recommended that users perform method correlation studies to estimate technology differences prior to switching to a new technology.

A negative result obtained with this product indicates that the target DNA is not detected in the DNA extracted from the sample; however, it cannot be excluded that the target DNA has a lower titer than the product detection limit (see [12 PERFORMANCE CHARACTERISTICS WITH ELITe InGenius, ELITe BeGenius and MyGenius PRO page 28](#)). In this case the result could be a false negative.

Results obtained with this product may sometimes be invalid due to failure of internal control. In this case the sample shall be retested, starting from extraction, which can lead to a delay in obtaining final results.

Possible polymorphisms, insertions or deletions within the region of the DNA targeted by the product primers and probes may impair detection and quantification of target DNA.

As with any other diagnostic medical device, the results obtained with this product must be interpreted in combination with all relevant clinical observations and laboratory results.

As with any other diagnostic medical device, there is a residual risk of obtaining invalid or erroneous results with this product. This residual risk cannot be eliminated or further reduced. In some cases, this residual risk could contribute to wrong decisions with potentially dangerous effects for the patient. However, this residual risk associated to the intended use of the product has been weighed against the potential benefits to the patient and it has been assessed acceptable.

18 TROUBLESHOOTING

ELITe InGenius and ELITE BeGenius

Table 51

Invalid Q-PCR Standard reaction, Standard curve or Positive Control reaction	
Possible Causes	Solutions
Instrument setting error.	Check the position of Q-PCR Mix, Q-PCR Standards and Positive Control. Check the volumes of Q-PCR Mix, Q-PCR Standards and Positive Control.
PCR Mix degradation.	Do not use the Q-PCR Mix for more than 5 independent sessions (3 hours each in the Inventory Area, Cool Block or in the Cooler Unit). Do not use the Q-PCR Mix for more than 3 consecutive sessions (7 hours in the Inventory Area Cool Block or in the Cooler Unit) Do not leave the Q-PCR Mix at room temperature for more than 30 minutes. Use a new tube of Q-PCR Mix.
Q-PCR Standards or Positive Control degradation.	Do not use the Q-PCR Standard for more than 4 independent sessions (2 hours each in the Extraction Area or in the Cooler Unit). Do not use the Positive Control for more than 4 independent sessions (3 hours each in the Extraction Area or in the Cooler Unit). Use new tubes of Q-PCR Standards or Positive Control.
Instrument error.	Contact ELITechGroup Technical Service.

Table 52

Invalid Negative Control reaction	
Possible Causes	Solutions
Instrument setting error.	Check the position of Q-PCR Mix and Negative Control. Check the volumes of Q-PCR Mix and Negative Control.
Contamination of the Negative Control.	Do not use the Negative Control for more than 1 session. Use a new aliquot of molecular biology grade water.
Contamination of the PCR Mix.	Use a new tube of Q-PCR Mix.
Contamination of the extraction area, Racks, Inventory Block or Cooler Unit.	Clean surfaces with aqueous detergents, wash lab coats, replace tubes and tips in use.
Instrument error.	Contact ELITechGroup Technical Service.

Table 53

Invalid Sample reaction	
Possible Causes	Solutions
Instrument setting error.	Check the position of Q-PCR Mix, Internal Control, and sample. Check the volumes of Q-PCR Mix, Internal Control and sample.
PCR Mix degradation.	Do not use the Q-PCR Mix for more than 5 independent sessions (3 hours each in the Inventory Area or in the Cooler Unit). Do not use the Q-PCR Mix for more than 3 consecutive sessions (7 hours in the Inventory Area Cool Block or in the Cooler Unit). Do not leave the Q-PCR Mix at room temperature for more than 30 minutes. Use a new tube of Q-PCR Mix.
Internal Control template degradation.	Use a new tube of Internal Control.
Inhibition due to interfering substances in the sample.	Repeat the amplification of eluted sample as is or with a 1:2 dilution in molecular biology grade water in a "PCR Only" session. Repeat the extraction of the sample with a 1:2 dilution in molecular biology grade water in an "Extract + PCR" session.
Instrument error.	Contact ELITechGroup Technical Service.

Table 54

Anomalous dissociation curve	
Possible causes	Solutions
Absence of a defined peak. Defined peak but T _m different from that of the other samples and that of the Standards or Positive Control.	Check for target Ct lower than 30. High quantity of amplification product at the end of the reaction may interfere with the melting curve analysis. Repeat the sample amplification to confirm the presence of target with a possible mutation. The target in the sample should be sequenced to confirm mutation.

Table 55

Error in Ct calculation	
Possible Causes	Solutions
Too high concentration of target in the sample or sample with anomalous fluorescence signal.	<p>If significant amplification is observed in PCR plot, select the track related to the sample and manually approve the result as positive.</p> <p>If no amplification is observed in PCR plot select the track related to the sample and manually approve the result as negative or leave it as invalid.</p> <p>If a Ct value is required:</p> <ul style="list-style-type: none"> - repeat the amplification of eluted sample with a 1:10 dilution in molecular biology grade water in a "PCR Only" session - repeat the extraction of the sample with a 1:10 dilution in molecular biology grade water in an "Extract + PCR" session.

Table 56

Abnormal high rate of positive results within the same session (reactions with similar late Ct values)	
Possible Causes	Solutions
Sample-to-sample contamination during preanalytical steps.	<p>Clean the micropipette with fresh 3% sodium hypochlorite solution (bleach) or DNA/RNA cleaner after pipetting each sample.</p> <p>Do not use Pasteur pipettes. The pipettes must be of the positive displacement type or used with aerosol filter tips.</p> <p>Introduce samples in the last positions of the instruments, as indicated by the GUI. Follow the loading sequence indicated by the software.</p>
Laboratory environmental contamination.	<p>Clean all surfaces in contact with the operator and samples (including the pipettes) with fresh 3% sodium hypochlorite solution (bleach) or DNA/RNA cleaner.</p> <p>Perform an U.V. decontamination cycle.</p> <p>Use a new tube of Q-PCR Mix and / or Internal Control</p>

MyGenius PRO**Table 57**

Invalid Q-PCR Standard reaction, Standard curve or Positive Control reaction	
Possible Causes	Solutions
Instrument setting error.	<p>Check the volumes of Q-PCR Mix, Q-PCR Standards and Positive Control.</p> <p>It is important that the tubes used on MyGenius PRO have not been previously used by other platforms or for setting sessions on open platforms because the instrument reads the label and automatically assigns the maximum number of replicates that can be run if that tube has never been loaded on MyGenius PRO. Therefore, if a non-new tube is used, the volume may not be sufficient.</p>
PCR Mix degradation.	<p>Do not use the Q-PCR Mix for more than 7 hours in the Carousel or for more than 3 hours in the Carousel for five times.</p> <p>Do not leave the Q-PCR Mix at room temperature for more than 30 minutes.</p> <p>Use a new tube of Q-PCR Mix.</p>

Table 57 (continued)

Invalid Q-PCR Standard reaction, Standard curve or Positive Control reaction	
Possible Causes	Solutions
Q-PCR Standards or Positive Control degradation.	Do not use the Q-PCR Standard for more than 4 independent sessions (2 hours each in the Carousel). Do not use the Positive Control for more than 4 independent sessions (3 hours each in the Carousel). Use new tubes of Q-PCR Standards or Positive Control.
Instrument error.	Contact ELITechGroup Technical Service.

Table 58

Invalid Negative Control reaction	
Possible Causes	Solutions
Instrument setting error.	Check the volumes of Q-PCR Mix and Negative Control. It is important that the tubes used on MyGenius PRO have not been previously used by other platforms or for setting sessions on open platforms because the instrument reads the label and automatically assigns the maximum number of replicates that can be run if that tube has never been loaded on MyGenius PRO. Therefore, if a non-new tube is used, the volume may not be sufficient.
Contamination of the Negative Control.	Do not use the Negative Control for more than 1 session and for more than 8 hours. Use a new tube of Negative Control.
Contamination of the PCR Mix.	Use a new tube of Q-PCR Mix.
Contamination of the extraction area, Racks, Inventory Block or Cooler Unit.	Clean surfaces with aqueous detergents, wash lab coats, replace tubes and tips in use.
Instrument error.	Contact ELITechGroup Technical Service.

Table 59

Invalid Sample reaction	
Possible Causes	Solutions
Instrument setting error.	Check the volumes of Q-PCR Mix, Internal Control and Samples. It is important that the tubes used on MyGenius PRO have not been previously used by other platforms or for setting sessions on open platforms because the instrument reads the label and automatically assigns the maximum number of replicates that can be run if that tube has never been loaded on MyGenius PRO. Therefore, if a non-new tube is used, the volume may not be sufficient.
PCR Mix degradation.	Do not use the Q-PCR Mix for more than 7 hours in the Carousel or for more than 3 hours in the Carousel for five times. Do not leave the Q-PCR Mix at room temperature for more than 30 minutes. Use a new tube of Q-PCR Mix.
Internal Control template degradation.	Do not use the Internal Control for more than 8 hours. Use a new tube of Internal Control.
Inhibition due to interfering substances in the sample.	Repeat the extraction and amplification of the sample. Repeat the extraction and amplification of the sample with a 1:2 dilution in molecular biology grade water.

Table 59 (continued)

Invalid Sample reaction	
Possible Causes	Solutions
Too high concentration of target in the sample or sample with anomalous fluorescence signal.	If a Ct value is required, repeat the extraction and amplification of the sample with a 1:10 dilution in molecular biology water.
Instrument error.	Contact ELITechGroup Technical Service.

Table 60

Anomalous dissociation curve	
Possible causes	Solutions
Absence of a defined peak. Defined peak but T _m different from that of the other samples and that of the Standards or Positive Control.	<p>Check for target Ct lower than 30.</p> <p>High quantity of amplification product at the end of the reaction may interfere with the melting curve analysis.</p> <p>Repeat the sample extraction and amplification to confirm the presence of target with a possible mutation.</p> <p>The target in the sample should be sequenced to confirm mutation</p>

Table 61

Error in Ct calculation	
Possible causes	Solutions
Too high concentration of target in the sample or sample with anomalous fluorescence signal.	If a significant amplification appears in the PCR plot, the result should be considered positive. If no amplification appears in the PCR plot, the result should be considered negative or invalid. If a Ct value is required, repeat the extraction and amplification of the sample with a 1:10 dilution in molecular biology water.

Table 62

Abnormal high rate of positive results within the same session (reactions with similar late Ct values)	
Possible Causes	Solutions
Sample-to-sample contamination during preanalytical steps.	<p>Clean the micropipette with fresh 3% sodium hypochlorite solution (bleach) or DNA/RNA cleaner after pipetting each sample.</p> <p>Do not use Pasteur pipettes. The pipettes must be of the positive displacement type or used with aerosol filter tips.</p> <p>Introduce samples as the last of the day.</p>
Laboratory environmental contamination.	<p>Clean all surfaces in contact with the operator and samples (including the pipettes) with fresh 3% sodium hypochlorite solution (bleach) or DNA/RNA cleaner.</p> <p>Perform an U.V. decontamination cycle.</p> <p>Use a new tube of Q-PCR Mix and / or Internal Control</p>

Open Platform

Table 63

Invalid Q-PCR Standard reaction, Standard curve or Positive Control reaction	
Possible Causes	Solutions
Incorrect dispensing into the microplate wells.	Check the volumes of PCR Mix, Q-PCR Standards and Positive Control dispensed in the Q-PCR microplate.
Q-PCR Mix degradation.	Do not freeze and thaw the PCR mix more than 5 times. Do not leave the Q-PCR Mix at room temperature for more than 30 minutes. Use a new tube of Q-PCR Mix.
Q-PCR Standards or Positive Control degradation.	Do not freeze and thaw the Q-PCR standard more than 4 times. Use new tubes of Q-PCR Standards or Positive Control.
Instrument setting error.	Check the position of PCR Mix, Q-PCR Standards and Positive Control on the instrument. Check the thermal cycle settings on the instrument.

Table 64

Invalid Negative Control reaction	
Possible Causes	Solutions
Instrument setting error.	Check the position of Q-PCR Mix and Negative Control. Check the volumes of Q-PCR Mix and Negative Control.
Microplate badly sealed.	Take care when sealing the Q-PCR microplate with the optical seal.
Contamination of the Negative Control.	Do not use the Negative Control for more than 1 session. Use a new aliquot of molecular biology grade water.
Contamination of the PCR Mix.	Use a new tube of Q-PCR Mix.
Contamination of the preparation area, racks and micropipette.	Clean surfaces and instruments with aqueous detergents, wash lab coats, replace test tubes and tips in use.

Table 65

Invalid Sample reaction	
Possible Causes	Solutions
Instrument setting error.	Check the position of Q-PCR Mix, Internal Control and sample. Check the volumes of Q-PCR Mix, Internal Control and sample.
PCR Mix degradation.	Do not freeze and thaw the PCR mix more than five times. Do not leave the Q-PCR Mix at room temperature for more than 30 minutes. Use a new tube of Q-PCR Mix
Internal Control template degradation.	Use a new tube of Internal Control.
Inhibition due to interfering substances in the sample.	Repeat the amplification of eluted sample with a 1:2 dilution in molecular biology grade water. Repeat the extraction of the sample with a 1:2 dilution in molecular biology grade water.

Table 66

Irregular or high background fluorescence in the reactions	
Possible causes	Solutions
Incorrect dispensing of sample.	Check the volumes of reagents and samples dispensed in the Q-PCR microplate.
Baseline setting error.	If the calculation range for the Baseline set from cycle 6 to cycle 15 is not proper to normalize the background, set the calculation range within cycles where the background fluorescence has already stabilized (check Results > Component) and the target fluorescence has not yet started to increase.

Table 67

Anomalous dissociation curve	
Possible causes	Solutions
Absence of a defined peak. Defined peak but different from that of the other samples and that of the Standards or Positive Control.	Check for target Ct lower than 30. High quantity of amplification product at the end of the reaction may interfere with the melting curve analysis. Repeat the sample amplification to confirm the presence of target with a possible mutation. The target in the sample should be sequenced to confirm mutation.

19 SYMBOLS



Catalogue Number.



Upper limit of temperature.



Batch code.



Use by (last day of month).



in vitro diagnostic medical device.



Fulfilling the requirements of the IVDR Regulation 2017/746/EC for *in vitro* diagnostic medical device. Certification released by TÜV SÜD Product Service GmbH, Germany.



Unique Device Identification



Contains sufficient for "N" tests.



Consult instructions for use.



Contents.



Keep away from sunlight.



Manufacturer.

20 NOTICE TO THE USERS

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and /or the patient is established. To inform ELITechGroup S. p. A., manufacturer of this device, please use the following mail address: egspa.vigilance@elitechgroup.com.

A “Summary of Safety and Performance” will be made available to the public via the European database on medical devices (Eudamed) when this informatic system will be functional. Before the notice of full functionality of Eudamed has been published, the “Summary of Safety and Performance” will be made available to the public upon request by email at emd.support@elitechgroup.com, without undue delay.

21 NOTICE TO PURCHASER: LIMITED LICENSE

This product contains reagents manufactured by Promega Corporation and supplied to ELITechGroup S.p.A. for use exclusively as components of ELITech branded diagnostic kits. The reagents may be used only within this kit. No additional rights or licenses—express or implied—are granted for any use, modification, or resale of the reagents outside this kit.

ELITe MGB® detection reagents and ELITech platforms (ELITe InGenius®, ELITe BeGenius®, MyGenius PRO®) are protected by granted patents and pending applications.

Any use of the reagents or related data outside the scope of this kit requires prior written authorization from ELITechGroup S.p.A.

ELITe MGB® detection reagents are covered by one or more of U. S. Patent numbers 7319022, 7348146, 7541454, 7671218, 7723038, 7767834, 8163910, 8969003, 9056887, 9085800, 9169256, 9328384, 10677728, 10738346, 10890529, and EP patent numbers 2689031, 2714939, 2736916, 2997161 as well as applications that are currently pending.

ELITe InGenius®, ELITe BeGenius® and MyGenius PRO® technologies are covered by patents and pending applications.

This limited license allows the person or entity to whom the product has been provided to use the product and data generated by the use of the product, solely for human diagnostics. Neither ELITechGroup S.p.A. nor its licensors grant any other licenses, expressed or implied for any other purposes.

MGB®, Eclipse Dark Quencher®, AquaPhluor®, ELITe MGB®, the ELITe MGB® logo, ELITe InGenius®, ELITe BeGenius® and the MyGenius PRO® (registration name ELIVERSE®) are registered trademarks of ELITechGroup within the European Union.
QIASymphony® is registered trademark of QIAGEN GmbH.
Ficoll® is registered trademark of GE Healthcare Bio-Sciences AB.

Appendix A EBV ELITE MGB Kit used in association with Genius series® platforms



CAUTION

This document is a simplified version of the official instruction for use. Please refer to the complete document before use: www.elitechgroup.com

Intended use

The product **EBV ELITE MGB® Kit** is an *in vitro* diagnostic medical device intended to be used by healthcare professionals as quantitative nucleic acids Real-Time PCR assay for the detection and quantification of the **DNA of Epstein-Barr virus (EBV)** extracted from clinical specimens.

The assay is validated in association with the **ELITE InGenius®** and **ELITE BeGenius®** instruments, automated and integrated systems for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid.

The assay is also validated in association with **MyGenius PRO®** (registration name ELIVERSE®) instrument, automated and integrated system for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA.

The assay is also validated in association with the **ELITE GALAXY**, automatic extraction and PCR set-up system and **7500 Fast Dx Real-Time PCR Instrument**, Real-Time PCR platform, using human specimens of whole blood collected in EDTA.

The product is intended for use as an aid in the diagnosis and monitoring of EBV infections in patients suspected of having or undergoing monitoring of EBV infections.

The results must be interpreted in combination with all relevant clinical observation and laboratory outcomes.




Amplified sequence

Sequence	Gene	Fluorophore	Channel
Target	EBNA-1 gene	FAM	EBV
Internal Control	promoter and 5' UTR region of the human beta-globin gene	AP525	IC

Validated matrix

- Whole blood collected in EDTA
- Plasma collected in EDTA
- CSF

Kit content and related products

EBV ELITe MGB Kit	EBV ELITe Standard	EBV- ELITe Positive Control
 X 4	 X 2	 X 2
Ready-to-use PCR Mix 4 tubes of 540 µL 96 reactions per kit 5 freeze-thaw cycles per tube	Ready-to-use 4 levels: 10 ⁵ , 10 ⁴ , 10 ³ , 10 ² 2 sets of 4 tubes of 200 µL 8 reactions per kit 4 freeze-thaw cycles per tube	Ready-to-use PC 2 tubes of 160 µL 8 reactions per kit 4 freeze-thaw cycles per tube

Maximum shelf-life: **24 months**Storage Temperature: **-20 °C**

Other products required not provided in the kit

<ul style="list-style-type: none"> • ELITe InGenius instrument: INT030. • ELITe BeGenius instrument: INT040. • ELITe InGenius SP 200: INT032SP200. • ELITe InGenius SP1000: INT033SP1000 	<ul style="list-style-type: none"> • CPE - Internal Control: CTCPE • ELITe InGenius and ELITe BeGenius Consumables (see ELITe InGenius and ELITe BeGenius Instruction for use)
--	--

ELITe InGenius and ELITe BeGenius protocol

<ul style="list-style-type: none"> › Sample volume › CPE volume › Total elution volume 	200 µL (InGenius and BeGenius) or 1000 µL (InGenius only) 10 µL 100 µL	<ul style="list-style-type: none"> › Eluate PCR input volume: › Q—PCR Mix volume: › Frequency of controls: › Frequency of calibration: 	20 µL 20 µL 15 days 30 days
---	--	--	--------------------------------------

ELITe InGenius and ELITe BeGenius Performances

Matrix	Limit of Detection		Diagnostic Sensitivity	Diagnostic Specificity	Linearity (IU/mL)		Conversion factor (IU/copies)
	IU/mL	copies/mL			IU / mL	copies/mL	
whole blood	104	36	98.5 %	97.3 %	104 - 1,000,001	36 — 344,828	2.9
Plasma (sample volume 200 µL)	124	65	98.9 %	98.3 %	124 — 10,000,000	65 — 5,263,158	1.9
Plasma (sample volume 1000 µL, InGenius only)	18	11			99 – 1,000,000	62 – 625,000	1.6
CSF	176	80	98.8 %	100 %	176 — 1,000,000	80 — 454,545	2.2

Sample preparation

This product is intended for use on the **ELITE InGenius** and **ELITE BeGenius** with the following clinical specimens identified according to laboratory guidelines, and collected, transported, and stored under the following conditions.

Sample type	Collection requirements	Transport/Storage conditions			
		+16 / +26 °C (room temperature)	+2 / +8 °C	-20 ± 10 °C	-70 ± 15 °C
Whole blood	EDTA	≤ 1 d	≤ 3 d	≤ 30 d	≤ 30 d
Plasma	EDTA	≤ 1 d	≤ 3 d	≤ 30 d	≤ 30 d
CSF	Minimal patient blood contamination	≤ 1 d	≤ 2 d	≤ 30 d	≤ 30 d

C EDTA, Ethylenediaminetetraacetic acid; CSF, cerebrospinal fluid, d, day.

ELITE InGenius Procedures

The user is guided step-by-step by the Graphic User Interface (GUI) of ELITE InGenius software to setup the run. All the steps: extraction, Real-Time PCR and result interpretation are automatically performed. Two operational modes are available: complete run (Extract + PCR) or PCR Only.

Before analysis

1. Switch on ELITE InGenius. Log in with username and password. Select the mode "Closed".	2. Verify calibrators: Q-PCR Standard in the "Calibration" menu. Verify controls: Positive Control and Negative Control in the "Controls" menu. Note: All must have been run, approved and not expired.	3. Thaw the PCR Mix and the CTRCPE tubes. Vortex gently. Spin down 5 sec.
---	--	---

Procedure 1 - Complete run: Extract + PCR (e.g., samples)

1. Select "Perform Run" on the touch screen	2. Verify the extraction volumes: Input: "200 µL" or "1000 µL, elution: "100 µL"	3. Scan the sample barcodes with hand-barcode reader or type the sample ID
4. Select the "Assay Protocol" of interest: EBV ELITE_WB_200_100 or EBV ELITE_PL_200_100 or EBV ELITE_PL_1000_100 or EBV ELITE_CSF_200_100	5. Select the method "Extract + PCR" and the sample position: Primary tube or Extraction Tube	6. Load the PCR Mix and the Internal Control in the Inventory Block
7. Load: PCR Cassette, Extraction cartridge, Elution tube, Tip Cassette, Extraction Tube racks and primary sample racks	8. Close the door. Start the run	9. View, approve and store the results

NOTE

If an Extract Only mode is needed, refer to the instrument user's manual for procedure.

Procedure 2: PCR Only (e.g., eluates, standards, controls)

1. Select "Perform Run" on the touch screen	2. Verify the extraction volumes: Input: "200 µL" or "1000 µL", elution: "100 µL"	3. Scan the sample barcodes with hand-barcode reader or type the sample ID
4. Select the "Assay protocol" of interest: EBV ELITE_PC and EBV ELITE_NC, or EBV ELITE_STD, or EBV ELITE_PL_200_100 or EBV ELITE_PL_1000_100, or EBV ELITE_WB_200_100, EBV ELITE_CSF_200_100.	5. Select the method "PCR Only" and the sample position "Elution Tube"	6. Load the PCR Mix in the Inventory Block
7. Load: PCR Cassette rack and Elution tube rack with the extracted nucleic acid	8. Close the door. Start the run	9. View, approve and store the results

ELITE BeGenius Procedures

The user is guided step-by-step by the Graphic User Interface (GUI) of ELITE BeGenius software to setup the run. All the steps, extraction, Real-Time PCR and result interpretation, are automatically performed. Two operational modes are available: complete run (Extract + PCR) or PCR Only.

Before analysis

1. Switch on ELITE BeGenius. Log in with username and password. Select the mode "Closed".	2. Verify calibrators: Q-PCR Standard in the "Calibration" menu. Verify controls: Positive Control and Negative Control in the "Controls" menu. Note: All must have been run, approved and not expired.	3. Thaw the PCR Mix and the CTRCPE tubes. Vortex gently. Spin down 5 sec.
---	--	---

Procedure 1 - Complete run: Extract + PCR (e.g., samples)

1. Select "Perform Run" on the touch screen and then click on the run mode «Extract + PCR»	2. Insert the Sample Rack with the barcoded samples in the Cooler Unit. The barcode scan is already active	3. Verify the extraction volumes: Input: "200 µL", Eluate: "100 µL"
4. Select the "Assay protocol" of interest EBV ELITE_Be_WB_200_100 or EBV ELITE_Be_PL_200_100 or EBV ELITE_Be_CSF_200_100 Note: If a second extraction is performed repeat steps from 2 to 4	5. Print the labels to barcode the empty elution tubes. Load the tubes in the Elution Rack and insert it in the Cooler Unit	6. Load the PCR Mix and the Internal Control in the Reagent/Elution Rack and insert it in the Cooler Unit
7. Load "PCR Rack" with "PCR Cassette" and the "Extraction Rack" with the "ELITE InGenius SP 200" extraction cartridges and the required extraction consumables	8. Close the door. Start the run	9. View, approve and store the results

NOTE

If an Extract Only mode is needed, refer to the instrument user's manual for procedure.

Procedure 2: PCR Only (e.g., eluates, standards, controls)

1. Select “Perform Run” on the touch screen	2. Verify the extraction volumes: Input: “200 µL”, elution: “100 µL”	3. Scan the sample barcodes with hand-barcode reader or type the sample ID
4. Select the “Assay protocol” of interest: EBV ELITe_Be_PC and EBV ELITe_Be_NC, or EBV ELITe_Be_STD.	5. Select the method “PCR Only” and the sample position “Elution Tube”	6. Load the PCR Mix in the Inventory Block
7. Load: PCR Cassette rack and the Elution tube rack with the extracted nucleic acid	8. Close the door. Start the run	9. View, approve and store the results

Appendix B EBV ELITE MGB Kit used in association with MyGenius PRO



CAUTION

This document is a simplified version of the official instruction for use. Please refer to the complete document before use: www.elitechgroup.com

Intended use

The product **EBV ELITE MGB® Kit** is an *in vitro* diagnostic medical device intended to be used by healthcare professionals as quantitative nucleic acids Real-Time PCR assay for the detection and quantification of the **DNA of Epstein-Barr virus (EBV)** extracted from clinical specimens.

The assay is validated in association with the **ELITE InGenius®** and **ELITE BeGenius®** instruments, automated and integrated systems for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid.

The assay is also validated in association with **MyGenius PRO®** (registration name ELIVERSE®) instrument, automated and integrated system for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA.

The assay is also validated in association with the **ELITE GALAXY**, automatic extraction and PCR set-up system and **7500 Fast Dx Real-Time PCR Instrument**, Real-Time PCR platform, using human specimens of whole blood collected in EDTA.

The product is intended for use as an aid in the diagnosis and monitoring of EBV infections in patients suspected of having or undergoing monitoring of EBV infections.

The results must be interpreted in combination with all relevant clinical observation and laboratory outcomes.




Amplified sequence

Sequence	Gene	Fluorophore	Channel
Target	EBNA-1 gene	FAM	EBV
Internal Control	promoter and 5' UTR region of the human beta-globin gene	AP525	IC

Validated matrix

- Whole blood collected in EDTA

Kit content and related products

EBV ELITe MGB Kit	EBV ELITe Standard	EBV- ELITe Positive Control
 X 4	 X 2	 X 2
Ready-to-use PCR Mix 4 tubes of 540 µL 96 reactions per kit 5 freeze-thaw cycles per tube	Ready-to-use 4 levels: 10 ⁵ , 10 ⁴ , 10 ³ , 10 ² 2 sets of 4 tubes of 200 µL 8 reactions per kit 4 freeze-thaw cycles per tube	Ready-to-use PC 2 tubes of 160 µL 8 reactions per kit 4 freeze-thaw cycles per tube

Maximum shelf-life: **24 months**Storage Temperature: **-20 °C**

Other products required not provided in the kit

<ul style="list-style-type: none"> • MyGenius PRO (EG SpA ref: INT050) • MyGenius PRO Software version BB-04 (or later) • ELITe InGenius SP 200 (EG SpA, ref. INT032SP200) 	<ul style="list-style-type: none"> • Negative Control (EG SpA, ref. CTRNEG) • Internal Control Maxi (EG SpA, ref. ICMAXI) • MyGenius PRO Consumables (see MyGenius PRO Instruction for use)
--	---

MyGenius PRO protocol

› Sample volume	200 µL	› Eluate PCR input volume	20 µL
› IC volume	10 µL	› Q—PCR Mix volume	20 µL
› Total elution volume	100 µL	› Frequency of controls	15 days
		› Frequency of calibration	60 days

MyGenius PRO Performances

Matrix	Limit of Detection		Diagnostic Sensitivity	Diagnostic Specificity	Linearity (IU/mL)		Conversion factor (IU/copies)
	IU/mL	copies/mL			IU / mL	copies/mL	
whole blood	104	36	98.5 %	97.3 %	104 - 1,000,001	36 — 344,828	2.9

Sample preparation

This product is intended for use on the **MyGenius PRO** with the following clinical specimens identified according to laboratory guidelines, and collected, transported, and stored under the following conditions.

Sample type	Collection requirements	Transport/Storage conditions			
		+16 / +26 °C (room temperature)	+2 / +8 °C	-20 ± 10 °C	-70 ± 15 °C
Whole blood	EDTA	≤ 1 d	≤ 3 d	≤ 30 d	≤ 30 d

C EDTA, Ethylenediaminetetraacetic acid; d, day.

MyGenius PRO Procedures

The user is guided step-by-step by the Graphic User Interface (GUI) of MyGenius PRO software to setup the run. All the steps: extraction, Real-Time PCR and result interpretation are automatically performed. Two operational modes are available: complete run (Extract + PCR) or PCR Only interpretation are automatically performed. Two operational modes are available: complete run (Extract + PCR) or PCR Only (only for Calibrators and Controls).

<p>1. Switch on MyGenius PRO. Log in in STAND-BY mode with username and password.</p>	<p>2. Load all consumables in the drawers and empty liquid waste tank and solid waste boxes if needed. Press Start button to perform Preparation: after preparation, the instrument will go in Operation mode.</p>	<p>3. Verify calibrators: Q-PCR Standard in the "Calibration" menu. Verify controls: Positive Control and Negative Control in the "Controls" menu Note: All must have been run, approved and not expired.</p>
<p>4. Thaw the PCR Mix and the IC MAXI tubes. Vortex gently. Spin down 5 sec</p>	<p>5. Load PCR mixes tubes in the reagent carousel by following GUI indications.</p>	<p>6. Load IC MAXI tube(s) in dedicated blue rack and load it in the Auto Sampler area.</p>
<p>7a. If the instrument is connected to the LIS, insert samples in Auto Sampler area using dedicated racks depending on tube diameters used. The extraction automatically will start. 7b. If the instrument is not connected to the LIS, in Sample list, press "Assign test", read the barcode of the samples with the external barcode reader, select Whole blood matrix and assign the Assay protocol "EBV ELITe_My_WB_IU_200_100" or "EBV ELITe_My_WB_cmL_200_100" insert samples in Auto Sampler area using dedicated racks depending on tube diameters used. The extraction automatically will start.</p>		

Appendix C EBV ELITe MGB Kit used in association with ABI 7500 Instrument



CAUTION

This document is a simplified version of the official instruction for use. Please refer to the complete document before use: www.elitechgroup.com

Intended use

The product **EBV ELITe MGB® Kit** is an *in vitro* diagnostic medical device intended to be used by healthcare professionals as quantitative nucleic acids Real-Time PCR assay for the detection and quantification of the **DNA of Epstein-Barr virus (EBV)** extracted from clinical specimens.

The assay is validated in association with the **ELITe InGenius®** and **ELITe BeGenius®** instruments, automated and integrated systems for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid.

The assay is also validated in association with **MyGenius PRO®** (registration name ELIVERSE®) instrument, automated and integrated system for extraction, Real-Time PCR and results interpretation, using human specimens of whole blood collected in EDTA.

The assay is also validated in association with the **ELITe GALAXY**, automatic extraction and PCR set-up system and **7500 Fast Dx Real-Time PCR Instrument**, Real-Time PCR platform, using human specimens of whole blood collected in EDTA.

The product is intended for use as an aid in the diagnosis and monitoring of EBV infections in patients suspected of having or undergoing monitoring of EBV infections.

The results must be interpreted in combination with all relevant clinical observation and laboratory outcomes.




Amplified sequence

Sequence	Gene	Fluorophore	Channel
Target	EBNA-1 gene	FAM	EBV
Internal Control	promoter and 5' UTR region of the human beta-globin gene	AP525	IC

Validated matrix

- Whole blood collected in EDTA

Kit content and related product

EBV ELITe MGB Kit	EBV ELITe Standard	EBV- ELITe Positive Control
 X 4	 X 2	 X 2
Ready-to-use PCR Mix 4 tubes of 540 µL 100 reactions per kit 5 freeze-thaw cycles per tube	Ready-to-use 4 levels: 10^5 , 10^4 , 10^3 , 10^2 2 sets of 4 tubes of 200 µL 8 freeze-thaw cycles per tube	Ready-to-use PC 2 tubes of 160 µL 12 reactions per kit (Galaxy) 8 freeze-thaw cycles per tube

Maximum shelf-life: **24 months**

Storage Temperature: **-20 °C**

Other products required not provided in the kit

<ul style="list-style-type: none"> • ELITe GALAXY: INT020 • ELITe GALAXY 300 extraction kit: INT021EX • ABI 7500 Fast Dx Real—Time PCR Instrument 	<ul style="list-style-type: none"> • CPE - Internal Control: CTRCPE • Molecular biology grade water
---	---

7500 Real-Time PCR Instrument Performances

Nucleic acid extraction method	Limit of Detection	Diagnostic Specificity	Diagnostic Sensitivity	Linearity (IU/mL)	Formula to quantity (gEq / mL)	Conversion factor gEq/mL to IU/mL
ELITe GALAXY	127 IU / mL 249 gEq / mL	100% (65/65)	100% (60/60)	178→ 17.8 x 10 ⁶	35 x quantity	0.51
QIASymphony SP/AS	—	98.3% (59/60)	100% (60/60)	-	24 x quantity	0.46

7500 Real-Time PCR Instrument Procedures

The procedure below summarized the main steps of the sample analysis with conventional PCR workflow: validated extraction systems, PCR instrument settings, PCR set-up and result interpretation.

Extraction - Validated systems

Extraction	Validated matrix	Sample volume processed	Min. sample volume	Total eluate volume	CPE Internal Control volume
ELITe Galaxy	WB	300 µL	400 µL	200 µL	10 µL

Amplification - Settings of 7500 Fast Dx

1. Switch on the thermal-cycler
2. Set "EBV" detector with "FAM" and quencher "none"
3. Set "Internal Control" detector with "VIC" and quencher "none"
4. Set passive fluorescence as "Cy5"
5. Set up the thermal profile as indicated. Fluorescence acquisition must be set during the hybridation step at 60°C

Stage	Temperature	Timing
Decontamination	50°C	2 min
Initial Denaturation	94°C	2 min
Amplification Detection 45 cycles	94°C	10 sec
	60°C	30 sec
	72°C	20 sec

The melt curve analysis is optional, refer to the complete IFU

Amplification - PCR Set-up (performed by ELITe Galaxy and QIASimphony)

To perform the PCR session set up:

1. thaw Q PCR-Mix and Positive Control / Q-PCR standard tubes
2. mix gently and spin-down
3. prepare the **Negative Control** (not provided)
4. prepare a **Q-PCR microplate**
5. The instrument automatically performs the PCR set-up dispensing in each well of the **Q-PCR microplate 20 µL of PCR Mix and 20 µL of extracted DNA / Q-PCR Standard / Controls**

After the PCR set-up performed by the instrument:

1. seal the **Q-PCR microplate** with an optical seal
2. transfer the **Q-PCR microplate** onto the **7500 Fast Dx Real-Time PCR Instrument** and start the PCR. Save the run file with a unique and recognizable name (e.g. "year-month-day-TARGET-EGSpA")

Amplification - Threshold for qualitative analysis

Instrument	EBV FAM	Internal Control VIC
7500 Fast Dx Real Time PCR	0.2	0.1

Interpretation

Qualitative results		
EBV Ct value	Internal Control Ct value	Interpretation
Determined	—	Positive
Undetermined	Ct ≤ 35	Negative
	Ct >35 or Undetermined	Invalid*

**Repeat the assay starting from the extraction*

Quantitative results
The EBV Ct value obtained for each sample and the standard curve generated are used to calculate the quantity of target DNA in the reaction
The sample quantification ranges from approximately 10 to 10 ⁶ copies/reaction.

ELITechGroup S.p.A.
C.so Svizzera, 185, 10149 Torino ITALY
Tel. +39-011 976 191
Fax +39-011 936 76 11
E. mail: emd.support@elitechgroup.com
WEB site: www.elitechgroup.com

