

## NOTICE of CHANGE dated 28/10/2022

### IMPORTANT COMMUNICATION FOR THE USERS OF PRODUCT:







# «HSV2 ELITe MGB® Kit» Ref. RTS032PLD

This new revision of the Instruction for Use (IFU) contains the following changes:

- Update for the use of the product for CSF matrix in association with «ELITe BeGenius®» instrument (REF INT040);
- Internal Control Ct cut-off update in association with «ELITe InGenius®» (REF INT030) and «ELITe BeGenius®» (REF INT040)

Composition, use and performance of the product remain unchanged.

## PLEASE NOTE


	LA REVISIONE DI QUESTO IFU E' COMPATIBILE ANCHE CON LA VERSIONE PRECEDENTE DEL KIT
	THE REVIEW OF THIS IFU IS ALSO COMPATIBLE WITH THE PREVIOUS VERSION OF THE KIT
	CET IFU MIS A JOUR ANNULE ET REMPLACE ET EST PARFAITEMENT COMPATIBLE AVEC LA VERSION PRECEDENTE DU KIT
	LA REVISIÓN DE ESTE IFU ES COMPATIBLE TAMBIÉN CON LA VERSIÓN ANTERIOR DEL KIT
	A REVISÃO DO ESTE IFU ÉTAMBÉM COMPATÍVEL COM A VERSÃO ANTERIOR DO KIT
	DIE REVIEW VON DIESER IFU IST KOMPATIBLE MIT DER VORIGE VERSION VON DEM TEST-KIT



## HSV2 ELITe MGB® Kit

reagent for DNA Real Time amplification

REF RTS032PLD



### TABLE OF CONTENTS

INTENDED USE	page 2
ASSAY PRINCIPLES	page 2
PRODUCT DESCRIPTION	page 3
MATERIALS PROVIDED IN THE PRODUCT	page 3
MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT	page 3
OTHER PRODUCTS REQUIRED	page 3
WARNINGS AND PRECAUTIONS	page 5
ELITe InGenius	page 6
SAMPLES AND CONTROLS	page 6
PROCEDURE	page 7
ELITe BeGenius	page 14
SAMPLES AND CONTROLS	page 14
PROCEDURE	page 15
PERFORMANCE CHARACTERISTICS ELITe InGenius and ELITe BeGenius	page 20
ABI 7500 Fast Dx Real-Time PCR Instrument ABI 7300 Real-Time System	page 29
SAMPLES AND CONTROLS	page 29
PROCEDURE	page 31
PERFORMANCE CHARACTERISTICS	page 39
Roche cobas z 480 analyzer	page 43
SAMPLES AND CONTROLS	page 43
PROCEDURE	page 44
PERFORMANCE CHARACTERISTICS	page 49
REFERENCES	page 51
PROCEDURE LIMITATIONS	page 51
TROUBLESHOOTING	page 52
SYMBOLS	page 54
NOTICE TO PURCHASER: LIMITED LICENSE	page 55

HSV2 ELITe MGB® Kit  
reagent for DNA Real Time amplification

REF RTS032PLD

### INTENDED USE

The «HSV2 ELITe MGB® Kit» product is a part of qualitative and quantitative nucleic acids amplification assay for the **detection and quantification of the DNA of type 2 Herpes Simplex human virus (HSV2)** in DNA samples extracted from cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

The product is intended for use in the diagnosis and monitoring of HSV2 infections, alongside clinical data of the patient and other laboratory tests outcomes.

### ASSAY PRINCIPLES


The assay consists of a real time amplification reaction with a programmable thermostat provided with a fluorescence detection optical system system (real time amplification thermal cycler).

In each well, two amplification reactions are performed starting from DNA extracted from the samples being tested: a specific reaction for a region of the **glycoprotein G (gpG)** of HSV2 and a specific reaction for a region of the **human beta Globin gene** (Internal Control of inhibition). The HSV2 specific probe with ELITe MGB® technology, labelled with FAM fluorophore, is activated when hybridizes with the specific product of the HSV2 amplification reaction. The Internal Control specific probe with ELITe MGB® technology, labelled with AP525 fluorophore (analogous to VIC), is activated when hybridizes with the specific product of the amplification reaction for the Internal Control. As the specific product of the amplification reaction increases, the fluorescence emission increases and is measured and recorded by the instrument. The processing of the data allows detecting the presence and the titre of HSV2 DNA in the starting sample.

At the end of the amplification session, dissociation curve (melting curve) analysis can be carried out in order to determine the dissociation temperature (melting temperature) and to confirm the presence of the correct target or to identify the presence of mutations.

The assay is validated with the systems described in this instruction for use.

In the following picture is synthetically showed the mechanism of activation and fluorescence emission of ELITe MGB® technology probe. Note that the probe is not hydrolyzed during the amplification cycle so as it can be utilized for the dissociation curve analysis.



## PRODUCT DESCRIPTION

The «**HSV2 ELITe MGB Kit**» product supplies the **ready to use** complete mixture "**HSV2 Q - PCR Mix**" for real time amplification in a stabilising solution, **aliquoted into four disposable test tubes**. Each tube contains **540 µL** of solution, sufficient for **24 tests** in association with «**ELITe InGenius®**» and «**ELITe BeGenius®**» systems and **25 tests** in association with other systems.

The primers and the HSV2 specific probe (stabilized by MGB® group, labelled by FAM fluorophore and quenched by a non fluorescent molecule) are specific for a region of the **gpG** of HSV2 (region US4).

The primers and the probe for the Internal Control (stabilized with MGB® group, labelled by AP525, analogous to VIC, fluorophore and quenched by a non fluorescent molecule) are specific for the **promoter and 5' UTR** region of the **human beta Globin gene**.

The reaction mixture provides buffer, magnesium chloride, triphosphate nucleotides, AP593 fluorophore (used instead of ROX or CY5) as passive reference for fluorescence normalisation, the enzyme Uracil N-glycosidase (UNG) to inactivate contamination by the amplification product, the "hot start" DNA polymerase enzyme.

The product is sufficient for **96 tests in association with «ELITe InGenius» and «ELITe BeGenius»** systems, including standards and controls.

The product is sufficient for **100 tests in association with other systems**, including standards and controls.

## MATERIALS PROVIDED IN THE PRODUCT

Component	Description	Quantity	Classification of hazards
HSV2 Q - PCR Mix	Complete reaction mixture	4 x 540 µL	-

## MATERIALS REQUIRED BUT NOT PROVIDED IN THE PRODUCT

- Laminar airflow hood.
- Disposable nitrile powder-free gloves or similar material.
- Vortex mixer.
- Bench microcentrifuge (12,000 - 14,000 RPM).
- Micropipettes and sterile tips with aerosol filter or sterile positive displacement tips (0.5-10 µL, 2-20 µL, 5-50 µL, 50-200 µL, 200-1000 µL).
- Molecular biology grade water.
- Programmable thermostat with optical fluorescence detection system 7300 Real Time PCR System or 7500 Fast Dx Real-Time PCR Instrument calibrated following manufacturer's instructions.
- Programmable thermostat with optical fluorescence detection system cobas z 480 analyzer, calibrated following manufacturer's instructions.

## OTHER PRODUCTS REQUIRED

The reagents for the extraction of DNA from the samples, the positive control of extraction, the positive control of the amplification, the known quantity DNA standards and the consumables **are not** included in this product.

For automatic sample analysis with the instrument «**ELITe InGenius**» (ELITechGroup S.p.A., ref. INT030) the following generic products are required: the extraction cartridges «**ELITe InGenius® SP 200**» (ELITechGroup S.p.A., ref. INT032SP200), the consumables for extraction and amplification of nucleic acids from biological samples «**ELITe InGenius® SP 200 Consumable Set**» (ELITechGroup S.p.A., ref. INT032CS), «**ELITe InGenius® Waste Box**» (ELITechGroup S.p.A., ref. F2102-000), «**ELITe InGenius® PCR Cassette**» (ELITechGroup S.p.A., ref. INT035PCR) and «**300 µL Filter Tips Axygen**» (Axygen BioScience Inc., CA, USA, ref. TF-350-L-R-S).

For automatic DNA extraction, Real Time amplification and interpretation of sample analysis to be analysed, the instrument «**ELITe InGenius**» (ELITechGroup S.p.A., ref. INT030) and the following specific Assay protocols (ELITechGroup S.p.A.), are required:

- for the calibrators «**HSV2 ELITe STD**»,
- for the positive control of amplification «**HSV2 ELITe PC**»,
- for negative control of amplification «**HSV2 ELITe NC**»,
- for samples analysis «**HVS2 ELITe\_WB\_200\_100**», «**HSV2 ELITe\_PL\_200\_100** » and «**HSV2 ELITe\_CSF\_200\_100**».

For automatic sample analysis with the instrument «**ELITe BeGenius**» (ELITechGroup S.p.A., ref. INT040) the following generic products are validated: the extraction cartridges «**ELITe InGenius® SP 200**» (ELITechGroup S.p.A., ref. INT032SP200), the consumables for extraction and amplification of nucleic acids from biological samples «**ELITe InGenius® SP 200 Consumable Set**» (ELITechGroup S.p.A., ref. INT032CS), «**ELITe InGenius® Waste Box**» (ELITechGroup S.p.A., ref. F2102-000), «**ELITe InGenius® PCR Cassette**» (ELITechGroup S.p.A., ref. INT035PCR) and «**1000 µL Filter Tips Tecan**» (Tecan, Switzerland, ref. 30180118).

For automatic DNA extraction, real time amplification and interpretation of sample analysis to be analysed, the instrument «**ELITe BeGenius**» (ELITechGroup S.p.A., ref. INT040) and the following specific Assay Protocols (ELITechGroup S.p.A.), are required:

- for the calibrators «**HSV2 ELITe Be STD**»,
- for the positive control of amplification «**HSV2 ELITe Be PC**»,
- for negative control of amplification «**HSV2 ELITe Be NC**»,
- for samples analysis «**HSV2 ELITe Be\_WB\_200\_100**», «**HSV2 ELITe Be\_PL\_200\_100**» and «**HSV2 ELITe Be\_CSF\_200\_100**».

For automatic DNA extraction from samples to be analyzed, it is validated the use of generic product «**ELITe STAR 200 Extraction kit**» (ELITechGroup S.p.A., ref. INT011EX), kit for extraction of DNA and RNA from non-cellular and cellular samples with the «**ELITe STAR**» instrument (ELITechGroup S.p.A., ref. INT010).

For automatic DNA extraction and preparation of microplates for amplification of samples to be analyzed, it is validated the use of generic product «**ELITe GALAXY 300 Extraction Kit**» (ELITechGroup S.p.A., ref. INT021EX), kit for extraction of DNA and RNA from non-cellular and cellular samples with the instrument «**ELITe GALAXY**» (ELITechGroup S.p.A., ref. INT020).

For automatic DNA extraction from samples to be analyzed, it is also validated the use of the generic products «**NucliSENS® easyMAG® Reagents**» (bioMérieux SA, ref. 280130, 280131, 280132, 280133, 280134, 280135), kits for extraction of nucleic acid from biological samples, with the instrument «**NucliSENS® easyMAG®**» (bioMérieux SA, ref. 200111).

For automatic DNA extraction from samples to be analyzed, the products «**QIAasymphony® DNA Mini Kit**» (QIAGEN GmbH, ref. 931236) and «**QIAasymphony® DSP Virus / Pathogen Midi kit**» (QIAGEN GmbH, ref. 937055), kits for extraction of nucleic acid from biological samples, with the instrument «**QIAasymphony® SP/AS**» (QIAGEN GmbH, ref. 9001297, 9001301) and related generic products are also validated.

For automatic DNA extraction from samples to be analyzed, the product «**MagNA Pure 24 Total NA Isolation Kit**» (Roche, ref. 07658036001), kit for extraction of nucleic acid from biological samples, with the instrument «**MagNA Pure 24 System**» (Roche, ref. 07290519001) is also validated.

As positive control of nucleic acids extraction from non-cellular samples and inhibition control, it is required the use of generic product «**CPE - Internal Control**» (ELITechGroup S.p.A., ref. CTRCPE), a stabilised solution containing two plasmid DNAs and the genomic RNA of MS2 phage.

When a 7300 Real-Time PCR System is used, it is required the use of generic product «**MicroAmp™ Optical 96-Well Reaction Plate**» (Life Technologies, ref. N8010560), microplates with 0.2 mL wells and adhesive sealing sheets for real time amplification.

When a 7500 Fast Dx Real-Time PCR Instrument is used, it is required the use of generic product: «**MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL**» (Life Technologies, ref. 4346906), microplates with 0.1 mL wells and adhesive sealing sheets for real time amplification.

When a cobas z 480 analyzer is used, it is required the use of generic product «**AD-plate 0.3ml**» (Roche, ref. 05232724001), microplates with 0.3 mL wells and adhesive sealing sheets for real time amplification.

If detection of HSV2 DNA is required for qualitative analysis, use the product «**HSV2 - ELITe Positive Control**» (ELITechGroup S.p.A., ref. CTR032PLD) or the product «**HSV2 - ELITe Positive Control RF**» (ELITechGroup S.p.A., ref. CTR032PLD-R), positive control of plasmid DNA.

If detection and quantification of HSV2 DNA is required for quantitative analysis, use the product «**HSV2 ELITe Standard**» (ELITeGroup S.p.A., ref. STD032PLD), four dilutions of known quantity plasmid DNA to obtain the standard curve.

## WARNINGS AND PRECAUTIONS

**This product is designed for *in-vitro* use.**

### General warnings and precautions

Handle and dispose of all biological samples as if they were able to transmit infective agents. Avoid direct contact with the biological samples. Avoid splashing or spraying. Materials that come into contact with the biological samples must be treated for at least 30 minutes with 3% sodium hypochlorite (bleach) or autoclaved for one hour at 121°C before disposal. Do not allow extraction reagents to contact sodium hypochlorite (bleach).

Handle and dispose of all reagents and all materials used to carry out the assay as if they were able to transmit infective agents. Avoid direct contact with the reagents. Avoid splashing or spraying. Waste must be handled and disposed of in compliance with adequate safety standards. Disposable combustible material must be incinerated. Liquid waste containing acids or bases must be neutralised before disposal.

Wear suitable protective clothes and gloves and protect eyes and face.

Never pipette solutions by mouth.

Do not eat, drink, smoke or apply cosmetic products in the work areas.

Carefully wash hands after handling samples and reagents.

Dispose of leftover reagents and waste in compliance with the regulations in force.

Carefully read all the instructions provided in the product before running the assay.

While running the assay, follow the instructions provided in the product.

Do not use the product after the indicated expiry date.

Only use the reagents provided in the product and those recommended by the manufacturer.

Do not use reagents from different batches.

Do not use reagents from other manufacturers.

### Warnings and precautions for molecular biology

Molecular biology procedures, such as nucleic acid extraction, amplification and detection, require qualified and trained staff to avoid the risk of erroneous results, especially due to the degradation of nucleic acids contained in the samples or sample contamination by amplification products.

When amplification session is manually setup, it is necessary to have available separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplification products. Never introduce an amplification product in the area designated for extraction / preparation of amplification reactions.

When amplification session is manually setup, it is necessary to have available lab coats, gloves and tools which are exclusively used for the extraction / preparation of the amplification reactions and for the amplification / detection of amplification products. Never transfer lab coats, gloves or tools from the area designated for the amplification / detection of amplification products to the area designated for the extraction / preparation of the amplification reactions.

The samples must be exclusively used for this type of analysis. Samples must be handled under a laminar airflow hood. Tubes containing different samples must never be opened at the same time. Pipettes used to handle samples must be exclusively used for this specific purpose. The pipettes must be of the positive displacement type or be used with aerosol filter tips. The tips used must be sterile, free from DNases and RNases, free from DNA and RNA.

The reagents must be handled under a laminar airflow hood. The pipettes used to handle the reagents must be exclusively used for this purpose. The pipettes must be of the positive displacement type or be used with aerosol filter tips. The tips used must be sterile, free from DNases and RNases, free from DNA and RNA.

Amplification products must be handled in such a way as to reduce as much as possible dispersion into the environment in order to avoid the possibility of contamination. The pipettes used to handle amplification products must be exclusively used for this purpose.

### Warnings and precautions specific for the components

The **HSV2 Q - PCR Mix** must be stored at -20°C in the dark.

The **HSV2 Q - PCR Mix** must be used within one month from the first opening.

The **HSV2 Q - PCR Mix** can be frozen and thawed for no more than **five times**: further freezing / thawing cycles may cause a loss of product performances.

## ELITe InGenius

## SAMPLES AND CONTROLS

### Samples

This product must be used with the following clinical samples:

#### Whole blood collected in EDTA

The whole blood samples for DNA extraction must be collected in EDTA and identified according to laboratory guidelines, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of three days, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from 200 µL of whole blood is carried out with the **ELITe InGenius** and with **ELITe InGenius Software** version 1.3 (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_WB\_200\_100**. This protocol processes 200 µL of sample, adds the **CPE** Internal Control at 10 µL / extraction and elute the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

#### Plasma collected in EDTA

The plasma samples for nucleic acids extraction must be collected in EDTA according to laboratory guidelines, transported at +2 / +8 °C and stored at +2 / +8 °C for a maximum of three days, otherwise they must be frozen and stored at -20 °C for a maximum of thirty days or at -70 °C for longer periods. It is recommended to split the samples into aliquots before freezing, in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from 200 µL of plasma is carried out with the **ELITe InGenius** and with **ELITe InGenius Software** version 1.3 (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_PL\_200\_100**. This protocol processes 200 µL of sample, adds the **CPE** Internal Control at 10 µL / extraction and elutes the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

#### Cerebrospinal fluid (CSF)

The CSF samples for nucleic acid extraction must be collected according to laboratory guidelines avoiding contamination by patient blood, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of four hours, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from CSF is carried out with the **ELITe InGenius** and with **ELITe InGenius Software** version 1.1 (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_CSF\_200\_100**. This protocol processes 200 µL of sample, add the **CPE** Internal Control at 10 µL / extraction and elute the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

### Other samples

At the moment there are no data available concerning product performance with other clinical samples as: suspensions of leucocytes, suspensions of granulocytes and amniotic fluid.



### Interfering substances

The sample must not contain heparin, in order to prevent the problem of inhibition and the possibility of frequent invalid results.

High quantity of human genomic DNA in the DNA extracted from the sample may inhibit the amplification reaction.

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

### Amplification calibrators and amplification controls

Before analysis of any sample, it is absolutely mandatory to generate and to approve the Calibration curve and the amplification controls for each lot of amplification reagent:

as calibrator set, use the four concentration levels of the **HSV2 ELITe Standard**, in association with protocol «**HSV2 ELITe STD**»  
as amplification Positive Control use the **HSV2 - ELITe Positive Control**, in association with protocol «**HSV2 ELITe PC**»  
as amplification Negative Control, use molecular grade water (not provided with this kit) in association with protocol «**HSV2 ELITe NC**».

**Note:** **ELITe InGenius** system requires approved and valid results of calibration curve and amplification controls for each lot of amplification reagent stored in its database.

The calibration curves, approved and stored in the database, will expire after **60 days**. At expiration date it is necessary to re-run the Q-PCR Standards in association with the amplification reagent lot.

The amplification control results, approved and stored in the database, will expire after **15 days**. At the expiration date it is necessary to re-run the Positive and Negative Controls in association with the amplification reagent lot.

Furthermore, the calibrators and amplification controls must be re-run when:

- a new lot of reagents is started,
- the results of Quality control analysis (see following paragraph) are out of specification,
- any major maintenance service is performed on the instrument.

### Quality controls

The planned validation of the extraction and amplification procedure is recommended. Tested samples or certified reference material can be used. External quality controls shall be used in accordance with local, state, federal accrediting organizations, as applicable.

## ELITe InGenius PROCEDURE

Using the «**HSV2 - ELITe MGB® Kit**» with the system **ELITe InGenius** consists of three steps:

- System readiness verification
- Set up of the session
- Review and approval of results

### Verification of the system readiness

Before starting the sample analysis session, referring to the instrument documentation, it is necessary to:

- switch on the **ELITe InGenius** and select the mode «**CLOSED**»;
- verify that the Calibrators (**HSV2 Q - PCR Standard**) have been run, approved and not expired (status). This can be checked under the "Calibration" menu in the Home page;
- verify that the amplification Controls (**HSV2 - Positive Control**, **HSV2 Negative Control**) have been run, approved and not expired (status). This can be checked under the "Control" menu in the Home page;
- choose the type of run and set up the run, following the instructions Graphical User Interface (GUI) for the session set up and using the Assay Protocols provided by ELITechGroup. These IVD protocols were specifically validated with **ELITe MGB** kits, matrices and **ELITe InGenius** instrument.

The Assay protocols available for sample testing with the product «**HSV2 ELITe MGB Kit**» are described in the table below.

Assay protocols for HSV2 ELITe MGB Kit			
Name	Matrix	Report unitage	Characteristics
<b>HSV2 ELITe MGB kit_WB_200_100</b>	Whole Blood	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Sonication: NO PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
<b>HSV2 ELITe MGB kit_PL_200_100</b>	Plasma	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Sonication: NO PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
<b>HSV2 ELITe MGB kit_CSF_200_100</b>	CSF	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Sonication: NO PCR Mix volume: 20 µL Sample PCR input volume: 20 µL

If the assay protocol of interest is not in the system, contact your local ELITechGroup Customer Service.

Protocols for qualitative analysis are available on request.

### Setup of the session

The **HSV2 ELITe MGB Kit** in association to the **ELITe InGenius** can be used in order to perform:

- A. Integrated run (Extract + PCR),
- B. Amplification run, (PCR only),
- C. Calibration run (PCR only),
- D. Amplification run for Positive and Negative Control run (PCR only),

All the parameters are included in the Assay Protocols available on the instrument and are loaded automatically when the Assay protocol is selected.

**Note:** the **ELITe InGenius** can be connected to the "Laboratory Information System" (LIS) which enables loading the session information. Refer to the instrument manual for more details.

The main steps for the setup of the four types of runs are described here below.

#### A. Integrated run

To set up the integrated run carry out the following steps as per the **SW Graphical User Interface**:

1. Thaw samples at room temperature (~+25 °C) and handle according to laboratory guidelines and to the "Samples and Controls" section.
2. Thaw HSV2 Q - PCR Mix tubes in sufficient number for the session. Each tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Protect the HSV2 Q - PCR Mix from light while thawing because this reagent is photosensitive.

3. Thaw the CPE tubes for the session. Each tube is sufficient for 12 extractions. Mix gently, spin down the content for 5 seconds.
4. Select "Perform Run" from the "Home".
5. Ensure that the Extraction Input Volume is 200 µL and the Extracted Elute Volume is 100 µL.
6. For each Track of interest fill in the "SampleID" (SID) by typing or by scanning the sample barcode.
7. Select the assay protocol to be used in the "Assay" column (i.e. HSV2 ELITe\_WB\_200\_100).
8. Ensure that the "Protocol" displayed is: "Extract + PCR".
9. Select the sample loading position in the "Sample Position" column:  
if a primary tube is used select "Primary Tube",  
if a secondary tube is used select "Sonication Tube".

Click "Next" to continue the setup.

10. Load CPE and HSV2 Q-PCR Mix on the Inventory Block selected by following the GUI instruction. Click "Next" button to continue the setup.
11. Load and check the Tip Racks in the Inventory Area selected by following the GUI instruction. Click "Next" button to continue the setup.
12. Load the "PCR Cassettes", the "ELITe InGenius SP 200" extraction cartridges, all the required consumables and the samples to be extracted in the positions specified in step 8, following the GUI instruction. Click "Next" to continue the setup.
13. Close the instrument door.
14. Press "Start" to start the run.

When the session is finished, the **ELITe InGenius** allows users to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Extracted Sample can be removed from the instrument, capped, identified and stored at -20 °C. Avoid the spilling of the extracted sample.

**Note:** At the end of the run the PCR Cassettes with the reaction products and other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

**Note:** At the end of the run the PCR mix can be kept on board in the refrigerated block up to 16 hours.

#### B. Amplification run

To set up the amplification run carry on the steps below following the GUI:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each tube is sufficient for 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Protect the HSV2 Q - PCR Mix from light while thawing because this reagent is photosensitive.

2. Select "Perform Run" from the "Home screen".
3. Ensure that the Extraction Input Volume is 200 µL and the Extracted Elute Volume is 100 µL.
4. For each Track of interest type the "SampleID" (SID) by typing or by scanning the sample barcode.
5. Select the assay protocol to be used in the "Assay" column (i.e. HSV2 ELITe\_WB\_200\_100).
6. Select "PCR Only" in the "Protocol" column.
7. Ensure the Eluted sample loading position in the "Sample Position" column is "ExtraTube (bottom row)". Click "Next" to continue the setup.
8. Load HSV2 Q-PCR Mix on the Inventory Block selected by following the GUI instruction. Click "Next" to continue the setup.
9. Load and check the Tip Racks in the Inventory Area selected by following the GUI instruction. Click "Next" to continue the setup.
10. Load the "PCR Cassettes" and the extracted Nucleic Acid samples following the GUI instruction. Click "Next" to continue the setup.
11. Close the instrument door.
12. Press "Start" to start the run.

When the session is finished, the **ELITe InGenius** allows users to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Extracted Sample can be removed from the instrument, capped and stored at -20 °C. Avoid the spilling of the Extracted Sample.

**Note:** At the end of the run, the PCR Cassettes and the consumables must be disposed of following all governmental and environmental regulations. Avoid spilling of the reaction products.

**Note:** At the end of the run the PCR mix can be kept on board in the refrigerated block up to 16 hours.

#### C. Calibration run

To set up the Calibration run carry on the steps below following the GUI:

1. Thaw HSV2 Q - PCR Mix tubes in sufficient number for the session. Each tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Protect the HSV2 Q - PCR Mix from light while thawing because this reagent is photosensitive.

2. Thaw HSV2 Q - PCR Standard tubes (Cal1: HSV2 Q - PCR Standards 10<sup>2</sup>, Cal2: HSV2 Q - PCR Standards 10<sup>3</sup>, Cal3: HSV2 Q - PCR Standards 10<sup>4</sup>, Cal4: HSV2 Q - PCR Standards 10<sup>5</sup>). Each tube is sufficient for 4 sessions. Mix gently, spin down the content for 5 seconds.
3. Select "Perform Run" from the "Home screen".
4. Ensure that the Extraction Input Volume is 200 µL and the Extracted Elute Volume is 100 µL.
5. Starting from the Track of interest, select the assay protocol to be used in the "Assay" column (HSV2 ELITe\_STD) and fill with the lot number and expiry date for the HSV2 Q - PCR Standard. Click "Next" button to continue the setup.
6. Load the HSV2 Q-PCR Mix on the Inventory Block selected by following the GUI instruction. Click "Next" to continue the setup.
7. Load and check the Tip Racks in the Inventory Area selected by following the GUI instruction. Click "Next" to continue the setup.
8. Load the Calibrator tubes and PCR Cassettes on board, following the GUI instruction. Click "Next" to continue the setup. Take care to load the PCR Standard fluids to the correct tracks as indicated in the GUI.
9. Close the instrument door.
10. Press "Start" to start the run.

After process completion, the **ELITe InGenius** allows to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Calibrators can be removed from the instrument, capped and stored at -20 °C.

**Note:** The Calibrators can be used for 4 separate sessions of 3 hours each.

**Note:** At the end of the run, the PCR Cassettes and other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

**Note:** At the end of the run the PCR mix can be kept on board in the refrigerated block up to 16 hours.

#### D. Amplification run for Positive Control and Negative Control

To set up the amplification Positive Control and Negative Control run, carry on the steps below following the GUI:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Protect the HSV2 Q - PCR Mix from light while thawing because this reagent is photosensitive.

2. Thaw the product HSV2 - Positive Control for Positive Control amplification. Each tube is sufficient for 4 sessions. Mix gently, spin down the content for 5 seconds.
3. Transfer at least 50 µL the molecular biology grade water for the sessions in one Elution tube, provided with the ELITe InGenius SP Consumable Set.
4. Select "Perform Run" from the "Home screen".
5. In the Track of interest, select the assay protocol to be used in the "Assay" column.
6. For the positive control, select HSV2 ELITe\_PC for the positive control and fill in the lot number and expiry date for the HSV2 Positive Control.
7. For the negative control, select HSV2 ELITe\_NC and fill in the lot number and expiry date for the molecular biology grade water.

8. Click "Next" to continue the setup.
9. Load HSV2 Q-PCR Mix on the Inventory Block selected by following the GUI instruction. Click "Next" to continue the setup.
10. Load and check the Tip Racks in the Inventory Area selected by following the GUI instruction. Click "Next" to continue the setup.
11. Load the amplification PCR cassette, the Positive Control and/or the Negative Control, following the GUI instruction. Click "Next" to continue the setup.
12. Close the instrument door.
13. Press "Start" to start the run.

After process completion, the **ELITe InGenius** allows to view, approve, store the results and to print and save the report.

**Note:** The Positive Control must be run as amplification control, to set up the Control Chart. Four (4) Positive Control values, from 4 different runs are requested to set up the chart. After that, the Positive control values are used for monitoring the amplification step. Refer to the user's manual of the instrument for more details.

**Note:** At the end of the run the remaining Positive Control can be removed from the instrument, capped and stored at -20 °C. Avoid the spilling of the Extracted Sample.

**Note:** The Positive Control can be used for 4 separate sessions of 3 hours each.

**Note:** At the end of the run, the PCR Cassettes and other consumables must be disposed of following all governmental and environmental regulations. Avoid spilling the reaction products.

**Note:** PCR mix can be kept on board in the refrigerated block up to 16 hours.

#### Review and approval of results

The ELITe InGenius monitors target and internal control fluorescence signals for each reaction and automatically applies the Assay Protocol parameters to generate PCR curves which are then interpreted into results.

At the end of the run, the "Results Display" screen is automatically shown. In this screen the sample / Calibrator / Control results and the information regarding the run are shown. From this screen is possible to approve the result, print or save the reports ("Sample Report" or "Track Report").

**Note:** The **ELITe InGenius** can be connected to the "Laboratory Information System" (LIS) which enables uploading the session results to the laboratory data center.. Refer to the instrument manual for more details.

**Note:** For detailed information refer to the **ELITe InGenius** instrument user manual.

The **ELITe InGenius** generates the results using the «**HSV2 ELITe MGB Kit**» through the following procedure:

- A. Validation of Calibration curve,
- B. Validation of amplification Positive Control and Negative Control results,

- C. Validation of sample results,
- D. Sample result reporting.

#### A. Validation of Calibration curve

The **ELITe InGenius software** interprets the PCR results for the specific HSV2 probe ("HSV2") and by the specific Internal Control probe ("IC") in the Calibrators amplification reaction with the HSV2 ELITe STD Assay Protocol parameters. The resulting Ct values are used to validate the system (reagents lot and instrument).

The Calibration curve, specific for the amplification reagent lot, is stored in the database after the approval of the "Administrator" or "Analyst" personnel by following the GUI instruction.

The Calibration curve, specific for the amplification reagent lot, will expire **after 60 days**.

Before analysing any sample, it is absolutely mandatory to generate and to approve the Calibration curve for the lot of amplification reagent used. The availability of Calibration curve results with "Approved" (Status) is shown in the "Calibration" window of the ELITe InGenius software.

**Note:** If the Calibration curve does not meet the acceptance criteria, the "Failed" message is shown on the "Calibration" screen. In this case, the results cannot be approved and the Calibrator amplification reactions must be repeated.

**Note:** When the Calibration Curve is run together with samples and its result is invalid, the entire session is invalid and the amplification of all samples must be repeated.

#### B. Validation of amplification Positive Control and Negative Control results

The **ELITe InGenius software** interprets the PCR results for the specific HSV2 probe ("HSV2") and by the specific Internal Control probe ("IC") in the Positive Control and Negative Control amplification reaction with the HSV2 ELITe\_PC and HSV2 ELITe\_NC Assay Protocols parameters. The resulting Ct values are used to validate the system (reagents lot and instrument).

The amplification Positive Control and Negative Control results, specific for the amplification reagent lot, are stored in the database (Controls) after the approval of the "Administrator" or "Analyst" personnel by following the GUI instruction.

The Positive Control and Negative Control results, specific for the PCR reagent lot, are recorded in the database (Controls). They can be viewed and approved by "Administrator" or "Analyst" users, following the GUI instructions.

The amplification Positive Control results, specific for the amplification reagent lot, will expire **after 15 days**.

Before analysing any sample and after approval of the Calibration curve, it is absolutely mandatory to generate and to approve an amplification Positive Control and Negative Control results for the lot of amplification reagent used. The availability of an amplification Positive Control and Negative Control results with "Approved" (Status) is shown in the "Controls" window of the ELITe InGenius software. If the amplification Positive Control and Negative Control results are missing, generate them as described above.

The **ELITe InGenius software** processes the Positive Control and Negative Control results and generates Control Charts. Four approved Positive Control and Negative Control results are used to set up the initial Control Chart. For subsequent controls, the results are analysed by the software to ensure the system performances are within the acceptance criteria, shown in the Control Chart plots. Refer to the instrument manual for more details.

**Note:** If the Positive Control or Negative Control result does not meet the acceptance criteria, the "Failed" message is shown on the "Controls" screen. In this case, it the results cannot be approved and the Positive Control or Negative Control runs must be repeated.

**Note:** If the Positive Control or Negative Control result is not valid and samples were included in the same run, the samples can be approved but their results are not validated. In this case, the failed Control(s) and samples must all be repeated.

#### C. Validation of Samples results

The ELITe InGenius software interprets the PCR results for the HSV2 probe (Channel "HSV2") and the Internal Control probe (Channel "IC") with the HSV2 ELITe\_WB\_200\_100, HSV2 ELITe\_PL\_200\_100 and HSV2 ELITe\_CSF\_200\_100 Assay Protocol parameters. The resulting HSV2 Ct values are converted to concentration.

Results are shown in "Result Display" module.

The Sample run can be approved when the three conditions reported in the table below are true.

1) Calibration curve	Status
HSV2 Q - PCR Standard	APPROVED
2) Positive Control	Status
HSV2 - Positive Control	APPROVED
3) Negative Control	Status
HSV2 - Negative Control	APPROVED

The sample results are automatically interpreted by the ELITe InGenius software using Assay Protocol parameters. The possible result messages are listed in the table below.

Result of Sample run	Interpretation
HSV2: DNA Detected, quantity equal to XXX copies / mL	<b>HSV2 DNA detected</b> within the measurement range of the assay, quantity as shown.
HSV2: DNA Detected, quantity below LLoQ copies / mL	<b>HSV2 DNA detected</b> below the lower limit of quantification of the assay
HSV2: DNA Detected, quantity beyond ULoQ copies / mL	<b>HSV2 DNA detected</b> beyond the upper limit of quantification of the assay
HSV2: DNA Not Detected or below LoD copies / mL	<b>HSV2 DNA was not detected</b> in the sample. The sample is negative for HSV2 DNA, or its concentration is below the Limit of Detection of the assay.
Invalid - Retest Sample	<b>Not valid assay result</b> caused by Internal Control failure (due to e.g., incorrect extraction, inhibitors carry-over). The test should be repeated.

Samples reported as "HSV2: DNA Detected, quantity below LLoQ" are not suitable for quantification. The concentration of HSV2 DNA detected in the sample is below the level at which it can be accurately quantified. If the sample was diluted before extraction or PCR, it can be retested without dilution.

Samples reported as "HSV2: DNA Detected, quantity beyond ULoQ" are not suitable for quantification. The concentration of HSV2 DNA detected in the sample is above the level at which it can be accurately quantified. The sample may be diluted before extraction or PCR and retested to yield results within the linear range of the assay.

Samples reported as "HSV2 DNA Not Detected or below LoD" are suitable for analysis but HSV2 DNA was not detected. In this case, the sample may be either negative for HSV2 DNA or the HSV2 DNA is present at a concentration below the limit of detection of the assay (see "Performance characteristics").

HSV2 DNA positive samples at a concentration below the LoD, if detected, are reported as "HSV2: DNA Detected, quantity below LLoQ" (see "Performance characteristics").

Samples reported as "Invalid - Retest Sample" are not suitable for result interpretation. In this case, the Internal Control DNA was not efficiently detected, which could be due to problems in the PCR or extraction step (degradation or loss of DNA during the extraction or inhibitors in the eluate), which may cause incorrect results.

If sufficient eluate volume remains, the eluate can be retested (as is or diluted) by an amplification run in "PCR Only" mode. If the second result is invalid, the sample must be retested starting from extraction of a new sample using "Extract + PCR" mode.

**Note:** The results obtained with this assay must be interpreted taking into consideration all the clinical data and the other laboratory test outcomes concerning the patient.

The Sample run results are stored in the database and, if valid, can be approved (Result Display) by "Administrator" or "Analyst" personnel by following the GUI instruction. From the Result Display" window it is possible to print and save the Sample run results as "Sample Report" and "Track Report".

#### D. Samples result reporting

The sample results are stored in the database and can be exported as "Sample Report" and "Track Report".

The "Sample Report" shows the details of a sample run sorted by Sample ID (SID).

The "Track Report" shows the details of a sample run track by selected track.

The "Sample Report" and "Track Report" can be printed and signed by authorized personnel.

## ELITe BeGenius

## SAMPLES AND CONTROLS

### Samples

This product must be used with the following clinical samples:

#### Whole blood collected in EDTA

The whole blood samples for DNA extraction must be collected in EDTA and identified according to laboratory guidelines, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of three days, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from 200 µL of whole blood is carried out with the **ELITe BeGenius** and with **ELITe BeGenius Software** version **2.0** (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_Be\_WB\_200\_100**. This protocol processes 200 µL of sample, adds the **CPE** Internal Control at 10 µL / extraction and elutes the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

#### Plasma collected in EDTA

The plasma samples for nucleic acids extraction must be collected in EDTA according to laboratory guidelines, transported at +2 / +8 °C and stored at +2 / +8 °C for a maximum of three days, otherwise they must be frozen and stored at -20 °C for a maximum of thirty days or at -70 °C for longer periods. It is recommended to split the samples into aliquots before freezing, in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from 200 µL of plasma is carried out with the **ELITe BeGenius** **ELITe** and with **ELITe BeGenius Software** version **2.0** (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_Be\_PL\_200\_100**. This protocol processes 200 µL of sample, adds the **CPE** Internal Control at 10 µL / extraction and elutes the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

#### Cerebrospinal fluid (CSF)

The CSF samples for nucleic acid extraction must be collected according to laboratory guidelines avoiding contamination by patient blood, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of four hours, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when the DNA extraction from CSF is carried out with the **ELITe BeGenius** and with **ELITe BeGenius Software** version **2.0.0** (or later equivalent versions), use the extraction protocol **HSV2 ELITe\_Be\_CSF\_200\_100**. This protocol processes 200 µL of sample, add the **CPE** Internal Control at 10 µL / extraction and elute the nucleic acids in 100 µL.

When the primary tube is used, the volume of the sample varies according to the type of the tube loaded. Refer to the instruction for use of the extraction kit for more information on how to set up and perform the extraction procedure.

#### Other samples:

At the moment there are no data available concerning product performance with other clinical samples such as: suspensions of leucocytes, suspensions of granulocytes and amniotic fluid.

#### Interfering substances



The sample must not contain heparin, in order to prevent the problem of inhibition and the possibility of frequent invalid results.

High quantity of human genomic DNA in the DNA extracted from the sample may inhibit the amplification reaction.

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

#### Amplification calibrators and amplification controls

Before analysis of any sample, it is absolutely mandatory to generate and to approve the Calibration curve and the amplification controls for each lot of amplification reagent:

as calibrator set, use the four concentration levels of the **HSV2 ELITe Standard**, in association with protocol «**HSV2 ELITe\_Be\_STD**»

as amplification Positive Control use the **HSV2 - ELITe Positive Control**, in association with protocol «**HSV2 ELITe\_Be\_PC**»

as amplification Negative Control, use molecular grade water (not provided with this kit) in association with protocol «**HSV2 ELITe\_Be\_NC**».

**Note:** **ELITe BeGenius** system requires approved and valid results of calibration curve and amplification controls for each lot of amplification reagent stored in its database.

The calibration curves, approved and stored in the database, will expire after **60 days**. At expiration date it is necessary to re-run the Q-PCR Standards in association with the amplification reagent lot.

The amplification control results, approved and stored in the database, will expire after **15 days**. At the expiration date it is necessary to re-run the Positive and Negative Controls in association with the amplification reagent lot.

Furthermore, the calibrators and amplification controls must be re-run when:

- a new lot of reagents is started,
- the results of Quality control analysis (see following paragraph) are out of specification,
- any major maintenance service is performed on the instrument.

#### Quality controls

The planned validation of the extraction and amplification procedure is recommended. Tested samples or certified reference material can be used. External quality controls shall be used in accordance with local, state, federal accrediting organizations, as applicable.

### ELITe BeGenius PROCEDURE

Using the «**HSV2 ELITe MGB Kit**» with the system **ELITe BeGenius** consists of three steps:

- System readiness verification
- Set up of the session
- Review and approval of results

#### Verification of the system readiness

Before starting the sample analysis session, referring to the instrument documentation, it is necessary to:

- switch on the **ELITe BeGenius** and select the mode «**CLOSED**»;
- verify that the Calibrators (**HSV2 Q-PCR Standard**) have been run, approved and not expired (status). This can be checked under the «Calibration» menu in the Home page;
- verify that the amplification Controls (**HSV2 - Positive Control**, **HSV2 Negative Control**) have been run, approved and not expired (status). This can be checked under the «Control» menu in the Home page;
- choose the type of run and set up the run, following the instructions Graphical User Interface (GUI) for the session set up and using the Assay Protocols provided by ELITechGroup. These IVD protocols were specifically validated with **ELITe MGB Kits**, matrices and **ELITe BeGenius** instrument.

The Assay protocols available for «**HSV2 ELITe MGB Kit**» are described in the table below.

Assay protocols for « <b>HSV2 ELITe MGB Kit</b> » and <b>ELITe BeGenius</b>			
Name	Matrix	Report unitage	Characteristics
<b>HSV2 ELITe_Be_PL_200_100</b>	Plasma	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Dilution Factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
<b>HSV2 ELITe_Be_WB_200_100</b>	Whole Blood	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Dilution Factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL
<b>HSV2 ELITe_Be_CSF_200_100</b>	CSF	copies/mL	Extraction Input Volume: 200 µL Extracted Elute Volume: 100 µL Internal Control: 10 µL Dilution Factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 20 µL

If the assay protocol of interest is not in the system, contact your local ELITechGroup Customer Service.

#### Setup of the session

The **HSV2 ELITe MGB Kit** can be used on **ELITe BeGenius** to perform:

- A. Sample run, (EXTR + PCR),
- B. Amplification run (PCR only),
- C. Calibration run (PCR only),
- D. Amplification run for Positive and Negative Control run (PCR only).

All the parameters needed for the session are included in the Assay protocol available on the instrument and are automatically recalled when the Assay protocol is selected.

**Note:** The **ELITe BeGenius** system can be linked to the «Location Information Server» (LIS) through which it is possible to load the work session information. Refer to the instrument user's manual for more details.

The main steps for the setup of the four types of runs are described here below.

#### A. Integrated run

To set up the integrated run, carry out the steps below following the **GUI**:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each new tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Thaw **HSV2 PCR Mix** in the dark because this reagent is sensitive to the light.

2. Thaw a sufficient number of CPE tubes for the session. Each new tube is sufficient for 12 extractions. Mix gently, spin down the content for 5 seconds.
3. Select «Perform Run» from the «Home screen».
4. Remove all the Racks from the «Cooler Unit» and place them on the preparation table.
5. Select the «run mode»: «Extract + PCR».
6. Load the samples into the cooling area starting from the L5 Sample Rack.
7. Insert the Rack into the «Cooler Unit». Click «Next» to continue the setup.

**Note:** If secondary tubes are loaded, flag «2 mL Tube». If secondary tubes are not barcoded, type manually the sample ID.

8. Check the Extraction Input Volume (200 µL) and the Extracted Elute Volume (100 µL).

9. Select the assay protocol to be used in the "Assay" column (i.e. HSV2 ELITe\_Be\_PL\_200\_100). Click "Next" to continue the setup.
10. If a second extraction has to be performed, repeat steps 6 to 9 using the L4 Sample Rack.
11. Load the barcoded eluate tubes into cooling area starting from L3 Elution Rack.

**Note:** Elution tubes can be labelled to improve traceability.

12. Insert the L3 Elution Rack into the "Cooler Unit". Click "Next" to continue the setup.
13. Repeat steps 11 and 12 using the L2 Reagent/Elution Rack.
14. Load CPE and HSV2 Q-PCR Mix into the into cooling area.
15. Insert the L1 Reagent Rack into the "Cooler Unit". Click "Next" to continue the setup.
16. Load and check the Tip Racks in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
17. Load the PCR Rack with "PCR Cassette" in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
18. Load the Extraction Rack with the "ELITe InGenius SP 200" extraction cartridges and the required extraction consumables by following the GUI instruction. Click "Next" to continue the setup.
19. Close the instrument door.
20. Press "Start" to start the run.

After process completion, the **ELITe BeGenius** allows the user to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Extracted Sample can be removed from the instrument, capped, identified and stored at -20 °C. Avoid the spilling of the Extracted Sample.

**Note:** At the end of the run the "PCR Cassette" with the reaction products and the consumables must be removed from the instrument and eliminated without producing environmental contaminations. Avoid the spilling of the reaction products.

**Note:** The PCR Mix can be used for 5 independent work sessions of 3 hours each or can be kept on board in the refrigerated block up to 3 consecutive work sessions of 3 hours each. Mix gently and spin down the content for 5 seconds before starting the next session.

## B. Amplification run

To set up the amplification run, with eluted samples, carry out the steps below following the GUI:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each new tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Thaw **HSV2 PCR Mix** in the dark because this reagent is sensitive to the light.

2. Select "Perform Run" from the "Home screen".
3. Remove Racks 1, 2 and 3 from the "Cooler Unit" and place them on the preparation table.
4. Select the "run mode": "PCR Only".
5. Load the samples into the cooling area starting from the L3 Elution Rack.
6. Insert the Rack into the "Cooler Unit". Click "Next" to continue the setup.
7. Even if extraction is not performed, check the Extraction Input Volume (200 µL) and the Extracted Elute Volume (100 µL).
8. Select the assay protocol to be used in the "Assay" column (e.g. HSV2 ELITe\_Be\_PL\_200\_100). Click "Next" to continue the setup.
9. Load HSV2 Q-PCR Mix into the cooling area.
10. Insert the Rack into the "Cooler Unit". Click "Next" to continue the setup.
11. Load and check the Tip Racks in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.

12. Load the PCR Rack with "PCR Cassette" in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
13. Close the instrument door.
14. Press "Start" to start the run.

After process completion, the **ELITe BeGenius** allows the user to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Extracted Sample can be removed from the instrument, capped, identified and stored at -20 °C. Avoid the spilling of the Extracted Sample.

**Note:** At the end of the run the "PCR Cassette" with the reaction products must be removed from the instrument and eliminated without producing environmental contaminations. Avoid the spilling of the reaction products.

**Note:** The PCR Mix can be used for 5 independent work sessions of 3 hours each or can be kept on board in the refrigerated block up to 3 consecutive work sessions of 3 hours each. Mix gently and spin down the content for 5 seconds before starting the next session.

## C. Calibration run

To set up the Calibration run, with the Q-PCR Standards, carry out the steps below following the GUI:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each new tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Thaw **HSV2 PCR Mix** in the dark because this reagent is sensitive to the light.

2. Thaw the HSV2 Q - PCR Standard tubes (Cal1: HSV2 Q-PCR Standards 10<sup>2</sup>, Cal2: HSV2 Q-PCR Standards 10<sup>3</sup>, Cal3: HSV2 Q-PCR Standards 10<sup>4</sup>, Cal4: HSV2 Q-PCR Standards 10<sup>5</sup>). Each tube is sufficient for 4 sessions. Mix gently, spin down the content for 5 seconds.
3. Select "Perform Run" from the "Home screen".
4. Remove Racks 1, 2 and 3 from the "Cooler Unit" and place them on the preparation table.
5. Select the "run mode": "PCR Only".
6. Load the Calibrator tubes into the L3 Elution Rack.
7. Insert the Rack into the "Cooler Unit". Click "Next" to continue the setup.
8. Even if extraction is not performed, check the Extraction Input Volume (200 µL) and the Extracted Elute Volume (100 µL).
9. Select the assay protocol to be used in the "Assay" column (HSV2 ELITe\_Be\_STD). Click "Next" button to continue the setup.
10. Load HSV2 Q-PCR Mix into the L2 Reagent/Elution Rack.
11. Insert the L2 Reagent/Elution Rack into the "Cooler Unit". Click "Next" to continue the setup.
12. Load and check the Tip Racks in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
13. Load the PCR Rack with "PCR Cassette" in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
14. Close the instrument door.
15. Press "Start" to start the run.

After process completion, the **ELITe BeGenius** allows the user to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Calibrators can be removed from the instrument, capped and stored at -20 °C. Avoid the spilling of the Q-PCR Standards.

**Note:** At the end of the run the "PCR Cassette" with the reaction products must be removed from the

instrument and disposed of without producing environmental contaminations. Avoid any spilling of the reaction products.

**Note:** The PCR Mix can be used for 5 independent work sessions of 3 hours each or can be kept on board in the refrigerated block up to 3 consecutive work sessions of 3 hours each. Mix gently and spin down the content for 5 seconds before starting the next session.

#### D. Amplification run for Positive Control and Negative Control

To set up the Positive Control and Negative Control run, carry out the steps below following the GUI:

1. Thaw a sufficient number of HSV2 Q - PCR Mix tubes for the session. Each new tube is sufficient for preparing 24 reactions in optimal reagent consumption conditions. Mix gently, spin down the content for 5 seconds.

**Note:** Thaw **HSV2 PCR Mix** in the dark because this reagent is sensitive to the light.

2. Thaw the product HSV2 - ELITe Positive Control, for Positive Control amplification. Each tube is sufficient for 4 sessions. Mix gently, spin down the content for 5 seconds.
3. Transfer at least 50 µL of the molecular biology grade water (as Negative Control) for the sessions in one Elution tube, provided with the ELITe InGenius SP Consumable Set.
4. Select "Perform Run" from the "Home screen".
5. Remove Racks 1, 2 and 3 from the "Cooler Unit" and place them on the preparation table.
6. Select the "run mode": "PCR Only".
7. Load the Positive Control and Negative Control tubes into the L3 Elution Rack.
8. Insert the Rack into the "Cooler Unit". Click "Next" to continue the setup.
9. Even if extraction is not performed, check the Extraction Input Volume (200 µL) and the Extracted Elute Volume (100 µL).
10. Select the assay protocol to be used in the "Assay" column (HSV2 ELITe\_Be\_PC and HSV2 ELITe\_Be\_NC). Click "Next" button to continue the setup.
11. Load HSV2 Q-PCR Mix into the L2 Reagent/Elution Rack.
12. Insert the L2 Reagent/Elution Rack into the "Cooler Unit". Click "Next" to continue the setup.
13. Load and check the Tip Racks in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
14. Load the PCR Rack with "PCR Cassette" in the Inventory Area by following the GUI instruction. Click "Next" to continue the setup.
15. Close the instrument door.
16. Press "Start" to start the run.

After process completion, the **ELITe BeGenius** allows the user to view, approve, store the results and to print and save the report.

**Note:** At the end of the run the remaining Positive Control can be removed from the instrument, capped and stored at -20 °C. Avoid the spilling of the Positive Controls.

**Note:** At the end of the run the "PCR Cassettes" with the reaction products must be removed from the instrument and disposed of without producing environmental contaminations. Avoid any spilling of the reaction products.

**Note:** The PCR Mix can be used for 5 independent work sessions of 3 hours each or can be kept on board in the refrigerated block up to 3 consecutive work sessions of 3 hours each. Mix gently and spin down the content for 5 seconds before starting the next session.

#### Review and approval of results

At the end of the run, the "Results Display" screen is automatically shown. In this screen the sample / Calibrator / Control results and the information regarding the run are shown. From this screen is possible to approve the result, print or save the reports ("Sample Report" or "Track Report").

**Note:** The **ELITe BeGenius** system can be linked to the "Location Information Server" (LIS) through which it

is possible send the work session results to the laboratory data center. Refer to the instrument user's manual for more details.

The **ELITe BeGenius** generates the results using the HSV2 ELITe MGB Kit through the following procedure:

- A. Validation of Calibration curve,
- B. Validation of amplification Positive Control and Negative Control results,
- C. Validation of sample results,
- D. Sample result reporting.

**Note:** Please, refer to the same **ELITe InGenius** chapters for the details.

### PERFORMANCE CHARACTERISTICS ELITe InGenius and ELITe BeGenius

#### Analytical sensitivity: limit of detection

The analytical sensitivity of this assay, as Limit of Detection (LoD) of the DNA amplification, allows detecting the presence of about 10 copies in 20 µL of DNA added to the amplification reaction.

The LoD of this assay was tested using plasmid DNA containing amplification product whose initial concentration was measured by spectrophotometer. The plasmid DNA was diluted to a titre of 10 copies / 20 µL in human genomic DNA at a titre of 500 ng / 20 µL. This sample was tested in 24 replicates carrying out the amplification by ELITechGroup S.p.A. products on two different instruments. The results are summed up in the following table.

Samples	N	positive	negative
10 copies plasmid DNA + 500 ng of human genomic DNA	24	24	0

The Limit of Detection (LoD) of HSV2 ELITe MGB Kit was verified in association with **Whole Blood, Plasma collected in EDTA** and **CSF** samples in association with **ELITe InGenius** and **ELITe BeGenius** systems (Extr + PCR mode).

#### For Whole Blood

The LoD of this assay was verified by testing 20 replicates of Whole blood sample spiked at 171 copies / mL on **ELITe InGenius** and **ELITe BeGenius** systems in "Extract + PCR" mode. The samples were spiked using the certified reference material "Heat Inactivated HSV Type 2 Culture Fluid" (ZeptoMetrix Corporation).

The LoD is confirmed if at least 18 out of 20 replicates give a positive result. The results are reported in the following tables.

Limit of Detection for Whole Blood samples and ELITe InGenius					
Sample	LoD	N	Valid	Positive	Negative
Whole blood collected in EDTA	171 copies / mL	20	20	20	0

Limit of Detection for Whole Blood samples and ELITe BeGenius					
Sample	LoD	N	Valid	Positive	Negative
Whole blood collected in EDTA	171 copies / mL	20	20	20	0

The LoD value for HSV2 target was confirmed at 171 copies / mL for Whole Blood collected in EDTA.

#### For Plasma

The LoD of this assay was verified by testing 20 replicates of Plasma sample spiked at 119 copies / mL on **ELITe InGenius** and **ELITe BeGenius** systems in "Extract + PCR" mode. The samples were spiked using the certified reference material "Heat Inactivated HSV Type 2 Culture Fluid" (ZeptoMetrix Corporation).

The LoD is confirmed if at least 18 out of 20 replicates give a positive result. The results are reported in the following tables.

Limit of Detection for Plasma samples and ELITe InGenius					
Sample	LoD	N	Valid	Positive	Negative
Plasma collected in EDTA	119 copies / mL	20	20	20	0

Limit of Detection for Plasma samples and ELITe BeGenius					
Sample	LoD	N	Valid	Positive	Negative
Plasma collected in EDTA	119 copies / mL	20	20	20	0

The LoD value for HSV2 target was confirmed at 119 copies / mL for Plasma collected in EDTA.

#### For Cerebrospinal Fluid (CSF)

The LoD of this assay was verified by testing 20 replicates of CSF sample spiked at 119 copies / mL on ELITe InGenius and ELITe BeGenius systems in "Extract + PCR" mode. The samples were spiked using the reference material Herpes Simplex Virus Type 2 (HSV-2) Culture Fluid Heat Inactivated, ZeptoMetrix.

The LoD is confirmed if at least 18 out of 20 replicates give a positive result.

The results are reported in the following tables.

Limit of Detection for CSF samples and ELITe InGenius					
Sample	LoD	N	Valid	Positive	Negative
Cerebrospinal Fluid	119 copies / mL	20	20	20	0

Limit of Detection for CSF samples and ELITe BeGenius					
Sample	LoD	N	Valid	Positive	Negative
Cerebrospinal Fluid	119 copies / mL	20	20	19	1

The LoD value for HSV2 target was confirmed at 119 copies / mL for Cerebrospinal Fluid.

#### Linear measuring range and Limits of quantification


The linear measuring range of HSV2 ELITe MGB Kit used in association with **Whole Blood, Plasma collected in EDTA** and **CSF** samples in association with **ELITe InGenius** and **ELITe BeGenius** was verified with a panel of HSV2 dilutions. The panel was prepared by diluting "Heat Inactivated HSV Type 2 Culture Fluid" (ZeptoMetrix Corporation) in HSV2 DNA - negative matrices.

##### For Whole Blood

The panel consisted of ten dilution points from about  $2.5 \times 10^7$  copies / mL to about 171 copies / mL. Each sample of the panel was tested in 3 replicates.

The analysis of the obtained data, performed by linear regression analysis, demonstrated that the assay in association with Whole Blood samples shows a linear response for all the dilutions with a Square Correlation Coefficient (R<sup>2</sup>) equal to 0.997 for **ELITe InGenius** and 0.997 for **ELITe BeGenius**.

The results are reported in the following graphs.



The Lower Limit of Quantification (LLOQ) was set at, the LoD concentration, that gives quantitative results precise (Standard Deviation equal to 0.1081 Log copies / mL for **ELITe InGenius** and 0.1900 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to -0.1221 Log copies / mL for **ELITe InGenius** and 0.1591 Log copies / mL for **ELITe BeGenius**): 171 copies / mL.


The Upper Limit of Quantification (ULOQ) was set at, the highest concentration tested, that gives quantitative results precise (Standard Deviation equal to 0.0472 Log copies / mL for **ELITe InGenius** and 0.1191 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to 0.2744 Log copies / mL for **ELITe InGenius** and 0.3029 Log copies / mL for **ELITe BeGenius**): 25,000,000 copies / mL.

The final results are summarized in the following table.

Linear measuring range for Whole Blood samples and ELITe InGenius and ELITe BeGenius		
Unit of measure	lower limit	upper limit
copies / mL	171	25,000,000

The results obtained by **ELITe InGenius** and **ELITe BeGenius** were analysed by orthogonal and linear regression in order to calculate the correlation between the methods.

The results are summed up in the following figure.



In this test, the orthogonal regression analysis generated a slope equal to 0.943 (95% CI: 0.919; 0.967) and an intercept equal 0.383 (95% CI: 0.276; 0.491). The linear regression analysis generated a R<sup>2</sup> of 0.996.


##### For Plasma

The panel consisted of ten dilution points from about  $2.5 \times 10^7$  copies / mL to about 119 copies / mL. Each sample of the panel was tested in 3 replicates.

The analysis of the obtained data, performed by linear regression analysis, demonstrated that the assay in association with Plasma samples shows a linear response for all the dilutions with a Square Correlation Coefficient (R<sup>2</sup>) equal to 0.992 for **ELITe InGenius** and 0.996 for **ELITe BeGenius**.



The results are reported in the following graphs.



The Lower Limit of Quantification (LLoQ) was set at the LoD concentration, that gives quantitative results precise (Standard Deviation equal to 0.2529 Log copies / mL for **ELITe InGenius** and 0.2931 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to 0.0503 Log copies / mL for **ELITe InGenius** and 0.0005 Log copies / mL for **ELITe BeGenius**): 119 copies / mL.


The Upper Limit of Quantification (ULoQ) was set at, the highest concentration tested, that gives quantitative results precise (Standard Deviation equal to 0.0234 Log copies / mL for **ELITe InGenius** and 0.0761 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to -0.0258 Log copies / mL for **ELITe InGenius** and -0.0646 Log copies / mL for **ELITe BeGenius**): 25,000,000 copies / mL.

The final results are summarized in the following table.

Linear measuring range for Plasma samples and ELITe InGenius and ELITe BeGenius		
Unit of measure	lower limit	upper limit
copies / mL	119	25,000,000

The results obtained by **ELITe InGenius** and **ELITe BeGenius** were analysed by orthogonal and linear regression in order to calculate the correlation between the methods.

The results are summed up in the following figure.




In this test, the orthogonal regression analysis generated a slope equal to 1.030 (95% CI: 1.005; 1.056) and an intercept equal to -0.137 (95% CI: -0.257; -0.017). The linear regression analysis generated a R2 of 0.996.

**For Cerebrospinal Fluid CSF:**

The panel consisted of ten dilution points from about 2.5 x 10<sup>7</sup> copies / mL to about 119 copies / mL. Each sample of the panel was tested in 4 replicates. Two replicates were excluded as outliers.

The analysis of the obtained data, performed by linear regression analysis, demonstrated that the assay in association with CSF samples shows a linear response for all the dilutions with a Square Correlation Coefficient (R2) equal to 0.997 for **ELITe InGenius** and 0.996 for **ELITe BeGenius**.

The results are reported in the following graphs.



The Lower Limit of Quantification (LLoQ) was set at the LoD concentration, that gives quantitative results precise (Standard Deviation equal to 0.0261 Log copies / mL for **ELITe InGenius** and 0.2349 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to -0.2938 Log copies / mL for **ELITe InGenius** and 0.1138 Log copies / mL for **ELITe BeGenius**): 119 copies / mL.


The Upper Limit of Quantification (ULoQ) was set at the highest concentration tested, that gives quantitative results precise (Standard Deviation equal to 0.0637 Log copies / mL for **ELITe InGenius** and 0.0277 Log copies / mL for **ELITe BeGenius**) and accurate (Bias equal to -0.0345 Log copies / mL for **ELITe InGenius** and 0.1823 Log copies / mL for **ELITe BeGenius**): 25,000,000 copies / mL.

The final results are summarized in the following table.

Linear measuring range for CSF samples and ELITe InGenius and ELITe BeGenius		
Unit of measure	lower limit	upper limit
copies / mL	119	25,000,000

The results obtained by **ELITe InGenius** and **ELITe BeGenius** were analysed by orthogonal and linear regression in order to calculate the correlation between the methods.

The results are summed up in the following figure.



In this test, the orthogonal regression analysis generated a slope equal to 1.004 (95% CI: 0.9822; 1.0267) and an intercept equal to 0.324 (95% CI: -0.4334; -0.2147). The linear regression analysis generated a R<sup>2</sup> of 0.995.

#### Repeatability

The Repeatability of results obtained by the product HSV2 ELiTe MGB Kit in association with the **ELiTe InGenius** and **ELiTe BeGenius** systems was tested by analysing a panel of Whole blood samples collected in EDTA. The panel included one negative sample and two samples spiked by HSV2 certified reference material "Heat Inactivated HSV Type 2 Culture Fluid, (Zeptomatrix) at concentration of 3 x LoD (about 513 cp / mL) and of 10 x LoD (about 1710 cp / mL).

The Intra – Session Repeatability on **ELiTe InGenius** was obtained through the analysis of panel samples in eight replicates, in two runs per day, with the same lot of product, with the same instrument, by the same operator, on the same day. Samples were processed in randomized positions.

The Inter – Session Repeatability on **ELiTe InGenius** was obtained through the analysis of panel samples in eight replicates, in two runs per day, with the same lot of product, with the same instrument, by the same operator, on two different days. Samples were processed in randomized positions.

The Ct values of the target and of Internal Control were used to calculate the %CV in order to evaluate the Repeatability as imprecision.

A summary of results is shown in the tables below.

Intra – Session Repeatability ELiTe InGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24 / 24	22.91	0.51	2.24
3 x LoD	8 / 8	35.65	0.87	2.45				
10 x LoD	8 / 8	33.96	0.20	0.58				

Inter – Session Repeatability ELiTe InGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 16	N.A.	N.A.	N.A.	48 / 48	23.40	0.82	3.49
3 x LoD	16 / 16	35.87	0.71	1.99				
10 x LoD	16 / 16	33.92	0.22	0.66				

In the Repeatability test on **ELiTe InGenius**, the assay detected the HSV2 target as expected and showed Ct values with %CV below 5% for HSV2 and for Internal Control.

The Intra – Session Repeatability on **ELiTe BeGenius** was obtained through the analysis of panel samples in eight replicates, in one run per day, with the same lot of product, with the same instrument, on the same day. Samples were processed in randomized positions.

The Inter – Session Repeatability on **ELiTe BeGenius** was obtained through the analysis of panel samples in eight replicates, in one run per day, with the same lot of product, with the same instrument, on two different days. Samples were processed in randomized positions.

The Ct values of the target and of Internal Control were used to calculate the %CV in order to evaluate the Repeatability as imprecision.

A summary of results is shown in the tables below.

Intra – Session Repeatability ELiTe BeGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24/24	27.16	0.87	3.19
3 x LoD	8 / 8	37.58	0.67	1.78				
10 x LoD	8 / 8	35.25	0.56	1.58				

Inter – Session Repeatability ELiTe BeGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 16	N.A.	N.A.	N.A.	48 / 48	26.87	0.94	3.49
3 x LoD	16 / 16	37.59	0.66	1.75				
10 x LoD	16 / 16	35.41	0.84	2.37				

In the Repeatability test on **ELiTe BeGenius**, the assay detected the HSV2 target as expected and showed Ct values with %CV below 5% for HSV2 and for Internal Control.

#### Reproducibility

The Reproducibility of results obtained by the product HSV2 ELiTe MGB Kit in association with the **ELiTe InGenius** and **ELiTe BeGenius** systems was tested by analysing a panel of Whole blood samples collected in EDTA. The panel included one negative sample and two samples spiked by HSV2 certified reference material "Heat Inactivated HSV Type 2 Culture Fluid, (Zeptomatrix) at concentration of 3 x LoD (about 513 cp / mL) and of 10 x LoD (about 1710 cp / mL).

The Inter – Instrument Reproducibility on **ELiTe InGenius** was obtained through the analysis of panel samples in eight replicates, in one run per day, in two days, with two different instruments by two different operators. Samples were processed in randomized positions on **ELiTe InGenius** system in "Extract + PCR" mode.

The Inter – Batch Reproducibility on **ELiTe InGenius** was obtained through the analysis of panel samples in eight replicates, in two runs per day, with two different lots and the same instrument. Samples were processed in randomized positions on **ELiTe InGenius** system in "Extract + PCR" mode.

The Ct values of the target and of Internal Control were used to calculate the %CV in order to evaluate the Reproducibility as imprecision.

A summary of results is shown in the table below.

Inter – Instrument Reproducibility ELiTe InGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24 / 24	23.91	0.35	1.45
3 x LoD	8 / 8	36.33	0.55	1.51				
10 x LoD	8 / 8	34.46	0.26	0.86				

Inter – Batch Reproducibility ELiTe InGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24 / 24	24.18	0.42	1.76
3 x LoD	8 / 8	36.45	0.46	1.26				
10 x LoD	8 / 8	34.66	0.26	0.76				

In the Reproducibility test on **ELITe InGenius**, the assay detected the HSV2 target as expected and showed Ct values with %CV below 5% for HSV2 and for Internal Control.

The Inter – Instrument Reproducibility on **ELITe BeGenius** was obtained through the analysis of panel samples in eight replicates, in one run per day, in two days, with two different instruments by two different operators. Samples were processed in randomized positions on **ELITe BeGenius** system in “Extract + PCR” mode.

The Inter – Batch Reproducibility on **ELITe BeGenius** was obtained through the analysis of panel samples in eight replicates, in two runs per day, with two different lots and the same instrument. Samples were processed in randomized positions on **ELITe BeGenius** system in “Extract + PCR” mode.

The Ct values of the target and of Internal Control were used to calculate the %CV in order to evaluate the Reproducibility as imprecision.

A summary of results is shown in the table below.

Inter – Instrument Repeatability ELITe BeGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24 / 24	28.32	0.41	1.46
3 x LoD	8 / 8	37.43	0.42	1.15				
10 x LoD	8 / 8	35.88	0.25	0.75				

Inter – Batch Repeatability ELITe BeGenius								
Sample	HSV2				Internal Control			
	Pos. / Rep.	Mean Ct	SD	% CV	Pos. / Rep.	Mean Ct	SD	% CV
Negative	0 / 8	N.A.	N.A.	N.A.	24 / 24	27.86	0.66	2.37
3 x LoD	8 / 8	36.90	0.58	1.59				
10 x LoD	8 / 8	35.17	0.49	1.39				

In the Reproducibility test on **ELITe BeGenius**, the assay detected the HSV2 target as expected and showed Ct values with %CV below 5% for HSV2 and for Internal Control.

#### Analytical sensitivity: reproducibility with certified reference material

The analytical sensitivity of the assay, as reproducibility of value of a calibrated reference material, was evaluated using as reference material the calibrated panel «AcroMetrix HSV2 Plasma Panel» (Life Technologies, US). Each sample of the panel was tested in 2 replicates carrying out the whole procedure of analysis, extraction, amplification, detection and result interpretation with **ELITe InGenius** and ELITechGroup S.p.A. products.

The results are reported in the following table.

Tests with calibrated reference materials and ELITe InGenius				
Sample	Nominal titre copies / mL	Nominal titre Log <sub>10</sub> copies / mL	Positive / Replicates	Mean results Log <sub>10</sub> copies / mL
HSV2 DNA 1E6	10 <sup>6</sup>	6.000	2/2	6.035
HSV2 DNA 1E5	10 <sup>5</sup>	5.000	2/2	5.067
HSV2 DNA 1E4	10 <sup>4</sup>	4.000	2/2	3.978
HSV2 DNA 1E3	10 <sup>3</sup>	3.000	2/2	3.190
HSV2 DNA 1E2	10 <sup>2</sup>	2.000	2/2	1.711

All samples were correctly detected as positive with a titre that was within the expected value  $\pm 0.5$  Log.

Further tests were carried out using as reference material QCMD 2014 Herpes Simplex Virus DNA EQA Panel (Qnostics Ltd, UK) a panel of HSV2 dilutions. Each sample of the panel was tested in 2 replicates carrying out the whole procedure of analysis, extraction, amplification, detection and result interpretation with **ELITe InGenius** and ELITechGroup S.p.A. products.

The results are reported in the following table.

Tests with calibrated reference materials and ELITe InGenius		
Sample	Sample Status	Positive / Replicates
HSV DNA14-01	HSV1 positive, HSV2 negative	0/2
HSV DNA14-02	HSV1 and HSV2 negative	0/2
HSV DNA14-03	HSV1 positive, HSV2 negative	0/2
HSV DNA14-04	HSV1 positive, HSV2 negative	0/2
HSV DNA14-05	HSV1 positive, HSV2 negative	0/2
HSV DNA14-06	HSV2 Frequently detected	2/2
HSV DNA14-07	HSV2 Detected	2/2
HSV DNA14-08	HSV2 Frequently detected	2/2
HSV DNA14-09	HSV1 positive, HSV2 negative	0/2
HSV DNA14-10	HSV1 and HSV2 negative	0/2

All samples were correctly detected.

#### Diagnostic sensitivity: confirmation of positive samples

The diagnostic sensitivity of the assay, as confirmation of positive clinical samples, was evaluated by analyzing some clinical samples of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid positive for HSV2 DNA in association with **ELITe InGenius**. As **ELITe BeGenius** has equivalent analytical performances to **ELITe InGenius**, the diagnostic performances of the assay performed on the two instruments are also considered equivalent. Therefore, the diagnostic sensitivity of the assay obtained in association with **ELITe InGenius** is also applicable to **ELITe BeGenius**.

The diagnostic sensitivity was evaluated using 23 samples of whole blood collected in EDTA negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV2MQP01-Medium sample, from “HSV2 Molecular Q Panel” (Qnostics Ltd, UK) at a titre of 750 copies / mL; 30 samples of plasma collected in EDTA negative for HSV2 DNA, that were spiked for HSV2 DNA adding “HSV2 ELITe-IQC High” (ELITech Group S.p.A.) at a titre of 750 copies / mL and 20 samples of CSF negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV2 ELITe-IQC High” (ELITech Group S.p.A.) at a titre of 750 copies / mL.

Each sample was tested carrying out the whole analysis procedure, extraction, amplification, detection and result interpretation with **ELITe InGenius** and ELITechGroup S.p.A. products. The results are summed up in the following table.

Samples	N	positive	negative
Whole blood collected in EDTA and positive for HSV2 DNA	23	22	1
Plasma collected in EDTA and spiked with HSV2 DNA	30	30	0
Cerebrospinal fluid and spiked with HSV2 DNA	20	20	0

All plasma and cerebrospinal fluid samples were valid and positive. All whole blood samples were valid; 22 samples were confirmed as positive. One sample tested discrepant negative possibly because of inhibition. The diagnostic sensitivity of the assay in this test was equal to 99%.

#### Diagnostic specificity: confirmation of negative samples

The diagnostic specificity of the assay, as confirmation of negative samples, was evaluated by analyzing some clinical samples of whole blood collected in EDTA, plasma collected in EDTA and cerebrospinal fluid negative for HSV2 DNA. in association with **ELITe InGenius**. As **ELITe BeGenius** has equivalent analytical performances to **ELITe InGenius**, the diagnostic performances of the assay performed on the two instruments are also considered equivalent. Therefore, the diagnostic specificity of the assay obtained in association with **ELITe InGenius** is also applicable to **ELITe BeGenius**.

The diagnostic specificity was evaluated using 34 whole blood samples collected in EDTA from healthy donors that were presumably negative for HSV2 DNA; 39 plasma samples collected in EDTA from healthy donors that were presumably negative for HSV2 DNA and 22 CSF samples from healthy donors that were presumably negative for HSV2 DNA. Each sample was tested carrying out the whole analysis procedure, extraction, amplification, detection and result interpretation with **ELITe InGenius** and ELITechGroup S.p.A. products.

The results are summed up in the following table.

Samples	N	positive	negative
Whole blood collected in EDTA and negative for HSV2 DNA	34	0	34
Plasma collected in EDTA and negative for HSV2 DNA	39	0	39
Cerebrospinal fluid and negative for HSV2 DNA	22	0	22

All whole blood, plasma and cerebrospinal fluid samples were valid and negative. The Internal Control Ct (IC Ct) cut-off value is set at 35. The diagnostic specificity of the assay in this test was equal to 100%.

**ABI 7500 Fast Dx Real-Time PCR Instrument**  
**ABI 7300 Real-Time System**

**SAMPLES AND CONTROLS**

**Samples**

This product must be used with **DNA extracted** from the following biological samples: cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

**Cerebrospinal fluid (CSF)**

The CSF samples for nucleic acid extraction must be collected according to laboratory guidelines avoiding contamination by patient blood, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of four hours, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods.

It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing.

When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when you carry out the DNA extraction from cerebrospinal fluid with **ELITe STAR** and with **software version 3.4.13** (or later equivalent versions) use the extraction protocol **UUNI\_E100S200\_ELI**, that uses 200 µL of sample and elutes the extract in 100 µL. The samples in primary tubes can be directly loaded on **ELITe STAR**. A minimum volume of 700 µL is always required for each sample. Add **200 µL** of **CPE** into Proteinase-Carrier tube as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from cerebrospinal fluid with the **ELITe GALAXY** with **software version 1.3.1** (or later equivalent versions) use the extraction protocol **xNA Extraction (Universal)**, that uses 300 µL of sample and elutes the extract in 200 µL. Samples in primary tubes can be directly loaded on «**ELITe GALAXY**». A minimum volume 400-650 µL, dependent on the tube class used, is always required for each sample. Add **10 µL / sample** of **CPE**. The **CPE** must be added to **IC + Carrier solution** as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from cerebrospinal fluid with the instrument «**NucliSENS® easyMAG®**», please follow the extraction protocol **Generic 2.0.1** and follow these directions: transfer **500 µL** of sample in the 8 well strip and run the extraction. After the 10 minutes incubation, add **5 µL** of **CPE** for the internal control before adding the **NucliSENS® easyMAG® Magnetic Silica** and proceed with the extraction. Elute the nucleic acids in **100 µL** of elution buffer.

**Whole blood collected in EDTA**

The whole blood samples for DNA extraction must be collected in EDTA according to laboratory guidelines, transported at +2 / +8 °C and stored at +2 / +8 °C for a maximum of three days, otherwise they must be frozen and stored at -20 °C for a maximum of thirty days or at -70 °C for longer periods.

It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing.

When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when you carry out the DNA extraction from whole blood with **ELITe STAR** and with **software version 3.4.13** (or later equivalent versions) use the extraction protocol **UUNI\_E100S200\_ELI**, that uses 200 µL of sample and elutes the extract in 100 µL. The samples in primary tubes can be directly loaded on «**ELITe STAR**». A minimum volume of 700 µL is always required for each sample. Add **200 µL** of **CPE** into Proteinase-Carrier tube as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from whole blood with the **ELITe GALAXY** with **software version 1.3.1** (or later equivalent versions) use the extraction protocol **xNA Extraction (Universal)**, that uses 300 µL of sample and elutes the extract in 200 µL. Samples in primary tubes can be directly loaded on «**ELITe GALAXY**». A minimum volume 400-650 µL, dependent on the tube class used, is always required for each sample. Add **10 µL / sample** of **CPE**. The **CPE** must be added to **IC + Carrier solution** as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from whole blood using «**EXTRABlood**» kit, please, follow the manual of Instructions for Use: start from **200 µL** of sample (no more than 2 millions of leucocytes), elute the DNA in **100 µL** of elution buffer.

**Plasma collected in EDTA**

The plasma samples for nucleic acid extraction must be collected in EDTA according to laboratory guidelines, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of three days, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods.

It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing.

When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when you carry out the DNA extraction from plasma with **ELITe STAR** and with **software version 3.4.13** (or later equivalent versions) use the extraction protocol **UUNI\_E100S200\_ELI**, that uses 200 µL of sample and elutes the extract in 100 µL. The samples in primary tubes can be directly loaded on «**ELITe STAR**». A minimum volume of 700 µL is always required for each sample. Add **200 µL** of **CPE** into Proteinase-Carrier tube as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from plasma with the **ELITe GALAXY** with **software version 1.3.1** (or later equivalent versions) use the extraction protocol **xNA Extraction (Universal)**, that uses 300 µL of sample and elutes the extract in 200 µL. Samples in primary tubes can be directly loaded on «**ELITe GALAXY**». A minimum volume 400-650 µL, dependent on the tube class used, is always required for each sample. Add **10 µL / sample** of **CPE**. The **CPE** must be added to **IC + Carrier solution** as indicated in the manual of the extraction kit. For the extraction procedure refer to the instruction for use manual of the extraction kit.

**Note:** when you carry out the DNA extraction from plasma with the instrument «**NucliSENS® easyMAG®**», please follow the extraction protocol **Generic 2.0.1** and follow these directions: transfer **500 µL** of sample in the 8 well strip and run the extraction. After the 10 minutes incubation, add **5 µL** of **CPE** for the internal control before adding the **NucliSENS® easyMAG® Magnetic Silica** and proceed with the extraction. Elute the nucleic acids in **100 µL** of elution buffer.

**Note:** when you carry out the DNA extraction from plasma with the instrument «**QIAsymphony® SP/AS**» and the kit «**QIAsymphony® DSP Virus / Pathogen Midi kit**» with **software version 3.5**, use the extraction protocol «**Virus Cell free 500\_V3\_DSP\_default IC**» and follow these directions: the instrument is able to use a primary tube, sample volume required for the extraction is **500 µL**, it's always requested a minimum dead volume of 100 µL. Prepare the solution containing AVE buffer and RNA carrier, according to the instruction manual of the extraction kit. Add **6 µL / sample** of **CPE** to the solution for each requested sample. Load on the instrument, in the "internal control" slot, the tubes containing the solution, as indicated in the instruction for use manual of the kit; indicate the position where eluates will be dispensed and specify the elution volume of **85 µL**. For details on the extraction procedure follow indications in the instruction for use manual of the kit.

**Interfering substances**

The DNA extracted from the sample must not contain heparin, haemoglobin, dextran, Ficoll®, ethanol or 2-propanol in order to prevent the problem of inhibition and the possibility of frequent invalid results.



High quantity of human genomic DNA in the DNA extracted from the sample may inhibit the amplification reaction.

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

#### Amplification controls

It is absolutely mandatory to validate each amplification with a negative control reaction and a positive control reaction.

For the negative control, use molecular biology grade water (not provided with this product) added to the reaction in place of the DNA extracted from the sample.

For the positive control, use the «**HSV2 - ELITe Positive Control**» product or the «**HSV2 ELITe Standard**» product.

#### Quality controls

It is recommended to validate the whole analysis procedure of each extraction and amplification session by testing Process Controls, i.e. a negative tested sample and a positive tested sample or a calibrated reference material.

### PROCEDURE

#### Setting of the real time amplification session

(To perform in the amplification / detection of amplification products area)

When **7300 Real-Time PCR System** instrument is used.

Before starting the session, referring to the instrument documentation, it is necessary to:

- switch on the real time thermal cycler, switch on the computer, run the dedicated software and open an "absolute quantification" session;
- set (Detector Manager) the "detector" for the HSV2 probe with the "reporter" = "FAM" and the "quencher" = "none" (non fluorescent) and call it "HSV2";
- set (Detector Manager) the "detector" for the Internal Control probe with the "reporter" = "VIC" (AP525 is analogous to VIC) and the "quencher" = "none" (non fluorescent) and call it "IC";
- for each well in use in the microplate, set (Well Inspector) the "detector" (type of fluorescence that is to be measured), the "passive reference" "ROX" (AP593 is used instead of ROX, normalisation of the measured fluorescence) and the type of reaction (sample, negative amplification control, positive amplification control or known quantity standard). Add this information to the **Work Sheet** enclosed at the end of this manual or print the microplate set up. The **Work Sheet** must be followed carefully during the transfer of the reaction mixture and samples into the wells.

**Note:** In order to determine the DNA titre in the starting sample, set up a series of reactions with the **Q - PCR Standards** (10<sup>5</sup> copies, 10<sup>4</sup> copies, 10<sup>3</sup> copies, 10<sup>2</sup> copies) to obtain the **Standard curve**.

See below, by way of example, how you can organise the quantitative analysis of 12 samples.

S1	S2	S3	S4	S5	S6	S7	S8	S9	S10	S11	S12
NC	10 <sup>2</sup>	10 <sup>3</sup>	10 <sup>4</sup>	10 <sup>5</sup>							

**Legend:** S1 - S12: Samples to be analysed; NC: Negative Control of amplification; 10<sup>2</sup>: 10<sup>2</sup> standard copies; 10<sup>3</sup>: 10<sup>3</sup> standard copies; 10<sup>4</sup>: 10<sup>4</sup> standard copies; 10<sup>5</sup>: 10<sup>5</sup> standard copies.

Referring to the instrument documentation, set on the dedicated software (Instrument > Thermal Cycler Protocol > Thermal Profile) the parameters of the **thermal cycle**:

- add to amplification stage the step (Add Step) of **extension at 72 °C**;

**Note:** the fluorescence acquisition (Instrument > Thermal Cycler Protocol > Settings > Data Collection) must be set during the step of hybridization at 60 °C.

- modify timing as indicated in the table "**Thermal cycle**";
- set the number cycles to **45**;
- set the volume for the software emulation of thermal transfer to reaction ("Sample volume") to **30 µL**;
- optional: add dissociation stage (Add Dissociation Stage) and set the temperature from **40 °C** to **80 °C**.

Thermal cycle		
Stage	Temperatures	Timing
Decontamination	50 °C	2 min.
Initial denaturation	94 °C	2 min.
Amplification and detection (45 cycles)	94 °C	10 sec.
	60 °C (data collection)	30 sec.
	72 °C	20 sec.
Dissociation (optional)	95 °C	15 sec.
	40 °C	30 sec.
	80 °C	15 sec.

When a **7500 Fast Dx Real-Time PCR Instrument** is used.

Before starting the session, referring to the instrument documentation, it is necessary to:

- switch on the real time thermal cycler, switch on the computer, run the dedicated software and open an "absolute quantification" session and set "Run mode: Fast 7500";
- set (Detector Manager) the "detector" for the HSV2 probe with the "reporter" = "FAM" and the "quencher" = "none" (non fluorescent) and call it "HSV2";
- set (Detector Manager) the "detector" for the internal control probe with the "reporter" = "VIC" (AP525 is similar to VIC) and the "quencher" = "none" (non fluorescent) and call it "IC";
- for each well in use in the microplate, set (Well Inspector) the "detector" (type of fluorescence that is to be measured), the "passive reference" = "CY5" (AP593 is used instead of CY5, normalisation of the measured fluorescence) and the type of reaction (sample, negative amplification control, positive amplification control or known quantity standard). Add this information to the **Work Sheet** enclosed at the end of this manual or print the microplate set up. The **Work Sheet** must be followed carefully during the transfer of the reaction mixture and samples into the wells.

**Note:** In order to determine the DNA titre in the starting sample, set up a series of reactions with the **Q - PCR Standards** (10<sup>5</sup> copies, 10<sup>4</sup> copies, 10<sup>3</sup> copies, 10<sup>2</sup> copies) to obtain the **Standard curve**.

The set up of the quantitative analysis of some samples is shown, by way of example, in the previous paragraph describing the procedure for the **7300 Real Time PCR System** instrument.

Referring to the instrument documentation, set on the dedicated software (Instrument > Thermal Cycler Protocol > Thermal Profile) the parameters of the **thermal cycle**:

- add to amplification stage the step (Add Step) of **extension at 72 °C**;

**Note:** the fluorescence acquisition (Instrument > Thermal Cycler Protocol > Settings > Data Collection) must be set during the step of hybridization at 60 °C.

- modify timing as indicated in the table "**Thermal cycle**";
- set the number cycles to **45**;
- set the volume for the software emulation of thermal transfer to reaction ("Sample volume") to **30 µL**;
- optional: add dissociation stage (Add Dissociation Stage) and set the temperature from **40 °C** to **80 °C**.

Thermal cycle		
Stage	Temperatures	Timing
Decontamination	50 °C	2 min.
Initial denaturation	94 °C	2 min.
Amplification and detection (45 cycles)	94 °C	10 sec.
	60 °C (data collection)	30 sec.
	72 °C	20 sec.
Dissociation (optional)	95 °C	15 sec.
	40 °C	1 min.
	80 °C	15 sec.
Dissociation (optional)	60 °C	15 sec.

#### Amplification set-up

(To be performed in extraction / preparation area of the amplification reaction)

Before starting the session, it is important to do the following:

- take and thaw the tubes containing the samples to be analysed. Mix gently, spin down the content for 5 seconds and keep them on ice;
- take and thaw the **HSV2 Q - PCR Mix** tubes required for the session, remembering that each tube is sufficient for preparing **25 reactions**. Mix gently, spin down the contents for 5 seconds and keep them on ice;
- take and thaw the **HSV2 - Positive Control** or the **HSV2 Q - PCR Standard** tubes. Mix them gently, centrifuge them for 5 seconds spinning down the contents and keep them on ice;
- take the **Amplification microplate** that will be used during the session, being careful to handle it with powder-free gloves and not to damage the wells.

1. Accurately pipet **20 µL** of **HSV2 Q - PCR Mix** on the bottom of the **Amplification microplate** wells, as previously established in the **Work Sheet**. Avoid creating bubbles.

**Note:** If not all the reaction mixture is used, store the remaining volume in the dark at -20°C for no longer than one month. Freeze and thaw the reaction mixture from a maximum of **5 times**.

2. Accurately pipet, by placing into the reaction mixture, **20 µL** of **DNA extract** from the first sample in the corresponding well of **Amplification microplate**, as previously established in the **Work Sheet**. Mix well the sample by pipetting the **extracted DNA** three times into the reaction mixture. Avoid creating bubbles. Proceed in the same way with the other samples of **extracted DNA**.
3. Accurately pipet, by placing into the reaction mixture, **20 µL** of **molecular biology grade water** (not provided with this product) in the well of **Amplification microplate** of the negative control of amplification, as previously established in the **Work Sheet**. Mix well the negative control by pipetting the **molecular biology grade water** three times into the reaction mixture. Avoid creating bubbles.
4. On the basis of the result required (qualitative or quantitative), one of these two options must be followed:


- When a **qualitative** result is required (detection of HSV2 DNA): accurately pipet, by placing into the reaction mixture, **20 µL** of **HSV2 - Positive Control** in the corresponding well of **Amplification microplate**, as previously established in the **Work Sheet**. Mix well the positive control by pipetting the **HSV2 - Positive Control** three times into the reaction mixture. Avoid creating bubbles.

- When a **quantitative** result is required (quantification of HSV2 DNA): accurately pipet, by placing into the reaction mixture, **20 µL** of **HSV2 Q - PCR Standard 10<sup>2</sup>** in the corresponding well of **Amplification microplate**, as previously established in the **Work Sheet**. Mix well the standard by pipetting the **HSV2 Q - PCR Standard 10<sup>2</sup>** three times into the reaction mixture. Avoid creating bubbles. Proceed in the same way with the other **HSV2 Q - PCR Standards (10<sup>3</sup>, 10<sup>4</sup>, 10<sup>5</sup>)**.

5. Accurately seal the **Amplification microplate** with the **Amplification Sealing Sheet**.
6. Transfer the **Amplification microplate** into the real time thermal cycler (in the amplification / detection of amplification products area) and start the thermal cycle for the amplification saving the session setting with an univocal and recognizable file name (e.g. "year-month-day-HSV2-EGSpA").

**Note:** At the end of the thermal cycle the **Amplification microplate** and the reaction products must be removed from the instrument and eliminated without producing environmental contaminations. In order to avoid the spilling of the reaction products, the **Amplification Sealing Sheet must not to be removed from the Amplification microplate**.

The following figure shows synthetically the preparation of the amplification reaction.



**Note:** if the preparation of the amplification is performed with the instrument «**QIA Symphony® SP/AS**», insert the microplate containing the extracts, the reagents and the amplification microplate in the dedicated slots, using the special adaptors, then follow indications in the instruction for use manual of the setup module and the steps required by the software.

**Note:** if the preparation of the amplification reaction is performed with the «**ELITe GALAXY**» instrument, load the elution microplate, the complete reaction mixture and the amplification microplate as indicated in the instrument user manual and following the steps required by the GUI.

### Qualitative analysis of the results

The recorded values of the fluorescence emitted by the specific HSV2 probe (FAM detector "HSV2") and by the specific Internal Control probe (VIC detector "IC") in the amplification reactions must be analysed by the instrument software.

Before starting the analysis, referring to the instrument documentation, it is necessary to:

- set manually (Results > Amplification plot > delta Rn vs Cycle) the calculation range for the **Baseline** (fluorescence background level) from cycle 6 to cycle 15;

**Note:** In the case of a positive sample with a high titre of HSV2 DNA, the FAM fluorescence of the HSV2 specific probe may begin to increase before the cycle 15. In this case the calculation range for the **Baseline** must be adapted from cycle 6 to the cycle in which the FAM fluorescence of the sample begins to increase, as detected by the instrument software (Results > Component).

When a **7300 Real-Time PCR System** instrument is used:

- set manually the **Threshold** for the FAM detector "HSV2" to **0.1**;
- set manually the **Threshold** for the VIC detector "IC" to **0.05**.

When a **7500 Fast Dx Real-Time PCR Instrument** is used:

- set manually the **Threshold** for the FAM detector "HSV2" to **0.2**;
- set manually the **Threshold** for the VIC detector "IC" to **0.1**.

The values of fluorescence emitted by the specific probes in the amplification reaction and the **Threshold** value of fluorescence allow to determine the **Threshold cycle (Ct)**, the cycle in which the fluorescence reached the **Threshold** value.

In the **Positive Control\*** amplification reaction, the **Ct** value of HSV2 (Results > Report) is used to validate the amplification and the detection as described in the following table:

Positive Control reaction detector FAM "HSV2"	Assay result	Amplification / Detection
<b>Ct ≤ 25</b>	<b>POSITIVE</b>	<b>CORRECT</b>

If the result of the **Positive control** amplification reaction is **Ct > 25** or **Ct Undetermined** for HSV2, the target DNA has not correctly detected. This means that problems occurred during the amplification or detection step (incorrect dispensation of the reaction mix or of the positive control, degradation of the reaction mix or of the positive control, incorrect setting of the position of the positive control, incorrect setting of the thermal cycle) which may lead to incorrect results. The session is not valid and needs to be repeated starting from the amplification step.

**\*Note:** When this product is used for the quantification of HSV2 DNA, the **Q - PCR Standard** reactions were set up instead of the **Positive Control** reaction. In this case, validate the amplification and the detection by referring to the amplification reaction of **Q - PCR Standard 10<sup>5</sup> (Ct ≤ 25)**.

In the **Negative control** amplification reaction, the **Ct** value of HSV2 (Results > Report) is used to validate the amplification and the detection as described in the following table:

Negative Control reaction detector FAM "HSV2"	Assay result	Amplification / Detection
<b>Ct Undetermined</b>	<b>NEGATIVE</b>	<b>CORRECT</b>

If the result of the **Negative control** amplification reaction is different from **Ct Undetermined (Undetermined)** for HSV2, the target DNA was detected. This means that problems have occurred during the amplification step (contamination), which may lead to incorrect results and false positives. The session is not valid and needs to be repeated starting from the amplification step.

In the amplification reaction of each **sample**, the **Ct** value of HSV2 is used to detect the target DNA while the **Ct** value of Internal Control is used to validate extraction, amplification and detection.

**Note:** Verify with the instrument software (Results > Amplification plot > delta Rn vs Cycle) that the **Ct** was determined by a fast and regular increase of the fluorescence values and not by peaks or an increase of the background (irregular or high background).

This product is able to detect a minimal quantity of about 10 copies of DNA of the gpG gene of HSV2 in the amplification reaction, corresponding to 10 genome Equivalents per reaction (limit of detection, see Performance Characteristics paragraph).

The results as **Ct** of the amplification reactions of each **sample** (Results > Report) are used as described in the following table:

Sample reaction		Sample suitability	Assay result	HSV2 DNA
detector FAM "HSV2"	detector VIC "IC"			
<b>Ct Undetermined</b>	<b>Ct &gt; 35 or Ct Undetermined</b>	<b>unsuitable</b>	<b>invalid</b>	<b>-</b>
	<b>Ct ≤ 35</b>	<b>suitable</b>	<b>valid, negative</b>	<b>NOT DETECTED</b>
<b>Ct Determined</b>	<b>Ct &gt; 35 or Ct Undetermined</b>	<b>suitable*</b>	<b>valid, positive</b>	<b>DETECTED</b>
	<b>Ct ≤ 35</b>	<b>suitable</b>	<b>valid, positive</b>	<b>DETECTED</b>

If the result of the amplification reaction of a sample is **Ct Undetermined** for HSV2 and **Ct > 35** or **Ct Undetermined** for the Internal Control, it means that it is impossible to detect efficiently the DNA for the Internal Control. In this case problems have occurred during the amplification step (inefficient or absent amplification) or during the extraction step (loss of DNA during the extraction or presence of inhibitors) which may lead to incorrect results and false negatives. The sample is not suitable, the assay, is invalid and it needs to be repeated starting from the extraction of a new sample.

If the result of the amplification reaction of a sample is **Ct Undetermined** for HSV2 and **Ct ≤ 35** for the Internal Control, it means that the HSV2 DNA is not detected in the DNA extracted from the sample; but it can not be excluded that the HSV2 DNA has a lower titre than the detection limit of the product (see Performance Characteristics). In this case the result could be a false negative.

The results obtained with this assay must be interpreted taking into consideration all the clinical data and the other laboratory tests outcomes about the patient.

**\*Note:** When in the amplification reaction of a sample the HSV2 DNA is detected, the Internal Control may result as **Ct > 35** or **Ct Undetermined**. In fact, the low efficiency amplification reaction for the Internal Control may be displaced by competition with the high efficiency amplification reaction for HSV2 DNA. In this case the sample is nevertheless suitable and the positive result of the assay is valid.

### Quantitative analysis of the results

After carrying out the procedure for qualitative analysis of the results it is possible to perform the quantitative analysis of the results of the positive samples.

In the amplification reactions of the four **Q - PCR standards**, the **Ct** values of HSV2 are used to calculate the **Standard Curve** (Results > Standard Curve) for the amplification session, to validate the amplification and the detection as described in the following table:

Standard Curve detector FAM "HSV2"	Acceptability range	Amplification / Detection
<b>Correlation coefficient (R2)</b>	<b>0.990 ≤ R2 ≤ 1.000</b>	<b>CORRECT</b>

If the **Correlation coefficient (R2)** value does not fall within the limits, this means that problems have occurred during the amplification or detection step (incorrect dispensation of the reaction mixture or of the standards, degradation of the reaction mixture or of the standards, incorrect setting of the position of the standards, incorrect setting of the thermal cycle) which may lead to incorrect results. The session is not valid and needs to be repeated starting from the amplification step.

The **Ct** values of HSV2 in the amplification reaction of each **sample** and the **Standard Curve** of the amplification session are used to calculate the **Quantity** of target DNA present in the amplification reactions of the samples.

This product is able to quantify from 1,000,000 to 10 copies of DNA of the gpG gene of HSV2 in the amplification reaction, corresponding to the genome Equivalents per reaction (linear measuring range, see Performance Characteristics), as described in the following table:

Sample result detector FAM "HSV2"	HSV2 genome Equivalents per reaction
Quantity > 1 x 10 <sup>6</sup>	MORE THAN 1,000,000
1 x 10 <sup>1</sup> ≤ Quantity ≤ 1 x 10 <sup>6</sup>	= Quantity
Quantity < 1 x 10 <sup>1</sup>	LESS THAN 10

The results (**Quantity**) of each **samples** (Results > Report) are used to calculate the genome Equivalents (**cp**) of HSV2 present in the extracted sample (**Nc**) according to this formula:

$$Nc \text{ (cp)} = \frac{Ve \times \text{Quantity}}{Vc \times Va \times Ep}$$

Where:

**Vc** is the quantity of the sample used in the extraction in rate to the required unit of measurement,  
**Ep** is the efficiency of the procedure, extraction and amplification, **expressed in decimal**;  
**Ve** is the total volume of the extraction product **expressed in µL**;  
**Va** is the volume of the extraction product used in the amplification reaction **expressed in µL**;  
**Quantity** is the result of the amplification reaction of the sample **expressed in cp per reaction**.

When «**ELITE STAR**» is used with whole blood, plasma collected in EDTA or cerebrospinal fluid samples collected in EDTA and the result **expressed in cp / mL** is required, the formula becomes:

Simplified formula for whole blood, plasma, cerebrospinal fluid and «ELITE STAR»
Nc (cp / mL) = 28 x Quantity

When «**ELITE GALAXY**» is used with whole blood, plasma collected in EDTA or cerebrospinal fluid samples collected in EDTA and the result **expressed in cp / mL** is required, the formula becomes:

Simplified formula for whole blood, plasma, cerebrospinal fluid and «ELITE GALAXY»
Nc (cp / mL) = 35 x Quantity

When «**EXTRAblood**» extraction kit is used with whole blood samples collected in EDTA and the result **expressed in cp / mL** is required, the formula becomes:

Simplified formula for whole blood and «EXTRAblood»
Nc (cp / mL) = 25 x Quantity

When «**NucliSENS® easyMAG®**» extraction system is used with plasma collected in EDTA or cerebrospinal fluid samples and the result **expressed in cp / mL** is required, the formula becomes:

Simplified formula for plasma, cerebrospinal fluid and «NucliSENS® easyMAG®»
Nc (cp / mL) = 10 x Quantity

When «**QIAasympy® SP/AS**» extraction system is used with plasma samples collected in EDTA and the result **expressed in cp / mL** is required, the formula becomes:

Simplified formula for plasma and «QIAasympy® SP/AS»
Nc (cp / mL) = 12 x Quantity

#### Calculation of the linear measuring range limits

The linear measuring range limits as cp / mL of the sample, when a particular extraction method is used, may be calculated from the linear measuring range of the amplification reaction according to the following formula:

$$\text{Lower limit (cp / mL)} = \frac{Ve \times 10 \text{ cp}}{Vc \times Va \times Ep}$$

$$\text{Upper limit (cp / mL)} = \frac{Ve \times 1,000,000 \text{ cp}}{Vc \times Va \times Ep}$$

When «**EXTRAblood**» extraction kit is used with whole blood samples collected in EDTA, the formula becomes:

Measuring range limits (cp / mL) with «EXTRAblood»
Lower limit (cp / mL) = 25 x 10 cp
Upper limit (cp / mL) = 25 x 1,000,000 cp
from 250 to 25,000,000 cp / mL

When «**ELITE STAR**» extraction system is used with whole blood, plasma collected in EDTA or cerebrospinal fluid samples, the formula becomes:

Measuring range limits (cp / mL) with «ELITE STAR»
Lower limit (cp / mL) = 28 x 10 cp
Upper limit (cp / mL) = 28 x 1,000,000 cp
from 280 to 28,000,000 cp / mL

When «**ELITE GALAXY**» extraction system is used with whole blood, plasma collected in EDTA or cerebrospinal fluid samples, the formula becomes:

Measuring range limits (cp / mL) with «ELITE GALAXY»
Lower limit (cp / mL) = 35 x 10 cp
Upper limit (cp / mL) = 35 x 1,000,000 cp
from 350 to 35,000,000 cp / mL

When «**NucliSENS® easyMAG®**» extraction system is used with plasma samples collected in EDTA or cerebrospinal fluid samples, the formula becomes:

Measuring range limits (cp / mL) with «NucliSENS® easyMAG®»
Lower limit (cp / mL) = 10 x 10 cp
Upper limit (cp / mL) = 10 x 1,000,000 cp
from 100 to 10,000,000 cp / mL

When «**QIAasympy® SP/AS**» extraction system is used with plasma samples collected in EDTA, the formula becomes:

Measuring range limits (cp / mL) with «QIAasympy® SP/AS»
Lower limit (cp / mL) = 12 x 10 cp
Upper limit (cp / mL) = 12 x 1,000,000 cp
from 120 to 12,000,000 cp / mL



**PERFORMANCE CHARACTERISTICS**

**Analytical sensitivity: limit of detection**

The analytical sensitivity of this assay allows detecting the presence of about 10 target DNA molecules in 20 µL of DNA added to the amplification reaction.

The analytical sensitivity of this assay, as detection limit, was tested using plasmidic DNA containing the amplification product whose initial concentration was measured by spectrophotometer. The plasmidic DNA was diluted to a titre of 10 copies / 20 µL in human genomic DNA at a titre of 500 ng / 20 µL. This sample was used in 50 replicates carrying out the amplification by ELITechGroup S.p.A. products.

The final results are summed up in the following table:

Samples	No.	positive	negative
10 copies plasmidic DNA + 500 ng of human genomic DNA	50	50	0

The analytical sensitivity of this assay used in association to whole blood samples and «**ELITe GALAXY**» was verified with a panel of HSV2 dilutions within the limiting concentration. The panel was prepared by diluting the HSV08-12 sample of the "QCMD 2008 Herpes Simplex Virus EQA Panel" (Qnostics, Ltd, UK) in HSV2 DNA - negative EDTA whole blood. The viral concentrations ranged from 10 cp / mL to 560 cp / mL. Each sample of the panel was tested in 12 replicates carrying out the whole analysis procedure, extraction and PCR Setup with «**ELITe GALAXY**» and amplification with ELITechGroup S.p.A. products. The statistical analysis was performed by the Probit regression. The limit of detection was calculated for the concentrations at which the probability of a positive result is the 95%.

The analytical sensitivity as cp/mL is reported below

Limit of Detection for whole blood samples and «ELITe GALAXY» (cp / mL)			
		95% confidence range	
		lower limit	upper limit
<b>95% positivity</b>	<b>171 cp / mL</b>	112 cp / mL	505 cp / mL

The analytical sensitivity of this assay used in association to plasma samples and «**ELITe GALAXY**» was verified with a panel of HSV2 dilutions within the limiting concentration. The panel was prepared by diluting the HSV08-12 sample of the "QCMD 2008 Herpes Simplex Virus EQA Panel" (Qnostics, Ltd, UK) in HSV2 DNA - negative EDTA plasma. The viral concentrations ranged from 10 cp / mL to 560 cp / mL. Each sample of the panel was tested in 12 replicates carrying out the whole analysis procedure, extraction and PCR Setup with «**ELITe GALAXY**» and amplification with ELITechGroup S.p.A. products.

The statistical analysis was performed by the Probit regression. The limit of detection was calculated for the concentrations at which the probability of a positive result is the 95%.

The analytical sensitivity as cp/mL is reported below

Limit of Detection for plasma samples and «ELITe GALAXY» (cp / mL)			
		95% confidence range	
		lower limit	upper limit
<b>95% positivity</b>	<b>119 cp / mL</b>	75 cp / mL	547 cp / mL

**Analytical sensitivity: linear measuring range**

The analytical sensitivity of this assay allows the quantification from 1,000,000 to 10 molecules of target DNA in the 20 µL of DNA added to the amplification reaction.

The analytical sensitivity of this assay, as linear measuring range, was determined using a panel of dilutions (1 log<sub>10</sub> between one dilution and the next) of a plasmidic DNA containing the amplification product whose initial concentration was measured by a spectrophotometer. The dilutions from 10<sup>7</sup> molecules per reaction to 10<sup>1</sup> molecules per reaction were tested in 9 replicates carrying out the amplification by the ELITechGroup S.p.A. products.

The analysis of the obtained data, performed by linear regression, demonstrated that the assay displays a linear response for all the panel points (linear correlation coefficient greater than 0.99).

The upper limit of the linear measuring range was set at 10<sup>6</sup> molecules per reaction corresponding to genome Equivalent per reaction, within one logarithm from the highest concentration Q - PCR Standard amplification standard (10<sup>5</sup> molecules / 20 µL).

The lower limit of the linear measuring range was set at 10 molecules per reaction corresponding to genome Equivalent per reaction, within one logarithm from the lowest concentration Q - PCR Standard amplification standard (10<sup>2</sup> molecules / 20 µL).

The final results are summed up in the following table:

Linear measuring range (cp / reaction)	
Upper limit	1,000,000 cp / reaction
Lower limit	10 cp / reaction

The linear measuring range limits as cp / mL referred to the used extraction kit are calculated at page 26.

**Analytical sensitivity: Precision and Accuracy**

The precision of the assay, as the variability of results obtained with several replicates of the same sample tested within the same session, allowed to obtain a mean percentage Coefficient of Variation (% CV) of about 18.6% of measured quantities, within the range from 10<sup>6</sup> molecules to 10 molecules in the 20 µL of DNA added to the amplification reaction.

The accuracy of the assay, as the difference between the mean of results obtained with several replicates of a sample within the same session and the theoretical concentration of the sample, allowed to obtain a mean percentage Inaccuracy (% Inacc.) of about 14.4% of measured quantities, within the range from 10<sup>6</sup> molecules to 10 molecules in the 20 µL of DNA added to the amplification reaction.

The precision and the accuracy were determined using data obtained for the study of the linear measuring range.

**Analytical sensitivity: reproducibility with calibrated reference material**

The analytical sensitivity of the assay, as reproducibility of results compared with results obtained using other assays in different laboratories, was checked by testing a proficiency panel.

The tests were carried out using as calibrated reference material a panel of dilutions of HSV2 within the concentration limit (QCMD 2009 Herpes Simplex Virus DNA EQA Panel, Qnostics Ltd, UK). Each sample was tested carrying out the whole analysis procedure: extraction with «**ELITe STAR**» and amplification with ELITechGroup S.p.A. products.

The results are reported in the following table.

Tests with calibrated reference materials and «ELITe STAR»				
Sample	Virus	Commercial assay consensus Log <sub>10</sub> virus conc.	Standard Deviation	Mean results Log <sub>10</sub> cp / mL
HSV09-01	HSV1	2.215	0.789	-
HSV09-02	HSV2	2.236	0.938	1.687
HSV09-03	HSV2	3.293	0.915	3.700
HSV09-04	HSV2	2.314	0.898	2.329
HSV09-05	Negative	NA	-	-
HSV09-06	HSV1	2.402	0.665	-
HSV09-07	HSV1	4.189	0.599	-
HSV09-08	HSV2	2.389	0.827	2.352
HSV09-09	Negative	NA	-	-
HSV09-10	HSV1	3.205	0.767	-

All samples were correctly detected. The quantitative results are within the range defined by the Consensus ± 1 Standard Deviation.

Further tests were carried out using as calibrated reference material a panel of dilutions of HSV2 within the concentration limit (QCMD 2012 Herpes Simplex Virus DNA EQA Panel, Qnostics Ltd, UK). Each sample was tested carrying out the whole analysis procedure: extraction with «**ELITe GALAXY**» and amplification with ELITechGroup S.p.A. products.

The results are reported in the following table in cp/mL.

Tests with calibrated reference materials and «ELITe GALAXY»				
Sample	Commercial assay consensus Log <sub>10</sub> virus conc.	Standard Deviation	Positive / Replicates	Mean results Log <sub>10</sub> cp / mL
HSV12-01	Negative, NA	-	0/2	-
HSV12-02	HSV1, 3.910	0.582	0/2	-
HSV12-03	HSV2, 1.948	0.305	1/2	2.128
HSV12-04	HSV1, 3.680	0.547	0/2	-
HSV12-05	HSV2, 1.352	0.629	0/2	-
HSV12-06	HSV1, 2.318	0.441	0/2	-
HSV12-07	HSV1, 2.014	0.296	0/2	-
HSV12-08	HSV2, 3.424	1.098	2/2	3.634
HSV12-09	Negative, NA	-	0/2	-
HSV12-10	HSV2, 3.417	1.042	2/2	3.861

All negative samples were correctly reported. The positive samples within the theoretical limit of detection of the system (350 copies/mL) were correctly detected within the range of the mean "Consensus" value of the commercial assay  $\pm$  1 Standard Deviation. One sample below the theoretical limit of detection of the system (22.5 copies/mL) was reported negative. Samples with titer below the limit of detection can be stochastically reported as positive or negative.

#### Diagnostic sensitivity: detection and quantification efficiency on different genotypes / subtypes

The diagnostic sensitivity of the assay, as detection and quantification efficiency on different genotypes / subtypes, was evaluated by comparison of sequences with nucleotide databases.

The analysis of the regions chosen for the hybridisation of the primers and of the fluorescent probe in the alignment of the sequences available in the database for the gpG gene of HSV2 showed conservation and absence of significant mutations.

#### Diagnostic sensitivity: confirmation of positive samples

The diagnostic sensitivity of the assay, as confirmation of positive clinical samples, was tested using some clinical samples whole blood collected in EDTA, tested positives for HSV2 DNA.

The diagnostic sensitivity was evaluated using as reference material 20 samples of whole blood collected in EDTA from presumably negative donors (Biological Sample Library Europe S.A.S., Lion, France), tested positives at low titre for HSV2 DNA with QCMD 2008 Herpes simplex virus EQA Panel (Qnostics Ltd, UK). Each sample of Whole blood was used to carry out the whole analysis, extraction with «EXTRAAblood» and amplification by ELITechGroup S.p.A. products.

The results are summed in the following table.

Samples	N	positive	negative
Whole blood collected in EDTA positive for HSV2 DNA	20	20	0

All spiked samples were correctly detected as positive for HSV2 DNA.  
The diagnostic sensitivity of the assay in this test was equal to 100%.

The diagnostic sensitivity was evaluated using 22 samples of cerebrospinal fluid negative for HSV2 DNA, that were spiked for HSV2 DNA adding IQC032RH sample, from ELITechGroup S.p.A., 30 samples of plasma negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV08-12 sample, from QCMD 2008 Herpes Simplex Virus EQA Panel (Qnostics Ltd, UK) and 30 whole blood samples negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV08-12, from QCMD 2008 Herpes Simplex Virus EQA Panel (Qnostics Ltd, UK). Each sample was used to carry out the whole analysis procedure: extraction with «ELITe STAR» and amplification with ELITechGroup S.p.A. products.

The results are summed up in the following table.

Samples	N	positive	negative
Cerebrospinal fluid spiked for HSV2 DNA	22	22	0
Plasma collected in EDTA spiked for HSV2 DNA	30	30	0
Whole blood collected in EDTA spiked for HSV2 DNA	30	29	0

One HSV2 positive sample resulted invalid.  
The diagnostic sensitivity of the assay in this test was equal to 100%.

The diagnostic sensitivity was evaluated using 21 samples of cerebrospinal fluid negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV2MQP01 - High sample, from Molecular Q Panel HSV2MQP01 (Qnostics Ltd, UK), 30 samples of plasma negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV08-12 sample, from QCMD 2008 Herpes Simplex Virus DNA EQA Panel, (Qnostics Ltd, UK) and 32 whole blood samples negative for HSV2 DNA, that were spiked for HSV2 DNA adding HSV08-12 sample, from QCMD 2008 Herpes Simplex Virus DNA EQA Panel, (Qnostics Ltd, UK) and HSV10-01 sample, from QCMD 2010 Herpes Simplex Human virus DNA EQA Panel, (Qnostics Ltd, UK). Each sample was used to carry out the whole analysis procedure: extraction with «ELITe GALAXY» and amplification with ELITechGroup S.p.A. products. The results are summed up in the following table.

Samples	N	positive	negative
Cerebrospinal fluid spiked for HSV2 DNA	21	21	0
Plasma collected in EDTA spiked for HSV2 DNA	30	30	0
Whole blood collected in EDTA spiked for HSV2 DNA	32	32	0

All spiked samples were correctly detected as positive for HSV2 DNA.  
The diagnostic sensitivity of the assay in this test was equal to 100 %.

#### Analytical specificity: absence of cross-reactivity potential interfering markers

The analytical specificity of the assay, as absence of cross-reactivity with other potential interference markers, was evaluated by comparison of sequences with nucleotide databases.

The analysis of the alignment of the sequences of the primers and of the fluorescent probe with the sequences available in databases for organisms other than HSV2, including the HSV1 and the VZV complete genomes, the human herpes virus that is most similar to HSV2, showed their specificity and the absence of significant homology.

The analytical specificity of the assay, as absence of cross-reactivity with other potential interference markers, was checked by testing a proficiency panel.

The analytical specificity was evaluated by using as calibrated reference material a panel including HSV1 positive and VZV positive samples (QCMD 2007 Herpes simplex virus EQA Panel, Qnostics Ltd, UK). Each sample was tested in duplicates carrying out the whole analysis procedure, extraction and amplification by ELITechGroup S.p.A. products.

The results are reported in the paragraph "Analytical sensitivity: reproducibility with calibrated reference material".

No cross-reactivity was detected with HSV1 and VZV positive samples.

#### Diagnostic specificity: confirmation of negative samples

The diagnostic specificity of the assay, as confirmation of negative samples, was tested using some HSV2 DNA negative clinical samples of whole blood collected in EDTA, tested negatives for HSV2 DNA.

The diagnostic specificity was evaluated using as reference material 24 samples of whole blood collected in EDTA from presumably negative donors for HSV2 DNA (Biological Sample Library Europe S.A.S., Lion, France). Each sample of Whole blood was used to carry out the whole analysis, extraction with «EXTRAAblood» and amplification by ELITechGroup S.p.A. products. The results are summed in the following table.

Samples	N	positive	negative
Whole blood collected in EDTA negative for HSV2 DNA	24	0	24

All samples were detected as negative for HSV2 DNA.

The diagnostic specificity of the assay in this test was equal to 100%.

The diagnostic specificity was evaluated using 24 cerebrospinal fluid samples that were negative for HSV2 DNA, 30 plasma samples collected in EDTA that were negative for HSV2 DNA and 30 whole blood samples collected in EDTA that were negative for HSV2 DNA (tested with a real time amplification CE IVD product). Each sample was used to carry out the whole analysis procedure: extraction «ELITe STAR» and amplification with ELITechGroup S.p.A. products.

The results are summed up in the following table.

Samples	N	positive	negative
Cerebrospinal fluid negative for HSV2 DNA	24	0	24
Plasma collected in EDTA negative for HSV2 DNA	30	0	30
Whole blood collected in EDTA negative for HSV2 DNA	30	0	30

All samples were detected as negative for HSV2 DNA.  
The diagnostic specificity of the assay in this test was equal to 100%.

The diagnostic specificity was evaluated using 22 cerebrospinal fluid samples that were negative for HSV2 DNA, 34 plasma samples collected in EDTA that were presumably negative for HSV2 DNA and 36 whole blood samples collected in EDTA that were presumably negative for HSV2 DNA. Each sample was used to carry out the whole analysis procedure: extraction with «**ELITe GALAXY**» and amplification with ELITechGroup S.p.A. products. The results are summed up in the following table.

Samples	N	positive	negative
Cerebrospinal fluid negative for HSV2 DNA	22	0	22
Plasma collected in EDTA presumably negative for HSV2 DNA	34	0	34
Whole blood collected in EDTA presumably negative for HSV2 DNA	36	0	36

All samples were correctly detected as negative for HSV2 DNA.  
The diagnostic specificity of the assay in this test was equal to 100 %.

**Note:** The complete data and results of the tests carried out to evaluate the product performance characteristics with matrices and instruments are recorded in Section 7 of the Product Technical File "HSV2 ELITe MGB Kit", FTP RTS032PLD.

## Roche cobas z 480 analyzer

## SAMPLES AND CONTROLS

### Samples

This product must be used with **DNA extracted** from the following clinical samples:

#### Whole blood collected in EDTA

The whole blood samples for DNA extraction must be collected in EDTA and identified according to laboratory guidelines, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of three days, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when you carry out DNA extraction from whole blood samples with the "**MagNA Pure 24 System**" instrument with **software version 1.0** (or equivalent later versions), use the "**Pathogen200**" extraction protocol and follow these instructions: dispense **350 µL** of sample into the MagNA Pure Tube 2.0 mL, load the tube into the instrument and begin the extraction. This protocol processes 200 µL of sample, adds **CPE** 20 µL / extraction and elutes the nucleic acids into 100 µL. The **CPE** must be diluted 1:2 in ultra-pure molecular biology grade water. For details of the extraction procedure, follow the instructions contained in the kit's User Manual carefully.

#### Plasma collected in EDTA

The plasma samples for nucleic acid extraction must be collected in EDTA according to laboratory guidelines, transported at +2° / +8°C and stored at +2° / +8°C for a maximum of three days, otherwise they must be frozen and stored at -20°C for a maximum of thirty days or at -70°C for longer periods. It is recommended to split the samples to be frozen into aliquots in order to prevent repeated cycles of freezing and thawing. When using frozen samples, thaw the samples just immediately before the extraction in order to avoid possible nucleic acid degradation.

**Note:** when you carry out DNA extraction from plasma samples with the "**MagNA Pure 24 System**" instrument with **software version 1.0** (or equivalent later versions), use the "**Pathogen200**" extraction protocol and follow these instructions: dispense **350 µL** of sample into the MagNA Pure Tube 2.0 mL, load

the tube into the instrument and begin extraction. This protocol processes 200 µL of sample, adds **CPE** 20 µL / extraction and elutes the nucleic acids into 100 µL. The **CPE** must be diluted 1:2 in ultra-pure molecular biology grade water. For details of the extraction procedure, follow the instructions contained in the kit's User Manual carefully.

### Interfering substances

The DNA extracted from the sample must not contain heparin, haemoglobin, dextran, Ficoll®, ethanol or 2-propanol in order to prevent inhibition problems and the possibility of frequent invalid results.

High quantity of human genomic DNA in the DNA extracted from the sample may inhibit the amplification reaction.

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

### Amplification controls

It is absolutely mandatory to validate each amplification session with a negative control reaction and a positive control reaction.

For the negative control, add ultra-pure molecular biology grade water (not included in the kit) to the reaction instead of the DNA extracted from the sample.

For the positive control, use the «**HSV2 - ELITe Positive Control**» or alternatively «**HSV2 - ELITe Positive Control RF**» product, or the «**HSV2 ELITe Standard**» product.

### Quality controls

It is recommended to validate the whole analysis procedure of each extraction and amplification session by testing Process Controls, i.e. a negative tested sample and a positive tested sample or a calibrated reference material.

## PROCEDURE

### Setting of the real time amplification session

(To perform in the amplification / detection of amplification products area)

Before starting the session, referring to the instrument documentation, it is necessary to:

- switch on the control computer and the real time thermal cycler. Open the dedicated software and in the main window, open an "New Experiment" session;
- set the reaction volume ("Reaction volume") to 40 µL;
- assign an identifier to each sample ("Sample editor");
- define the reaction's Thermal Cycle according to the following table:

Thermal Cycle		
Stage	Temperatures	Periods
Decontamination	50°C	2 mins.
Initial denaturation	94°C	2 mins.
Amplification and detection (45 cycles)	94°C	10 sec.
	60°C (fluorescence acquisition)	30 sec.
	72°C	20 sec.
Dissociation (optional)	95 °C	15 sec.
	40 °C	30 sec.
	80 °C	15 sec.

**Note:** fluorescence acquisition occurs individually, set Ramp Rate (°C/sec) to 4.4 °C/sec.


- select the signal detection channels: "detector" for the HSV2 sensor with "channel FAM 465-510" and "detector" for the IC internal control sensor with "channel VIC 540-580";

Fill in the **Work Plan** attached at the end of this User Manual, transcribing this information or printing the microplate's layout. This **Work Plan** must be followed carefully when transferring the reaction mixture and samples into the wells.

**Note:** to determine the concentration of DNA in the source sample, you must perform a series of reactions

with **Q - PCR Standard** ( $10^5$  copies,  $10^4$  copies,  $10^3$  copies,  $10^2$  copies) to obtain the **Standard Curve**.

See below, by way of example, how you can organise the quantitative analysis of 12 samples.



**Legend:** C1 - C12: Samples to be analyzed; NC: Negative amplification control;  $10^2$ : Standard  $10^2$  copies;  $10^3$ : Standard  $10^3$  copies;  $10^4$ : Standard  $10^4$  copies;  $10^5$ : Standard  $10^5$  copies.

#### Amplification set-up

(To be performed in the extraction / preparation of the amplification reaction area)

Before starting the session, it is necessary to:

- retrieve and thaw the test tubes containing the samples to be analyzed. Shake the tubes gently and then place them in the centrifuge for 5 seconds to send the contents to the bottom and then keep them on ice;
- retrieve and thaw the test tubes containing **HSV2 Q - PCR Mix** required for the session, remembering that the contents of each tube is enough to perform **25 reactions**. Shake the tubes gently and then place them in the centrifuge for 5 seconds to send the contents to the bottom and then keep them on ice;
- retrieve and thaw the test tubes containing **HSV2 - Positive Control** or alternatively **HSV2 - ELITe Positive Control RF** or the test tubes containing **HSV2 Q - PCR Standard**. Shake the tubes gently and then place them in the centrifuge for 5 seconds to send the contents to the bottom and then keep them on ice;
- retrieve the **AD-plate** to be used in the session, making sure you handle it wearing dust-free gloves and do not damage the wells.

- Without creating any bubbles and depositing it precisely on the bottom, transfer **20  $\mu$ L** of reaction mixture **HSV2 Q - PCR Mix** into the wells on the **AD-plate** as previously established in the **Work Plan**.

**Note:** If not using all the reaction mixture, store any remaining mixture at  $-20^\circ\text{C}$  for a maximum of one month. Freeze and thaw the reaction mixture a maximum of **5 TIMES**.

- Depositing it precisely into the reaction mixture, transfer **20  $\mu$ L** of **extracted DNA** from the first sample in the corresponding well on the **AD-plate** as previously established in the **Work Plan**. Mix the sample well by pipetting the **extracted DNA** three times into the reaction mixture. Be sure not to create any bubbles. Proceed in the same manner with all the other **extracted DNA**.

- Depositing it precisely into the reaction mixture, transfer **20  $\mu$ L** of **ultra-pure molecular biology grade water** (not supplied with the product) into the well on the **AD-plate** containing the negative amplification control as previously established in the **Work Plan**. Mix the negative control well by pipetting the **ultra-pure molecular biology grade water** three times into the reaction mixture. Be sure not to create any bubbles.

- On the basis of the result required (qualitative or quantitative), one of these two options must be followed:


- When a **qualitative** result is required (detection of HSV2 DNA): accurately pipet, by placing into the reaction mixture, **20  $\mu$ L** of **HSV2 - Positive Control** or alternatively «**HSV2 - ELITe Positive Control RF**» in the corresponding well of **Amplification microplate**, as previously established in the **Work Sheet**. Mix well the positive control by pipetting the **HSV2 - Positive Control** three times into the reaction mixture. Avoid creating bubbles.
- When a **quantitative** result is required (quantification of HSV2 DNA): accurately pipet, by placing into

the reaction mixture, **20  $\mu$ L** of **HSV2 Q - PCR Standard  $10^2$**  in the corresponding well of **Amplification microplate**, as previously established in the **Work Sheet**. Mix well the standard by pipetting the **HSV2 Q - PCR Standard  $10^2$**  three times into the reaction mixture. Avoid creating bubbles. Proceed in the same way with the other **HSV2 Q - PCR Standards ( $10^3$ ,  $10^4$ ,  $10^5$ )**.

- Carefully seal the **AD-plate** using the **Sealing Film**.
- Transfer the **AD-plate** into the real-time Thermal Cycler in the amplification/detection of amplification products area and start the amplification thermal cycle, saving the session settings under a unique and recognizable identifier (e.g. "year-month-day-HSV2-EGSpA").

**Note:** At the end of the thermal cycle, the **AD-plate** and reaction products must be removed from the instrument and disposed of in a way that does not cause environmental pollution. **Never remove the Sealing Film from the Amplification microplate** to avoid any leakage of the reaction products.

The following figure shows synthetically the preparation of the amplification reaction..



#### Qualitative results analysis

The emitted fluorescence values recorded by the HSV2 detector and Internal Control (IC) detector during the amplification reactions must be analyzed by the instrument's software.

Select the menu "Analysis" and choose "Absolute Quant/Fit Points" (2 points)

Select the group of samples to be analyzed

In accordance with the instrument's documentation, before starting the analysis you must:

- manually enter the calculation range (Background button) for the **Background Fluorescence Level** from cycle 2 to cycle 6.

For **Plasma** samples

- manually set the **Threshold** and **Noiseband** for the FAM "HSV2" detector to **0.55**;
- manually set the **Threshold** and **Noiseband** for the VIC "IC" detector to **1.2**

For **Whole blood** samples

- manually set the **Threshold** and **Noiseband** for the FAM "HSV2" detector to **0.80**;
- manually set the **Threshold** and **Noiseband** for the VIC "IC" detector to **1.5**



The fluorescence values emitted by the specific detectors in the amplification reaction and the **Threshold** and **Noiseband** fluorescence values are used to determine the **Threshold Cycle (Ct)**, i.e. the cycle in which the fluorescence **Threshold** is reached.

The **Ct** values for HSV2 in the amplification reactions of the four **Q - PCR Standard** are used to calculate the **Standard Curve** (Results > Standard Curve) of that amplification session and to validate the amplification and detection as shown in the following table:

Reaction Q - PCR Standard 10 <sup>5</sup> "HSV2" detector	Assay result	Amplification / Detection
Ct ≤ 25	POSITIVE	CORRECT

If the result of the **Positive control** amplification reaction is **Ct > 25** or **Ct Undetermined** for HSV2, the target DNA was not correctly detected. This means that problems occurred during the amplification or detection step (incorrect dispensation of the reaction mix or of the positive control, degradation of the reaction mix or of the positive control, incorrect setting of the position of the positive control, incorrect setting of the thermal cycle) which may lead to incorrect results. The session is not valid and needs to be repeated starting from the amplification step.

\* **Note:** When this product is used for the quantification of HSV2 DNA, the **Q - PCR Standard** reactions were set up instead of the **Positive Control** reaction. In this case, validate the amplification and the detection by referring to the amplification reaction of **Q - PCR Standard 10<sup>5</sup> (Ct ≤ 25)**.

During the **Negative Control** amplification reaction, the value of **Ct** for HSV2 (Analysis window) is used to validate amplification and detection as shown in the following table:

Negative Control Reaction "HSV2" detector	Assay result	Amplification / Detection
Ct Undetermined	NEGATIVE	CORRECT

If the result of the **Negative Control** amplification reaction is other than **Ct Undetermined** for HSV2, the presence of the DNA target has been detected. Problems occurred during the amplification stage (contamination) which can lead to incorrect results and false positives. The session is invalid and must be repeated from the amplification stage.

During the amplification reactions for each **sample**, the value of **Ct** for HSV2 is used to detect the presence of the DNA target, whilst the value of **Ct** for the Internal Control is used to validate the extraction, amplification and detection.

**Note:** Check using the instrument's software (Analysis window) that the **Ct** is determined by a rapid and regular increase in fluorescence values and not by peaks or an increase of the background signal (irregular or noisy background).

Results like **Ct** from each **sample's** amplification reactions (Analysis window) are used as shown in the following table:

Sample reaction		Sample suitability	Assay result	HSV2 DNA
"HSV2" detector	"IC" detector			
Ct Undetermined	Ct > 35 or Ct Undetermined	not suitable	invalid	-
	Ct ≤ 35	suitable	valid, negative	NOT DETECTED
Ct Determined	Ct > 35 or Ct Undetermined	suitable	valid, positive	DETECTED
	Ct ≤ 35	suitable	valid, positive	DETECTED

If the result of a sample's amplification reaction is **Ct Undetermined** for HSV2 and **Ct > 35** or **Ct Undetermined** for the Internal Control, it was not possible to detect the Internal Control DNA efficiently. In this case, problems occurred during the amplification stage (inefficient or null amplification) or in the extraction stage (degraded sample DNA, sample with insufficient number of cells, loss of DNA during extraction or presence of inhibitors in the extracted DNA) which can cause incorrect results and false

negatives. The sample is not suitable, the assay is not valid and must be repeated starting from the extraction of a new sample.

If the result of a sample's amplification reaction is **Ct Undetermined** for HSV2 and **Ct ≤ 35** for the Internal Control, the HSV2 DNA was not detected in the DNA extracted from the sample but it cannot be excluded that the HSV2 DNA is present at a concentration lower than the product's limit of detection (see Performance Characteristics). In this case, the result would constitute a false negative.

The results obtained with this assay must be interpreted by considering all the clinical data and the results of other laboratory tests connected to the patient.

**Note:** When HSV2 DNA is detected during the amplification reaction of a sample, amplification of the Internal Control can produce a result of **Ct > 35** or **Ct Undetermined**. In fact, the low-efficiency Internal Control amplification reaction can be eliminated from the competition with the high-efficiency HSV2 reaction. In this case, the sample is then suitable and the positive assay result is valid.

#### Quantitative results analysis

After having performed the qualitative analysis procedure, you can carry out the quantitative analysis of the results relating to the positive sample.

If the result of the amplification reaction for the **Q - PCR Standard 10<sup>5</sup>** is **Ct > 25** or **Ct Undetermined** or if the **Ct** values of the four **Q - PCR standards** don't fit regularly the standard curve the DNA target was not correctly detected. Problems occurred during the amplification or detection stage (incorrect dispensing of the reaction mixture or standards, degradation of the reaction mixture or standards, incorrect setting of the standard positions, incorrect setting of the thermal cycle) which can cause incorrect results. The session is invalid and must be repeated from the amplification stage.

The **Ct** values for HSV2 in the amplification reactions of each **sample** and the **Standard Curve** (**Standard Curve** button) from the amplification session are used to calculate the **Quantity** of DNA target present in the amplification reactions relating to the samples.

This product is able to quantify from 1,000,000 down to around 10 copies per reaction, from 25,000,000 to 250 copies per mL for whole blood and plasma using the **MagNA Pure 24** extraction system (see Performance Characteristics), as shown in the following table:

Sample result FAM "HSV2" detector	HSV2 copies per reaction
Quantity > 1 x 10 <sup>6</sup>	GREATER THAN 1,000,000
1.0 x 10 <sup>1</sup> ≤ Quantity ≤ 1 x 10 <sup>6</sup>	= Quantity
Quantity < 1.0 x 10 <sup>1</sup>	LESS THAN 10

The results (**Quantity**) relating to each **sample** (Analysis window) are used to calculate the **copies** of HSV2 present in the source sample (**Nc**) according to this formula:

$$Nc = \frac{V_e \times \text{Quantity}}{V_c \times V_a \times E_p}$$

Where:

**Vc** is the quantity of sample used in the extraction in relation to the required unit of measure;  
**Ep** is the efficiency of the procedure, extraction and amplification, **expressed in decimals**,  
**Ve** is the total volume obtained from the extraction **expressed in µL**;  
**Va** is the volume of extraction product used in the amplification reaction **expressed in µL**;  
**Quantity** is the result of the amplification reaction relating to the sample **expressed in copies per reaction**.

When using samples of whole blood collected in EDTA or plasma collected in EDTA and the **MagNA Pure 24** extraction system and the result is to be **expressed in copies / mL**, the formula becomes:

<b>Simplified formula for whole blood and plasma and MagNA Pure 24</b>
<b>Nc (copies / mL) = 25 x Quantity</b>

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity: limit of detection

The analytical sensitivity of this assay, as limit of detection, allows for the detection of around 10 copies in 20 µL of DNA added to the amplification reaction.

The analytical sensitivity of this assay, as limit of detection, has been tested using a plasmid DNA containing the amplification product whose initial concentration was measured using a spectrophotometer. The plasmid DNA was diluted to a concentration of 10 copies / 20 µL in 150,000 copies of pBETAGLOBIN / 20 µL. This sample was used in 27 replicates to carry out amplification using ELITechGroup S.p.A. products. The final results are summarized in the following table.

Samples	N	positives	negatives
10 copies of plasmid DNA + 150,000 copies of Beta-globin	27	27	0

### Analytical sensitivity: linear measuring range

The analytical sensitivity of this assay, as linear measuring range, allows for the quantification from around 25,000,000 to 25 copies in 20 µL of DNA added to the amplification reaction.

The analytical sensitivity of this assay was evaluated using a panel of dilutions (1 Log<sub>10</sub> between one dilution and the next) of plasmid DNA containing the amplification product, whose initial concentration was measured using a spectrophotometer. The points of the panel from 10<sup>7</sup> molecules per reaction to 10<sup>1</sup> molecules per reaction were used in 9 replicates to carry out amplification using ELITechGroup S.p.A. products. Analysis of the obtained data, performed using linear regression, showed that the assay has a linear response for all panel points (linear correlation coefficient greater than 0.99).

The lower limit of the linear measuring range was set at around 10 copies / reaction within one logarithm from the lowest concentration of Q - PCR Standard amplification standard (10<sup>2</sup> copies / 20 µL).

The upper limit of the linear measuring range was set at 10<sup>6</sup> copies / reaction within one logarithm from the highest concentration of Q - PCR Standard amplification standard (10<sup>5</sup> copies / 20 µL).

The results are shown in the following table.

Linear measuring range using MagNA Pure 24		
	Lower limit	Upper limit
copies / mL	250	25,000,000
copies / reaction	10	1,000,000

Conversions from copies / mL to copies / reaction and vice versa were calculated as shown on page 39.

### Analytical sensitivity: Precision and Accuracy

The precision of this assay, in terms of the variability of the results obtained in the same amplification session using different replicates of a sample, allowed to obtain a mean Variation Coefficient percentage (CV%) of the values of Ct lower than 1% in the range from 10<sup>6</sup> molecules to 10<sup>1</sup> molecules in 20 µL of DNA added to the amplification reaction.

The precision of this assay, in terms of the variability of the results obtained in the same amplification session using different replicates of a sample, allowed to obtain a mean Variation Coefficient percentage (CV%) of the measured quantities of around 9% in the range from 10<sup>6</sup> molecules to 10<sup>1</sup> molecules in 20 µL of DNA added to the amplification reaction.

The accuracy of this assay, in terms of the difference between the mean of the results obtained in the same amplification session using different replicates of a sample and the sample's theoretical concentration value, allowed to obtain a mean Inaccuracy percentage of the measured logarithmic quantity of around 6.4% in the range from 10<sup>6</sup> molecules to 10<sup>1</sup> molecules in 20 µL of DNA added to the amplification reaction.

Precision and accuracy were determined using the data obtained during the experiments assessing the linear measuring range.

### Analytical sensitivity: reproducibility with certified reference material

The analytical sensitivity of the assay, as reproducibility of value of a calibrated reference material, was evaluated using as reference material the «HSV2 Molecular 'Q' Panel» (Qnostics Ltd, UK). Each sample of the panel was tested in 2 replicates carrying out the whole analysis procedure: extraction using the **MagNA Pure 24** automatic extraction system and amplification using ELITechGroup S.p.A. products.

The results are reported in the following table.

Tests with calibrated reference materials and «MagNA Pure 24»	
Sample	Positives / Replicates
HSV2MQP01-High	2/2
HSV2MQP01-Medium	2/2
HSV2MQP01-Low	2/2
HSV2MQP01-Negative	0/2

All samples were correctly detected.

The analytical sensitivity of the assay, as reproducibility of value of a calibrated reference material, was evaluated using as reference material the QCMD 2017 Herpes Simplex Virus DNA Panel (HSV DNA17S, Qnostics Ltd, UK) a panel of HSV2 dilutions. Each sample of the panel was tested in 2 replicates carrying out the whole analysis procedure: extraction using the **MagNA Pure 24** automatic extraction system and amplification using ELITechGroup S.p.A. products.

The results are reported in the following table.

Tests with calibrated reference materials and «MagNA Pure 24»		
Sample	Sample Status	Positive / Replicates
HSV DNA17S -01	HSV1 positive, HSV2 negative	0/2
HSV DNA17S -02	HSV1 positive, HSV2 negative	0/2
HSV DNA17S -03	HSV1 positive, HSV2 negative	0/2
HSV DNA17S -04	HSV1 positive, HSV2 negative	0/2
HSV DNA17S -05	HSV1 positive, HSV2 negative	0/2
HSV DNA17S -06	HSV2 Frequently detected	2/2
HSV DNA17S -07	HSV2 Frequently detected	2/2
HSV DNA17S -08	HSV2 Frequently detected	2/2
HSV DNA17S -09	HSV2 Detected	2/2
HSV DNA17S -10	HSV2 Frequently detected	2/2

All samples were correctly detected.

### Diagnostic sensitivity: confirmation of positive samples

The diagnostic sensitivity was evaluated using as reference material 30 samples of whole blood collected in EDTA, negative for HSV2 DNA, which were spiked for HSV2 DNA adding HSV2MQP01-High sample (Qnostics Ltd, UK), and 30 samples of plasma collected in EDTA, negative for HSV2 DNA, which were spiked for HSV2 DNA adding HSV2MQP01-High sample (Qnostics Ltd, UK).

Each sample was used carrying out the whole analysis procedure: extraction using the **MagNA Pure 24** automatic extraction system and amplification using ELITechGroup S.p.A. products. The results are summarized in the following table.

Samples	N	positives	negatives
Whole blood collected in EDTA spiked for HSV2 DNA	30	30	0
Plasma collected in EDTA spiked for HSV2 DNA	30	30	0

All samples were valid at first test and confirmed positive for HSV2 DNA.

The total diagnostic sensitivity of the assay was 100%.

#### Diagnostic specificity: confirmation of negative samples

The diagnostic specificity was evaluated using as reference material 40 samples of whole blood collected in EDTA presumably negative for HSV2 DNA and 34 samples of plasma collected in EDTA presumably negative for HSV2 DNA.

Each sample was used carrying out the whole analysis procedure: extraction using the **MagNA Pure 24** automatic extraction system and amplification using ELITechGroup S.p.A. products. The results are summarized in the following table.

Samples	N	positives	negatives
Whole blood collected in EDTA presumably negative for HSV2 DNA	40	0	40
Plasma collected in EDTA presumably negative for HSV2 DNA	34	0	34

All samples were valid at first test and confirmed negative for HSV2 DNA.

The total diagnostic specificity of the assay was 100%.

**NOTE:** The complete data and results of the tests carried out to evaluate the product performance characteristics with matrices and instruments are recorded in Section 7 of the Product Technical File "HSV2 ELITe MGB Kit", FTP RTS032PLD.

#### REFERENCES

- E. Aurelius et al. (1993) *J. Med. Virology* **39**: 179 - 186  
E. A. Lukhtanov et al. (2007) *Nucleic Acids Res.* **35**: e30

#### PROCEDURE LIMITATIONS

Use this product only with DNA extracted from the following clinical samples: cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

Do not use DNA extracted from heparinized samples with this product: heparin inhibits the amplification reaction of nucleic acids and causes invalid results.

Do not use extracted DNA that is contaminated with haemoglobin, dextran, Ficoll®, ethanol or 2-propanol with this product: these substances inhibit the amplification reaction of nucleic acids and may cause invalid results.

Do not use with this product extracted DNA containing high quantity of human genomic DNA that may inhibit the amplification reaction of nucleic acids.

There are no data available concerning product performance with DNA extracted from the following clinical samples: suspensions of leucocytes, suspensions of granulocytes and amniotic fluid.

Use this product only with the validated instruments and associated clinical samples indicated in the section "Samples and Controls".

There are no data available concerning inhibition caused by antiviral, antibiotic, chemotherapeutic or immunosuppressant drugs.

The results obtained with this product depend on an adequate identification, collection, transport storage and processing of the samples. To avoid incorrect results, it is therefore necessary to take care during these steps and to carefully follow the instructions for use provided with the products for nucleic acid extraction.

Owing to its high analytical sensitivity, the real time amplification method used in this product is sensitive to cross-contaminations from the HSV2 positive samples, the positive controls and the same amplification products. Cross-contaminations cause false positive results. The product format is able to limit cross-contaminations. However, the cross-contaminations can be avoided only by good laboratory practices and following these instructions for use.

This product must be handled by qualified personnel trained in the processing of potentially infective biological samples and chemical preparations classified as dangerous to prevent accidents with potentially serious consequences for the user and other persons.

This product requires the use of work clothes and areas that are suitable for the processing of potentially infective biological samples and chemical preparations classified as dangerous to prevent accidents with potentially serious consequences for the user and other persons.

This product must be handled by qualified personnel trained in molecular biology techniques, such as extraction, amplification and detection of nucleic acids, to avoid incorrect results.

It is necessary to have separate areas for the extraction/preparation of amplification reactions and for the amplification / detection of amplification products to prevent false positive results.

This product requires the use of special clothing and instruments for extraction/preparation of amplification reactions and for amplification / detection of amplification products to avoid false positive results.

Due to inherent differences between technologies, it is recommended that users perform method correlation studies to estimate technology differences prior to switching to a new technology.

A negative result obtained with this product means that the HSV2 DNA is not detected in the DNA extracted from the sample; but it can not be excluded that the HSV2 DNA has a lower titre than the product detection limit (see Performance Characteristics). In this case the result could be a false negative.

Results obtained with this product may sometimes be invalid due to failed internal control and require retesting, starting from extraction, that can lead to a delay in obtaining final results.

Possible polymorphisms within the region of the viral genome covered by the product primers and probes may impair detection and quantification of HSV2 DNA.

As with any other diagnostic medical device, the results obtained with this product must be interpreted taking into consideration all the clinical data and other laboratory tests done on the patient.

As with any other diagnostic medical device, there is a residual risk of invalid, false positive and false negative results obtained with this product. This residual risk can not be eliminated or further reduced. In some cases, as the prenatal or emergency diagnosis, this residual risk could contribute to wrong decisions with potentially dangerous effects for the patient.

#### TROUBLESHOOTING

##### Target DNA not detected in the Positive Control or Q - PCR Standard reactions or invalid correlation coefficient of the Standard curve

Possible Causes	Solutions
Incorrect dispensing into the microplate wells.	Take care when dispensing reactions into the microplate wells and comply with the work sheet. Check the volumes of reaction mixture dispensed. Check the volumes of positive control or standard dispensed.
Incorrect session setup on ELITe InGenius.	Check the position of reaction mixture, positive control or standards. Check the volumes of reaction mixture, positive control or standards.
Probe degradation.	Use a new aliquot of reaction mixture.
Positive control or standard degradation.	Use a new aliquot of positive control or standard.
Instrument setting error.	Check the position settings for the positive control or standard reactions on the instrument. Check the thermal cycle settings on the instrument.
Instrument error.	Contact ELITechGroup Technical Service.

##### Target DNA detected in the Negative control reaction

Possible Causes	Solutions
Incorrect dispensing into the microplate wells.	Avoid spilling the contents of the sample test tube. Always change tips between one sample and another. Take care when dispensing samples, negative controls, positive controls and standards into the microplate wells and comply with the work sheet.
Incorrect session setup on ELITe InGenius	Check the position of reaction mixture or negative control. Check the volumes of reaction mixture or negative control.
Error while setting the instrument	Check the position settings of the samples, negative controls, positive controls and standards on the instrument
Microplate badly sealed.	Take care when sealing the microplate.
Contamination of the molecular biology grade water.	Use a new aliquot of water.
Contamination of the reaction mixture.	Use a new aliquot of reaction mixture.

Contamination of the extraction / preparation area for amplification reactions.	Clean surfaces and instruments with aqueous detergents, wash lab coats, replace test tubes and tips in use.
Instrument error.	Contact ELITechGroup Technical Service.








Target and Internal Control DNA not detected in the sample reactions	
Possible Causes	Solutions
Incorrect dispensing into the microplate wells.	Avoid spilling the contents of the sample test tube. Always change tips between one sample and another. Take care when dispensing samples into the microplate wells and comply with the work sheet.
Incorrect session setup on ELITe InGenius.	Check the position of reaction mixture or samples. Check the volumes of reaction mixture or samples.
Internal Control degradation.	Use new aliquots of Internal Control.
Inhibition due to sample interfering substances.	Repeat the amplification with a 1:2 dilution in molecular biology grade water of eluted sample in a "PCR only" session. Repeat the extraction and amplification of sample.
Incorrect reagent storage.	Verify that reaction mix was not exposed to room temperature for more than 30 minutes.
Problems during extraction.	Verify quality and concentration of extracted DNA.
Instrument error.	Contact ELITechGroup Technical Service.

Irregular or high levels of background fluorescence in the reactions	
Possible causes	Solutions
Incorrect dispensing of sample.	Take care, by pipetting three times, when mixing samples, negative controls and positive controls or standards into the reaction mixture. Avoid creating bubbles.
Baseline setting error.	Set the baseline calculation range within cycles where the background fluorescence has already stabilized (check the "Results", "Component" data) and the signal fluorescence has not yet started to increase, e.g. from cycle 6 to cycle 15. Use the automatic baseline calculation by setting the "Auto Baseline" option.

Anomalous dissociation curve	
Possible causes	Solutions
Absence of a defined peak. Defined peak but different from that of the other samples and of the standards or positive control.	Check for detector FAM Ct lower than 30. High quantity of amplification product at the end of the reaction may interfere with the melting curve analysis. Repeat the sample amplification to confirm the presence of target DNA with a possible mutation. The target DNA of the sample should be sequenced to confirm mutation.

With ELITe InGenius: Error 30103	
Possible Causes	Solutions
Too high concentration of target in the sample.	If significant amplification is observed in PCR plot: - repeat the amplification of eluted sample in molecular biology grade water, in a "PCR only" session or - repeat the extraction with a dilution of the primary sample in molecular biology grade water, in a "Extract + PCR" session.

## SYMBOLS

<b>REF</b>	Catalogue Number.
	Upper limit of temperature.
<b>LOT</b>	Batch code.
	Use by (last day of month).
<b>IVD</b>	<i>in vitro</i> diagnostic medical device.
	Fulfilling the requirements of the European Directive 98/79/EC for <i>in vitro</i> diagnostic medical device.
	Contains sufficient for "N" tests.
	Attention, consult instructions for use.
<b>CONT</b>	Contents.
	Keep away from sunlight.
	Manufacturer.



**NOTICE TO PURCHASER: LIMITED LICENSE**

This product contains reagents manufactured by Life Technologies Corporation and are sold under licensing arrangements between ELITechGroup S.p.A. and its Affiliates and Life Technologies Corporation. The purchase price of this product includes limited, nontransferable rights to use only this amount of the product solely for activities of the purchaser which are directly related to human diagnostics. For information on purchasing a license to this product for purposes other than those stated above, contact Licensing Department, Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008. Phone: +1(760)603-7200. Fax: +1(760)602-6500. Email: outlicensing@thermofisher.com.

ELITe® MGB detection reagents are covered by one or more of U.S. Patent numbers 6,127,121, 6,485,906, 6,660,845, 6,699,975, 6,727,356, 6,790,945, 6,949,367, 6,972,328, 7,045,610, 7,319,022, 7,368,549, 7,381,818, 7,662,942, 7,671,218, 7,715,989, 7,723,038, 7,759,126, 7,767,834, 7,897,736, 8,008,522, 8,067,177, 8,163,910, 8,389,745, 8,969,003, 8,980,855, 9,056,887, 9,085,800, 9,169,256 and EP patent numbers 1068358, 1144429, 1232157, 1261616, 1430147, 1781675, 1789587, 1975256, 2714939 as well as applications that are currently pending.

This limited license permits the person or legal entity to which this product has been provided to use the product, and the data generated by use of the product, only for human diagnostics. Neither ELITechGroup S.p.A. nor its licensors grants any other licenses, expressed or implied for any other purposes.

ELITe MGB®, the ELITe MGB® logo, ELITe InGenius® and ELITe BeGenius® are registered as trademarks within the European Union.

«NucISSENS® easyMAG®» are registered trademarks of bioMérieux SA.

«QIASymphony®» is a registered trademark of QIAGEN GmbH.

Ficoll® is a registered trademark of GE Healthcare.

# HSV2 ELITE MGB® kit used with Genius series platforms

Ref: RTS032PLD



This document is a simplified version of the official instruction for use. Please refer to the complete document before use: [www.elitechgroup.com](http://www.elitechgroup.com)  
This document is available only in English.

## A. Intended use

The «**HSV2 ELITE MGB® Kit**» product is part of a qualitative and quantitative nucleic acids amplification assay for the **detection and quantification of the DNA of type 2 Herpes Simplex human virus (HSV2)** in DNA samples extracted from cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

The product is intended for use in the diagnosis and monitoring of HSV2 infections alongside clinical data of the patient and other laboratory tests outcomes. The assay is CE-IVD validated in combination with the instruments **ELITE InGenius®** and **ELITE BeGenius®**.

## B. Amplified sequence


Target	Gene	Fluorophore
HSV2	Glicoprotein G (gpG)	FAM
Internal Control	Human beta globin gene	AP525

## C. Validated matrix

- › **Whole Blood** EDTA, **Plasma** EDTA, **CSF**

## D. Kit content

**HSV2 Q-PCR Mix**



**X 4**

Ready-to-use PCR Master Mix  
4 tubes of 540 µL  
96 reactions per kit  
5 freeze-thaw cycles per tube

- › Maximum shelf-life: 24 months
- › Storage temperature: - 20°C

## E. Material required not provided in the kit

- |                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>› <b>ELITE InGenius®</b> instrument: INT030</li> <li>› <b>ELITE BeGenius®</b> instrument: INT040</li> <li>› <b>ELITE InGenius SP200</b> extraction cartridge: INT032SP200</li> <li>› <b>ELITE InGenius PCR Cassette</b> amplification cartridge: INT035PCR</li> <li>› <b>ELITE InGenius SP200 Consumable Set</b> consumable for extraction: INT032CS</li> </ul> | <ul style="list-style-type: none"> <li>› <b>HSV2 - ELITE Standard:</b> STD032PLD</li> <li>› <b>HSV2 - ELITE Positive Control:</b> CTR032PLD</li> <li>› <b>CPE – Internal Control:</b> CTRCPE</li> <li>› <b>ELITE InGenius Waste Box:</b> F2102-000</li> <li>› <b>300 µL Filter Tips Axygen:</b> TF-350-L-R-S</li> <li>› <b>1000 µL Filter Tips Tecan :</b> 30180118</li> </ul> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## F. Protocol

- |                                                                                                                                                                                                                                        |                                                                                                                                                                                |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>› Sample volume 200 µL</li> <li>› CPE Internal Control volume 10 µL</li> <li>› Total eluate volume 100 µL</li> <li>› PCR eluate input volume 20 µL</li> <li>› Q-PCR Mix volume 20 µL</li> </ul> | <ul style="list-style-type: none"> <li>› Unit of quantitative result Copies/mL</li> <li>› Frequency of controls 15 days</li> <li>› Frequency of calibration 60 days</li> </ul> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## G. Performance ELITE InGenius® and ELITE BeGenius®

Matrix	Limit of Detection	Linearity Range	Diagnostic Sensitivity	Diagnostic Specificity
Whole Blood	171 cp / mL	171 – 25,000,000	96% 22/23*	100% 34/34*
Plasma	119 cp /mL	119 – 25,000,000	100% 30/30*	100% 39/39*
CSF	119 cp /mL	119 – 25,000,000	100% 20/20*	100% 22/22*

\*confirmed samples/ tested samples










## H. Procedures ELITE InGenius®

The user is guided step-by-step by the ELITE InGenius software to prepare the run. All the steps: extraction, amplification and result interpretation are automatically performed. Three operational mode are available: complete run, or extraction only, or PCR only.

### Before analysis

1. Switch on ELITE InGenius Identification with username and password Select the mode "Closed"	2. Verify calibrators: HSV2 Q-PCR Standard in the "Calibration menu" Verify controls: HSV2 positive and negative controls in the "Control menu" <i>N.B:</i> Both have been run, approved and not expired	3. Thaw the Q-PCR-Mix and the Internal Control tubes Vortex gently Spin down 5 sec
---------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------

### Procedure 1 - Complete run: Extraction + PCR

1. Select "Perform Run" on the touch screen 	2. Verify the extraction volume: Input: "200 µL", eluate: "100 µL" 	3. Scan the sample barcodes with hand-held barcode reader or type the sample ID 
4. Select the "Assay protocol" of interest 	5. Select the sample position: Primary tube or extraction tube 	6. Load the Q-PCR Mix and the Internal Control in the inventory block 
7. Load: PCR cassette, Extraction cartridge, Elution tube, Tip, extraction tube and primary sample racks 	8. Close the door Start the run 	9. View, approve and store the results 

### Procedure 2 - PCR only

1 to 4 : Follow the Complete Run procedure described above	5. Select the protocol "PCR only" and set the sample position "Extra tube"	6. Load the extracted nucleic acid tubes in the Elution tubes rack
7. Load the PCR cassette rack Load the Q-PCR Mix in the inventory block	8. Close the door Start the run	9. View, approve and store the results

## Procedure 3 - Extraction only

<b>1 to 4 :</b> Follow the Complete Run procedure described above	<b>5.</b> Select the protocol “Extraction Only” and set the sample position : Primary tube or Secondary tube	<b>6.</b> Load the Internal Control in the inventory block
<b>7.</b> Load: Extraction cartridge, Elution tube, Tip cassette, extraction tube and primary sample racks	<b>8.</b> Close the door Start the run	<b>9.</b> Archive the eluate sample










## L. Procedures ELITE BeGenius®

The user is guided step-by-step by the ELITE BeGenius software to prepare the run. All the steps: extraction, amplification and result interpretation are automatically performed. Three operational mode are available: complete run, or extraction only, or PCR only.

### Before analysis

<b>1.</b> Switch on ELITE BeGenius Identification with username and password Select the mode “Closed”	<b>2.</b> Verify calibrators: HSV2 Q-PCR standard in the “Calibration menu” Verify controls: HSV2 pos. and neg. controls in the “Control menu” <i>NB:</i> Both have been run, approved and not expired	<b>3.</b> Thaw the HSV2 Q- PCR-Mix and the CPE Internal Control tubes Vortex gently Spin down 5 sec
----------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------

### Procedure 1 - Complete run: Extraction + PCR

<b>1.</b> Select “Perform Run” on the touch screen and then click on the run mode «Extraction and PCR» 	<b>2.</b> Insert the Sample Rack with the barcoded samples in the cooling area. The barcode scan is already active 	<b>3.</b> Verify the extraction volumes: Input: “200 µL”, Eluate: “100 µL” 
<b>4.</b> Select the “Assay protocol” of interest 	<b>5.</b> Print the labels to barcode the empty elution tubes. Load the tubes in the Elution Rack and insert it in the cooling area 	<b>6.</b> Load the Q-PCR-Mix and the CPE Internal Control in Reagent Rack and insert it in the cooling area 
<b>7.</b> Load: Filter Tips, Extraction rack, and PCR rack 	<b>8.</b> Close the door. Start the run 	<b>9.</b> View, approve and store the results 

**Note:** if a second extraction is performed repeat steps from 2 to 4



### Procedure 2 - PCR only

<b>1.</b> Select “Perform Run” on the touch screen and the click on the run mode «PCR Only»	<b>2.</b> Load the extracted nucleic acid barcoded tubes in the Elution Rack and insert it in the cooling area	<b>3.</b> Select the “Assay protocol” of interest
<b>4.</b> Load the Q-PCR-Mix in Reagent Rack and insert it in the cooling area Load filter tips and the PCR rack	<b>5.</b> Close the door. Start the run	<b>6.</b> View, approve and store the results

### Procedure 3 - Extraction only

<b>1 to 4 :</b> Follow the Complete Run procedure described above	<b>5.</b> Select the protocol “Extraction Only” in the Assay Protocol selection screen.	<b>6.</b> Load the CPE Internal Control in the Elution Rack and insert it in the cooling area
<b>7.</b> Load : Filter Tips and the Extraction Rack	<b>8.</b> Close the door Start the run	<b>9.</b> Archive the eluate sample

# HSV2 ELITE MGB® Kit used with ABI PCR instrument

Ref: RTS032PLD



This document is a simplified version of the official instruction for use. Please refer to the complete document before use: [www.elitechgroup.com](http://www.elitechgroup.com)  
This document is available only in English.

## A. Intended use

The «**HSV2 ELITE MGB® Kit**» product is part of a qualitative and quantitative nucleic acids amplification assay for the **detection and quantification of the DNA of type 2 Herpes Simplex human virus (HSV2)** in DNA samples extracted from cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

The product is intended for use in the diagnosis and monitoring of HSV2 infections alongside clinical data of the patient and other laboratory tests outcomes. The assay is CE-IVD validated in combination with **ABI PCR thermal cyclers** (Thermo-Fisher) and the following extraction systems: **ELITE STAR** (ELITechGroup), **ELITE GALAXY** (ELITechGroup), **easyMAG** (BioMérieux) or **QIAasymphony** (Qiagen).

## B. Amplified sequence

Target	Gene	Fluorophore
<b>HSV2</b>	Glicoprotein G (gpG)	FAM
<b>Internal Control</b>	Human beta globin gene	AP525

## C. Validated matrix

› Whole blood EDTA

› Plasma EDTA

› Cerebrospinal fluid

## D. Kit content

### HSV2 Q-PCR Mix



X 4

Ready-to-use PCR Master Mix  
4 tubes of 540 µL  
100 reactions per kit  
5 freeze-thaw cycles per tube

- › Maximum shelf-life: 24 months
- › Storage Temperature: - 20°C

## E. Material required not provided in the kit

- › **7500 Fast Dx and 7300 PCR Instrument**
- › **ELITE STAR**: INT010
- › **ELITE STAR 200 extraction kit**: INT011EX
- › **ELITE GALAXY**: INT020
- › **ELITE GALAXY 300 extraction kit**: INT021EX
- ›

- › **HSV2 ELITE Standard**: STD032PLD
- › **HSV2 - ELITE Positive Control**: CTR032PLD
- › **CPE - Internal Control**: CTCRPE
- › **easyMAG** - Generic protocol 2.0.1
- › **QIAasymphony** - DNA Mini kit or DSP Virus/Pathogen Midi kit
- › **Molecular biology grade water**

## F. Performance

System	Matrix	Limit of Detection	Diagnostic Sensitivity	Diagnostic Specificity
ELITE STAR - ABI	Whole blood	<b>10 gEq/reaction</b>	<b>100%</b> (30/30)*	<b>100%</b> (30/30)*
	Plasma	<b>10 gEq/reaction</b>	<b>100%</b> (30/30)*	<b>100%</b> (30/30)
	CSF	<b>10 gEq/reaction</b>	<b>100%</b> (22/22)*	<b>100%</b> (24/24)*
ELITE GALAXY - ABI	Whole blood	<b>171 gEq/mL</b>	<b>100%</b> (32/32)*	<b>100%</b> (36/36)*
	Plasma	<b>119 gEq/mL</b>	<b>100%</b> (30/30)*	<b>100%</b> (34/34)*
	CSF	<b>10 gEq/reaction</b>	<b>100%</b> (21/21)*	<b>100%</b> (22/22)*

\*confirmed samples/tested samples

## G. Procedure

The procedure below summarized the main steps of the sample analysis with conventional PCR workflow: validated extraction systems, PCR instrument settings, PCR set-up and result interpretation.

### Extraction - Validated systems

Extraction	Validated matrix	Sample volume processed	Min. sample volume	Total eluate volume	CPE Internal Control volume
ELITE Star	WB, Plasma, CSF	200 µL	700 µL	100 µL	200 µL
ELITE Galaxy	WB, Plasma	300 µL	400 µL	200 µL	10 µL
EasyMAG	CSF, Plasma	500 µL	-	100 µL	5 µL
QIAasympphony	Plasma	500 µL	700 µL	85 µL	10 µL

### Amplification - Settings of 7500 Fast Dx and 7300 PCR instruments


1. Switch on the thermal-cycler
2. Set "HSV2" detector with "FAM" and quencher "none"
3. Set "Internal Control" detector with "VIC" and quencher "none"
4. Set passive fluorescence as "Cy5" with 7500 Fast Dx and as "ROX" with 7300 instrument
5. Set up the thermal profile as indicated. Fluorescence acquisition must be set during hybridization step at 60°C

Stage	Temperature	Timing
Decontamination	50°C	2 min
Denaturation	94°C	2 min
Amplification and detection	94°C	10 sec
	60°C	30 sec
45 cycles	72°C	20 sec

*The melt curve analysis is optional, refer to the complete IFU*

### Amplification - PCR Set -up

1. Thaw HSV2 Q-PCR Mix and Q-PCR standard tubes
2. Mix gently and spin-down
3. Pipet **20 µL** of Q-PCR-Mix in all microplate wells in use
4. Add, **20 µL** of extracted DNA in sample wells, **20 µL** of molecular grade water in Negative Control well, and **20 µL** of the 4 Q-PCR standards in standard curve wells. Each one has to be mixed by pipetting 3 times into the reaction mixture
5. Seal the microplate with the amplification sealing sheet
6. Transfer the microplate in the thermocycler and start



### Amplification – Baseline and Threshold for qualitative analysis

Instrument	Baseline	HSV2 FAM	Internal Control VIC
7500 Fast Dx Real Time PCR	6 - 15	0.2	0.1
7300 Real Time PCR	6 - 15	0.1	0.05

### Interpretation - Qualitative results

HSV2 Ct value	Internal Control Ct value	Interpretation
Determined	–	Positive
Undetermined	Ct ≤ 35	Negative
	Ct >35 or Undetermined	Invalid*

*\*Repeat the assay starting from the extraction*

### Interpretation - Quantitative results

The HSV2 Ct value obtained for each sample and the standard curve generated are used to calculate the quantity of target DNA in the reaction.

The sample quantification ranges from approximately 10 to 10<sup>6</sup> cp/reaction or approximately from 100 to 10<sup>7</sup> cp/mL.

# HSV2 ELITE MGB® Kit used with Cobas-Z 480 PCR instruments

Ref.: RTS032PLD



This document is a simplified version of the official instruction for use. Please refer to the complete document before use: [www.elitechgroup.com](http://www.elitechgroup.com)  
This document is available only in English.

## A. Intended use

The «**HSV2 ELITE MGB® Kit**» product is part of a qualitative and quantitative nucleic acids amplification assay for **the detection and quantification of the DNA of type 2 Herpes Simplex human virus (HSV2)** in DNA samples extracted from cerebrospinal fluid (CSF), whole blood collected in EDTA, plasma collected in EDTA.

The product is intended for use in the diagnosis and monitoring of HSV2 infections alongside clinical data of the patient and other laboratory tests outcomes. The assay is CE-IVD validated in combination with **Cobas – Z 480 analyzer (Roche)** and the following extraction systems: **MagNA Pure 24 System**.

## B. Amplified sequence

Target	Gene	Fluorophore
HSV2	Glicoprotein G (gpG)	FAM
Internal Control	human beta globin gene	AP525

## C. Validated matrix

- › Whole blood EDTA
- › Plasma EDTA

## D. Kit content

### HSV2 Q-PCR Mix



X 4

Ready-to-use PCR Master Mix  
4 tubes of 540 µL  
100 reactions per kit  
5 freeze-thaw cycles per tube

- › Maximum shelf-life: 24 months
- › Storage Temperature: - 20°C

## E. Material required not provided in the kit

- › Cobas – Z 480 analyzer PCR Instrument
- › MagNA Pure 24 System, software 1.0
- › HSV2 - ELITE Positive Control: CTR032PLD
- › HSV2 ELITE Standard: STD032PLD
- › CPE Internal Control: CTRCPE

## F. Performance

System	Matrix	Limit of Detection	Diagnostic Sensitivity	Diagnostic Specificity
MagNA Pure 24	Whole blood	10 cp/reaction	100% (30/30)*	100% (40/40)*
	Plasma	10 cp/reaction	100% (30/30)*	100% (34/34)*

\*confirmed samples/tested samples



## G. Procedure

The procedure below summarized the main steps of the sample analysis with conventional PCR workflow: validated extraction systems, PCR instrument settings, PCR set-up and result interpretation.

### Extraction - Validated systems

Extraction	Validated matrix	Sample volume processed	Min. sample volume	Total eluate volume	CPE Internal Control volume
MagNA Pure 24	WB, Plasma	200 µL	350 µL	100 µL	20 µL

### Amplification - Settings of Cobas-Z 480 PCR instruments


1. Switch on the thermal-cycler
2. Set "HSV2" detector with "FAM" and quencher "465 - 510"
3. Set "Internal Control" detector with "VIC" and quencher "540 -580"

Stage	Temperature	Timing
Decontamination	50°C	2 min
Denaturation	94°C	2 min
Amplification and detection	94°C	10 sec
	60°C	30 sec
45 cycles	72°C	20 sec

*The melt curve analysis is optional, refer to the complete IFU*

### Amplification - PCR Set-up

1. Thaw HSV2 Q-PCR Mix and Q-PCR standard tubes
2. Mix gently and spin-down
3. Pipet **20 µL** of Q-PCR Mix in all microplate wells in use
4. Add, **20 µL** of extracted DNA in sample wells, **20 µL** of molecular grade water in Negative Control well, and **20 µL** of the 4 Q-PCR standards in standard curve wells  
Each one has to be mixed by pipetting 3 times into the reaction mixture
5. Seal the microplate with the amplification sealing sheet
6. Transfer the microplate in the thermocycler and start



### Amplification – Background and Threshold for qualitative analysis\*

Instrument	Matrix	Background	HSV2 FAM	Internal Control VIC
Cobas-Z 480 PCR instruments	Plasma	2 - 6	0.55	1.2
Cobas-Z 480 PCR instruments	WB	2 - 6	0.8	1.5

*\*manually set the Threshold and Noiseband*

### Interpretation - Qualitative results

HSV-1 Ct value	Internal Control Ct value	Interpretation
Determined	-	Positive
Undetermined	Ct ≤ 35	Negative
	Ct >35 or Undetermined	Invalid*

*\*Repeat the assay starting from the extraction*

### Interpretation - Quantitative results

The HSV2 Ct value obtained for each sample and the standard curve generated are used to calculate the quantity of target DNA in the reaction. The sample quantification ranges from approximately 10 to 10<sup>6</sup> copies/reaction or approximately from 100 to 10<sup>7</sup> copies/mL.

# WORK SHEET

[illegible]