



**Analyte Specific Reagent**  
**Analytical and performance characteristics are not established**

## **MGB Alert<sup>®</sup> HSV 1 & 2 Probe Mix ASR**

This package insert applies to PN 800426

**SPIN TUBES  
PRIOR TO  
OPENING**

 ELITechGroup<sup>®</sup> MDx LLC  
21720 23<sup>rd</sup> Drive SE, Suite 150  
Bothell, WA 98021 USA | Phone: 1.800.453.2725

### **Product Description**

MGB Alert<sup>®</sup> HSV 1 & 2 Probe Mix is an analyte specific reagent that contains a chemically synthesized DNA oligonucleotide designed to detect nucleic acid sequences from HSV 1 and 2 (Herpes Simplex Virus 1 and 2). The mix consists of one probe at 0.2 nmol per vial in a 1X TE buffer solution. The probe dye is labeled with 6-carboxyfluorescein and has an excitation wavelength maxima of 492 nm and an emission wavelength maxima of 516 nm.

### **Storage and Storage Stability**

The reagent contained in this package should be protected from light and stored at  $\leq -10^{\circ}\text{C}$ . The reagent must be used by the expiration date printed on the vial label. Product use is limited to five freeze-thaws.

#### **Technical Assistance Center**

For Technical Assistance, please call toll free 1.800.453.2725 or  
E-mail: MDx@elitechgroup.com

### **Limited Product Warranty**

ELITechGroup MDx LLC (EGI MDx) warrants that this product will meet the specifications stated above. If any component of this product does not conform to these specifications, EGI MDx will at its sole discretion, as its sole and exclusive liability and as the users' sole and exclusive remedy, replace the product at no charge or refund the cost of the product; provided that notice of non-conformance is given to EGI MDx, within sixty (60) days of receipt of the product.

This warranty limits EGI MDx' liability to the replacement of this product or refund of the cost of the product. NO OTHER WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT, ARE PROVIDED BY EGI MDx. EGI MDx shall have no liability for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use this product and its components.

In no event shall EGI MDx be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or use or the failure of EGI MDx products to perform in accordance with the stated specifications.

Some components of nucleic acid analysis, such as specific methods and compositions for manipulating or visualizing nucleic acids for analysis, may be covered by one or more patents owned by other parties. Similarly, nucleic acids containing specific nucleotides sequences may be patented. Making, using, offering for sale, or selling such components or nucleic acids may require one or more licenses. Nothing in this document should be construed as an authorization or implicit license to make, use or sell any so covered component or nucleic acid under any such patents.

### **Licensing and legal notice**

MGB Alert detection reagents are covered by one or more of U.S. Patents Nos. 6,127,121, 6,485,906, 6,660,845, 6,699,975, 6,727,356, 6,790,945, 6,949,367, 6,972,328, 7,045,610, 7,319,022, 7,368,549, 7,381,818, 7,662,942, 7,671,218, 7,715,989, 7,723,038, 7,759,126, 7,767,834, 7,897,736, 8,008,522, 8,067,177, 8,163,910, 8,389,745, 8,969,003, 8,980,855, 9,056,887, 9,085,800 and 9,169,256 as well as applications that are currently pending.

This limited license permits the person or legal entity to which this product has been provided to use the product, and the data generated by use of the product, only: (a) for internal ASR assay validation purposes involving hybridization-based analysis of nucleic acids; and (b) for clinical diagnostics as an Analyte Specific Reagent, as that term is defined and restricted by the U.S. Food and Drug Administration in 21 CFR 864.4020 and 21 CFR 809.30. Neither EGI MDx nor any of its licensors grants any other licenses, whether express or implied, for any other purposes. No license under the PCR process is conveyed through purchase of this product. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department, Roche Diagnostics, Grenzacherstrasse 124, CH-4070, Basel, Switzerland.

#### **Trademarks:**

MGB Alert<sup>®</sup> is a registered trademark of ELITechGroup, Inc.