ELI.H.A RF

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ELI.H.A RF

Quantitative detection of Rheumatoid Factor in serum samples by indirect haemagglutination

220 tests
(Ref. 04115)

8000180-en-2011-11

1 - AIM
ELI.H.A RF enables the quantitative determination of IgM Rheumatoid Factor in serum samples (4, 5) by indirect haemagglutination (10, 11). Each kit allows 220 tests to be carried out or 20 reactions of 11 dilutions.

2 - INTRODUCTION
The Rheumatoid Factor is a high diagnostic interest biological marker for adult rheumatoid arthritis (5).

3 - PRINCIPLE
ELI.H.A RF is based on the Rheumatoid Factor specific haemagglutination properties, used in reactions of the Waaler-Rose type (1, 2, 8). The reagent **R1** is composed of sheep red blood cells sensitized by a gamma-globulins fraction of rabbit anti-sheep red blood cells serum. The presence of Rheumatoid Factor in serum samples is revealed by a haemagglutination of the reagent **R1**: a reddish-brown film can be observed in the well. In the absence of Rheumatoid Factor, these red blood cells deposit, forming a button in well bottom. The reagent **R2**, composed of unsensitized sheep red blood cells, ensures the reaction specificity and allows the elimination of interferences due to natural anti-sheep agglutinins (Forsman heteroantibodies, infectious mononucleosis antibodies and other macroglobulins). The reaction is carried out in a V-microplate. Handling is simple and fast, with results within 2 hours.

4 - REAGENTS AND MATERIAL

Description	Quantity
R1: Vial of 4.4 mL of active reagent	1
R2: Vial of 1 mL of control reagent	1
BUF: Vial of 55 mL of phosphate buffer pH 7.2	1
R3: Vial of 2 mL of adsorbent	1
CONTROL +: Vial of 0.5 mL of titrated positive control	1
CONTROL -: Vial of 0.5 mL of negative control	1
MICROPLATE: Microplate with a V-bottom	3
DROPPER: Special dropper	2

5 - PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only and must be handled by authorized personnel.
- Tests are for single use only.
- All the reagents, except the **BUF** reagent, contain raw materials of animal origin and must be handled with caution.
- Patient samples are potentially infectious. They must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- The **R1** reagent contains raw materials of human origin which has been screened for and found not to contain anti-HIV antibodies, anti-HCV antibodies and HBs Ag. Nonetheless it must be handled as a potentially infectious product.
- The **CONTROL** reagents contain sodium azide (<0.1%).
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Prior to use, allow the serum and the reagents to reach room temperature.
- Carefully shake the **R1** and **R2** reagents before use.
- When dispensing the **R1** and **R2** reagents, make sure that the dropper is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery volumes.

6 – SAMPLE COLLECTION AND TREATMENT
Use fresh serum or serum preserved at -20°C, and not showing any sign of haemolysis, cloudiness or of contamination.
Avoid repeated freezing and defrosting.
Do not decompartment the serum.

7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS
The reagents are ready-to-use.
All the reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box. Do not freeze.

8 - MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
- Contaminated waste containers;
- Centrifuge;
- Haemolysis tubes.

9 - METHOD
Allow the reagents to reach room temperature before use.

9.1 – Sample preparation
Carry out a 1:20 dilution of the serum to be tested:

- 50 µL of serum;
- 950 µL of **BUF** reagent.

9.2 – Realization of the test on a microplate

- Using a multichannel micropipette, add 50 µL of **BUF** reagent to 12 wells of the microplate.

- Using a micropipette, add 50 µL of diluted serum to the 1st well. Mix the serum with the **BUF** reagent and carry out a serial dilution, preferably using a microdiluter, by transferring 50 µL from the 1st well into the 2nd, then 50 µL from the 2nd to the 3rd, and so on until the 10th well is reached. 50 µL from the 10th well is then discarded. In this way, dilutions from 1:40 to 1:20480 are obtained.

- Add 50 µL of diluted serum to the 11th well. Mix the serum with the **BUF** reagent and then discard 50 µL. This dilution (1:80) is the serum control, whose role is to detect the natural anti-sheep agglutinins that could be present in certain serum samples.

- Carefully shake the **R1** and **R2** reagents.
 - Add 1 drop of **R1** reagent to the first 10 wells.
 - Add 1 drop of **R2** reagent to the 11th well (serum control).
 - Add 1 drop of **R1** reagent to the 12th well (reagent control) whose role is to control the validity of the **BUF** and **R1** reagent.

Note: Only carry out one reagent control for each series of tests.

- Very carefully, shake the contents of the wells manually, by tapping laterally the side of the microplate that has been posed flat on the bench.

- Now leave the plate to rest, away from any sources of vibration.

- The plate can be read after 2 hours.


9.3 - Adsorption of the natural anti-sheep agglutinins in the event of agglutination of the serum control

- Carefully shake the **R3** reagent.
- In a tube, add and mix:
 - 0.1 mL of serum;
 - 0.3 mL of **R3** reagent.
- Incubate at room temperature for 60 minutes.
- Centrifuge at 2000 rpm for 15 minutes.
- Collect the supernatant; the serum is now at a 1:4 dilution.
- Carry out a 1:5 dilution of the supernatant in **BUF** reagent to obtain an adsorbed stock dilution (1:20).
- Follow the steps described in "Realization of the test on a microplate", but replace the stock dilution by the adsorbed stock dilution.

10 - READING AND INTERPRETATION OF RESULTS

Negative reaction: Absence of haemagglutination.
Presence of a button at the bottom of the well.

Positive reaction: Presence of haemagglutination.
Presence of a cloudy red/brown deposit coating the well.

The titer is given by the last dilution showing haemagglutination.
The titer in IU/mL is equal to the inverse of this limit dilution multiplied by the sensitivity threshold printed on the box.

Example: If a serum is positive until the 1:1280 dilution, and if the sensitivity threshold is 0.2 IU/mL, the titer of the serum will be: 1280 x 0.2 = 256 IU/mL.

11 – INTERNAL QUALITY CONTROL
The **CONTROL +** and **CONTROL -** reagents must be treated like test serums. The titer of the **CONTROL +** reagent must be the same as the titer printed on the vial label ± one dilution. There must not be any haemagglutination of the **CONTROL -**. If haemagglutination is present then the test is not valid.

12 – CAUSES OF ERROR AND TEST LIMITS

- Poor conservation of the serum.
- Poor conservation of the reagents after opening.
- Only use the droppers provided in the kit.
- Do not interchange the droppers between the **R1** and **R2** reagents.
- In the case of a positive reaction in the first 10 wells, carry out a further serial dilution in order to determine the titer limit of haemagglutination.
- The serum control must give a negative reaction (button of red blood cells). In the event of haemagglutination of this control, it will be necessary to renew the test after having eliminated the natural anti-sheep agglutinins from the serum by adsorption.
- The reagent control must give a negative reaction (button of red blood cells). In the event of haemagglutination of this control, the **ELI.H.A RF** cannot be used.
- Certain serums, whose antibody concentration is very high, can give rise to a zone phenomenon (with disappearance of the clouding) in the initial dilutions, which disappears in the subsequent dilutions.
- The Rheumatoid Factor is a high diagnostic interest biological marker for adult rheumatoid arthritis. Nevertheless, the Rheumatoid Factor is not rheumatoid arthritis specific and may also be present, at low levels, with healthy people (more than 70 years old) or in case of connectivities (Gougerot Sjögren syndrome, disseminated lupus erythematosus, mixed connectivities...).
- On the other hand, low titers may be observed at the beginning of the disease or in treated patients.
- It is generally admitted that the level of Rheumatoid Factor is significant when it is higher than about 12 IU/mL (6), although the probability is higher for a more elevated titer (usually 30 IU/mL or more) (7). Each batch sensitivity threshold, expressed in IU/mL, (determined by W.H.O. Standard) is printed on the box (6).
- In addition, the use of **WAALER-ROSE Bicolor** is recommended to perform the screening of the Rheumatoid Factor in 2 minutes by a standardized Waaler-Rose slide reaction.
- The quality of the reagents makes it possible to carry out the reaction in the evening and to read the test the following morning, provided that the microplate is not moved in any way and is protected from any sources of vibration
- In all cases, it is necessary that the clinical, epidemiological and biological data are taken fully into consideration before establishing the final diagnosis.

13 - PERFORMANCE
The results of the **ELI.H.A RF** evaluations show a sensitivity of 97.6% and a specificity of 97.4%.

14 - WASTE ELIMINATION
Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.
If the **BUF** reagent is spilled, clean the work area with absorbent paper and rinse with water.
If a serum or another reagent is spilled on the work area, clean using bleach and absorbent paper.

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