STREPTO Color

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STREPTO Color

Haemagglutination slide test for the detection of antibodies to streptococcal extracellular antigens in serum samples and plasma

50 tests
(Ref. 04145)

8000330-en-2011-12

1 – AIM

STREPTO Color is a haemagglutination slide test for the rapid qualitative and semi-quantitative detection of antibodies to the Group A Streptococcus antigens. Each kit, Ref. 04145, allows 50 tests to be carried out.

2 – INTRODUCTION

STREPTO Color offers a simple method for rapid detection of Group A Streptococcus, such may develop in streptococcal pharyngitis, rheumatic fever, pyoderma, glomerulonephritis, and other related condition (4)

3 – PRINCIPLE

The **REAG** reagent, brown-green colour, is composed of sheep red blood cells sensitized by exoenzyme of Group A Streptococcus (streptolysin, streptokinase, hyaluronidase, streptodornase, ANAdase), in a green dye.

The presence of anti-exoenzyme of Group A Streptococcus antibodies is revealed by haemagglutination on a light green background, visible to the naked eye. Otherwise, no agglutination is observed and the suspension remains homogeneous and an uniform brown-green colour.

Handling is simple and fast. The results are obtained in 2 minutes.

4 – REAGENTS AND MATERIAL

Description – Kit of 50 tests (Ref. 04145)	Quantity
REAG: vial containing 1.0 mL of sensitized red blood cells	1
CONTROL +: vial containing 0.3 mL of positive control	1
CONTROL -: vial containing 0.3 mL of negative control	1
TEST CARD: disposable reaction cards	7
STICK: disposable stirrers	50

5 – PRECAUTIONS

- The reagents are intended solely for in vitro diagnostic use and must be handled by authorized personnel.
- Tests are for a single use only.
- **REAG** reagent contains raw materials of animal origin. The samples are potentially infectious. They must be handled with caution, in observance of hygiene rules and current regulations for this type of product in the country of use.
- The reagents **CONTROL +** and **CONTROL -** contain raw materials of human origin which has been screened for and found not to contain anti-HIV antibodies, anti-HCV antibodies and HBs Ag. Nonetheless it must be handled as a potentially infectious product.
- The product contains sodium azide (< 0.1 %).
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Allow reagents to reach room temperature.
- Carefully shake the **REAG** reagent before use.

6 – SAMPLE COLLECTION

Use fresh or inactivated serum, or non-heparinized plasma. The serum or plasma can be stored during 24 hours at 2-8°C. If the test is not performed within 24 hours after the sample collection, they must be frozen at -20°C.

It is recommended to prepare aliquots in order to avoid repeated freezing and defrosting.

Do not use a serum or a plasma showing any sign of haemolysis, cloudiness or contamination.

7 – STABILITY, STORAGE AND PREPARATION OF REAGENTS

All reagents are ready-to-use.

All reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box. **Do not freeze.**

8 – MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
- Haemolysis tubes;
- Physiological water;
- Contaminated waste container.

9 – METHOD

Allow the reagents to reach room temperature before use.

QUALITATIVE METHOD

1 – Sample preparation (1:100)

Carry out a 1:100 dilution of the sample to be tested:

- 10 µL of serum or plasma;
- 990 µL of physiological water.

2 – Realization of the test on a slide

- Using a micropipette, add 20 µL of diluted sample to one disposable slide circle.
- After ensuring that the **REAG** reagent is well mixed, add 20 µL using a micropipette.
- Mix the 2 drops with a disposable stirrer and spread them to cover the entire surface of the circle.
- Apply **a slow rocking motion from top to bottom (10 per minute)** for 2 minutes and then put the slide on a flat area.
- Observe within 10 seconds the possible appearance of brown agglutinates on a light green background.

SEMI-QUANTITATIVE METHOD

In case of a positive result on 1:100 diluted sample, it is possible to evaluate the amount of streptococcal antibodies by testing increasing dilutions in physiological water, until a negative reaction is obtained:

a. Prepare a serie of sample dilutions in physiological water, as indicated below:

Tube 1	Dilution 1:100	10 µL of sample + 990 µL of physiological water
Tube 2	Dilution 1:200	0.5 mL of tube 1 + 0.5 mL of physiological water
Tube 3	Dilution 1:400...	0.5 mL of tube 2 + 0.5 mL of physiological water

b. Carry out a slide test for each dilution following the procedure described in §“QUALITATIVE METHOD”.

10 – READING

Positive reaction: Formation of brown agglutinates visible to the naked eye on more or less light green background.

Negative reaction: No agglutination. The suspension remains homogeneous and brown-green.

Warning: A slight granulation, without a light green background, should not be confounded with a positive reaction.

The titer in IU/mL is the higher dilution showing a clear agglutination, visible to the naked eye, with a light green background (normal values are lower than 1:100).

11 – INTERNAL QUALITY CONTROL

The **CONTROL** reagents are ready-to-use and **have to be used pure**. They allow to validate the test. An agglutination must be observed with the **CONTROL+** reagent and no agglutination must be observed with the **CONTROL-** reagent. If not, then the test is not valid.

12 – INTERPRATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	Presence of at least one of five anti-exoenzyme of Group A Streptococcus antibodies (SD, ASLO, ASK, ANAdase and ASH)
NEGATIVE	Absence of the five anti-exoenzyme of Group A Streptococcus antibodies

13 – CAUSES OF ERROR AND TEST LIMITS

- It is recognized that a single determination of anti-exoenzymes Group A Streptococcus antibodies is not as significant as serial titrations performed for up to 6 weeks after diagnosis. Sequential determinations can give the trend of the patient's antibody production to allow for clinical conclusions with regard to the progress of the disease and treatment.
- In all cases, it is necessary that the clinical, epidemiologic and biological data as well as other test results are fully taken into consideration before establishing the final diagnosis.

14 – PERFORMANCE

STREPTO Color enables the detection of anti-exoenzymes Group A Streptococcus antibodies, which appear frequently (more than 90% of cases) in the serum of patients with streptococcal infections (especially rheumatic fever, glomerulonephritis, bacterial endocarditis, scarlet fever...).

The performance studies of **STREPTO Color** demonstrated a sensitivity and specificity above 95%.

15 – WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use. If the reagent or the serum is spilled, clean using bleach and absorbent paper.

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