ELITex Bicolor Mono

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ELITex Bicolor Mono

Latex slide agglutination test for the detection of infectious mononucleosis

50 tests
(Ref. 04155)

100 tests
(Ref. 04156)

8000080-en-2022-05

1 - AIM

ELITex Bicolor Mono is a latex slide agglutination test for the rapid detection of infectious mononucleosis heterophile antibodies which appear at the beginning of the disease.

Each kit, Ref. 04155, allows 50 tests to be carried out.
Each kit, Ref. 04156, allows 100 tests to be carried out.

2 - INTRODUCTION

Infectious mononucleosis is a disease caused by Epstein-Barr virus, relatively frequent, usually harmless and often causing a high asthenia. From the earliest days of disease, heterophile antibodies of infectious mononucleosis are present in the serum in sufficient quantity to allow diagnosis. They disappear gradually but can persist up to 3 months (9-11).

These heterophile antibodies can agglutinate sheep, horse, beef and rabbit red blood cells (4).

3 - PRINCIPLE

ELITex Bicolor Mono is a direct agglutination slide test. **TEST LATEX** reagent consists of latex particles sensitized by erythrocytes stroma specially treated, making unnecessary the adsorption of other heterophile antibodies (2) and ensuring the reaction specificity (3).

The presence of heterophile antibodies in sera is revealed by a red agglutination on a more or less intense blue background, visible to the naked eye. Otherwise, no agglutination is observed and the suspension remains homogeneous and a uniform purple colour.

Handling is simple and fast. It is performed on pure serum and the results are obtained in 2 minutes.

4 - REAGENTS AND MATERIAL

Description - Kit of 50 tests (Ref. 04155)	Quantity
TEST LATEX: dispenser vial containing 1.5 mL of sensitized latex	1
CONTROL +: dispenser vial containing 0.5 mL of positive control	1
CONTROL -: dispenser vial containing 0.5 mL of negative control	1
TEST CARD: disposable reaction cards	7
STICK: Disposable stirrers	50

Description - Kit of 100 tests (Ref. 04156)	Quantity
TEST LATEX: dispenser vial containing 3 mL of sensitized latex	1
CONTROL +: dispenser vial containing 1 mL of positive control	1
CONTROL -: dispenser vial containing 1 mL of negative control	1
TEST CARD: disposable reaction cards	13
STICK: Disposable stirrers	100

5 - PRECAUTIONS

- The reagents are intended solely for *in vitro* use and must be handled by authorized personnel.
- Tests are for a single use only.
- The **TEST LATEX** and **CONTROL-** reagents contain raw materials of animal origin. The samples are potentially infectious. They must be handled with caution,

in observance of hygiene rules and current regulations for this type of product in the country of use.

- The **CONTROL+** reagent contains raw materials of human origin which has been screened for and found not to contain anti-HIV antibodies, anti-HCV antibodies and HBs Ag. Nonetheless it must be handled as a potentially infectious product.
- The reagents contain sodium azide (< 0.1 %).
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Allow the reagents to reach room temperature.
- Carefully shake the **TEST LATEX** reagent before use.
- When dispensing reagents, make sure that the dispenser vial is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery volumes. As a precautionary measure, wipe the dispenser vial tip after use.

6 - SAMPLE COLLECTION

Use fresh serum.

The serum can be stored during 24 hours at 2-8°C. If the test is not performed within 24 hours after the serum collection, the serum must be frozen at -20°C. It is recommended to prepare aliquots in order to avoid repeated freezing and defrosting. Do not deplete the serum. Do not use a serum showing any signs of haemolysis, cloudiness or contamination.

7 - CONSERVATION AND PREPARATION OF REAGENTS

All reagents are ready-to-use.

Reagents conserved at 2-8°C, in their original state, are stable up to the expiry date indicated on the box. **Do not freeze.**

8 - MATERIAL REQUIRED BUT NOT SUPPLIED

Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
Physiological water;- Contaminated waste container.

9 - METHOD

Allow the reagents to reach room temperature before use.

QUALITATIVE METHOD

- Using a micropipette, add 50 µL of serum to one of the **circle** of the disposable slide.
- After ensuring that the **TEST LATEX** reagent is well mixed, add 1 drop of reagent using the dispenser vial.
- Mix the 2 drops **using a disposable stirrer** to cover the entire surface reaction. **Apply a slow rotation for 2 minutes.** A red agglutination may appear on a more or less intense blue background.

SEMI-QUANTITATIVE METHOD

In case of a positive result, it is possible to evaluate the amount of heterophile antibodies to determine its evolution during the disease, by testing increasing dilutions.

- Prepare serial dilutions from 1:2 to 1:2 using physiological water.
- Carry out a slide test for each dilution following the procedure described in § «QUALITATIVE METHOD».

10 - READING

Positive reaction: **Formation of red agglutinates** on a more or less blue intense background.

Negative reaction: **No agglutination.** The suspension remains homogeneous and purple.

For the semi-quantitative method, serum titer is comprised between the inverse of the last positive dilution and the inverse of the first negative dilution.

11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	Presence of infectious mononucleosis heterophile antibodies
NEGATIVE	Absence of infectious mononucleosis heterophile antibodies

12 - INTERNAL QUALITY CONTROL

The **CONTROL** reagent are ready-to-use and should be used pure. They allow to validate the test. An agglutination must be observed with the **CONTROL +** and no agglutination must be observed with the **CONTROL-**. If not, then the test is not valid.

13 - CAUSES OF ERROR AND TEST LIMITS

- The presence of a mononucleosis syndrome (7) associated to a positive result with **ELITEX Bicolor Mono** allows, in most cases, to give a rapid and early diagnosis of infectious mononucleosis. Some forms of mononucleosis without heterophile antibodies have yet been reported by different authors, particularly among young children and in some cases with moderate clinical manifestations. In this event, it is recommended to associate an anti-EBV antibodies detection (5, 8).
- In all cases, it is necessary that the clinical, epidemiologic and biological data are taken fully into consideration before establishing the final diagnosis.

14 - PERFORMANCE

Evaluation has demonstrated that the test **ELITex Bicolor Mono** has a sensitivity of 98.4% and a specificity of 99.7%.

15 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use. If the reagent or the serum is spilled, clean using bleach and absorbent paper.

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Changes from the previous version are highlighted in gray