# ELITex Bicolor RF

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# ELITex Bicolor RF

Latex slide agglutination test for  
the detection of Rheumatoid Factor in serum samples

**50 tests**  
(Ref. 04165)

**100 tests**  
(Ref. 04166)

8000090-en-2022-05

## 1 - AIM

**ELITex Bicolor RF** is a latex slide agglutination test for the rapid detection of Rheumatoid Factor in serum samples.

Each kit, Ref. 04165, allows 50 tests to be carried out.

Each kit, Ref. 04166, allows 100 tests to be carried out.

## 2 - INTRODUCTION

The Rheumatoid Factor is a high diagnostic interest biological marker for adult rheumatoid arthritis (1-4-5).

## 3 - PRINCIPLE

**ELITex Bicolor RF** is an agglutination slide test. It consists of latex particles allowing the rapid qualitative or semi-quantitative detection of IgM Rheumatoid Factor in sera from most of patients with rheumatoid arthritis. The reactif **TEST LATEX** is constituted by latex particles sensitized by human gamma-globulins, which agglutinate specifically Rheumatoid factor (FII-type reaction). The presence of Rheumatoid Factor in sera is revealed by a red agglutination on a more or less intense blue background, visible to the naked eye. Otherwise, no agglutination is observed and the suspension remains homogeneous with an uniform purple colour.

Handling is simple and fast. It is performed on pure serum and the results are obtained in 2 minutes.

## 4 - REAGENTS AND MATERIAL

Description - Kit of 50 tests (Ref. 04165)	Quantity
<b>TEST LATEX:</b> dispenser vial containing 2.5 mL of sensitized latex	1
<b>CONTROL +:</b> vial containing 0.2 mL of positive control	1
<b>CONTROL -:</b> vial containing 0.2 mL of negative control	1
<b>TEST CARD:</b> disposable reaction cards	8
<b>STICK:</b> Disposable stirrers	50

Description - Kit of 100 tests (Ref. 04166)	Quantity
<b>TEST LATEX:</b> dispenser vials containing 2.5 mL of sensitized latex	2
<b>CONTROL +:</b> vial containing 0.2 mL of positive control	1
<b>CONTROL -:</b> vial containing 0.2 mL of negative control	1
<b>TEST CARD:</b> disposable reaction cards	13
<b>STICK:</b> Disposable stirrers	100

## 5 - PRECAUTIONS

- The reagents are intended solely for *in vitro diagnostic* use and must be handled by authorized personnel.
- Tests are for a single use only.
- Reagents contain raw materials of animal origin. The samples are potentially infectious. They must be handled with caution, in observance of hygiene rules and current regulations for this type of product in the country of use.
- The **TEST LATEX** reagent contains raw materials of human origin which has been screened for and found not to contain anti-HIV antibodies, anti-HCV antibodies and HBs Ag. Nonetheless it must be handled as a potentially infectious product.
- Reagents contain sodium azide (< 0.1 %)
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Allow reagents to reach room temperature.
- Carefully shake the **TEST LATEX** reagent before use.

When dispensing the **TEST LATEX** reagent, make sure that the dispenser vial is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery

volumes. Wipe the tip of the **TEST LATEX** reagent dispenser vial in order to obtain calibrated drops.

## 6 - SAMPLE COLLECTION

Use fresh serum.

The serum can be stored during 24 hours at 2-8°C. If the test is not performed within 24 hours after the sample collection, the serum must be frozen at - 20°C. It is recommended to prepare aliquots in order to avoid repeated freezing and defrosting.

Do not decomplex the serum.

Do not use a serum showing any sign of haemolysis, cloudiness or contamination.

## 7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

All reagents are ready-to-use

All reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box.

**Do not freeze.**

## 8- MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
- Haemolysis tubes;
- Physiological water;
- Contaminated waste container.

## 9 - METHOD

**Allow the reagents to reach room temperature before use.**

### QUALITATIVE METHOD

a. Using a micropipette, add 10 µL of serum to one of the **circle** of the disposable slide.

b. After ensuring that the **TEST LATEX** reagent is well mixed, add 1 drop of reagent using the dispenser vial.

c. Mix the 2 drops using a disposable stirrer to cover the entire reaction area.

Apply a slow rotation for 2 minutes, a red agglutination may appear on a more or less blue background.

### SEMI-QUANTITATIVE METHOD

In case of a positive result on pure serum, it is possible to evaluate the amount of Rheumatoid Factor by testing increasing dilutions in physiological water, until a negative reaction is obtained:

a. Prepare a serie of sample dilutions in physiological water, as indicated below:

Tube 1	Dilution 1:2	100µL of sample + 100 µL of physiological water
Tube 2	Dilution 1:4	100µL of tube 1 + 100 µL of physiological water
Tube 3	Dilution 1:8	100µL of tube 2 + 100 µL of physiological water
Tube 4	Dilution 1:16...	100µL of tube 3 + 100 µL of physiological water

b. Carry out a slide test for each dilution following the procedure described in § « QUALITATIVE METHOD »

## 10 – READING

**Positive reaction :** Formation of red agglutinates on more or less blue intense background.

**Negative reaction :** No agglutination. The suspension remains homogeneous and purple.

**The titer IU/mL is equal to the inverse of the higher dilution showing a clear agglutination multiplied by the sensitivity threshold printed on the box** (determined by WHO standard).

Example: if a serum is positive until 1:8 dilution, and if the sensitivity threshold is 8 IU/mL, the titer of the serum will be :  $8 \times 8 = 64$  IU/mL.

Warning: A slight granulation in the purple background should not be confounded with a positive reaction.

## 11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	Presence of Rheumatoid Factor
NEGATIVE	Absence of Rheumatoid Factor

## 12 - INTERNAL QUALITY CONTROL

The **CONTROL** reagents are ready-to-use and have to be used pure. They allow to validate the test. An agglutination must be observed with the **CONTROL+** reagent and no agglutination must be observed with the **CONTROL-** reagent. If not, then the test is not valid.

## 13 - CAUSES OF ERROR AND TEST LIMITS

- The Rheumatoid Factor is a high diagnostic interest biological marker for adult rheumatoid arthritis. When associating a Waaler-Rose reaction (**WAALER-ROSE Bicolor**), it has been observed for 80% of the patients after a one year evolution period (2). Nevertheless, the Rheumatoid Factor is not specific and may also be present, at low levels, in healthy people (more than 70 years old) or in case of connectivities (Gougerot Sjögren syndrome, disseminated lupus erythematosus, mixed connectivities...). On the other hand, low titers may be observed at the beginning of the disease in treated patients.

- It is generally admitted that the level of Rheumatoid Factor is significant when higher than about 12 IU/mL (3), although the probability is higher for a more elevated titer (usually 30 IU/mL) or more. Each batch sensitivity threshold, expressed in IU/mL (determined by W.H.O Standard) is printed on the box.

- In all cases, it is necessary that the clinical, epidemiologic and biological data as well as other test results are fully taken into consideration before establishing the final diagnosis.

## 14 – PERFORMANCE

Evaluation has demonstrated that the test **ELITex Bicolor RF** has a sensitivity of 98.75 % and a specificity of 98.37 %.

## 15 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

If the reagent or the serum is spilled, clean using bleach and absorbent paper.

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