ELI.H.A Distoma

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ELI.H.A Distoma

Serodiagnosis of distomatosis by indirect haemagglutination

102 tests
(Ref: 66606)

8000150-en-2012-09

1 - AIM
ELI.H.A Distoma enables the quantitative determination of anti-*Fasciola hepatica* serum antibodies by indirect haemagglutination.

Each kit allows 102 tests to be carried out or 17 reactions of 6 dilutions.

2 - INTRODUCTION

The distomatosis are parasitic diseases caused by trematodes. The location of the adult worms in the host gives rise to the following classifications:

- hepato-biliary distomatosis
- intestinal distomatosis
- pulmonary distomatosis

Of the hepato-biliary distomatosis, infection with *Fasciola hepatica* is characterized by a progressive life-cycle in the human definitive host, which includes a transperitoneal and transhepatic larval migratory phases. Diagnosis of *F. hepatica* is primarily immunological during the invasive phase of infection and relies upon the detection of antibodies (phase where there is no faecal egg transmission). Immunological diagnosis is also of great interest during the early disease state phase, because the eggs are only emitted in small quantities.

3 - PRINCIPLE
ELI.H.A Distoma is based on the indirect haemagglutination principle. The sensitized red blood cells consist of sheep red blood cells covered with a *Fasciola hepatica* antigen. The presence of anti-*Fasciola hepatica* serum antibodies results in agglutination of the sensitized red blood cells resulting in a cloudy red/brown deposit coating the well. In the absence of specific antibodies, the red blood cells form a ring-like deposit at the bottom of the well. The non-sensitized red blood cells ensure the specificity of the reaction making it possible to eliminate any interference from the natural anti-sheep agglutinins (Forsman heteroantibodies, erythrocyte mononuclear antibodies, ...).

4 - REAGENTS AND MATERIAL

Handling is simple and fast, with results within 2 hours. The reaction is carried out in a U-microplate.

Description	Quantity
R1 : Vial of 2 mL of sensitized red blood cells	1
R2 : Vial of 1 mL of non-sensitized red blood cells	1
BUF : Vial of 55 mL of phosphate buffer pH 7.2	1
R3 : Vial of 2 mL of adsorbent	1
CONTROL + : Vial of 0.2 mL of titrated positive control	1
CONTROL - : Vial of 0.2 mL of negative control	1
MICROPLATE : Microplate with a U-bottom	2
DROPPER : Special dropper	2

5 - PRECAUTIONS

The reagents are intended for *in vitro* diagnostic use only and must be handled by authorized personnel.

- Tests are for a single use only.
- All the reagents, except the **BUF** reagent, contain raw materials of animal origin and must be handled with caution.
- Patient samples are potentially infectious. They must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- Do not use reagents from different batch numbers.
- Prior to use, allow the serum and the reagents to reach room temperature.
- Carefully shake the **R1** and **R2** reagents before use.
- When dispensing the **R1** and **R2** reagents, make sure that the dropper is perfectly vertical.

6 - SAMPLE COLLECTION AND TREATMENT
Use fresh serum or serum preserved at -20°C, and not showing any sign of haemolysis, cloudiness or contamination.
Avoid repeated freezing and defrosting.
Do not decontaminate the serum.

7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

The reagents are ready-to-use.
All the reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box. Do not freeze.

8 - MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
- Contaminated waste containers;
- Centrifuge;
- Haemolysis tubes.

9 - METHOD

Allow the reagents to reach room temperature before use.

9.1 - Sample preparation

Carry out a 1:40 dilution of the serum to be tested:

- 50 µL of serum;
- 1.95 mL of **BUF** reagent.

9.2 - Realization of the test on a microplate

Using a multichannel micropipette, add 50 µL of **BUF** reagent to 8 wells of the microplate.

- Using micropipette, add 50 µL of diluted serum to the 1st well.
- Mix the serum with the **BUF** reagent and then discard 50 µL.
- In this way, dilutions from 1:80 to 1:2560 are obtained.
- Carefully shake the **R1** and **R2** reagents.
- Add 1 drop of **R1** reagent to the first 6 wells.
- Add 1 drop of **R2** reagent to the 7th well (serum control).
- Add 1 drop of **R1** reagent to the 8th well (reagent control) whose role is to control the validity of the **BUF** and **R1** reagent.

Very carefully, shake the contents of the wells:

- either manually, by tapping laterally the side of the microplate that has been posed flat on the bench;
- or by using a vibrating plate shaker for microtiter plates (for example at 1300 rpm for 10 seconds). Do not use an orbital shaker.

Now leave the plate to rest, away from any sources of vibration.

The plate can be read after 2 hours.

9.3 - Adsorption of the natural anti-sheep agglutinins in the event of agglutination of the serum control

- Carefully shake the **R3** reagent.
- In a tube, add and mix:
- 0.1 mL of serum;
- 0.3 mL of **R3** reagent.
- Incubate at room temperature for 60 minutes.
- Centrifuge at 2000 rpm for 15 minutes.
- Collect the supernatant; the serum is now at a 1:4 dilution.
- Carry out a 1:10 dilution of the supernatant in **BUF** reagent to obtain an adsorbed stock dilution (1:40).
- Follow the steps described in "Realization of the test on a microplate", but replace the stock dilution by the adsorbed stock dilution.

10 - READING

Negative reaction:
Absence of haemagglutination.

Presence of a more or less large ring at the bottom of the well.

Positive reaction:
Presence of haemagglutination.

Presence of a cloudy red/brown deposit coating the well, sometimes there is the presence of a fine peripheral border.

Example: Serum positive at a dilution of 1:1280

11 - INTERPRETATION OF RESULTS

Non significant reaction.

Probable absence of distomatosis.

Renew the test 2 to 3 weeks later and also carry out an electrosyneresis or an immunoelectrophoresis test.

Doubtful reaction.

Renew the test 2 to 3 weeks later and also carry out an electrosyneresis or an immunoelectrophoresis test.

Titre ≥ 1:320:

Significant reaction in favour of progressive distomatosis.

12 - INTERNAL QUALITY CONTROL

The **CONTROL +** and **CONTROL -** reagents must be treated like test serums. The titer of the **CONTROL +** reagent must be the same as the titer printed on the vial label ± one dilution. There must not be any haemagglutination of the **CONTROL -**. If haemagglutination is present then the test is not valid.

13 - CAUSES OF ERROR AND TEST LIMITS

- Poor conservation of the reagents after opening.
- Only use the droppers provided in the kit.
- Do not interchange the droppers between the **R1** and **R2** reagents.
- In the case of a positive reaction in the first 6 wells, carry out a further serial dilution in order to determine the titer limit of haemagglutination.
- The serum control must give a negative reaction (ring). In the event of haemagglutination of this control, it will be necessary to renew the test after having eliminated the natural anti-sheep agglutinins from the serum by adsorption.
- The reagent control must give a negative reaction (ring). In the event of haemagglutination of this control, the ELI.H.A Distoma cannot be used.
- Certain serums, whose antibody concentration is very high, can give rise to a zone phenomenon (with disappearance of the clouding) in the initial dilutions, which which disappears in the subsequent dilutions.
- The quality of the reagents makes it possible to carry out the reaction in the evening and to read the test the following morning, provided that the microplate is not moved in any way and is protected from any sources of vibration
- In all cases, it is necessary that the clinical, epidemiologic and biological data are taken fully into consideration before establishing the final diagnosis.

14 - PERFORMANCE

ELI.H.A Distoma consists of red blood cells sensitized with a *Fasciola hepatica* antigen, that ensures the specificity and sensitivity of the reaction.

Evaluation has demonstrated that the test has a sensitivity of 96.1% and a specificity of 96.6%.

15 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

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