ELITex Bicolor Amoeba

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☎ : 04 94 88 55 00
Fax : 04 94 88 55 22

<p>ELITex Bicolor Amoeba</p> <p>Latex slide agglutination test for the detection of amoebic antibodies in serum</p> <p>25 tests (Ref: 66608)</p>
--

8000040-en-2012-05 

ELITex Bicolor Amoeba is a slide agglutination test for the detection of amoebic antibodies in serum. Each kit allows 25 tests to be carried out.

2 - INTRODUCTION
Amoebiasis is a parasitic infection caused by a protozoon specific to man: *Entamoeba histolytica*. It is the only pathogenic amoeba of man. Serodiagnosis enables the confirmation of hepatic and pulmonary amoeba infections.

3 - PRINCIPLE
The **TEST LATEX** reagent, brown in colour, is made up of red dyed latex particles in suspension with a green counter dye. These particles are sensitized with an *Entamoeba histolytica* total mixed antigen.
The presence of amoebic antibodies in serum is characterized by a red agglutination, that could progressively spread towards the periphery, on a green background.

In absence of specific antibodies, the suspension remains homogeneous and of an uniform brown colour.

4 - REAGENTS AND MATERIAL

Description	Quantity
TEST LATEX: Vial of 0.5 mL of sensitized latex	1
CONTROL +: Vial of 0.5 mL of positive control	1
CONTROL -: Vial of 0.5 mL of negative control	1
DIL: Vial of 5 mL of special diluent (dispensing drops of ± 40µL)	1
TEST CARD: Disposable reaction card	4
STICK: Disposable mixing sticks	25
DROPPER: Special dropper	1

5 - PRECAUTIONS
- The reagents are intended for *in vitro diagnostic* use only and must be handled by authorized personnel.
- Tests are for a single use only.
- The **TEST LATEX** and **CONTROL** reagents contain raw materials of animal origin and must be handled with caution.
- The reagents contain sodium azide (> 0.1 %).
- The serum dilution must be done with the **DIL** reagent provided in the kit.
- Patient samples are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Prior to use, allow the reagents to reach room temperature.
- Carefully shake the **TEST LATEX** reagent before use. Do not turn the vial upside down.
- When dispensing the **TEST LATEX** reagent, make sure that the dropper is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery volumes.

6 - SAMPLE COLLECTION

Use fresh serum.
The serum can be stored during 24 hours at 2-8°C. If the test is not performed within 24 hours after the sample collection, the serum must be frozen at - 20°C. It is recommended to prepare aliquots in order to avoid repeated freezing and defrosting.
Do not decompartment the serum.
Do not use a serum showing any signs of haemolysis, cloudiness or of contamination.

7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

All reagents are ready-to-use. All reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box. **Do not freeze.**

8 - MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured;
- Pasteur pipette or wire loop;
- Contaminated waste container.

9 - METHOD

Allow the reagents to reach room temperature before use.

QUALITATIVE METHOD

Preparation of 1:5 dilution of test serum

- 2 drops of **DIL** reagent;
- 20 µL of test serum.

- After ensuring that the **TEST LATEX** reagent is well mixed, add 1 drop of reagent with the dropper provided, or with the aid of a micropipette, add 20 µL.
- Using a micropipette, add 20 µL of diluted serum to one of the circles of the slide.

SEMI-QUANTITATIVE METHOD

- Using the provided **DIL** reagent, prepare serial dilutions from 1:2 to 1:2, starting from the 1:5 initial dilution.
- Carry out a slide test for each dilution following the procedure described in § QUALITATIVE METHOD.

10 - READING

For the qualitative method:
Positive reaction : Formation of a red agglutination (which can progressively spread towards the periphery) on a green background.

Negative reaction : Absence of agglutination. The suspension remains homogeneous and brown.

For the semi-quantitative method:
Serum titer is comprised between the inverse of the last positive dilution and the inverse of the first negative dilution.

11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	Presence of amoebic antibodies
NEGATIVE	Absence of amoebic antibodies
	Probable absence of amoeba tissue infection

The **CONTROL** reagents are ready-to-use and have to be used pure. They validate the test. An agglutination must be observed with the **CONTROL+** reagent and no agglutination must be observed with the **CONTROL-** reagent. If not, then the test is not valid.

13 - CAUSES OF ERROR AND TEST LIMITS

- Poor conservation of the serum.
- Poor conservation of the reagents after opening. Do not use a horizontal movement agitator (Kline type agitator).
- A positive reaction is an immunologic argument in favor of amoebiasis. A control must be carried out with the indirect haemagglutination, electrosynere, and indirect immunofluorescence techniques.
- In all cases, it is necessary that the clinical, epidemiologic and biological data are taken fully into consideration before establishing the final diagnosis.

14 - PERFORMANCE

ELITex Bicolor Amoeba consists of latex particles sensitized with an *Entamoeba histolytica* mixed total antigen, that ensures the specificity and sensitivity of the reaction. Evaluation has demonstrated that the test **ELITex Bicolor Amoeba** has a sensitivity of 98 % and a specificity of 96 %.

15 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.
If the **DIL** reagent is spilled, clean the work area with absorbent paper and rinse with water. If serum or another reagent is spilled, clean using bleach and absorbent paper.

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