

DSQ Alert™ Mycoplasma/Ureaplasma v2.0 RUO Monoreagent

For Research Use Only. Not for use in diagnostic procedures.



ELITechGroup MDx LLC
21720 23rd Dr SE, Suite 150
Bothell, WA 98021 USA
Telephone: 1-800-453-2725



M400791



**SPIN TUBES
PRIOR TO
OPENING**

Intended Use

The **DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent** is intended for use in a nucleic acid amplification test, to detect and distinguish DNA from *Mycoplasma* and *Ureaplasma* species in a nucleic acid sample. This product is intended for use with a real-time PCR system with appropriate optical specifications and melt curve analysis capability.

Assay Principle

The **DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent** is a multiplex real-time PCR reagent designed with DSQ hydrolysis probe chemistry, the next generation of MGB hydrolysis probes, to detect and distinguish DNA from ***Mycoplasma genitalium*, *Mycoplasma hominis*, *Ureaplasma urealiticum* and *Ureaplasma parvum* species**. To use this product effectively, thermal cycler parameters must include PCR thermal cycling with five color fluorescence detection. The reagent contains primers and a probe specific to the GAPDH gene for *Mycoplasma* species, and the urease gene for *Ureaplasma* species, labeled with a fluorophore and a duplex stabilizing quencher (DSQ). The reagent also contains a primer set and probe specific to an internal control (IC, sold separately).

The DSQ hydrolysis probe chemistry in this product is unique. During each cycle of PCR, the primers and probe anneal to their target template, if present, and DNA is synthesized from the primers by a polymerase. During synthesis, the polymerase encounters the probe annealed to the template downstream of the primer, and the exonuclease activity of the polymerase hydrolyzes the probe, releasing the fluorophore from the proximity of the DSQ and allowing fluorescence emission. The PCR cycles result in exponential amplification of the target DNA and fluorescence levels.

Product Description

The **DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent** is a ready-to-use mix of primer and probe sets specific to the DNA of the target pathogens, and to a synthetic sequence that serves as an internal control to monitor assay performance. (The IC DNA template is sold separately, see below.) Probes are labeled with FAM or an **AquaPhluor® (AP) fluorophore** (Table 1), and a DSQ.

Table 1. DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent components description. The number in the AP fluorophore name indicates its peak excitation wavelength.

Target template	Probe type and chemistry	Probe fluorophore	Analogous fluorophore (for optical channel selection)	Additional probe labels
<i>M. genitalium</i> GAPDH gene species-specific region	DSQ hydrolysis	FAM	FAM	DSQ
<i>M. hominis</i> GAPDH gene species-specific region	DSQ hydrolysis	AP530	VIC, JOE, HEX	DSQ
<i>U. parvum</i> urease gene species-specific region	DSQ hydrolysis	AP593	ROX, Texas Red	DSQ
<i>U. urealyticum</i> urease gene species-specific region	DSQ hydrolysis	AP639	Cy5, Quasar 670, Alexa Fluor 647	DSQ
Internal control IC1	DSQ hydrolysis	AP690	Cy5.5, Quasar 705	DSQ

DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent is provided at a volume of 1,000 µL per tube, 1 tube per kit, and contains all the necessary components for PCR except the template. The 1.250X monoreagent concentration is relative to the optimal final concentration of the primers and probes in the PCR.

Recommended Materials Not Provided

Table 2. Additional materials recommended for real-time PCR not provided in the DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent.

Material	Use	Vendor	Part Number
Internal Control IC1 DNA	Internal control DNA template to monitor nucleic acid extraction and PCR performance	ELITechGroup	M800735
Molecular biology grade water	Reaction mix preparation, negative controls	NA	NA
Positive controls	Positive control DNA for each target genotype if available	NA	NA

Recommended Reaction Setup

For optimal performance, protect all reagents from light, store at ≤-10°C while not in use, and limit the number of freeze-thaw cycles.

The following is an example of how to set up a real-time PCR using the DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent for 50 µL reactions. Preparation of the reaction mix should be done in an area separate from preparation and addition of samples and controls.

Table 3. Example recipe for real-time PCR reaction mix.

Reagent	Stock concentration	Volume per reaction (µL)
DSQ Alert RUO Monoreagent	1.250X	40.0
Sample/control template	--	10.0

1. Prepare reaction mix as above (Table 3) or adjust volumes per reaction based on final reaction volume, multiplying the volumes per reaction by the number of samples + controls being run and an appropriate overage to add the needed dead volume.
2. Array 40 µL of the monoreagent into the wells of an optical plate or tubes.
3. Prepare positive and negative controls as appropriate.
4. Pipette 10 µL of sample or control into the appropriate well or tube containing reaction mix.
5. Seal the plate with optical adhesive film or cap PCR tubes.
6. Load the plate/tubes onto the real-time PCR instrument and program the thermal cycling as below (Table 4). Start the run.

Table 4. Recommended thermal cycling conditions. Adjustments may be required to optimize the PCR for various real-time PCR instruments. Refer to the instrument manual to set up the real-time PCR.

Stage		Temperature	Time
UNG activation*	Hold	50°C	10 min
Polymerase activation	Hold	95°C	5 min
PCR (45 cycles)	Denaturation	95°C	10 sec
	Annealing/Extension**	63°C	25 sec

* The UNG activation step is optional and recommended when using a PCR master mix with UNG.

** Read fluorescence at the annealing/extension stage of PCR.

Data Analysis Guidelines

Analysis of results from the DSQ Alert Mycoplasma/Ureaplasma v2.0 RUO Monoreagent should be performed from the PCR stage. Amplification of the FAM fluorescence signal during PCR is indicative of the presence of *M. genitalium* DNA in the nucleic acid sample. Amplification of the AP530 fluorescence signal during PCR is indicative of the presence of *M. hominis* DNA in the nucleic acid sample. Amplification of the AP593 fluorescence signal during PCR is indicative of the presence of *U. parvum* DNA in the nucleic acid sample. Amplification of the AP639 fluorescence signal during PCR is indicative of the presence of *U. urealyticum* DNA in the nucleic acid sample. Amplification of the internal control AP690 signal indicates the PCR performed as expected. Amplification of the internal control AP690 signal may or may not be observed in samples that test positive for *Mycoplasma/Ureaplasma* DNA but must be observed in samples that test negative for *Mycoplasma/Ureaplasma* DNA to ensure the PCR performed as expected.

Warnings and Precautions

- This product is for Research Use Only, and not for use in diagnostic procedures.
- Use of this product requires personnel trained in molecular biology techniques.
- This product shall be protected from light and stored at ≤-10°C while not in use.
- This product shall not be used after its expiration date.

- This product shall be used in accordance with local, state, and federal regulations or accreditation requirements.
- Disposal of all waste material shall be done in accordance with local, state, and federal regulations or accreditation requirements.

Technical Support

For technical support, call or email the ELITechGroup MDx (EG MDx) Technical Support Center: 1.800.453.2725 or mdx@elitechgroup.com, or contact your EG MDx Field Applications Specialist.

Legal Notices

Limited Product Warranty

EG MDx warrants that this product will meet the specifications stated above. If any component of this product does not conform to these specifications, EG MDx will, at its sole discretion, as its sole and exclusive liability and as the users' sole and exclusive remedy, replace the product at no charge or refund the cost of the product; provided that notice of non-conformance is given to EG MDx within sixty (60) days of receipt of the product.

This warranty limits EG MDx's liability to the replacement of this product or refund of the cost of the product. NO OTHER WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT, ARE PROVIDED BY ELITECHGROUP MDx LLC. EG MDx shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use or the inability to use this product and its components.

In no event shall EG MDx be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue, or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or use or the failure of EG MDx products to perform in accordance with the stated specifications.

Licensing and Legal Notice

Some components of nucleic acid analysis, such as specific methods and compositions for manipulating or visualizing nucleic acids for analysis, may be covered by one or more patents owned by other parties. Similarly, nucleic acids containing specific nucleotide sequences may be patented. Making, using, offering for sale, or selling such components or nucleic acids may require one or more licenses. Nothing in this document should be construed as an authorization or implicit license to make, use, or sell any so covered component or nucleic acid under any such patents.

DSQ Alert detection reagent are covered by one or more of U.S. Patents Numbers 6972339, 7319022, 7348146, 7381818, 7541454, 7582739, 7601851, 7671218, 7718374, 7723038, 7759126, 7767834, 7851606, 8008522, 8067177, 8163910, 8389745, 8569516, 8969003, 9056887, 9085800, 9169256, 9328384, 10677728, 10738346, 10890529, and 11320376 as well as applications that are currently pending.

This limited license permits the person or legal entity to which this product has been provided to use the product, and the data generated by use of the product, only for internal RUO assay validation purposes involving hybridization-based analysis of nucleic acids as defined and restricted by the U.S. Food and Drug Administration in 21 CFR 864.4020 and 21 CFR 809.30. Neither EG MDx nor any of its licensors grants any other licenses, whether express or implied, for any other purposes.

Although patents covering the basic polymerase chain reaction (PCR) have expired, patents covering the use of certain enzymes and other uses of the PCR process owned by Hoffman-LaRoche and others remain in effect and may require a license. Purchase of this product does not include or provide a license with respect to these

patents. EG MDx does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is recommended for persons that either have the license to perform PCR or are not required to obtain a license. No license under the patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of this product. Nothing herein is to be construed as recommending any practice or any products in violation of any patent or in violation of any law or regulation.








Trademarks

ELITechGroup is a trademark of ELITechGroup SAS.

DSQ Alert and AquaPhluor are trademarks of ELITechGroup BV.

Symbols

The following symbols are used within ELITechGroup MDx DSQ Alert labeling.

	Catalog number		Upper limit of temperature
	Lot or Batch Code		Expiration Date YYYY-MM-DD
	Manufacturer		Keep away from sunlight
	Contains sufficient for <N> tests		