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## Transport of urogenital mycoplasmas

### UMMT Revolution 50 tests (REF 00061)

CP 0386 EN-2014-02

For *in vitro* diagnostic use only, for professional use only

#### 1 - INTENDED USE

The UMMT Revolution kit is a complement of UMMT (3 mL) transport medium for MYCOFAST Revolution (REF 00060) or MYCOFAST Screening Revolution (REF 00063) kits.

The UMMT (3 mL) medium allows the transport of *Ureaplasma urealyticum* (U.u.) and of *Mycoplasma hominis* (M.h.) present in pathological patient samples, as well as the inoculation of MYCOFAST Screening Revolution kits.

#### 2 - PRINCIPLE

The mycoplasmas are fragile bacteria lacking a cell wall and are thus sensitive to variations of osmotic pressure. The samples must be discharged into an isotonic medium containing sodium chloride, such as the UMMT transport medium. Moreover, mycoplasmas are sensitive to the products of catabolism: a strong concentration of catabolic products can induce auto-intoxication. Inoculation in a liquid medium produces a dilution effect.

#### 3 - REAGENT

UMMT (3 mL): Vial of liquid medium containing Mycoplasma broth, antibiotics and a preservative, pH: 6.0 ± 0.1.

#### 4 - PRECAUTIONS

- The reagents are intended solely for *in vitro* use and must be handled by authorised personnel.
- The patient samples and inoculated reagents are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- Reagents containing raw materials of animal origin must be handled with caution.
- Do not use reagents after the expiry date.
- Do not use reagents that have been damaged or that have been poorly conserved before use.

#### 5 - SAMPLE COLLECTION AND HANDLING

##### Sample collection

Endocervical and vaginal sample collection: Use only a Dacron or rayon swab or a cyobrush to collect samples.

The cervix should be carefully cleaned with a swab, to remove secretions, before collecting the sample with a new swab. As mycoplasmas adhere strongly to mucous cells, the mucous lining should be vigorously swabbed to obtain a rich specimen.

##### Urinal sample collection:

Urethral sample collection: Clean the meatus and swab or scrape the area to obtain cells. Sperm/urine: Collect sperm or first micturition in a sterile tube or bottle.

Gastric secretions: Collect gastric secretions from the neonate by aspiration with a catheter and transfer to a sterile bottle.

##### Inoculation and transport of the specimen

Inoculate a ml of UMMT (3 mL) medium with a swab or 300 µL of liquid sample. Storage in UMMT medium. For storage at 2-8 °C, first add two drops of "MYCOPLASMA Stabilizer".

#### 6 - PREPARATION AND STORAGE OF REAGENTS

Before opening, the UMMT medium may be kept for 20 hours at room temperature (18-25 °C) or 56 hours at 2-8 °C.

For storage during 3 days at -20 °C, first add two drops of "MYCOPLASMA Stabilizer". The UMMT (3 mL) may be stored at 2-8 °C, in their original packaging until the expiry date shown on the kit.

#### 7 - MATERIAL REQUIRED BUT NOT PROVIDED

- Sample collection (Dacron or rayon swabs, cyobrushes, sterile containers for liquid samples)
- MYCOPLASMA Stabilizer (REF 00064)
- Pipettes and tubes
- Freezer at -20 °C
- Waste container for contaminated waste

#### 8 - METHOD

Refer to the usage instructions of the MYCOFAST Revolution or MYCOFAST Screening Revolution kits.

#### 9 - PERFORMANCE

Refer to the performance data of the MYCOFAST Revolution or MYCOFAST Screening Revolution kits.

#### 10 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

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The changes from the previous version are highlighted in grey.

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