UMMt AMIES *Revolution*

Transport des mycoplasmes urogénitaux
UMMt AMIES (2.6 mL) 50 tests (REF 00083)

CPB 0412_FR-2018-03

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Transport of urogenital mycoplasmas
UMMt AMIES (2.6 mL) 50 tests (REF 00083)

CPB 0412_EN-2018-03

For *in vitro* diagnostic use only, for professional use only

1 - INTENDED USE

The UMMt AMIES *Revolution* kit is a complement of UMMt AMIES (2.6mL) medium for the MYCOFAST *Revolution* 2 AMIES kit (REF 00081) or for the COMPLÉMENT MYCOFAST *Revolution* 2 kit (REF 00082) or the MYCOFAST Screening *Revolution* kit (REF 00063).

The UMMt AMIES (2.6mL) medium allows the transport of *Ureaplasma urealyticum* (Uu) and *Mycoplasma hominis* (Mh) present in samples prepared with AMIES transport medium or in a universal transport medium for viruses, chlamydiae, mycoplasmas and ureaplasmas.

2 - PRINCIPLE

The mycoplasmas are fragile bacteria lacking a cell wall and are thus sensitive to variations of osmotic pressure(1, 3). The samples must be discharged into an isotonic medium containing sodium chloride, such as the UMMt transport medium. Moreover, mycoplasmas are sensitive to the products of catabolism. A strong concentration of catabolic products can induce auto-intoxication.

Inoculation in a liquid medium produces a dilution effect.

3 - REAGENT

UMMt AMIES (2.6 mL): Vial of liquid medium containing Mycoplasma broth, antibiotics and a preservative. pH: 6.0 ± 0.1.

4 - PRECAUTIONS

- The reagents are intended solely for *in vitro* use and must be handled by authorised personnel.
- The patient samples and inoculated reagents are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- Reagents containing raw materials of animal origin must be handled with caution.
- Do not use reagents after the expiry date.
- Do not use reagents that have been damaged or that have been poorly conserved before use.

5 - SAMPLE COLLECTION AND HANDLING

Sample collection

Endocervical and vaginal sample collection:

Only use the swab provided with the transport medium.

The cervix should be carefully cleaned with a swab, to remove secretions, before collecting the sample with a new swab. As mycoplasmas adhere strongly to mucous cells, the mucous lining should be vigorously swabbed to obtain a rich specimen.

Urethral sample collection:

Only use the swab provided with the transport medium.

Clean the meatus and swab or scrape the area to obtain cells.

Inoculation and transport of the specimens

Inoculate the UMMt AMIES (2.6mL) medium with 300µL of AMIES transport medium or a universal transport medium for viruses, chlamydiae, mycoplasmas and ureaplasmas.

Storage in UMMt AMIES (2.6mL) medium

The inoculated UMMt AMIES (2.6mL) medium may be kept for 20 hours at room temperature (18-25 °C) or 56 hours at 2-8 °C.

For storage during 3 days at -20 °C, first add two drops of "MYCOPLASMA Stabilizer".

6 - PREPARATION AND STORAGE OF REAGENTS

The UMMt AMIES (2.6mL) medium is stable when stored at 2-8 °C, in its original packaging until the expiry date shown on the kit.

Before opening, the UMMt AMIES (2.6mL) medium may be stored temporarily (3 months) at room temperature (18-25 °C).

The UMMt AMIES medium is ready-to-use.

7 - MATERIAL REQUIRED BUT NOT PROVIDED

- Sample collection materials (swabs provided with the AMIES transport medium or a universal transport medium for viruses, chlamydiae, mycoplasmas and ureaplasmas), pipettes and tips
- MYCOPLASMA Stabilizer (REF 00064)
- Pipettes and tips
- Freezer at -20 °C
- Incubator at 37 °C ± 1 °C
- Waste container for contaminated waste

8 - METHOD

Refer to the usage instructions of the MYCOFAST *Revolution* 2 AMIES or MYCOFAST Screening *Revolution* kits.

9 - PERFORMANCE

Refer to the performance data of the MYCOFAST *Revolution* 2 AMIES or MYCOFAST Screening *Revolution* kits.

10 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

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