# ELITex Bicolor albi-dubli

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Fax : 04 94 88 55 22

# ELITex Bicolor albi-dubli

Latex slide agglutination test for the identification of *Candida* of the *C. albicans* / *C. dubliniensis* group

**60 tests**  
(Ref. 44500)

**360 tests**  
(Ref. 44501)

8000030-en-2014-06

## 1 - AIM

ELITex Bicolor albi-dubli is a slide agglutination test for the rapid identification of *Candida* strains (Group *C. albicans* / *C. dubliniensis*), directly from isolated colonies.

Each kit, Ref. 44500, allows 60 tests to be carried out.

Each kit, Ref. 44501, allows 360 tests to be carried out.

## 2 - INTRODUCTION

Yeast of the genus *Candida* can be responsible for cutaneous candidiasis, mucosal candidiasis, and candidemias or invasive candidiasis.

*Candidas* are usually commensal yeast of the digestive and urogenital mucosae. They become pathogenic only when favourable conditions arise in the host. Among the factors that favour *Candida* infection, there are intrinsic physiological factors or pathologies (advanced age, pregnancy, diabetes, immunological deficiencies and malignant diseases). As for extrinsic factors, these are primarily iatrogenic in nature. The prevalence of candidiasis has increased considerably over the last twenty years, as a consequence of the emergence of pathologies such as AIDS, the generalization of antibiotic treatments and oral contraceptives, the development of therapeutic immunosuppressive agents, parenteral alimentation, and the multiplication of aggressive investigative methods and surgical procedures. For example, the candidemias account for approximately 10% of nosocomial infections; a percentage that can reach as high as 20% according to certain studies. Moreover, the prognosis for the candidiasis is poor, with mortality in affected patients oscillating between 38 and 50%.

*C. albicans* is the species that is most often isolated. It accounts for 60 to 80% of clinical isolates.

*C. albicans* identification can be accomplished at the same time as isolation on chromogenic media or by using phenotypical characteristics such as the capacity to produce germinative tubes. Biochemical tests can also be used, such as the carbon auxanogram (study of the assimilation of monosaccharides as a source of carbon and energy) and the zymogram (study of the use of monosaccharides in anaerobiosis). Their disadvantage is that there is a long 24 to 48 hours wait before the tests can be read. *C. albicans* is the most pathogenic *Candida* species. Its rapid identification is thus essential to the choice of therapeutic regimen.

ELITex Bicolor albi-dubli is a rapid test that enables the identification of *Candida* of the *C. albicans* and *C. dubliniensis* group. This test uses latex particles sensitized with a monoclonal antibody specific for an antigen that is common to these two species. Differentiation between these two species is possible via the use of the ELITex Bicolor *dubliniensis* test.

## 3 - PRINCIPLE

ELITex Bicolor albi-dubli test is carried out with the aid of 2 reagents:

- a **TEST LATEX** reagent (brown), made up of red dyed latex particles in suspension with a green counter dye. These particles are sensitized with a monoclonal antibody specific for a ***C. albicans* / *C. dubliniensis*** antigen that is present primarily in the yeast cell wall.

- a **R-DIS** reagent, containing enzymes that are able to dissociate the cells and to uncover the intraparietal antigen that is recognized by the monoclonal antibody.

When one adds the **TEST LATEX** reagent to colonies of ***C. albicans*** or ***C. dubliniensis*** previously placed in suspension with the **R-DIS** reagent, agitation of the suspension results in coagglutination between the blastospores, carrying the antigen that is recognized by the monoclonal antibody, and the sensitized latex particles.

A positive reaction results in a red agglutination on a more or less intense green background. With yeast colonies other than *C. albicans* and *C. dubliniensis*, no agglutination is observed and the suspension remains homogeneous and of a uniform brown colour.

Handling is simple and fast, with results within 5 minutes.

## 4 - REAGENTS AND MATERIAL

Description - Kit of 60 tests (Ref. 44500)	Quantity
<b>TEST LATEX:</b> vial of 1 mL of sensitized latex	1
<b>R-DIS:</b> Vial of lyophilized dissociating reagent, to be reconstituted with 0.45 mL of distilled water	3
<b>TEST CARD:</b> Disposable reaction card	8
<b>STICK:</b> Disposable mixing sticks	60
<b>DROPPER:</b> Dropper	1

Description - Kit of 360 tests (Ref. 44501)	Quantity
<b>TEST LATEX:</b> vial of 3 mL of sensitized latex	2
<b>R-DIS:</b> Vial of lyophilized dissociating reagent, to be reconstituted with 2.5 mL of distilled water	3
<b>TEST CARD:</b> Disposable reaction card	50
<b>STICK:</b> Disposable mixing sticks	360
<b>DROPPER:</b> Dropper	1

## 5 - PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only and must be handled by authorized personnel.
- Tests are for a single use only.
- The **TEST LATEX** reagent contains raw materials of animal origin and must be handled with caution.
- Patient samples are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- The reagents contain sodium azide (< 0.1 %).
- Do not use reagents after the expiry date.
- Do not use reagents from different batch numbers.
- Prior to use, allow the reagents to reach room temperature.
- Carefully shake the **TEST LATEX** reagent before use.
- When dispensing the **TEST LATEX** reagent, make sure that the dropper is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery volumes.

## 6 - SAMPLE COLLECTION

The test can be carried out:

- either directly from a 24 to 48 h primary culture (Sabouraud agar, blood agar...);
- or after re-seeding on a Sabouraud agar (24 to 48 h culture).

## 7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

All reagents stored at 2-8°C, in their original packaging, are stable until the expiry date indicated on the box. Do not freeze.

After reconstitution, the **R-DIS** reagent can be preserved for 45 days if it is stored at 2-8°C.

Do not expose the reagents to a strong light.

## 8 - MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured
- Distilled water - Pasteur pipette or wire loop - Contaminated waste container

## 9 - METHOD

Allow the reagents to reach room temperature before use.

### 1. Reconstitution of the R-DIS reagent

- Remove the external capsule then raise the stopper slightly without removing it completely so that air can enter the vial. Then, remove the stopper whilst taking care not to spill any of the lyophilizate. Add the exact volume (as indicated in paragraph 4) of distilled water to the vial.
- Close the vial and shake the contents by turning the vial upside down in order to dissolve the lyophilizate.

### 2. Realization of the test on a slide

- For each culture to be tested, add 20 µL of **R-DIS** reagent into separate circle on the slide.
- Using a Pasteur pipette or a wire loop, take **3 - 4 colonies**.
- Transfer the colonies** in a drop of **R-DIS** reagent and **spread them** to cover the entire surface of the circle.
- After ensuring that the **TEST LATEX** reagent is well mixed, add with the aid of the dropper provided, 1 drop of reagent to the circle.
- Using a disposable stick, mix and spread them to cover the entire surface of the circle until a homogeneous suspension is obtained.
- Manually, gently rock the slide for **5 minutes**. A red agglutination may appear on a green background.

## 10 - READING

**Positive reaction:** Formation of a red agglutination on a more or less intense green background.

**Negative reaction:** Absence of agglutination. The suspension remains homogeneous and brown.

## 11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
<b>POSITIVE</b>	The tested strain belongs to <b><i>Candida albicans</i></b> species or <b><i>Candida dubliniensis</i></b> species. Differentiation of these two species is possible with the <b>ELITex Bicolor <i>dubliniensis</i></b> test.
<b>NEGATIVE</b>	The tested strain is <b>not <i>Candida albicans</i></b> . The tested strain is <b>not <i>Candida dubliniensis</i></b> .

## 12 - CAUSES OF ERROR AND TEST LIMITS

- Certain yeast strains such as *C. parapsilosis*, whose culture is difficult to dissociate, cause the formation of white or red aggregates, but with the absence of a green background (red or brown background). The reaction is thus negative.

- In all cases, it is necessary that the clinical, epidemiologic and biological data are taken fully into consideration before establishing the final diagnosis.

## 13 - PERFORMANCE

ELITex Bicolor albi-dubli consists of latex particles sensitized with monoclonal antibodies from the "LIB-3H8" clone, developed by "Micro Seccion Departamental 1" of the faculty of pharmacy at the University of Valencia in collaboration with the "Laboratorio Integrado de Bioingeniería (LIB) of Universidad Politécnica de Valencia", Spain. It is this monoclonal antibody that ensures the sensitivity and specificity of the reaction.

## A. Sensitivity compared to the Blastèse test (germ tube test)

Protocol:

192 identified strains or isolates of ***C. albicans*** or ***C. dubliniensis*** (ID.32C system, BioMérieux) were tested with the **ELITex Bicolor albi-dubli** and the Blastèse (germ tube test) tests (Culture 24h at 37°C on Sabouraud + incubation of the blastospores for 3h at 37°C on human serum + reading of the tests under an inverted microscope).

Results:	Positive	%	Negative	%
Blastèse test	176	91.66	16	8.33
<b>ELITex Bicolor albi-dubli</b>	192	100	0	
ID.32C	192		0	

Conclusion:

**ELITex Bicolor albi-dubli** has a better sensitivity than the Blastèse test (germ tube test).

## B. Sensitivity/specificity compared to the ID.32C identification system (14)

Protocol:

2404 yeast isolates, of which 1747 were ***Candida*** of the ***C. albicans*/*C. dubliniensis*** group and 657 other yeast that were not ***C. albicans*** and ***C. dubliniensis*** (identified by the ID.32C system, BioMérieux) were tested with **ELITex Bicolor albi-dubli**, after culture on Sabouraud agar at 37°C for 24 or 48 hours.

3 or 4 colonies were taken and tested with **ELITex Bicolor albi-dubli**.

Results:	IDENTIFICATION ID.32C	
	<i>C. albicans</i> or <i>C. dubliniensis</i> (1747)	Non <i>C. albicans</i> and non <i>C. dubliniensis</i> yeast (657)
<b>ELITex Bicolor albi-dubli</b>	1744 (99.8%)	4*
	3	653 (99.4%)

\* 4 yeast were identified with the ID.32C as being:

- *Candida parapsilosis*
- *Candida glabrata*
- *Candida kefyr*
- *Candida guilliermondii*

Conclusion:

The results indicate:

- a very good correlation (99.7%) between **ELITex Bicolor albi-dubli** and the ID.32C identification system;
- a sensitivity and a specificity greater than 99% for the **ELITex Bicolor albi-dubli** product.

## 14 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use. If the **TEST LATEX** reagent is spilled or the work area is contaminated by colonies, clean using bleach and absorbent paper.

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