ELITex Bicolor dubliniensis

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ELITex Bicolor dubliniensis

Slide agglutination test for the rapid identification of *Candida dubliniensis*

100 tests
(Ref. 44502)

8000020-en-2012-01

1 - AIM

ELITex Bicolor *dubliniensis* is a slide agglutination test for the rapid identification of *Candida dubliniensis*, directly from isolated colonies. Each kit allows 100 tests to be carried out.

2 - INTRODUCTION

Yeast of the genus *Candida* can be responsible for cutaneous candidiasis, mucosal candidiasis, and candidemias or invasive candidiasis.

Candidas are usually commensal yeast of the digestive and urogenital mucosae. They become pathogenic only when favourable conditions arise in the host. Among the factors that favour *Candida* infection, there are intrinsic physiological factors or pathologies (advanced age, pregnancy, diabetes, immunological deficiencies and malignant diseases). As for extrinsic factors, these are primarily iatrogenic in nature. The prevalence of candidiasis has increased considerably over the last twenty years, as a consequence of the emergence of pathologies such as AIDS, the generalization of antibiotic treatments and oral contraceptives, the development of therapeutic immunosuppressive agents, parenteral alimentation, and the multiplication of aggressive investigative methods and surgical procedures. For example, the candidemias account for approximately 10% of nosocomial infections; a percentage that can reach 20% according to certain studies. Moreover, the prognosis for the candidiasis is poor, with mortality in affected patients oscillating between 38 and 50%.

C. albicans is the species that is most often isolated. It accounts for 60 to 80% of clinical isolates. However, other species such as *C. dubliniensis* (a species very close to *C. albicans*), *C. glabrata*, *C. tropicalis*, *C. krusei*, *C. lusitanae*, *C. parapsilosis*, and *C. lipolytica*, are increasingly being reported as mycotic organisms.

Since 1995, a new *Candida* has been identified: *C. dubliniensis*. This species, close to *C. albicans*, shares with *C. albicans* many phenotypic characteristics, such as capacity to produce germ tubes and chlamydo-spores.

These two species can easily be differentiated with molecular biology techniques: however it is time-consuming and requires special equipments.

ELITex Bicolor *dubliniensis* is a rapid test that enables the identification of *C. dubliniensis*, using latex particles sensitized with a specific *C. dubliniensis* monoclonal antibody. This test allows the differentiation between *C. albicans* and *C. dubliniensis*.

3 - PRINCIPLE

ELITex Bicolor *dubliniensis* test is based on the coagglutination of *C. dubliniensis* blastospores with blue latex particles in suspension in a red dye. The latex particles are sensitized with a monoclonal antibody which specifically detects a *C. dubliniensis* surface antigen.

The dissociation of *C. dubliniensis* colonies in the TEST LATEX reagent leads to an agglutination between blastospores and the latex particles. This positive reaction results in a blue agglutination (that could progressively spread towards the periphery) on a red background. The agglutination is visible to the naked eye.

With non-*C. dubliniensis* yeast colonies, and particularly with *C. albicans* colonies, there is no agglutination. The suspension remains homogeneous and purple in colour.

Handling is simple and fast, with results within 5 minutes.

4 - REAGENTS AND MATERIAL

Description	Quantity
TEST LATEX: 2 mL vial of sensitized latex	1
TEST CARD: Disposable reaction slide	13
STICK: Disposable mixing stick	100

5 - PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only and must be handled by authorized personnel. Tests are for single use only.
- The TEST LATEX reagent contains raw materials of animal origin and must be handled with caution.
- The TEST LATEX reagent contains sodium azide (< 0.1 %).
- Patient samples are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- Do not use TEST LATEX reagent after the expiry date.
- Prior to use, allow the TEST LATEX reagent to reach room temperature.

- Carefully shake the TEST LATEX reagent before use.

6 - SAMPLE COLLECTION

The test can be carried out:

- either directly from a 24, 48 to 72 h primary culture (Sabouraud agar, Candiselect®, CHROMagar Candida®, Candida ID®, blood agar...);
- or after re-seeding on a Sabouraud agar (24, 48 to 72 h culture).

7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

The reagent is ready to use.

The reagent stored at 2-8 °C, in its original packaging, is stable until the expiry date indicated on the box. Do not freeze.

Do not expose the TEST LATEX reagent to strong light.

8 - MATERIAL REQUIRED BUT NOT SUPPLIED

- Automatic pipette(s) with a pipetting volume adapted to the volume that will be measured.
- Pasteur pipette or wire loop
- Contaminated waste container

9 - METHOD

Allow the reagent to reach room temperature before use.

- For each culture to be tested add 20 µL of homogenized TEST LATEX reagent into separate circle on the slide.
- Using a Pasteur pipette or a wire loop or a stick, take 2-3 colonies from a 24, 48 to 72 hours culture.
- Transfer the colonies in the drop of TEST LATEX reagent and spread them to cover the entire surface of the circle until a homogeneous suspension is obtained.
- Manually, gently rock the slide for 3 to 5 minutes. A blue agglutination on a pink or red or purple background may appear.

10 - READING

Positive reaction: Formation of a blue agglutination on a pink or red or purple background, visible to the naked eye, which can spread towards the periphery forming a blue edge around a central pink or purple background.

Negative reaction: Absence of agglutination. The suspension remains homogeneous and purple.

11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	The tested strain belongs to <i>Candida dubliniensis</i> species.
NEGATIVE	The tested strain is not <i>Candida dubliniensis</i> ⇒ If the strain was previously identified as belonging to the group <i>Candida albicans</i> / <i>Candida dubliniensis</i> group (chromogenic media, ELITex Bicolor albi-dubli...), a negative result allows confirming the strain is <i>Candida albicans</i> .

12 - CAUSES OF ERROR AND TEST LIMITS

- To detect a possible association *Candida albicans* / *Candida dubliniensis*, it is recommended to realize from one culture several tests.

- Certain non-*C. dubliniensis* yeast strains (for example *C. parapsilosis*), whose culture is difficult to dissociate, cause the formation of white aggregates which cannot be confused with blue agglutination. The reaction is thus negative.

- In all cases, it is necessary that the clinical, epidemiologic and biological data are fully taken into consideration before establishing the final diagnosis.

13 - PERFORMANCE

ELITex Bicolor *dubliniensis* consists of latex particles sensitized with monoclonal antibodies that ensure the specificity and sensitivity of the reaction. A study, carried out on 200 yeast isolates, compared of ELITex Bicolor *dubliniensis* performances to molecular biology technique and to the ID 32C system (BIOMERIEUX). The results indicate sensitivity of 98 % and specificity of 100 %. The high sensitivity of the test allows performing the reaction with 1 to 10 strains without zone phenomenon.

14 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use. If the TEST LATEX reagent is spilled or the work area contaminated by colonies, clean using bleach and absorbent paper.

15 - BIBLIOGRAPHY

- ADOU-BRYN, K., DOUCHET, C., FERRER, A., GRIMAUD, L., ROBERT, R. AND RICHARD-LENOBLE, D., 2003. Morphological features of *Candida dubliniensis* on a modified Pal's medium. Preliminary study. *J. Mycol. Méd.*, 13 : 99-103.

- AL-MOSAID, A., SULLIVAN, D., COLEMAN, D.C., 2003. Differentiation of *Candida dubliniensis* from *Candida albicans* on Pal's agar. *J. Clin. Microbiol.*, 41: 4787-4789.
- AL-MOSAID, A., SULLIVAN, D., SALKIN, I. F., SHANLEY, D., COLEMAN, D. C., 2001. Differentiation of *Candida dubliniensis* from *Candida albicans* on Staib agar and caffeic acid-ferrocitrate agar. *J. Clin. Microbiol.*, 39 : 323-327.
- BIKANDI, J., SAN MILLÁN, R., MORAGUES, M. D., CEBAS, G., CLARKE, M., COLEMAN, D. C., SULLIVAN, D. J., QUINDÓS, G., PONTÓN, J., 1998. Rapid identification of *Candida dubliniensis* by indirect immunofluorescence based on differential localization of antigens on *C. dubliniensis* blastospores and *Candida albicans* germ tubes. *J. Clin. Microbiol.*, 36 : 2428-2433.
- COLEMAN, D., SULLIVAN, D., BENNET, G. P., MORAN, G., BARRY, H., SHALEY, D., 1997a. Candidiasis, the emergence of a novel species, *Candida dubliniensis*. *AIDS*, 11 : 557-567.
- COLEMAN, D., SULLIVAN, D., HAYENS, M., HENMAN, M., SHANLEY, D., BENNET, D., 1997b. Molecular and phenotypic analysis of *Candida dubliniensis*, a recently identified species linked with oral candidiasis in HIV-infected and AIDS patients. *Oral Dis*, 3 : 96-101.
- DONNELLY, S. M., SULLIVAN, D. J., SHANLEY, D. B., COLEMAN, D. C., 1999. Phylogenetic analysis and rapid identification of *Candida dubliniensis* based on analysis of *ACT1* intron and exon sequences. *Microbiology*, 145 : 1871-1882.
- GUTIERREZ, J., MORALES, P., GONZALES, M. A., QUINDOS, G., 2002. *Candida dubliniensis*, a new fungal pathogen : review. *J. Basic Microbiol.*, 42 : 207-227.
- JOLY, S., PUJOL, C., RYSZ, M., VARGAS, K., SOLL, D. R., 1999. Development and characterization of complex DNA fingerprinting probes for the infectious yeast *Candida dubliniensis*. *J. Clin. Microbiol.*, 37 : 1035-1044.
- KURZAI, O., HEINZ, W. J., SULLIVAN, D. J., COLEMAN, D. C., FROSCH, M., MUHLSCHEGEL, F. A., 1999. Rapid PCR test for discriminating between *Candida albicans* and *Candida dubliniensis* isolates using primers derived from the pH-regulated PHR1 and PHR2 genes of *Candida albicans*. *J. Clin. Microbiol.*, 37 : 1587-1590.
- MAROT-LEBLOND, A., GRIMAUD, L., NAIL, S., BOUTERIGE, S., APAIRE-MARCHAIS, V., SULLIVAN, D. J., ROBERT, R., 2000. New monoclonal antibody specific for *Candida albicans* germ tube. *J. Clin. Microbiol.*, 38 : 61-67.
- MEYER, W., MASZEWSKA, K., SORRELL, T. C., 2001. PCR fingerprinting, a convenient molecular tool to distinguish between *Candida dubliniensis* and *Candida albicans*. *Med. Mycol.*, 39 : 185-193.
- ODDS, F. C., NUFFEL, L. V., DAMS G., 1998. Prevalence of *Candida dubliniensis* isolates in a yeast stock collection. *J. Clin. Microbiol.*, 36 : 2869-2873.
- PARK, S., WONG, M., MARRAS, S. A., CROSS, E. W., KIEHN, T. E., CHATURVEDI, V., TYAGI, S., PERLIN, D. S., 2000. Rapid identification of *Candida dubliniensis* using a species-specific molecular beacon. *J. Clin. Microbiol.*, 38 : 2829-2836.
- PELTROCHE-LLACSAHUANGA, H., SCHMIDT, S., SEIBOLD, M., LUTTICKEN, R., HAASE, G., 2000. Differentiation between *C. dubliniensis* and *C. albicans* by fatty acid methyl ester analysis using gas liquid chromatography. *J. Clin. Microbiol.*, 38 : 3696-3704.
- PINJON, E., SULLIVAN, D., SALKIN, I., SHANLEY, D., COLEMAN, D. C., 1998. Simple inexpensive, reliable method for differentiation of *Candida dubliniensis* from *Candida albicans*. *J. Clin. Microbiol.*, 6 : 2093-2095.
- STAIB, P., MORSCHHÄUSER, J., 1999. Chlamydo-spore formation on Staib agar as species-specific characteristic of *Candida dubliniensis*. *Mycoses*, 42 : 521-524.
- SULLIVAN, D., COLEMAN, D., 1998. *Candida dubliniensis*, characteristics and identification. *J. Clin. Microbiol.*, 36 : 329-334.
- SULLIVAN, D. J., MORAN, G., DONNELLY, S., GEE, S., PINJON, E., MCCARTAN, B., SHANLEY, D. B., COLEMAN, D. C., 1999. *Candida dubliniensis*, An update. *Rev. Iberoam. Mycol.*, 16 : 72-76.
- SULLIVAN, D. J., WESTERNENG, T. J., HAYNES, K. A., BENNET, D. E., COLEMAN, D. C., 1995. *Candida dubliniensis* sp. nov., phenotypic and molecular characterization of a novel species associated with oral candidosis in HIV-infected individuals. *Microbiology*, 141 : 1507-1521.
- WILLIAMS, D. W., COULTER, W. A., WILSON, M. J., POTTS, A. J., LEWIS, M. A. O., 2001. Identification of *Candida dubliniensis*, based on ribosomal DNA sequence analysis. *Br. J. Biomed. Sci.*, 58 : 11-16.
- MAROT-LEBLOND, A., BEUCHER, B., DAVID, S., NAIL-BILLAUD, S., ROBERT, R., 2006. Development and evaluation of a rapid latex agglutination test using a monoclonal antibody to identify *Candida dubliniensis* colonies. *J. Clin. Microbiol.*, 44 : 138-142.
- SAHAND, I.H., MORAGUES, M.D., ROBERT, R., QUINDOS, G., PONTON, J., 2006. Evaluation of Bichro-dubli distinguish *Candida dubliniensis* from *Candida albicans*. *Diagnostic Microbiology and Infectious Disease*, 55 : 165-167.
- CHRYSSANTHOU, E., FERNANDEZ, V., PETRINI, B., 2007. Performance of commercial latex agglutination tests for the differentiation of *Candida dubliniensis* and *Candida albicans* in routine diagnostics. *APMIS*, 115 : 1281-1284.

ELITech MICROBIO

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83870 SIGNES

FRANCE

☎ : 04 94 88 55 00

Fax : 04 94 88 55 22