# ELITex *krusei*

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## ELITex krusei

Latex slide agglutination test for the identification of *Candida krusei*

**30 tests** (Réf. 44504)

**1 - AIM** ELITex *krusei* is a slide agglutination test for the rapid identification of *Candida krusei*, directly from colonies.

Each kit allows 30 tests to be carried out.

## 2 -INTRODUCTION

Yeast of the genus *Candida* can be responsible for cutaneous candidiasis, mucosal candidiasis, and candidemias or invasive candidases.

*Candidas* are usually commensal yeast of the digestive and urogenital mucosae. They become pathogenic only when favourable conditions arise in the host. Among the factors that favour *Candida* infection, there are intrinsic physiological factors or pathologies (advanced age, pregnancy, diabetes, immunological deficiencies and malignant diseases). As for extrinsic factors, these are primarily iatrogenic in nature. The prevalence of candidases have increased considerably over the last twenty years, as a consequence of the emergence of pathologies such as AIDS, the generalization of antibiotic treatments and oral contraceptives, the development of therapeutic immunosuppressive agents, parenteral alimentation, and the multiplication of aggressive investigative methods and surgical procedures. For example, the candidemias account for approximately 10% of nosocomial infections; a percentage that can reach 20% according to certain studies. Moreover, the prognosis for the candidases is poor, with mortality in affected patients oscillating between 38 and 50%.

*C. krusei* is the species that is most often isolated. It accounts for 60 to 80% of clinical isolates. However, other species such as *C. dubliniensis* (a species very close to *C. albicans*), *C. glabrata*, *C. tropicalis*, *C. krusei*, *C. lusitanae*, *C. parapsilosis*, and *C. lipolytica*, are increasingly being reported as mycotic organisms.

It is also found in nature. This yeast is seldom isolated from clinical samples such as haemocultures (2%). However, the mortality rate attributed to *C. krusei* is high (40%), with primary resistance to certain azoles being a problem (fluconazole, itraconazole, ...). The isolation and rapid identification of *C. krusei* is thus essential to the choice of therapeutic regimen. Traditionally, biochemical tests are used in the identification of *C. krusei*. The disadvantage of the carbon auxanogram (study of the assimilation of monosaccharides as a source of carbon and energy) and the zyogram (study of the use of monosaccharides in anaerobiosis), is that there is a long 24 to 48 hour wait before the tests can be read.

**3 -PRINCIPLE** ELITex *krusei* is a rapid identification test of *C. krusei* that uses latex particles sensitized with a specific *C. krusei* monoclonal antibody.

ELITex *krusei* test is based on the coagglutination of a *C. krusei* surface antigen with red latex particles that are sensitized with a monoclonal antibody. The dissociation of *C. krusei* colonies in the ELITex *krusei* reagent leads to an agglutination between blastospores and the latex particles. This positive reaction results in a red agglutination that is visible to the naked eye. With non-*C. krusei* yeast colonies, there is no agglutination and the suspension

## 4 - REAGENTS AND MATERIAL

Remains homogeneous and red. Handling is simple and fast, with results within 5 minutes.

Description	Quantity
TEST LATEX: Dispenser vial of 1,2 mL of sensitized latex	1
TEST CARD: Disposable reaction card	4

## 5 - PRECAUTIONS

The reagent is intended for *in vitro* diagnostic use only and must be handled by authorized personnel.

Tests are for a single use only.

The reagent contains raw materials of animal origin and must be handled with caution.

Patient samples are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.

The reagent contains sodium azide (< 0,1 %).

Do not use the reagent after the expiry date.

Prior to use, allow the reagent to reach room temperature.

Carefully shake the reagent before use.

When dispensing the reagent, make sure that the dispenser vial is perfectly vertical. Check for the absence of air bubbles in the drops to ensure constant delivery volumes. Wipe the tip of the dispenser vial before use in order to obtain well gauged drops.

## 6 - SAMPLE COLLECTION

The test can be carried out:

- either directly from a 24 to 48 h primary culture (Sabouraud agar, Candidchrom albicans , chromagar Candida , Albicans ID ...);
- or after re-seeding on a Sabouraud agar (24 to 48 h culture).

Note: It is not recommended that the test be carried out on isolated blood (fresh or cooked) agar colonies. It should be noted that agglutinations can be weak or slow to appear.

## 7 - STABILITY, STORAGE AND PREPARATION OF REAGENTS

The reagent is ready-to-use.

The reagent stored at 2-8°C, in its original packaging, is stable until the expiry date indicated on the box.

Do not freeze.

**Do not expose** the reagent to strong light.

## 8 - MATERIAL REQUIRED BUT NOT SUPPLIED

Pasteur pipette or wire loop

Contaminated waste container

## 9 - METHOD

Allow the reagent to reach room temperature before use.

a. For each culture to be tested, add 1 drop of homogenized **TEST LATEX**

reagent into separate circle on the slide.

b. Using a Pasteur pipette or a wire loop, take **2-3 colonies**.

c. Transfer the colonies in the drop of **TEST LATEX** reagent and spread them to cover the entire surface of the circle until a homogeneous suspension is obtained.

d. Manually, gently rock the slide for **5 minutes**. A red agglutination may appear.

## 10 - READING

**Positive reaction:** Formation of a red agglutination visible to the naked eye.

**Negative reaction:** Absence of agglutination.

The suspension remains homogeneous and red.

## 11 - INTERPRETATION OF RESULTS

RESULT	INTERPRETATION
POSITIVE	The test strain belongs to <i>Candida krusei</i> species
NEGATIVE	The test strain is not <i>Candida krusei</i>

## 12 - CAUSES OF ERROR AND TEST LIMITS

- Certain non-*C. krusei* yeast strains (for example *C. parapsilosis*), whose culture is difficult to dissociate, cause the formation of white aggregates which cannot be confused with red agglutination. The reaction is thus negative.

- In all cases, it is necessary that the clinical, epidemiologic and biological data are taken fully into consideration before establishing the final diagnosis.

## 13 - PERFORMANCE

ELITex *krusei* consists of latex particles sensitized with monoclonal antibodies from the 6B3 clone (developed by the SR2B - Avuille - FRANCE). It is these monoclonal antibodies that ensure the specificity and sensitivity of the reaction.

A study, carried out on 884 yeast isolates, compared ELITex *krusei* performances to the ID 32C system (BIOMERIEUX). The results indicate sensitivity of 98,8 % and specificity of 100 %.

## 14 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

If the **TEST LATEX** reagent is spilled or the work area is contaminated by colonies, clean using bleach and absorbent paper.

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