ELIchrom *glabrata*

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15 - FÖRTECKNING

- Abi-said D., Anaïsse E., Uzun O., Raad L., Pinzcowski H., and Vartivarian S. 1997. The epidemiology of hematogenous candidiasis caused by different *Candida* species. *Clin. Infect. Dis.* 24:1122-1128.
- Bouchara J., P., Declerck P., Cimon B., Planchenault C., de Gentile L., Chabasse D. 1996. Routine use of CHROMagar *Candida* medium for presumptive identification of *Candida* yeast species and detection of mixed fungal populations. *Clin Microbiol Infect* 2:202-208
- Fenn J. P., Billeldeux E., Segal H., Skodack-Jones L., Padilla PE., Bale M., Carroll K. 1999. Comparison of four methodologies for rapid and cost-effective identification of *Candida glabrata*. *J Clin Microbiol.* 37:3387-9.
- Fidel P. L., Vazquez J. A., Sobel J. D. 1999. *Candida glabrata*: review of epidemiology, pathogenesis, and clinical disease with comparison to *C. albicans*. *Clin Microbiol Rev.* 12:80-96.
- Freydière A. M., Guinet R., and Boiron P. 2001. Yeast identification in the clinical microbiology laboratory: phenotypical methods. *Med. Mycol.* 39:9-33.
- Lopez J., Dalle F., Mantelin P., Mouiroux P., Nierlich A. C., Pacot A., Cuisenier B., Vagner O., and Bonnin A. 2001. Rapid identification of *Candida glabrata* based on trehalose and sucrose assimilation using Rosco diagnostic tablets. *J. Clin. Microbiol.* 39:1172-1174.
- Parant F., Freydière A. M., Gille Y., Boiron P., and Odds F. C. 2001. A 'one minute' trehalase detection test for the identification of *Candida glabrata*. *J. Mycol. Med.* 11:26-31.
- Peltroche-Llacsahuanga H., Schnitzler N., Lütticken R., Haase G. 1999. Rapid identification of *Candida glabrata* by using a Dipstick to detect trehalase-generated glucose. *J Clin Microbiol.* 37:202-5.
- Pfaller M. A., Diekema D. J., Jones R. N., sader H., S., Fluit A., C., Hollis R., J., Messer S., A. 2001. International surveillance of blood stream infections due to *Candida* species: frequency of occurrence and *in vitro* susceptibilities to fluconazole, ravuconazole, and voriconazole of isolates collected from 1997 through 1999 in the SENTRY antimicrobial surveillance program. *J Clin. Microbiol.* 39:3254-9.
- Pfaller M. A., Messer S. A., Hollis R. J., Jones R. N., Doern G. V., Brandt M. E., and Hajjeh R. A. 1999. Trends in species distribution and susceptibility to fluconazole among blood stream isolates of *Candida* species in the United States. *Diag. Microbiol. Infect. Dis.* 33:217-222.
- Senet J-M., Robert R. 1995. Physiopathology of candidosis. *J. Mycol. Med.* 5:145-166.
- Freydière A-M., Robert R., Ploton C., Marot-Leblond A., Monereau F., Vandenesch F. 2003 Aug. Rapid identification of *Candida glabrata* with a new commercial test, GLABRATA RTT. *J Clin Microbiol.* 41(8) : 3861-3.
- Freydière A-M., Perry J-D, Faure O., Willinger B., Tortorano A-M., Nicholson A., Peman J., Verweij P-E. 2004 Oct. Routine use of a commercial test, GLABRATA RTT, for rapid identification of *Candida glabrata* in six laboratories. *J Clin Microbiol.* 42(10) : 4870-2.
- Peman J., Aparisi N., Garcia-Esteban C., Gombardo M., 2004 Jun. Rapid identification of *Candida glabrata* using a new commercial kit *Rev lberoam Micol.* 21(2): 82-4.
- Willinger B., Wein S., Hirschl A-M., Rotter M-L., Manafi M. 2005 Jan. Comparison of a new commercial test, GLABRATA RTT, with a dipstick test for rapid identification of *Candida glabrata*. *J Clin Microbiol.* 43(1) : 499-501.

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