SEROLOGY

Détection quantitative des anticorps anti-*Ureaplasma urealyticum* et anti-*Mycoplasma hominis*

14 tests (REF 00014)
120 tests (REF 00097)

FR-2010-09

1 - BUT

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SEROLOGY

Quantitative detection of antibodies against *Ureaplasma urealyticum* and *Mycoplasma hominis*

14 tests (REF 00014)
120 tests (REF 00097)

EN-2010-09

1 - INTENDED USE

The SEROLOGY kit enables the quantitative detection in serum of the antibodies against the principal urogenital mycoplasma species, *Ureaplasma urealyticum* (U.u.) and *Mycoplasma hominis* (M.h.). Each kit contains 14 or 120 tests.

2 - INTRODUCTION

Mycoplasmas are fragile germs lacking a cell wall. They are demanding in their requirements for growth promoters, endowed with a strong capacity of adherence and most often present as commensal organisms. The species isolated from the urogenital tract, *Ureaplasma urealyticum* and *Mycoplasma hominis* are the most frequently encountered. U.u. and M.h. are sexually transmitted and can be pathogenic in nature.

Given their lack of immunogenicity, serological testing of mycoplasmas does not replace direct diagnosis from patient specimens. The serologic test for urogenital mycoplasmas is especially used to diagnose a deep genital infection, like salpingitis or prostatitis. The technique of metabolic inhibition using lyophilized reagents and reference strains is the most appropriate technique for use by non-specialized laboratories (1).

3 - PRINCIPLE

The quantitative detection of anti-U.u. and/or anti-M.h. antibodies is achieved by a metabolic inhibition technique (2). The reaction is carried out in the wells of a tray containing in the *U. urealyticum* rows, lincomycin and in the *M. hominis* rows, erythromycin. The first 8 wells, of each side of the tray, are inoculated with UMM culture medium containing a pH indicator and seeded beforehand with a known quantity of U.u. and M.h. strains.

The serum to be tested is serially diluted in the tray. The presence of antibodies in the serum causes an inhibition of mycoplasma multiplication and is visualised by a yellow coloration of the medium. However, in the absence of antibodies the mycoplasmas multiply, degrading urea and/or arginine present in the medium, which results in an orangey, red or fuschia coloration of the medium.

For each row of U.u. and M.h. tests there is a growth control (well 1, labelled T). These are carried out in order to verify the viability of the U.u. and M.h. strains, as well as the correct operation of the test.

4 - REAGENTS

Description	Quantity #00014	#00097
UMMt: Vial of 2 mL of liquid mycoplasma base medium	14	120
UMMMyo: Vial of lyophilized mycoplasma growth medium (+ 2 mL UMMt)	14	120
U.u./M.h.: Vial containing lyophilized <i>U. urealyticum</i> and <i>M. hominis</i> strains	14	120
SEROLOGY: 20-well tray wrapped individually in an aluminium sachet	14	120
Positive control: Vial containing a lyophilized solution of known titre for verification of the reagents	1	5

Reconstituted medium (UMMt + UMMMyo) - Composition

Mycoplasma broth 20 g/L, Foal serum 200 mL, Yeast extract 5,8 g/L, Cystine 0,3 g/L, Arginine 9 g/L, Urea 3,6 g/L, Phenol red 0,04 g/L, Antibiotics 10 mL. pH : 6.1 ± 0.1

U.u. / M.h. strains

Lyophilized reference strains: *U. urealyticum* serotype 8 and *M. hominis* PG 21 strain.

SEROLOGY tray (§9.1)

- *U. urealyticum* row: Series of first 8 wells containing a *M. hominis* specific lincomycin inhibitor.

- *M. hominis* row: Series of first 8 wells containing an *U. urealyticum* specific erythromycin inhibitor.

Positive control

Mycoplasma medium with a mixture of specific antibiotics.

5 - PRECAUTIONS

● The reagents in this kit are intended solely for *in vitro* use and must be handled by authorised personnel.

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- Samples and inoculated reagents are potentially infectious and must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.
- The UMMt and UMMMyo media, as well as the positive control contain raw materials of animal origin and must be handled with caution.
- The U.u./M.h. vials contain group 2 bacteria (*M. hominis*) and must be handled with caution. Risk of contamination by a biological agent.
- Do not use reagents after the expiry date.
- Do not use reagents that have been damaged or that have been poorly preserved before use.
- After opening the kit, only remove from the refrigerator the reagents required. The U.u./M.h. strains are sensitive to heat and to variations in temperature.
- A positive result with the SEROLOGY method should not be used in isolation to reach a clinical diagnosis. A doctor, taking into account the biological results and the clinical symptoms should carry this out.

6 - SAMPLE COLLECTION AND HANDLING

- Collect the blood sample in a plain, non-additive, tube.
- Centrifuge for 10 minutes at x2500g.
- Collect the serum.

The serum can be stored for 3 hours at 2 to 8 °C or frozen at -20 °C. Specimens from the same patient taken at an interval of 15 days must be tested simultaneously.

7 - PREPARATION AND STORAGE OF REAGENTS

- The reagents, stored at 2 to 8 °C in their original packaging, are stable until the expiry date shown on the kit.
- The UMMMyo (+ 2 mL UMMt) should be regenerated with the UMMt medium (2 mL).
- The vial of U.u. and M.h. strains should be reconstituted with 1 ml of sterile distilled water. Mix well by gently inverting the vial. Do not conserve after use.
- The positive control vial should be reconstituted with 1 mL of sterile distilled water. Mix well. Do not conserve the vial after use.

8 - MATERIAL REQUIRED BUT NOT PROVIDED

- Sterile pipettes
- Sterile distilled water
- Paraffin oil
- Incubator calibrated to 37 °C
- Waste container for contaminated waste


9 - METHOD

Allow the reagents to reach room temperature (18 to 25 °C) before use.

9.1. Preparation and distribution of the suspension

- Inoculate the regenerated UMMMyo medium with 400 µL of the germ suspension (U.u./M.h. strains).

- Label each tray, lift up the adhesive film, then add 100 µL of prepared reagent to each of the 8 wells of the U.u. row and of the M.h. row. The last 2 wells of each row (empty wells) don't need to be filled.

9.2. Distribution of the serum to be tested (or of the positive control)

- Add 100 µL of serum to the second well of one of the two rows of wells. Serial dilute 100 µL of liquid from well 2 into well 3 and continue in this way for each well through to well 8 (1/128). Discard the extra 100 µL from well 8.
- Repeat the same operation for the second row of wells.
- Before using the kit to test patient serum, the wells of the two rows of the SEROLOGY tray should be inoculated with the positive control.
- Add 2 drops of paraffin oil to each well of the tray
- Cover the tray with the adhesive film.

9.3. Incubation of the tray

- Incubate the SEROLOGY tray for 24 hours (48 hours if necessary) (§10.1) at 37 °C.
- The reagent positive control tray should only be incubated for 24 hours.

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10 - READING AND INTERPRETATION

10.1. Validation

Verify that the contents of the wells are limp.

Do not read the tray unless the growth controls (T wells) are positive:

U.u.: Colour change from yellow to red or fuschia.

M.h.: Colour change from yellow to orangey-red or fuschia.

The incubation times of the two strains are independent and depend upon the colour change of the medium in the T wells of each row.

10.2. Reading of the tray

The presence of antibodies is determined by the absence of a colour change of the medium. In the absence of antibodies the medium becomes orangey-red or fuschia coloured. The strongest inhibited dilution (the last yellow well) represents the titre of antibodies.

10.3. Interpretation

Interpretation must be based upon the testing of two specimens taken at an interval of 15 days and tested simultaneously. Only a variation in antibody titre of these specimens, greater or equal to 2 dilution's, is any indication of a progressive infection.

11 - QUALITY CONTROL

11.1. Reagent positive control

Verification of the reagents must be carried out upon reception of the kit using the positive control vial in the place of serum. Reading of the results is carried out after incubation for 24 hours. The expected titre for the two rows of the tray is specified on the label.

Note: The titre of the kit's positive control, as specified on the label, can vary by one dilution.

11.2. Test positive control

The growth of the strain in the first well (T) of each U.u. and M.h. row validates the test (§10.1).

12 - LIMITS OF THE METHOD

● The patients whose samples are to be tested should not be under antibiotic treatment. The presence of antibiotics in the serum can inhibit the multiplication of the mycoplasmas and lead to misleadingly elevated levels.

● Interpretation must be based upon the testing of two specimens taken at an interval of 15 days and tested simultaneously.

● The serological testing of mycoplasma must not under any circumstances be

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