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2 - **Clinical and Laboratory Standards Institute (CLSI). 2011** Methods for Antimicrobial Susceptibility Testing for Human Mycoplasmas; Approved Guideline. CLSI Document M43-A. Vol.31 - N°19.

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6 - **WAITES KEN B, DONNAM. CRABB, and LYNN B. DUFFY. 2008.** Comparative In Vitro Activities of the Investigational Fluoroquinolone DC-159a and Other Antimicrobial Agents against Human Mycoplasmas and Ureaplasmas. ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, Vol. 52, No. 10, 3776–3778.

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**Urogenital mycoplasma diagnosis**  
**MYCOFAST® Screening Revolution**  
 Screening and differentiation  
 50 tests (REF 00063)

**COMPLEMENT MYCOFAST® Revolution 2**  
 Enumeration, identification and susceptibility testing  
 25 tests (REF 00082)

**UMMt AMIES Revolution**  
 50 tests (REF 00083)

CPB 0396-3\_EN-2018-03

For in vitro diagnostic use only, for professional use only

### 1 - INTENDED USE

MYCOFAST Screening *Revolution* (REF 00063) has been designed for the screening and differentiation of *Ureaplasma urealyticum* / *Ureaplasma parvum* (Uu) and *Mycoplasma hominis* (Mh) in various clinical specimens prepared with AMIES transport medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas. This kit should be used in association with the media contained in the UMMt AMIES *Revolution* kit (REF 00083). In the case of positive screening, analysis can be completed with trays contained in the COMPLEMENT MYCOFAST *Revolution* kit (REF 00062) allowing the enumeration and identification of Uu and/or Mh as well as the antimicrobial susceptibility testing according to the recommendations of the CLSI (Clinical and Laboratory Standards Institute) (2).

### 2 - INTRODUCTION

Mycoplasmas that include several species that have been identified in humans, all belong to the mollicutes class. They differ from other bacteria in their lack of a cell wall and hence a natural resistance to  $\beta$ -lactams, as well as by the presence of a membrane rich in sterol obtained through their adhesion to eukaryotic cells. Since mycoplasmas are relatively fragile, they will only grow in acellular culture in the presence of various growth factors and at an optimal temperature of 37°C (4). Most human mycoplasmas are commensal. *U. urealyticum* and *M. hominis* are the most commonly encountered species that have been isolated from the urogenital tract. *U. urealyticum* species are divided into two biovars: *U. urealyticum* and *U. parvum* (Uu.). Uu. and Mh. can be pathogenic. They are responsible for male genital infections (non-gonococcal urethritis, epididymitis, prostatitis, infertility); female genital infections (bacterial vaginosis, endometritis, salpingitis); fertility problems (chorioamnionitis, post-partum endometritis, preterm birth, spontaneous abortion), neonatal problems (low birth weight, respiratory and neurological infections, bacteremias, abscesses); extragenital infections (septic arthritis, reactive arthritis, other infection loci) (1). The diagnosis of mycoplasma infections depends upon the determination of the pathological threshold, followed by enumeration. The resistance of Uu./Mh. to certain drugs necessitates antimicrobial susceptibility testing (5, 6). The drugs tested and the interpretation criteria are adapted for the treatment of infections caused by mycoplasmas encountered in the urogenital tract or in extragenital sites (2).

### 3 - PRINCIPLE

MYCOFAST Screening *Revolution* is a liquid method based on the ability of Uu. and Mh. to metabolize urea and arginine respectively. The mycoplasma growth results in a colour change of the medium, containing phenol red indicator, from yellow-orange to red. This colour change is due to liberation of ammonia resulting in an alkaline pH of the medium.

Mycoplasma growth thus viewed enables:

- detection and differentiation; then, if a positive result is obtained:
- the enumeration of mycoplasma based on the rate of urea or arginine hydrolysis, which is proportional to the number of germs contained in the sample.
- the identification based on the sensitivity or otherwise of the germ to three antimicrobial agents.
- the Uu and Mh susceptibility testing to antimicrobial agents.

### 4 - REAGENTS

Description	Amount		
	réf 00083	réf 00082	réf 00063
<b>UMMt AMIES</b> : Vial of 2.6 mL mycoplasma broth with antimicrobial agents and preservative solution. pH : 6.0 ± 0.1	50		
<b>MYCOFAST SCREENING Revolution</b> : Divisible tray of 10 wells for 5 tests, individually packed in an aluminium sachet with an integrated desiccant			10
<b>Labels</b> : Sheet of 5 divisible labels			10
<b>S.Mh.</b> : Mycoplasma hominis growth activator (4.5 mL)			1
<b>MYCOFAST Revolution 2</b> : Tray of 24 wells for 1 test, packed in an aluminium sachet with an integrated desiccant		25	
<b>Closing System</b> : Protective translucent plastic lid for MYCOFAST Revolution tray		25	

#### MYCOFAST Screening Revolution tray:

Tray consisting of 5 rows of 2 wells: a *Ureaplasma urealyticum* (Uu) well containing lincomycin and urea and a *Mycoplasma hominis* (Mh) well containing erythromycin and arginine.

#### The MYCOFAST Revolution 2 tray:

Tray, in each of the 24 wells, contains dehydrated mycoplasma culture medium (foal serum, yeast extract, cysteine, arginine, urea, phenol red, antibiotics, pH: 6.1 ± 0.1) and includes 2 separate parts:

- the part intended for enumeration and susceptibility testing of the Uu species (wells identified on the label as black).
- the part intended for enumeration and susceptibility testing of the Mh species (wells identified on the label as red).

#### Diagnosis of Uu species (black part of tray label):

- Wells 1/2/3** : Identification and enumeration of Uu at  $10^3$ ,  $10^4$  and  $\geq 10^5$  CCU/mL (buffered solution and lincomycin inhibiting Mh growth).  
**Wells 4/5** : Evaluation of Uu susceptibility to Levofloxacin (LVX) at 2 / 4  $\mu$ g/mL  
**Wells 6/7** : Evaluation of Uu susceptibility to Moxifloxacin (MXF) at 2 / 4  $\mu$ g/mL  
**Wells 8/9** : Evaluation of Uu susceptibility to Erythromycin (ERY) at 8 / 16  $\mu$ g/mL  
**Wells 10/11** : Evaluation of Uu susceptibility to Tetracycline (TET) at 1 / 2  $\mu$ g/mL  
**Wells 12/13** : Evaluation of Uu susceptibility to Doxycycline (DOX) at 1 / 2  $\mu$ g/mL

#### Diagnosis of Mh species (red part of tray label):

- Wells 14** : Identification and Enumeration of Mh at  $\geq 10^4$  CCU/mL  
**Wells 15/16** : Evaluation of Mh susceptibility to Doxycycline (DOX) at 4 / 8  $\mu$ g/mL  
**Wells 17/18** : Evaluation of Mh susceptibility to Levofloxacin (LVX) at 1 / 2  $\mu$ g/mL  
**Wells 19/20** : Evaluation of Mh susceptibility to Moxifloxacin (MXF) at 0.25 / 0.5  $\mu$ g/mL  
**Wells 21/22** : Evaluation of Mh susceptibility to Clindamycin (CLI) at 0.25 / 0.5  $\mu$ g/mL  
**Wells 23/24** : Evaluation of Mh susceptibility to Tetracycline (TET) at 4 / 8  $\mu$ g/mL

		15	16	17	18	19	20	21	22	23	24		
		DOX		LVX		MXF		CLI		TET			
		4	8	1	2	0.25	0.5	0.25	0.5	4	8		
14	Mh $\geq 10^4$	<b>MYCOFAST® Revolution 2</b>										2	13
	Uu $10^3$											1	
		Uu $10^4$	Uu $\geq 10^5$	2	4	2	4	8	16	1	2		
		2	3	4	5	6	7	8	9	10	11		

### 5 - PRECAUTIONS

The reagents are intended solely for in vitro use and must be handled by authorized personnel.

The patient samples and inoculated reagents are potentially infectious; they must be handled with caution, in observance of hygiene rules and the current regulations for this type of product in the country of use.

Reagents containing raw materials of animal origin must be handled with caution.

Do not use reagents after the expiry date.

Do not use reagents that have been damaged or that have been poorly conserved before use.

A positive result with the MYCOFAST method indicates colonization by urogenital mycoplasmas, but cannot alone be used to make a clinical diagnosis. This must be made by a doctor and is a function of the biological results and clinical signs.

### 6 - SAMPLE COLLECTION AND HANDLING

#### 6.1 Sample collection

##### Cervicovaginal sample collection

Only use the swab provided with the transport medium.

The cervix should be carefully cleaned with a swab, to remove secretions, before collecting the sample with a new swab. As mycoplasmas adhere strongly to mucous cells, the mucous lining should be vigorously scrubbed to obtain a rich specimen.

##### Urethral sample collection :

Only use the swab provided with the transport medium.

Clean the meatus and swab or scrape the area to obtain cells.

#### T6.2. Transport of samples

Transport in AMIES medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas:

Refer to the manufacturer's operating instructions

##### Transport in UMMt AMIES medium:

Place 300  $\mu$ L of the inoculated transport medium in a vial of UMMt AMIES medium.

#### 6.3 Conservation of the samples

Conservation in AMIES medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas

Refer to the manufacturer's operating instructions

##### Conservation in UMMt AMIES medium

The inoculated UMMt AMIES medium may be kept for 20 hours at room

temperature (18-25°C) or 56 hours at 2-8°C. For storage during 3 days at -20°C, first add 2 drops of "MYCOPLASMA Stabilizer".

## 7 - PREPARATION AND STORAGE OF REAGENTS

All the reagents are ready-to-use. The vials may be stored at 2-8 °C, in their original packaging until the expiry date shown on the kit. Should only one or two, three, or four rows of (Uu) (Mh) wells be used, the remaining MYCOFAST Screening *RevolutioN* tray may be stored for 4 weeks at 2-8 °C in its original packaging and hermetically resealed. The UMMt medium may be stored temporarily (3 month) at room temperature but is more stable at 2 - 8°C. The S.Mh supplement is stable for 3 months after opening. Do not freeze the reagents contained in the kit.

## 8 - MATERIAL REQUIRED BUT NOT PROVIDED

Sample collection materials (swabs provided with the AMIES transport medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas)  
MYCOPLASMA Stabilizer (REF 00064)  
Pipettes and tips  
Waste container for contaminated waste  
Mineral oil  
Incubator calibrated at 37°C ± 1°C

## 9 - METHOD

**Allow the reagents to reach room temperature (20-30 minutes).**

### 9.1 SCREENING - MYCOFAST Screening *RevolutioN* tray

- Prepare as many rows of wells as samples to be tested. - If required separate one or several rows of (Uu)/(Mh) wells with the aid of the marks found on the tray.

#### 9.1.1 Inoculation of the UMMt AMIES vial

If the sample has been transported in UMMt AMIES medium inoculated with 300µl of AMIES medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas; **then pass directly to step 9.1.2.**

If the sample has been transported in AMIES medium or a universal transport medium for viruses, chlamydias, mycoplasmas and ureaplasmas, then transfer 300 µl of this transport medium into a vial of UMMt AMIES.

#### 9.1.2 Inoculation of Uu/Mh wells

- Distribute successively:  
(Uu) well: 100 µL of seeded UMMt medium.  
(Mh) well: 100 µL of seeded UMMt medium.  
50 µL of Mh supplement.  
- Add 2 drops of mineral oil to the two wells.  
- Cover the wells with the divisible labels and label the sample in order to identify it.  
- **Store excess seeded UMMt medium at 2-8°C** in order to continue the analysis in case of positive screening.

#### 9.1.3 Incubation of Uu/Mh wells

Incubate the wells of the tray for 24 hours at 37°C ± 1°C. Tray incubation can be extended for up to 48 hours only in the case of liquid samples that are negative after 24 hours.

#### 9.1.4 Reading and interpretation of Uu/Mh wells

- Check that the 2 (Uu) (Mh) wells are limpid. A cloudy appearance in a well indicates bacterial contamination. In this case repeat the analysis. Observe the colour change of the medium in the Uu and Mh wells:  
Uu wells are orangey or red: Presence of *Ureaplasma urealyticum*  
Mh wells are orangey or red: Presence of *Mycoplasma hominis*  
Uu / Mh wells are yellow: Absence of mycoplasma

**In the case of positive screening continue diagnosis with the MYCOFAST *RevolutioN* 2 tray.**

## 9.2 ENUMERATION, IDENTIFICATION AND SUSCEPTIBILITY TESTING

### 9.2.1 Inoculation of the MYCOFAST *RevolutioN* 2 tray

Remove the adhesive film by pulling on the tab and add the following to the wells of each row:  
Wells 1-24 100 µL of inoculated UMMt AMIES medium  
Wells 1-24 2 drops of mineral oil  
Cover the seeded tray with the "closing system".  
Label the sample.  
**Store excess UMMt AMIES medium at 2-8°C** for at least 48 hours for possible verification.

### 9.2.2 Incubation of the tray

Incubate the tray at 37°C ± 1°C for 24 hours.  
For Uu and Mh enumeration, read the results in 24 hours. Tray incubation can be extended for up to 48 hours only in the case of liquid samples that are negative after 24 hours.

### 9.2.3 Reading and interpretation

Check that all the wells in the row are limpid. A cloudy appearance in a well indicates bacterial contamination. In this case repeat the analysis. The results are read by the colour obtained in the different wells. Urogenital Mycoplasma growth is indicated when the medium turns red (alkaline). The medium remains yellow when no growth of urogenital mycoplasma occurs. An orangey coloration should be considered as a positive test (rate limit).  
For the interpretation of the results refer to the results sheet.

### Enumeration (wells 1, 2, 3 and 14)

Mark the wells that have turned orange or red and interpret:  
1 Uu value 10<sup>3</sup> CCU/mL  
1 and 2 Uu value 10<sup>4</sup> CCU/mL  
1, 2 and 3 Uu value ≥10<sup>5</sup> CCU/mL  
14 Mh value ≥ 10<sup>4</sup> CCU/mL

The pathological role of mycoplasmas in urogenital infections is subject to interpretation according to specific recommendations (1,3,7). The pathological thresholds usually quoted for *U. urealyticum* are: ≥10<sup>4</sup> CCU/mL for a urethral specimen or endotracheal specimen, ≥10<sup>3</sup> CCU/mL in a first urine stream or sperm (although a new local recommendation mentions a threshold ≥ 10<sup>4</sup> CCU / mL for semen (7)). The presence of *M. hominis* at a threshold ≥10<sup>4</sup>CCU/mL in a cervicovaginal specimen is abnormal (1, 3).

### Susceptibility testing (wells 4 to 13 and 15 to 24)

The red colour change of the medium in the wells containing an antibiotic indicates the presence of bacterial growth and hence resistance to the antibiotic concentration being tested. The yellow colour of the medium indicates the absence of bacterial growth and hence susceptibility to the antibiotic concentration being tested. The strains are characterized as being Susceptible or resistant to the antibiotics according to the following criteria defined by the CLSI (2):

### Table of MIC (µg/mL) interpretative criteria

Antibiotic		Uu		Mh		Comments
Class	Drug	S	R	S	R	
Quinolones	Levofloxacin	≤2	≥4	≤1	≥2	/
	Moxifloxacin	≤2	≥4	≤0.25	≥0.5	/

Lincosamides	Clindamycin	/	/	≤0.25	≥0.5	/
Tetracyclines	Tetracycline	≤1	≥2	≤4	≥8	/
	Doxycycline	≤1	≥2	≤4	≥8	/
Macrolides	Erythromycin	≤8	≥16	/	/	Organisms susceptibles to erythromycin will also be susceptible to azythromycin

Help with interpretation:

### Susceptibility testing for Uu

Antibiotic	LVX			MXF			ERY			TET			DOX		
Concentration (µg/mL)	2	4	int <sup>†</sup>	2	4	int <sup>†</sup>	8	16	int <sup>†</sup>	1	2	int <sup>†</sup>	1	2	int <sup>†</sup>
Profile	-	-	S	-	-	S	-	-	S	-	-	S	-	-	S
	+	-	R	+	-	R	+	-	R	+	-	R	+	-	R
	+	+	R	+	+	R	+	+	R	+	+	R	+	+	R

int<sup>†</sup> = interpretation

### Susceptibility testing for Mh

Antibiotic	LVX			MXF			CLI			TET			DOX		
Concentration (µg/mL)	1	2	int <sup>†</sup>	0.25	0.5	int <sup>†</sup>	0.25	0.5	int <sup>†</sup>	4	8	int <sup>†</sup>	4	8	int <sup>†</sup>
Profile	-	-	S	-	-	S	-	-	S	-	-	S	-	-	S
	+	-	R	+	-	R	+	-	R	+	-	R	+	-	R
	+	+	R	+	+	R	+	+	R	+	+	R	+	+	R

int<sup>†</sup> = interpretation

The strain is said to be susceptible when its growth is inhibited by the higher and lower critical concentrations of the antibiotic. The strain is said to be resistant when its growth is inhibited by the higher critical concentration of the antibiotic, but not the lower critical concentration or when its growth is not inhibited by either the higher or the lower critical concentrations of the antibiotic.  
*M. hominis* strains are innately resistant to macrolides (14 -15 carbon atoms), including erythromycin. In some patient populations, tetracycline resistance is as high as 45% for Uu and 39.6% for Mh (2). Uu/Mh quinolone (5, 6) and clindamycin resistance have been described but the prevalence is not known.

## 10 - PARTICULAR CASES

For high Uu and Mh levels, there is a red turn of the wells concerned by the germ. It is recommended that the sample be diluted in order to obtain more specific results. In this case, proceed as follows:  
Inoculate a new UMMt 3mL vial with 300 µL of the original UMMt medium stored at 2-8°C (see § 9.1).  
Inoculate a new tray with the new inoculated UMMt medium.  
Take the dilution (1:10) into account in the interpretation of the enumeration results.  
If necessary, confirm the presence of mycoplasmas on an A7 agar plate by re-isolating from the original UMMt medium stored at 2-8°C (§ 9.1). Anon-constant incubation temperature or <36°C (frequent opening and temperature heterogeneity of the incubator) can slow down the mycoplasmas growth kinetics.

## 11 - QUALITY CONTROL

Quality control can be carried out with the lyophilized *U. urealyticum* or *M. hominis* strains of the MYCOPLASMA CONTROL kit (REF 00900) or from a lyophilized reference strain (*U. urealyticum* ATCC 27815 or *M. hominis* ATCC 23114) previously calibrated at 10<sup>4-5</sup> CCU/mL.

Inoculate the MYCOFAST *Revolution* 2 tray and perform the test as indicated in these instructions (§ 9 and 10).

Expected results (ATCC): MYCOFAST *Revolution* 2

	Uu 10 <sup>3</sup>	Uu 10 <sup>4</sup>	Uu ≥10 <sup>5</sup>	Mh ≥10 <sup>4</sup>	LVX	MXF	ERY	TET	DOX	CLI
Strain Uu ATCC 27815	+	+	+/-	-	S	S	S	S/R	S	NI*
Strain Mh ATCC 23114	-	-	-	+	S/R	S	NI*	S	S	S

NI\* (Not interpretable)

## 12 - LIMITATIONS OF THE PROCEDURE

### 12.1 - Screening:

The MYCOFAST Screening *Revolution* tray allows detection at a threshold of  $\leq 10^3$  UCC/mL and does not enable numeration. The numeration obtained with the MYCOFAST *Revolution* tray can appear negative following positive screening.

### 12.2 Enumeration, identification and susceptibility testing:

Some bacteria that are present in quantities of  $> 10^{6-7}$  CFU/mL and contain urease may cause all the wells in the tray to change colour. The presence of these can be verified by re-isolating on chocolate agar from the original UMMt medium stored at 2-8°C (§ 9.2).

An alkaline sample pH (pH  $\geq 8$ ) may lead the medium to change colour. Should this occur, dilute the sample (1:10) in fresh UMMt AMIES and interpret the results taking the dilution into account.

A sample with an acidic pH (pH  $\leq 5$ ) can slow down the appearance of the colour change.

A sample containing blood may cause a colour change in the wells of the MYCOFAST *Revolution* 2 tray and could be interpreted as a positive result. In this case dilute the sample (1:10) in another UMMt medium and interpret the results, taking into account the dilution.

A sample with a low mycoplasma load ( $< 10^3$  CCU/mL) may lead to a random colour change in the different wells of the tray. As for all germ detection methods, the quality of the sample can influence the test result. A negative test does not therefore necessarily indicate the absence of infection.

## 13 - PERFORMANCES

### 13.1 Screening and differentiation

#### MYCOFAST Screening *Revolution* tray

A comparative study was carried out using dry swabs provided with their transport medium (Transwab etc) on clinical vaginal samples (n= 87 for both Uu and Mh species).

The results obtained with MYCOFAST Screening *Revolution* were compared with those obtained with the liquid micro-dilution method.

- For the screening of both Uu and Mh species, there was 97.7% concordance.

For Uu at  $10^3$  CCU/mL, 2 samples that were positive with the liquid micro-dilution method appeared negative with the MYCOFAST Screening *Revolution*. However, it is important to underline that the  $10^3$  CCU/mL concentration corresponds to an infra-pathological threshold usually quoted for Uu. For Uu and Mh, the global agreement for the supra-pathological threshold is 100%.

- Differentiation of the mycoplasma species from all tested samples were correctly identified in the corresponding Uu and Mh wells of the MYCO-

FAST Screening *Revolution* tray.

### 13.2 Identification - Enumeration

#### Direct method MYCOFAST *Revolution* 2 tray

% of overall agreement		Uu	Mh	Uu/Mh
Isolated strains (threshold $\leq 10^3$ CCU/mL) (see § 14.1.1 and note 1)	24h reading	88.9	NA*	NA*
	read as in § 9.3	86.7	NA*	NA*
Isolated strains (threshold $\geq 10^4$ CCU/mL) (see § 14.1.1 and note 1)	24h reading	91	96.4	93.7
	read as in § 9.3	82.1	92.5	87.3
vaginal clinical samples (see § 14.1.2)	read as in § 9.3	88.2	100	94

NA \*: not applicable

#### 13.2.1 Isolated strains

A comparative study was carried out using 26 isolated strains (ATCC strains and collection strains) tested separately (Uu or Mh) on 3 transport medium references (Sigma Transwab and Sigma VCM from Medical Wire and ESwab Collection Kit from BD) at several concentrations (a total of 279 tests).

The results obtained are compared with those obtained using the micro-dilution liquid method.

For interpretation with a pathological threshold set at  $10^3$  CCU/mL and a reading of the result at 24h, the overall agreement for Uu is 88.9% (we listed 14 false positives: 13 at  $10^2$  CCU/mL and one below  $10^2$  CCU/mL and 17 false negatives: 11 at  $10^3$  CCU/mL - 5 f at  $10^4$  CCU/mL and one at  $10^5$  CCU/mL enumeration method).

For interpretation with a pathological threshold set at  $10^4$  CCU/mL and a reading of the result at 24h; the global concordance for Uu is 91% (we have listed 18 false positives: at  $10^2$  -  $10^3$  CCU/mL and 7 false negatives: at  $10^4$  -  $10^5$  CCU/mL with the micro-dilution enumeration method).

The overall agreement for Mh is 96.4% for a reading of the result at 24h (we recorded one false positive at  $10^3$  CCU/mL, and 9 false negatives at  $10^4$  CCU/mL with the micro-dilution enumeration method).

The global agreement Uu + Mh with a reading of the result at 24h is 93.7%.

For interpretation with a pathological threshold set at  $10^3$  CCU/mL and a reading of the result according to the protocol described in the leaflet (§ 9.3); the overall agreement for Uu is 86.7% (35 false positives were recorded: 34 at  $10^2$  CCU/mL and one below  $10^2$  CCU/mL and 2 false negatives: at  $10^3$  CCU/mL with the micro-dilution enumeration method). For interpretation with a pathological threshold set at  $10^4$  CCU/mL and a reading of the result according to the protocol described in the leaflet (§ 9.2); the overall agreement for Uu is 82.1% (we listed 49 false positives: at rate at  $10^2$  -  $10^3$  CCU/mL and 1 false negative: at  $10^4$  CCU/mL with the micro-dilution enumeration method).

The overall agreement for Mh is 92.5% for a reading of the result according to the protocol described in the leaflet (§ 9.2) (we have listed 21 false positives at a rate of  $10^2$  -  $10^3$  CCU/mL with the micro-dilution enumeration method).

The global agreement Uu + Mh with a reading of the result according to the protocol described in the notice (§ 9.2) is 87.3%.

#### 13.2.2 Clinical samples

An initial comparative study was performed using vaginal clinical specimens (n =59) carried out on 3 types of transport medium and their associated swabs (Sigma Transwab and Sigma VCM from Medical Wire and ESwab Collection Kit from BD). The results obtained with MYCO-

FAST *Revolution* 2 AMIES were compared with the liquid micro-dilution enumeration method.

The overall agreement for Uu is 88.2% (we identified 3 false negative at  $10^4$  -  $10^5$  -  $10^6$  CCU/mL and 4 false positive at  $< 10^2$  -  $10^2$  and  $10^3$  CCU/mL with micro-dilution numeration method).

The overall agreement for Mh is 100%

The overall agreement for Uu and Mh is 94%.

### 13.3 Susceptibility testing

A comparative study was carried out in a national reference laboratory between the method for determining the minimum inhibitory concentrations (MIC) in liquid medium and the MYCOFAST *Revolution* 2.

The tested strains (7 *U. urealyticum*, 11 *U. parvum* and 16 *M. hominis*) were reference strains, wild-type clinical strains or strains with acquired resistance. Each strain was tested at  $10^3$ ,  $10^4$  and  $10^5$  CCU/mL dilutions in UMMt 3 mL.

For  $10^4$  and  $10^5$  CCU/mL rates, results were read and interpreted after 24 hours of incubation.

For  $10^3$  CCU/mL rate, results were read and interpreted after 48 hours incubation in case of negative test in 24 hours.

The results of both methods were interpreted as susceptible (S) or resistant (R) according to CLSI recommendations.

The overall agreement for *Ureaplasma urealyticum* / *Ureaplasma parvum* is: 95.5%.

The overall agreement for *Mycoplasma hominis* for rates at  $10^4$ - $10^5$  CCU/mL is: 100%.

Agreement	<i>Ureaplasma urealyticum</i> / <i>parvum</i> (n=40)					<i>Mycoplasma hominis</i> (n=28)				
	TET	DOX	MXF	LVX	ERY	TET	DOX	MXF	LVX	CLI
	34	38	40	39	40	28	28	28	28	28
ME	5 <sup>a</sup>	0	0	0	0	0	0	0	0	0
VME	1 <sup>b</sup>	2 <sup>c</sup>	0	1 <sup>d</sup>	0	0	0	0	0	0

ME: Major Error, VME : Very Major Error

<sup>a</sup>: 1 discrepancy at  $10^3$  CCU/mL (MIC of reference 0.5 µg/mL), 4 discrepancies at  $10^5$  CCU/mL (MIC of reference 0.5 - 1 and 8 µg/mL).

<sup>b</sup>: 1 discrepancy at  $10^5$  CCU/mL (MIC of reference 8 µg/mL).

<sup>c</sup>: 1 discrepancy at  $10^3$  CCU/mL (MIC of reference 8 µg/mL); 1

discrepancy at  $10^5$  CCU/mL (MIC of reference 2 µg/mL)

<sup>d</sup>: 1 discrepancy at  $10^5$  CCU/mL (MIC of reference 4 µg/mL).

## 14 - WASTE ELIMINATION

Waste should be disposed of in accordance with the hygiene rules and current regulations for this kind of product in the country of use.

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2 - Clinical and Laboratory Standards Institute (CLSI). 2011 Methods for Antimicrobial Susceptibility Testing for Human Mycoplasmas; Approved Guideline. CLSI Document M43-A. Vol.31 - N°19.

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of Infectious Diseases, 4th ed., vol. 2, Churchill Livingstone, New York.

**5 - WAITES KEN B. , BRENDAKATZ AND ROBERT L. SCHELONKA. 2005.** Mycoplasmas and Ureaplasmas as Neonatal Pathogens. Clin. Microbiol. Rev. Vol.18 -N°4 -757-789.

**6 - WAITES KEN B, DONNAM. CRABB, and LYNN B. DUFFY. 2008.** Comparative In Vitro Activities of the Investigational Fluoroquinolone DC-159a and Other Antimicrobial Agents against Human Mycoplasmas and Ureaplasmas. ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, Vol. 52, No. 10, 3776–3778.

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MYCOFAST® is a trademark of ELITech MICROBIO

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