BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFÍA/BIBLIOGRAFIA

1.- Schreiber,W.E., *Iron and Porphyrin Metabolism*, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc. eds.), (2010), 755 and appendix.
2.- Higgins, T., Beutler, E. and Doumas, B.T., *Hemoglobin, Iron and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed, Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
3.- Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 634.
4.- Siek, G., *et al.*, Clin. Chem., (2002), **48**, 161.
5.- Yamanishi, H. *et al.*, Clin. Chem, (2003), **49**, 175.
6.- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7.- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
8.- Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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IMPORTANT NOTE
- See § PROCEDURE: Manual entry required

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TIBC (µg/dL) = transferrine (mg/dL) x 1.4
TIBC (µmol/L) = transferrine (g/L) x 25

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English - EN

INTENDED USE

ELITech Clinical Systems Direct TIBC is an *in vitro* diagnostic reagent intended for the quantitative determination of Total Iron Binding Capacity in human serum samples on analyzers

This *in vitro* diagnostic device is for professional use only.

CLINICAL SIGNIFICANCE (1-3)

In the body, 65 to 70 % of total iron is circulating in erythrocyte hemoglobin, 20 to 25% is stored in cells mainly as an iron-ferritin complex, less than 0.1% is transported by transferrin, the remaining iron being associated with other proteins, including myoglobin. Total iron-binding capacity (TIBC) is the measure of the maximum concentration of iron that the serum proteins, mostly transferrin, can bind. TIBC is increased in iron deficiency, during pregnancy or by hormonal contraception. A decrease is associated with iron overload (hemochromatosis), inflammation, nephrotic syndrome, malnutrition or some malignancies. TIBC is used to assess the amount of iron that can be bound by serum proteins (mainly Transferrin) and to investigate iron deficiency or iron overload.

LIMITATION OF USE

The quantitative assay of TIBC alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE (4-5)

Chromazurol B – End Point.

The sample is added to reagent R1 which contains an acidic buffer, an iron-binding dye and ferric chloride, the concentration of which is in excess compared to the transferrin iron-binding capacity. Under acidic condition, transferrin-bound iron is released and binds to the dye. The colored complex then represents both the serum iron and excess iron from reagent R1. Addition of reagent R2 increases the pH, restoring the strong affinity of iron for serum transferrin. Separation of iron from the dye results in a decrease in absorbance, directly proportional to the total iron-binding capacity (TIBC) of the sample. Since transferrin is the only serum protein with actual specificity for iron, there is a direct relationship between TIBC and serum transferrin:

TIBC (µg/dL) = transferrin (mg/dL) x 1.4
TIBC (µmol/L) = transferrin (g/L) x 25

The accuracy of this conversion may however vary with the nonspecific influence of other serum proteins or transferrin's functional defects.

COMPOSITION

Reagent : R1

Acetate buffer (pH 4.5)		100 mmol/L
Chromazurol B	≥	108 µmol/L
Cetrimide		734 µmol/L
Ferric chloride		15.9 µmol/L

Reagent: R2

Buffer Sodium bicarbonate	338 mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS
- Consult Safety Data Sheet (SDS) for a proper handling.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.

STABILITY
Store at 2-8 °C and protect from light. Do not freeze.
Do not use after expiration dates indicated on the vial labels.
On board stability:
The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION
The reagents are ready to use.

PRODUCT DETERIORATION
- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES
Specimen (1,3)
- Serum
- Using any other specimen type should be validated by the laboratory.

Warnings and precautions
Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability (1,3)
- 4 days at room temperature
- 7 days at 2-8°C
- 1 month at -20°C

REFERENCE VALUES (1)

Serum	µg/dL	µmol/L
Adult :	250 - 450	44.8 - 80.6

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

PROCEDURE
For ELITech Clinical Systems Selectra Analyzers applications are available on request.
Wavelength 620 nm
Temperature: 37 °C
Read against reagent blank

	BLANK	CALIBRATION	TEST
Reagent R1	220 µL	220 µL	220 µL
Distilled water	18 µL	-	-
Calibrator	-	18 µL	-
Sample	-	-	18 µL

Mix, incubate 4 minutes and 42 seconds then add:

Reagent R2	66 µL

Mix and read the absorbance (A1) after 6 seconds following R2 addition, and again 390 seconds following R2 addition (A2).

- With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- In the application, the offset must be manually set to: -50 µg/dL (-9.0 µmol / L).

CALCULATION
$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Calibrator}}$ x n n = calibrator concentration

Conversion factor : µg/dL x 0.179 = µmol/L

CALIBRATION
ELICAL 2 is traceable to ERM-DA470k/IFCC reference material (transferrin concentration (mg/dL) converted into TIBC (µg/dL) with factor 1.4).

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL
It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.
Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT
Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES
Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range
70 - 790 µg/dL (12.5 - 141.4 µmol/L)
Do not report results outside the measuring range.

Limit of Detection (LoD) and Limit of Quantification (LoQ)
LoD = 41 µg/dL (7.3 µmol/L)
LoQ = 70 µg/dL (12.5 µmol/L)

- Precision
Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented in the following table.

		Mean		Within-run	Total
	n	µg/dL	µmol/L	CV (%)	
Level 1	80	157	28.1	2.8	6.1
Level 2	80	194	34.7	1.7	5.3
Level 3	80	356	63.7	1.4	3.3

- Correlation
A comparative study has been performed between Direct TIBC reagent on a Selectra ProM analyzer and a commercially available Transferrin reagent on 100 human serum samples. For the comparison, transferrin values (mg/dL) were transformed into TIBC values (µg/dL) through the 1.4 conversion factor (see §Method and §Bibliography, reference 5).
The sample concentrations ranged from 45 to 846 µg/dL (8.1 - 151.4 µmol/L).
The results are as follows :
Correlation coefficient : (r) = 0.989
Linear regression: y = 1.025x - 8 µg/dL (1.4 µmol/L).

- Limitations/Interferences
- Studies have been performed to determine the level of interference from different compounds.
The following TIBC levels were tested: 125 µg/dL and 450 µg/dL.
No significant interference is defined by a recovery ≤±25 µg/dL of initial value at TIBC concentration of 125 µg/dL and ≤± 10% of initial value at TIBC concentration of 450 µg/dL.

Unconjugated bilirubin: No significant interference up to 22.5 mg/dL (385 µmol/L).
Conjugated bilirubin: No significant interference up to 29.5 mg/dL (505 µmol/L)
Hemoglobin: No significant interference up to 500 mg/dL.
Turbidity: No significant interference up to 600 mg/dL (6.8 mmol/L) triglyceride equivalent
Ascorbic acid: No significant interference up to 9.9 mg/dL
Acetaminophen: No significant interference up to 30 mg/dL
Copper: No significant interference up to 500 µg/dL (78.7 µmol/L)
Iron: No significant interference up to 170 µg/dL (30.4 µmol/L)
Acetylsalicylic acid: No significant interference up to 65 mg/dL. Acetylsalicylic acid can interfere when present in concentrations above therapeutic doses.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.(6)

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.(7-8)

- On board stability/Calibration frequency
On Board Stability: 21 days
Calibration frequency: 7 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

DECLARATION OF SERIOUS INCIDENT
Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements.
By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

TECHNICAL ASSISTANCE
Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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Entre em contato com o seu distribuidor local ou com a ELITech Clinical Systems SAS.