


SEROLOGY, PREGNANCY & OBSTETRIC

SEROLOGY

REFERENCE	PRODUCT	KIT SIZE UNIT
ASOL-0100	ASO Latex	100 Tests
LXCR-0112	CRP Latex	100 Tests
LXRF-0112	RF Latex	100 Tests
LXWR-0112	WAALER ROSE	100 Tests
RPRL-0100	RPRL VDRL Carbon	100 Tests
TPHA-0004	TPHA Microplates	5 Units
TPHA-0100	TPHA 100 Tests	100 Tests + 5 Microplates



PREGNANCY & OBSTETRIC

REFERENCE	PRODUCT	KIT SIZE UNIT	METHODOLOGY	SPECIFICATIONS
RTHC-0040	PREGTEST PLUS HCG	30 Tests	Immuno-chromatography	10 mIU/mL in serum, 20 mIU/mL in urine
PRMT-0050	PMR TEST	50 Tests	Sensitive pH indicator, Nitrazine yellow	Premature Membrane Rupture
PRMT-0100	PMR TEST	100 Tests	Sensitive pH indicator, Nitrazine yellow	Premature Membrane Rupture



About ELITech Clinical Systems...

ELITech Clinical Systems is a global leader in BenchTop Chemistry Systems, with over 15,000 Chemistry Systems delivered. Our mission is to enable better medical decisions by bringing high-value diagnostic solutions and services to laboratories that are closer to the patient (proximity laboratories). ELITech Clinical Systems leadership originates from our ability to continuously provide differentiated products and services focused to meet the needs of proximity laboratories. ELITech Clinical Systems markets the Selectra Pro and Flexor EL analysers along with ready-to-use reagents, designed to deliver trouble free results laboratories can trust.

Delivering Trouble Free Results Laboratories can Trust

Robust, Always gets the job done	Trouble Free Operations	Right Results, On Time
<ul style="list-style-type: none"> ✓ Broad menu of chemistries with proven methodologies ✓ True workhorse when combined with Selectra Pro and Flexor EL systems ✓ Support you can count on 	<ul style="list-style-type: none"> ✓ Simplifies lab operations and reduces errors ✓ Liquid stable ready-to-use assays ✓ Optimum sized packaging further enhances efficiency 	<ul style="list-style-type: none"> ✓ Error reducing capabilities with automated programming of assays and 2D bar coded reagents ✓ Increased on-board stability ✓ Proven methodologies, referenced and traceable to industry standards

ELITech Clinical Systems test menu features a broad range of chemistries with liquid stable ready-to-use reagents and proven methodologies that are referenced to industry standards. The enhanced reagent supply provides flexibility in reagent containers and increased on-board stability. Optimum sized packaging further enhances efficiency. Together with the Selectra Pro and Flexor EL analysers, ELITech Clinical Systems test menu provides an ideal chemistry solution for laboratories that are closer to the patient. ELITech also develops and manufactures serology, pregnancy and obstetric kits to meet the needs of proximity laboratories. Over the past decades, the ELITech Group Companies have established a solid distribution network that operates world-wide. Thousands of customers have experienced the quality, convenience, reliability of ELITech products. Local support is provided by dedicated and well trained sales and service organizations.

BIOCHEMISTRY Reagents SEROLOGY, PREGNANCY & OBSTETRIC Rapid Tests



Part nr: ECS Gen-300-2010/00EN

DISTRIBUTOR CONTACT INFORMATION:

WORLDWIDE OFFICES:

Belgium T: + 32 9 2820 531
Brazil T: + 55 27 3025 1415
France T: + 33 4 90 17 54 50
Italy T: + 39 02 48 40 35 42

Middle East & Africa T: + 971 4 375 2744
Serbia T: + 381 11 2467119
Switzerland T: + 41 26 663 86 60
The Netherlands T: + 31 485 511 675

UK T: + 44 1268 493377
United States T: + 1 609 216 7360
855-ELITTECH (855-354-8324)

Headquarters T: + 33 1 41 45 07 13
www.elitechgroup.com
info@elitechgroup.com

