



► Scalable Solutions

ELITechGroup, offers a comprehensive range of benchtop chemistry systems for the physician office laboratory and small hospitals. Our benchtop chemistry systems are designed to meet the needs of laboratories running from 15,000 to 200,000 tests per year. All systems are supported by technical service and remote diagnostics.

Leading Family of Benchtop Chemistry Systems




INSTRUMENT	PATIENTS PER DAY	TESTS PER HOUR	ANNUAL TEST VOLUME
Selectra Pro S	5 - 15	140	15,000 - 60,000
Selectra Pro M	10 - 40	180	40,000 - 100,000
Selectra Pro S	20 - 80	490	50,000 - 200,000



Selectra Pro S

Fully-featured, compact system

- 5-15 patients per day
- 30 refrigerated reagent positions
- 25 position sample rotor



Selectra Pro M

Fully-featured, mid-volume system

- 10-40 patients per day
- 32 refrigerated reagent positions
- 50 position sample rotor



ENVOY500+

Higher throughput in benchtop chemistry

- 20-80 patients per day
- 40 refrigerated reagent positions
- 52 position sample wheel
- Clot detection sensor
- Reusable glass cuvettes



STANDARD FEATURES FOR ALL ANALYZERS

- Touch screen user interface
- Ready-to-use, liquid-stable reagents
- Positive sample identification, host query¹
- Barcoded reagents
- Liquid level sensing
- Programmable result checks, auto-dilutions
- Remote diagnostics

¹Selectra systems

► Comprehensive Chemistry Systems Service and Support

- Technical Support, Monday – Friday, 8AM – 7PM EST.
- Remote Diagnostics and on-site service
- Full one-year warranty; post-warranty options

► ELITechGroup

ELITech is a privately held group of worldwide manufacturers and distributors of in vitro diagnostic equipment and reagents. By bringing together IVD specialty companies that offer innovative products and solutions, ELITechGroup has become a major contributor in advancing clinical diagnostics to laboratories in the proximity market, those operating closer to the patient. The ELITechGroup manufactures and distributes diagnostic products for clinical chemistry, microbiology, immunology, and molecular biology through direct sales and a distribution network encompassing more than 100 countries. For more information about the ELITechGroup, visit our website: www.elitechgroup.com.

© 2021. ELITechGroup, Inc. All rights reserved.
 Envoy is a registered trademark of ELITechGroup.
 Windows is a registered trademark of the Microsoft Corporation.
 L7103, Rev. G



ENVOY500+ is available in the USA only. To explore ELITechGroup products available in your country, please contact your local sales representative.

ELITechGroup North America,
 370 West 1700 South Logan, UT 84321 USA
 Email: infoUS@elitechgroup.com
 Phone: +1-800-345-2822