



Molecular Diagnostics

**MGB**

# HCV ELITe MGB® Kit

Real-Time PCR  
Assay for HCV  
Monitoring

**IDEAL SOLUTION**  
for all small and  
medium labs

**EASE OF USE** on  
ELITe InGenius® system

**OPTIMAL SENSITIVITY**  
for HCV patient  
monitoring

**COMPLETE PANEL** in  
combination with HBV  
and HIV1 ELITe MGB® Kits

CE IVD  
0344

*In Vitro* Diagnostic Medical Device  
Not available in all countries. Not  
available in the United States.



**ELITechGroup**  
EMPOWERING IVD



"Early diagnosis can prevent health problems that may result from infection and prevent transmission of the virus. WHO recommends testing people who may be at increased risk of infection. Quantitative HCV RNA testing is recommended prior to initiating antiviral therapy and at the end of treatment after 12 or 24 weeks. WHO Global goal is to stop HCV by 2030 increasing diagnosis and treatment."

World Health Organization

**HCV ELITe MGB® Kit** is a quantitative assay detecting the RNA of Hepatitis C virus (HCV), genotypes 1, 2, 3, 4, 5 and 6 providing a new solution for the management of HCV- infected individuals undergoing antiviral therapy.

**HCV ELITe MGB® Kit** in combination with **ELITe InGenius®** shows optimal sensitivity and quantification of HCV RNA.



### PERFORMANCES

**LoD: 26 IU/ml in plasma and serum samples**  
**LoD: 11 copies/ml in plasma and serum samples**  
**Diagnostic Specificity: 100%**

### ALL YOU NEED IN ONE BOX

The Kit contains ready-to-use reagents, standards and controls.

### RESULTS IN LESS THAN 3 HOURS

Fast results support better disease management and better patient care



#### Ordering information

Reference	Product	Quantity
RTK601ING	<b>HCV ELITe MGB® Kit</b>	96 test

ELITechGroup Molecular Diagnostics, a major innovator in the molecular diagnostic market, offers a comprehensive real-time PCR range of products for diagnosis and monitoring of viral and microbiological infections.



ELITechGroup is certified by DEKRA Certification BV and applies a quality management system in compliance with ISO 9001:2015 and EN ISO 13485:2016

Limited license : Please visit <http://www.elitechgroup.com/corporate/elitemgb-legalnotice> for complete licensing and warranty information. ELITe MGB® is a registered trademark of ELITechGroup. Products not available in the US.



#### WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.



**Manufacturer:** ELITechGroup S.p.A.


#### Administrative and operational site

Corso Svizzera 185  
10149 Turin (TO) - Italy

#### Registered office

Corso Italia, 22  
20122 Milan (MI) - Italy

Headquarters	+33 1 41 45 07 10		
Australia	+61 1800 815 098	New Zealand	+64 800 555 611
Belgium, Luxembourg	+32 9 282 05 31	Serbia	+381 11 2467119
Brazil	+55 27 3025 1415	The Netherlands	+31 313 430 500
France	+33 4 83 36 10 82	UK	+44 1442 869 320
Italy	+39 011 97 61 91	United States	+1 800 453 2725



**ELITechGroup**  
EMPOWERING IVD  
[www.elitechgroup.com](http://www.elitechgroup.com)  
[info@elitechgroup.com](mailto:info@elitechgroup.com)