



Microbiology

MYCOFAST® REVOLUTION 2

Urogenital mycoplasma diagnosis
DIRECT METHOD



FULFILLING YOUR NEEDS IN THE MOST DEMANDING SITUATIONS

Detection, Enumeration, Identification and Antimicrobial Susceptibility Testing

- ◆ Complete diagnosis in an All-in-1 product
- ◆ Compatible for use with liquid-Amies collection and transport media
- ◆ Optimized management of co-infections
- ◆ Antibiotics compliant with worldwide recommendations
- ◆ The right pathologic threshold tested for each specimen type
- ◆ Fast results for clinicians
- ◆ Visual, colorimetric reading

MYCOFAST® REVOLUTION 2



► BENEFITS

◆ Complete diagnosis in 24h*

- ◆ Detection, Enumeration
- ◆ Identification
- ◆ Antimicrobial Susceptibility Testing

◆ Cutting-edge performance

Sensitivity

- ◆ Detection and enumeration of *U. urealyticum* at 10^3 , 10^4 and $\geq 10^5$ CCU/mL
- ◆ Detection and enumeration of *M. hominis* at $\geq 10^4$ CCU/mL
- ◆ The right pathologic threshold tested for each specimen type

◆ Specificity

- ◆ Selective growth media for each species tested separately

Reliability

- ◆ 5 antibiotics (LVX/MXF/E/CM/TE) tested in compliance with CLSI guidelines (Clinical & Laboratory Standards Institute)
- ◆ Doxycyclin added to the panel

Fast results

- ◆ Integrated growth activator for *M. hominis*

◆ Ease of use


- ◆ Ready-to-use reagents
- ◆ Each well pre-filled with growth medium
- ◆ Visual colorimetric reading
- ◆ No instrument required

◆ Flexibility and Safety

- ◆ Optimized transport medium: Specimens preserved up to 72 hours before analysis
- ◆ Clip-on closing system for the tray

► METHODOLOGY

1 Inoculation of the transport medium



* Reading of all samples in 24h.
Only read liquid samples having negative results within 24h.

► PRODUCT RANGE

REF.	NAME	QTY	PURPOSE
00080	MYCOFAST® Revolution N 2	25 tests	Complete kit
00081	MYCOFAST® Revolution N 2 AMIES	25 tests	Complete kit
00061	UMMt RevolutionN (optional)	50 vials	Additional transport medium
00083	UMMt AMIES RevolutionN (optional)	50 vials	Additional transport medium (for Amies specimens)
00082	COMPLEMENT MYCOFAST® Revolution N 2 (optional)	25 tests	Additional test trays
00064	MYCOPLASMA Stabilizer (optional)	100 tests	Additive for sample storage > 56 hours
00900	MYCOPLASMA Control (optional)	12 tests	U.u. and M.h. control strains

OPTIONAL



WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.

Headquarters

Australia
BELUX
Brazil
France
Italy

+33 1 41 45 07 10

1800 815 098
+32 92 820 531
+55 27 3025 1415
+33 4 83 36 10 82
+39 011 67 61 913

New Zealand
Serbia
The Netherlands
UK
United States

0800 555 611
+381 11 2467119
+31 313 430 500
+44 1442 869320
+1 435 0752 6011



: ELITech Microbio (France)



ELITechGroup

EMPOWERING IVD

www.elitechgroup.com
info@elitechgroup.com