

# HDL LDL CALIBRATOR

REF HLCA-0041





FV-HLCA-v4 (11/2021)

HDL	LDL value for reagent lots <b>before 21-0894</b> . For newer lots refer to page 2
36 mg/dL 0.93 mmol/L	62 mg/dL 1.61 mmol/L
Traceability : reference method approved by the CDC (Centers for Disease Control and Prevention)	

LOT 21-4019

2023-09



[www.elitechgroup.com](http://www.elitechgroup.com)

**LDL Reagent lots before 21-0894**

selectra  
**PRO**



HDL LDL Calibrator 0  
S25 FV-HLCA  
21-4019 2023-09

Place for future  
**Selectra MACH Barcode**



# HDL LDL CALIBRATOR

REF HLCA-0041





FV-HLCA-v4 (11/2021)

<b>HDL</b>	<b>LDL</b> value for reagent lots from 21-0894. For former lots refer to page 1
36 mg/dL 0.93 mmol/L	67 mg/dL 1.73 mmol/L
<b>Traceability</b> : reference method approved by the CDC (Centers for Disease Control and Prevention)	

LOT 21-4019


2023-09



[www.elitechgroup.com](http://www.elitechgroup.com)

LDL Reagent lots from 21-0894

selectra  
**PRO**



HDL LDL Calibrator 0  
S25 FV-HLCA  
21-4019 2023-09

Place for future  
**Selectra MACH Barcode**

