

# RHEUMATOID FACTOR CALIBRATOR

REF RHFA-M130



FV-RHFA-v1 (07/2021)

200 IU/mL


LOT







22-0160

Traceability : Rheumatoid Arthritis Serum, 1<sup>st</sup> British Standard




2024-08

Preparation 

Rhumatoid Factor Cal	μL	-	20	40	80	150	300
NaCl 9 g/L	μL						
Dilution Factor		0	1/16	1/8	1/4	1/2	1

Precautions



 [www.elitechgroup.com](http://www.elitechgroup.com)

Selectra MACH

Place du code-barre 2D

