

CRP IP CALIBRATOR SET

REF ICRP-0043



PITV-CRPCalS-v13 (06/2021)

| CAL | 1 | 2 | 3 | 4 | 5 |
|-------------------------------|------|------|------|-------|-------|
| mg/L | 7.6 | 15.2 | 30.3 | 119.9 | 200.0 |
| mg/dL | 0.76 | 1.52 | 3.03 | 11.99 | 20.00 |
| Traceability : ERM-DA474/IFCC | | | | | |

LOT 22-4003



2024-10


Preparation



Place for future

Selectra MACH Barcode

Precautions



www.elitechgroup.com

