

RHEUMATOLOGY CONTROL II

REF IRCT-0047





PITV-IRCT2-v4 (06/2021)

Selectra PRO		
ASLO-0250	ANTI-STREPTOLYSIN O	305 (213 - 397) IU/mL
IRFA-0230	RHEUMATOID FACTOR	52 (36 - 68) IU/mL


LOT 22-4016

2024-09



www.elitechgroup.com

selectra
PRO



RHEUMATOLOGY CTRL II 0
C85 PITV-IRCT2
22-4016 2024-09



RHEUMATOLOGY CONTROL II

REF IRCT-0047



PITV-IRCT2-v4 (06/2021)

Selectra MACH


ASLO-0250	ANTI-STREPTOLYSIN O	305 (213 - 397) IU/mL
ICRP-M230	CRP IP	52.7 (36.9 - 68.5) mg/L
RHFA-M130	RHEUMATOID FACTOR	52 (36 - 68) IU/mL

LOT 22-4016



2024-09


Preparation



Precautions



10
days
días
jours
tage



2°C



3
months
meses
mois
monat



-20°C



www.elitechgroup.com

Place for future

Selectra MACH Barcode

