

# CK-MB CONTROL

REF CKMB-0900



FVCE-CKMB-v6 (06/2021)

Selectra PRO	Selectra MACH
144 (115 - 173) U/L	144 (115 - 173) U/L

LOT 22-4017

2024-07

**Preparation**

3 mL → 30 min → 15°C - 25°C → ✓

**Precautions**


5 days / 5 jours / 5 tage  
2°C → 8°C

⊗



[www.elitechgroup.com](http://www.elitechgroup.com)

selectra  
**PRO**



CK-MB Control 0  
C30 FVCE-CKMB  
22-4017 2024-07

Selectra MACH

Place du code-barre 2D



ELITECH Clinical Systems SAS - Zone Industrielle - 61500 SEES FRANCE