

# VITAMIN D CONTROL SET

REF VITD-0049



PITV-VITDCo-v3 (07/2021)

	Control L	Control H
ng/mL	25.7 (19.3 - 32.1)	50.4 (37.8 - 63)
nmol/L	64.3 (48.3 - 80.3)	126.0 (94.5 - 157.5)


LOT 22-4020

2024-01

Preparation




Precautions



[www.elitechgroup.com](http://www.elitechgroup.com)

Control L  
selectra  
**PRO**



VITD CONTROL L 0  
C81 PITV-VITDCo  
22-4020 2024-01

Control H  
selectra  
**PRO**



VITD CONTROL H 0  
C86 PITV-VITDCo  
22-4020 2024-01



ELITech Clinical Systems SAS - Zone Industrielle - 61500 SEES FRANCE