

VITAMIN D CALIBRATOR SET

REF VITD-0043



PITV-VITDCa-v3 (06/2021)

CAL	1	2	3	4	5
ng/mL	0.0	31.7	66.4	103.5	156.5
nmol/L	0.0	79.3	166.0	258.8	391.3
Traceability : NIST SRM972a					

LOT 22-4021




2024-02

Preparation



selectra
PRO




VITD CALIBRATOR 0

S83 PITV-VITDCa

22-4021 2024-02

Precautions



www.elitechgroup.com



ELITech Clinical Systems SAS - Zone Industrielle - 61500 SEES FRANCE