

INTENDED PURPOSE

The ISE Tubing Kit is used for performing 4-month maintenance of the ISE Module on the Selectra Pro Series and EL Series Analyzers.

SUMMARY AND PRINCIPLE

The ISE Tubing Kit contains disposable tubing for the Selectra Pro Series and EL Series Analyzers. The components of the kit are replaced every 4 months per instructions in the user manual.

PROCEDURE

MATERIALS PROVIDED

Quantity	Part No.	Description
1	6003-450	Peristaltic Pump Tube
1	6003-451	Reference Tube
1	6003-452	Y-Connector Tube Assembly
1	6003-453	Cap Assembly: Reference Bottle Cap, Straws, and Tubing

MATERIALS REQUIRED BUT NOT PROVIDED

- Selectra Pro Series or EL Series Analyzer
- General Laboratory Equipment

INSTALLATION AND USE

1. Consult the Dry Electrode ISE / Maintenance chapter in the user manual for installation instructions and documentation.
2. The diagram below should be consulted in addition to reading the manual.
3. Verify and check that the electrode order is as shown in the diagram.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

SYMBOLS

Symbols used on the labeling are defined in the ISO 15223-1 standard.

- ☞ Modification from previous version.

