

# Accu-Sed® Plus

## Normal / Abnormal ESR Control, 4x8.5 mL



REF DS-71006

Kit Base Lot Number  
Kit numéro de lot de base  
Kit basen lot-nummer  
Kit bas lotnummer

LOT 2123500

Normal Control Contrôle normal Normal kontroll Normal kontroll	Abnormal Control Contrôle anormal Unormal kontroll Onormal kontroll
LOT 212010	LOT 213110
2022-06-07	2022-06-07

### INSTRUMENT ID CODES AND BARCODES CODES D'IDENTIFICATION ET CODE-BARRES DES INSTRUMENTS INSTRUMENT ID-KODER AND BARKODER ID-KODER OCH STRECKKODER FÖR INSTRUMENT

Instrument Equipement	Normal QC ID Code / Barcode Codes d'identification et code-barres QC Normal Normal kvalitetskontroll ID-kode/streckkode Normal QC ID kod/ streckkod		Abnormal QC ID Code / Barcode Codes d'identification et code-barres QC Anormal Unormal kvalitetskontroll ID-kode/streckkode Onormal QC ID kod / streckkod	
Excyte® 20	212106026		213270090	
Excyte® 40	212106026		213273093	
Monitor Family	212106015		213268042	
MixRate Family	212106015		213259097	
BD Sedi Family	212106026		213249052	
Greiner Bio-One Sed-Rate Screener Family	212107016		213267096	

\* ESR values corrected on temperature of 18°C, in accordance with Manley Table.

\* ESR valeurs corrigées de la température de 18°C, conformément au tableau Manley.

\* SR-verdiene korrigeret på temperaturer på 18 °C i samsvar med Manley-tabellen.

\* ESR-värden korrigerade för temperatur på 18 °C, i enlighet med tabellen Manley Table

# Accu-Sed® Plus

## Normal / Abnormal ESR Control, 4x8.5 mL



REF DS-71006

**ASSAY TABLE**  
**TABLEAU DE TESTS**  
**ANALYSETABELL**  
**ANALYSTABELL**

Normal Control Contrôle normal Normal kontroll Normal kontroll		Abnormal Control Contrôle anormal Unormal kontroll Onormal kontroll	
LOT	212010	LOT	213110
	2022-06-07		2022-06-07

Diluted Westergren Methods	Units	Mean	Range	Mean	Range
Méthodes de Westergren diluées	Unités	Valeur moyenne	Valeur prévue	Valeur moyenne	Valeur prévue
Utvannet Westergren metoder	Enheter	Gj.snitt	Område	Gj.snitt	Område
Utspádda Westergren Metoder	Enheter	Medeltal	Intervall	Medeltal	Intervall
Excyte® 10/M/Mini	mm/hr	7	3 – 11	67	49 – 85
Excyte® 20	mm/hr	6	2 – 10	70	52 – 88
Excyte® 40	mm/hr	6	2 – 10	73	55 – 91
Monitor Family	mm/hr	6	0 – 12	68	50 – 86
MixRate Family	mm/hr	6	0 – 12	59	41 – 77
BD Sedi Family	mm/hr	6	2 – 10	49	39 – 59
Greiner Bio-One Sed-Rate Screener Family	mm/hr	7	1 – 13	67	49 – 85
Polymedco Sediplast™	mm/hr	4	0 – 8	41	26 – 56
Dispette™ 2	mm/hr	6	1 – 11	50	35 – 65

### REFERENCES / RÉFÉRENCES / REFERANSER / REFERENSER

1. U. S. Department of Health and Human Services: Recommendation for Prevention of HIV Transmission in Health Care Settings. MMW Report, Aug 21, 1987, Vol. 36, No. 2S.

Accu-Sed and Excyte are registered trademarks of ELITechGroup.  
Sediplast is a trademark of LP Italiana SPA.  
Dispette is a trademark of Ulster Medical Products.

1. U. S. Department of Health and Human Services (Département de la Santé et des Services sociaux des États-Unis): Recommendation for Prevention of HIV Transmission in Health Care Settings. MMW Report, Aug 21, 1987, Vol. 36, No. 2S.

Accu-Sed et Excyte sont des marques déposées d'ELITechGroup.  
Sediplast est une marque de LP Italiana SPA.  
Dispette est une marque de Ulster Medical Products.

1. U. S. Department of Health and Human Services (avdeling for helse om omsorgstjenester): Recommendation for Prevention of HIV Transmission in Health Care Settings. MMW Report, Aug 21, 1987, Vol. 36, No. 2S.

Accu-Sed og Excyte er registrerte varemerker som tilhører ELITechGroup.  
Sediplast er et varemerke for LP Italiana SPA.  
Dispette er et varemerke for Ulster Medical Products

1. U. S. Department of Health and Human Services (Amerikanska ministeriet för hälso- och människovård): Recommendation for Prevention of HIV Transmission in Health Care Settings. MMW Report, Aug 21, 1987, Vol. 36, No. 2S.

Accu-Sed och Excyte är registrerade varumärken för ELITechGroup.  
Sediplast är ett varumärke som tillhör LP Italiana SPA.  
Dispette är ett varumärke som tillhör Ulster Medical Products