


ELICAL 2

REF : CALI-05550/1550


LOT : 03-3180

 2023-12

FOR ENVOY ANALYZERS

Reagent REF	REAGENT NAME	VALUE UNITS	VALUE UNITS	METHOD	TRACEABILITY
BITV-0850	TOTAL BILIRUBIN ENVOY	4.35 mg/dL	74.3 μmol/L	<i>Diazo (Blanked)</i>	<i>NIST SRM 916a</i>
BIDV-0850	DIRECT BILIRUBIN ENVOY	2.57 mg/dL	43.9 μmol/L	<i>Diazo (Blanked)</i>	<i>Jendrassik-Grof Manual method</i>
FEFE-0850	IRON ENVOY	191 μg/dL	34.2 μmol/L	<i>Ferene</i>	<i>NIST SRM 937</i>
MAGX-0850	MAGNESIUM ENVOY	2.72 mg/dL	1.12 mmol/L	<i>Xylidyl Blue</i>	<i>Atomic absorption</i>

SRM: Standard Reference Material

 Modification from previous version