




Real-Time PCR Solution for BCR-ABL diagnosis
and transcription monitoring

BCR-ABL P210 & P190 **ELITe MGB®** Kits



ONE STEP assay design

ACCURATE quantification of BCR-ABL transcripts
and correlation to the International Scale

BETTER performance with industry-leading MGB technology

IMPROVED patient monitoring during therapy

Molecular monitoring by quantitative Real-time PCR of BCR-ABL fusion gene expression is today the standard of care according to international guidelines.

BCR-ABL P210 ELITe MGB® and BCR-ABL P190 ELITe MGB® are CE-IVD one-step reverse transcription & quantitative Real-Time PCR assays, developed for BCR-ABL transcripts expression monitoring. The unique test design, associated with MGB technology, enables the laboratory to improve disease monitoring during therapy and result reproducibility across multiple operators and different laboratories.

Improved workflow

- Simplify reaction set-up with one-step reagents
- Drastically reduce manual steps & hands-on-time
- Improve time-to-result: 90 minutes from extracted RNA

Standardized results

- Direct correlation to the International Scale with a P210 RNA material reference
- Compliance with latest international recommendations including EAC guidelines
- Standard curves calibrated with European Reference Material

Accurate quantification

- Excellent clinical performance
 - BCR-ABL P210: 100% sensitivity and 100% specificity
 - BCR-ABL P190: 96% sensitivity and 100% specificity
- Extended measurement range
 - BCR-ABL P210 value from 0,0018% to 61%
 - BCR-ABL P190 value from 0,0032% to 55%

Comprehensive onco-hematology line

- One-step assays for BCR-ABL P210 (b2a2, b3a2) & BCR-ABL P190 (e1a2) diagnosis & expression monitoring
- Unique secondary P210 RNA reference material directly calibrated against the W.H.O. standard*
- Large range of nested PCR assays for diverse onco-hematologic marker identification and diagnosis

Ordering information

Reference	Product	Quantity
RTSG07PLD210	BCR-ABL P210 ELITe MGB Kit	50 x 2 reactions
STDG07PLD210	BCR-ABL P210 ELITe Standard	4 x 2 quantification
SPG07-210	PHILADELPHIA P210 RNA Reference	5 sessions
RTSG07PLD190	BCR-ABL P190 ELITe MGB Kit	50 x 2 reactions
STDG07PLD190	BCR-ABL P190 ELITe Standard	4 x 2 quantification

ELITechGroup Molecular Diagnostics, a major innovator in the molecular diagnostic market, offers a comprehensive Real-Time PCR range of products for diagnosis and patient monitoring in microbiology, onco-hematology and genetic disorders.

Limited License : please visit <http://www.elitechgroup.com/corporate/elitemgb-legalnotice> for complete licensing and warranty information

ELITe MGB® is a registered trademark of ELITechGroup

Product not available in the US - Please contact your sales representative for product availability in your country



ELITechGroup is certified by DEKRA and applies a quality management system in compliance with ISO 9001:2008 and ISO 13485:2012

*"1st W.H.O International Genetic Reference Panel for quantitation of BCR-ABL translocation by RQ-PCR" (NIBSC, UK, code: 09/138)

WORLDWIDE OFFICES

Headquarters
Australia T: +33 1 41 45 07 10
Belgium & Luxembourg T: +32 9 282 05 31
Brazil T: +55 27 3025 1415
France T: +33 4 83 36 10 82
Italy T: +39 02 48 40 35 42
New Zealand T: +64 800 555 611

Serbia T: +381 11 2467119
Switzerland T: +41 26 663 86 60
The Netherlands T: +31 313 430 500
UK T: +44 1442 869320
United Arab Emirates T: +971 4 375 2744
United States T: +1 800 453 2725

 **ELITechGroup**
SOLUTIONS
tailored to your needs
www.elitechgroup.com
info@elitechgroup.com