




Real-Time PCR Assay for *Toxoplasma g.* detection

Toxoplasma g. **ELITe MGB[®]** Kit

MGB



EXCLUSIVE monoreagent & ready-to-use
Toxoplasma gondii Real-Time PCR test

UNIQUE dual test design

BETTER performance with industry-leading MGB technology

IMPROVED detection of *Toxoplasma gondii* during pregnancy
and the follow-up of immunocompromised patient

Molecular diagnosis by a real-time PCR test is considered as the most relevant method to identify a congenital transmission from mother to fetus with an amniotic fluid sample as well as to establish an early diagnosis of toxoplasmosis in immunocompromised patients.

Toxoplasma g. ELITe MGB®* is the first mono-reagent and ready-to-use CE-IVD real-time PCR kit for **Toxoplasma** detection. The dual test design which targets a highly conserved and repeated sequence plus an internal control, combined with the industry-leading MGB technology, guarantee the most sensitive, specific and accurate diagnosis.

Easy to use test

- First CE-IVD mono-reagent and ready-to-use Real-Time PCR test *Toxoplasma gondii* detection
- Validated with amniotic fluid and whole blood sample
- Validated with automated and manual extraction

Accurate result

- Excellent analytical sensitivity : up to 1,91 *Toxoplasma gondii* / extraction
- Extended dynamic range: 10 to 10 Million copies of the highly conserved and repeated (RE) sequence
- Monitoring of the entire process with an integrated internal control
- Industry-leading MGB Technology

Comprehensive range

- Broadest CE-IVD Real-time PCR immunocompromised panel
- Universal cycling conditions to analyze several pathogens in the same amplification run
- Fully manual and automated CE-IVD validation from extraction to result analysis

Ordering information

Reference	Product	Quantity
RTST01PLD	Toxoplasma g. ELITe MGB® Kit	100 tests
CTRT01PLD	Toxoplasma g. RE - ELITe Positive Control	12 Amplification runs

ELITechGroup Molecular Diagnostics, a major innovator in the molecular diagnostic market, offers a comprehensive real-time PCR range of products for diagnosis and monitoring of viral and microbiological infections.



Limited License

Please visit <http://www.elitechgroup.com/corporate/elitemgb-legalnotice> for complete licensing and warranty information

*ELITe MGB® is a registered trademark of Elitech Holding B.V.

Product not available in the US - Please contact your sales representative for product availability in your country

WORLDWIDE OFFICES

Headquarters
 Australia T: +33 1 41 45 07 13
 Benelux T: +61 2 9894 6988
 Brazil T: +31 313 430 574
 France T: +55 27 3025 1415
 Italy T: +33 4 83 36 10 82
 Middle East & Africa T: +39 02 48 40 35 42
 T: +971 4 375 2744

New Zealand T: +64 800 555 611
 Serbia T: +381 11 2467119
 Switzerland T: +41 26 663 86 60
 The Netherlands T: +31 313 430574
 UK T: +44 1442 869320
 United States T: +1 609 216 7361

 **ELITechGroup**
 SOLUTIONS
tailored to your needs
www.elitechgroup.com
info@elitechgroup.com